

BioReference®

Directory of Services 2022-2023 EDITION



BioReference®

DIRECTORY OF SERVICES

2022-2023 EDITION

BIOREFERENCE HEALTH, LLC

DIRECTORY OF SERVICES 2022-2023 EDITION

TABLE OF CONTENTS

WELCOME	A5
MESSAGE FROM DR. JAMES WEISBERGER, M.D., CHIEF MEDICAL OFFICER	A6
COMPANY OVERVIEW	A7
MEDICAL, SCIENTIFIC, LABORATORY OPERATIONS, AND SERVICE STAFF	A9
LICENSURE AND REGULATORY INFORMATION	A12
PATIENT SERVICE CENTERS	A13
PATIENT SERVICE CENTER OFFERINGS	A14
PATIENT SERVICE CENTER LOCATIONS AND SCHEDULING	A15
SUPPLIES	A16
IMPORTANT REGULATORY NOTICE	A17
COMMONLY USED COLLECTION SUPPLIES	A18
SPECIMEN PREPARATION	A23
SPECIMEN PREPARATION AND LABELING	A24
PROPER PHLEBOTOMY TECHNIQUES	A24
SPECIMEN SUBMISSION CRITERIA	A24
DRUGS OF ABUSE TESTING	A26
INSTRUCTIONS FOR CLEAN CATCH URINE COLLECTIONS, URINE CULTURES, AND FROZEN SPECIMENS	A29
CYTOLOGY INSTRUCTIONS	A30
ANATOMIC PATHOLOGY (BIOPSIES AND SURGICAL SPECIMENS)	A34
COLLECTION AND TRANSPORT OF SPECIMENS	A34
SPECIMEN QUALITY AND REJECTION	A39
CYTOLOGY SPECIMEN QUALITY AND REJECTION	A41
CONNECTIVITY SOLUTIONS	A42
INSIGHTDX®: ONLINE TEST ORDERING AND RESULTS	A43
ORDER ENTRY	A44
COLLECTION GUIDANCE FOR SPECIMEN REQUIREMENTS	A45
RESULTS REVIEW	A46
LABORATORY ANALYTICS	A47
PATIENT PORTAL	A48
SCARLET HEALTH®	A49

DIRECTORY OF SERVICES 2022-2023 EDITION

TABLE OF CONTENTS

TEST REQUISITIONS AND RESULT REPORTS	A51
TEST REQUISITION OVERVIEW	A52
TEST ADDITIONS AFTER SUBMISSION	A53
RESULTS REPORTS	A54
CLINICAL REPORT	A54
CLARITEST CORE REPORT	A55
PAP, HPV HIGH RISK & CT/GC REPORT	A56
INHERIGEN® CARRIER SCREENING REPORT	A57
AEROBIC VAGINITIS REPORT	A58
BILLING AND PAYER POLICIES	A59
BILLING POLICIES AND INSURANCE COVERAGE	A60
PAYER POLICY	A62
ALLERGY EVALUATIONS	A63
ALLERGY EVALUATION	A64
REGION 1 RESPIRATORY PROFILE (NORTHEAST)	A65
REGION 2 RESPIRATORY PROFILE (MID-ATLANTIC)	A66
REGION 3 RESPIRATORY PROFILE (SOUTHEAST)	A67
REGION 4 RESPIRATORY PROFILE (SOUTH FLORIDA)	A68
REGION 10 RESPIRATORY PROFILE (WEST SOUTH CENTRAL)	A69
REGION 12 RESPIRATORY PROFILE (SOUTH CA, AZ DESERT)	A70
REGION 13 RESPIRATORY PROFILE (SOUTH CA COAST)	A71
REGION 14 RESPIRATORY PROFILE (CENTRAL CA)	A72
REGION 17 RESPIRATORY PROFILE (NORTH CA AND PACIFIC NORTHWEST)	A73
HOW TO USE THE TEST DIRECTORY	A74
HOW TO USE THE TEST DIRECTORY	A75
AMA DISEASE-SPECIFIC PROFILES	A76
AMA DISEASE-SPECIFIC PROFILES (APPROVED FOR GOVERNMENT HEALTH PLANS)	A77
AMA DISEASE-SPECIFIC PROFILES (APPROVED FOR GOVERNMENT HEALTH PLANS)	A78
LABORATORY TESTS IN ALPHABETICAL ORDER	B1
INDEX BY TEST CODE	C1

BioReference®

Welcome

MESSAGE FROM DR. JAMES WEISBERGER, M.D., CHIEF MEDICAL OFFICER


To our valued customers,

The COVID-19 pandemic has challenged the healthcare industry in unprecedented ways. It has also demonstrated the courage, resilience, and compassion of healthcare providers like you, and we would like to thank you for continuing to choose BioReference as your laboratory services provider.

BioReference has always placed a priority on the health and safety of the customers, patients, and communities we serve. We have worked diligently to deliver testing solutions and technologies that you need to care for your patients. With innovative, new technologies, you will be able to better anticipate patient needs, improve efficiencies, enable end-to-end virtual visits, and provide more personalized care.

The Directory of Services will show you our capabilities as one of the largest full service specialty laboratories in the country. Please do not hesitate to call on your Account Representative if there is additional information that you need that is not provided in this directory. On behalf of the medical, scientific, and technical staff, as well as the numerous service and support personnel at BioReference, we hope that you are pleased with our service and we look forward to working with you this year and in the future.

Sincerely,

A handwritten signature in black ink, appearing to read 'James Weisberger', with a stylized, flowing script.

James Weisberger, M.D.
Chief Medical Officer

COMPANY OVERVIEW

TRUSTED SERVICE. SCIENTIFIC EXCELLENCE.

For more than four decades, BioReference has built a legacy based on scientific excellence, innovation and world-class service in laboratory testing solutions. Amidst a tidal wave of change in healthcare in recent years, BioReference has evolved, adding significant services and solutions aimed at addressing the needs of today's customers and patients. Laboratory testing remains the cornerstone to the services we provide. As we challenge the limits of specialty diagnostics, we are making strategic investments to continue drive innovation and cultivate a unique customer experience, as well as expanding our reach to match the dynamic needs of an ever-changing healthcare system.

BioReference offers scientific expertise and laboratory innovation in oncology, urology, and women's health. Our organization provides credible and tailored solutions for a wide spectrum of customers and patients, including medical practices small and large, hospitals and health systems, correctional institutions, government agencies, educational systems, sport leagues, travel and leisure industries, and retail markets. In addition to an extensive test menu with 99% of tests performed in-house, BioReference's differentiated offerings include large-scale health screening programs, on-demand mobile phlebotomy, and transformative business solutions that optimize laboratory testing and management.

BioReference is part of OPKO Health, Inc., a multinational biopharmaceutical and diagnostics company that seeks to establish industry-leading positions in large, rapidly growing markets. Headquartered in Elmwood Park, New Jersey, BioReference operates 10 regional laboratory facilities around the country, is in-network with the largest health plans in the United States, and serves over 19 million patients annually.

BioReference looks at providing laboratory services through the eyes of the customer. Our specialized brands and products are best-in-class, and designed with service, science, and innovation in mind.

<https://www.bioreference.com/about/our-brands/>



GenPath® Oncology offers a comprehensive portfolio to support the continuum of care of cancer patients – from risk assessment to treatment planning. Backed by a team of specialized pathologists, GenPath Oncology delivers cutting-edge solutions that meet the unique needs of oncologists and pathologists, ranging from routine clinical and special coagulation to complex genomic testing for tumor sequencing and hereditary cancer syndromes.



GenPath® Urology provides cutting-edge uropathology, including the proprietary prostate cancer biomarker test, The 4Kscore Test®. As a follow-up to a suspicious PSA test or DRE, The 4Kscore Test helps support the biopsy decision-making process by determining a patient-specific probability for finding aggressive prostate cancer on biopsy.

GenPath®

WOMEN'S HEALTH

GenPath® Women's Health provides end-to-end laboratory solutions for all women for every stage in their life. With an evidence-based portfolio designed for OBGYNs, MFMs and women's healthcare providers, GenPath Women's Health is pleased to offer testing for cervical and vaginal health, reproductive health, and hereditary cancer screening.

SCARLET®

Scarlet® is an on-demand, mobile laboratory specimen collection service that gives patients the option to have provider-ordered lab test specimens collected when and where they choose – from the home or workplace.* Safe, easy, convenient, and backed by BioReference, Scarlet helps empower patients to take control over their healthcare.

**In select regions of the United States.*

MEDICAL, SCIENTIFIC, LABORATORY OPERATIONS, AND SERVICE STAFF

SENIOR LEADERSHIP TEAM	
Craig Allen	President and Chief Executive Officer
James Weisberger, MD, FACP, FASCP	Chief Medical Officer
Ellen Beausang	Chief Commercial Officer
Richard Schwabacher	Chief Business Officer
Michael Schmid	Chief Financial Officer
Jane Pine Wood	Chief Legal Officer
Robert Rossi	Chief Compliance and Privacy Officer
Vincent D'Orazio	Senior Vice President, Operations
Greg Cahill	Senior Vice President, Human Resources
Merve Ozkus	Chief Information Officer

MEDICAL AND SCIENTIFIC STAFF	
James Weisberger, MD, FACP, FASCP	Chief Medical Officer; Senior Vice President; Laboratory Director - Elmwood Park, NJ and Melbourne, FL
Tamera Paczos, MD, MS	Associate Chief Medical Officer, Vice President, Laboratory Director, Mineola, NY
Wojciech Gorczyca, MD, PhD	Chief of Flow Cytometry and Hematopathology Services
Luoquan Wang, MD	Chief of Anatomic Pathology and Cytopathology Services
James W. Sharp, MD, MBA	Vice President, Bioinformatics
Rebecca Khan, MD, FACOG	Senior Director of Women's Health, Medical Affairs
Ayala Aviram-Goldring, PhD	Director, Constitutional and Prenatal Genetics
Jeffrey Gilbert, MD	Senior Director, Sexually Transmitted Infections
Dehui Ku, PhD	Director, Special Coagulation
Bodhraj Acharya, PhD, FAACC	Assistant Director, Special Coagulation
Priscilla Fernandes, PhD, FACMG	Director, Reproductive Genetics
Karim Ouachi, MD, FACMG	Director, Cancer Genetics
Hashim Othman, PhD	Director, Special Chemistry and Toxicology
Wojciech Swat, PhD	Director, Cancer Genomics and Molecular Pathology
Katherine Berezowski, MD	Laboratory Director - Clarksburg, MD, and Gaithersburg, MD
Andrew Grainger, MD	Laboratory Director - Columbus, OH
Rajendra Persad, MBBS	Laboratory Director, Core Lab and Cytology - Houston, TX
Yan Fan, MD	Laboratory Director - Park Avenue, New York, NY
Blythe Gorman, MD	Laboratory Director - Houston, TX
Anila Guruji, MD	Laboratory Director - Burbank, CA
Joseph Weisensel, DO	Laboratory Director - Hackensack, NJ

LABORATORY AND SERVICE OPERATIONS	
Amy Frank	Vice President, Clinical Laboratory Operations
Tanuj Pasricha	Vice President, Logistics and Supply Chain Operations
Karen Muthusammy	Vice President, Customer Solutions
Ada Gazzillo	Vice President, Pre-Analytic Operations
Victorine Ndjatou	Senior Director, Microbiology
Susan Aveta	Vice President, Phlebotomy
Frank Buccini	Vice President, Technology
Jamie Schwartzmeyer	Vice President, Quality Systems
Adria Marallo	Senior Director, Pre-Analytical Technical Excellence

HOW TO CONTACT US	
Main Number	Telephone: (800)-229-5227 • Fax: (201)-345-7048
Español	Telephone: (888)-729-1201 • Fax: (201)-345-7048
Customer Service Clinical Laboratory Supplies	Telephone: (800)-229-5227 Option 1 • Fax: (201)-345-7048 Email: CustomerService@bioreference.com
Customer Service Women's Health Laboratory Supplies	Telephone: (800)-633-4522 Option 1 • Fax: (201)-345-7152 Email: WomensHealthCustomerService@genpathdiagnostics.com
Customer Service Oncology Laboratory Supplies	Telephone: (800)-627-1479 Option 1 • Fax: (201)-345-7166 Email: GenPathCustomerService@genpathdiagnostics.com
Billing	Telephone: (833)-4MYLABS (833-469-5227) • Fax: (201)-703-7130 Email: BillingCS@bioreference.com
Specimen Pick-Up Clinical	Telephone: (800)-229-5227 Option 3
Specimen Pick-Up Women's Health	Telephone: (800)-633-4522 Option 3
Specimen Pick-up Oncology	Telephone: (800)-627-1479 Option 3
Sales/Marketing/Account Setup	Telephone: (800)-229-5227 Option 4 Email: marketing@bioreference.com
Employee Directory By Name	Telephone: (800)-229-5227 Option 6
Human Resources	Telephone: (800)-229-5227 Ext. 8855
Compliance	Telephone: (800)-229-5227 Ext. 8222 Email: ComplianceDepartment@bioreference.com
Tech Support	Telephone: (800)-229-5227 Ext. 8462
HEDIS	Telephone: (800)-229-5227 Option 5

LICENSURE AND REGULATORY INFORMATION

BioReference Health, LLC, located in Elmwood Park, New Jersey, currently holds licenses for testing in the following states, as well as a CLIA Certificate issued by the Centers for Medicare and Medicaid Services (CMS). For a full list of licensure, please visit: <https://www.bioreference.com/about/licensure/>

CLIA CERTIFICATE & STATE LICENSE/PERMIT	ID NUMBER
CLIA Certificate of Accreditation	31D0652945
CAP Certificate of Accreditation	1237201
Medicare PTAN	301910
New York State	PFI 3130
New Jersey State License	0000283
Maryland State Permit	482
Pennsylvania State Permit	22757A
Rhode Island State License	LCO00305
California State License	CDS00800242
West Virginia State License	HIV Screening
Vermont State Certificate	Drug Screening
FDA License (Donor Services)	FEI:3003652672

BioReference®

Patient Service Centers

PATIENT SERVICE CENTER OFFERINGS

BioReference has a network of Patient Service Centers across the United States, providing convenient access and options for patients and healthcare providers.

BioReference Patient Service Centers are designed with patient comfort in mind. From newly remodeled waiting rooms to courteous and professional phlebotomists, BioReference is dedicated to delivering a consistently positive experience to every patient.

BioReference is contracted as an in-network provider with most national and regional insurance carriers. Patients having difficulty meeting their financial obligations or suffering from financial hardship may contact the laboratory to discuss payment options and determine eligibility for financial assistance by: calling 800-229-5227 (**option 2**) or emailing BillingCS@bioreference.com.

All patients must have a laboratory test order from a licensed healthcare provider at the time of the appointment. Please visit www.bioreference.com for a complete list of locations and hours of operation.

BioReference realizes every patient is unique. The following service offerings are available to your patients:



Pediatric Draws



Coagulation Draws



Tuberculosis Screening




Chain-of-Custody Testing



STI Screening



Genetic Testing



*COVID-19 Antibody Testing

* Patients with orders for COVID-19 testing seeking evaluation of active disease infection will not be collected for molecular diagnostic testing via nasal-swab at BioReference Patient Service Centers. Nasal-swabbed specimens should be collected at healthcare provider offices, hospitals or other clinic settings.

PATIENT SERVICE CENTER LOCATIONS AND SCHEDULING


Finding a location is easier than ever with our online search tool.

Simply visit: <https://www.bioreference.com/patients/locations-patient-services/find-a-location/>.

Then choose Patients, Find a Location, then enter your zip code, and hit Search.

APPOINTMENTS

Patient Service Centers Mobile Collection Mobile Collection (COVID-19)



an OPKO Health Company

For appointments at BioReference Patient Service Centers Only. All patients must have a laboratory order from a healthcare provider at the time of the appointment.

[CLICK HERE TO SCHEDULE AN APPOINTMENT](#)

At this time, mylabology operates on a walk-in basis only, with the exception of COVID-19 antibody testing.

LOCATIONS


Zip Code Category (Optional)

FIND A LOCATION

Patients may also schedule an appointment by selecting the red button labeled, **“CLICK HERE TO SCHEDULE AN APPOINTMENT”**.

This will open a new window where the patient will select **“Get Started”**. On the next page, the patient enters the requested information on the left-hand side of the page, and then selects **“Find nearest locations”**. Location options will populate in the middle of the screen, and the patient may then select the location of their choice by choosing **“Select”**.

On the following page, the patient will choose the date and time they prefer, and then select **“Complete Booking”**. The patient will receive a text and email confirmation of their appointment, which also provides helpful information about preparing for the appointment and accessing results.



Contact Information

First name:

Middle name:

Last name:

Date of birth:

Sex at birth:

Email address:

Confirm email address:

Phone number:

I authorize BioReference to text me updates about my appointment. Messaging and data rates may apply.

Home address:

Apt/Suite/Other:

Select location

Find locations near:

- PSC ELMWOOD PARK - JJ1003**
0.1 miles away - Show on map
491-B EDWARD H.ROSS DR., Elmwood Park, New Jersey 07407
[Calendar icon] [Sun] [Mon] [Tue] [Select]
- PSC CLIFTON, BROAD ST - JJ565**
3.65 miles away - Show on map
716 BROAD STREET, Clifton, New Jersey 07013
[Calendar icon] [Sun] [Mon] [Tue] [Select]
- PSC CLIFTON, CLIFTON AVE - JJ366**
3.85 miles away - Show on map
1033 CLIFTON AVE #213, Clifton, New Jersey 07013
[Calendar icon] [Sun] [Mon] [Tue] [Select]
- PSC WAYNE - J7260**
6.54 miles away - Show on map
510 HAMBURG TRKE #104, Wayne, New Jersey 07470
[Calendar icon] [Sun] [Mon] [Tue] [Select]
- PSC WESTWOOD - J7306**
7.71 miles away - Show on map
270 OLD HOOK ROAD, Westwood, New Jersey 07675
[Calendar icon] [Sun] [Mon] [Tue] [Select]
- PSC HAMILTON PLACE, W. 137TH ST. - CX400**
9.76 miles away - Show on map
17 HAMILTON PLACE, New York, New York 10031
[Calendar icon] [Sun] [Mon] [Tue] [Select]
- PSC UNIVERSITY AVE - CX576**
10.49 miles away - Show on map
1550 UNIVERSITY AVE, Bronx, New York 10453
[Calendar icon] [Sun] [Mon] [Tue] [Select]
- PSC HARLEM, W. 127TH ST. - CX682**
10.5 miles away - Show on map
154 WEST 127TH ST, New York, New York 10027
[Calendar icon] [Sun] [Mon] [Tue] [Select]
- PSC BURNSIDE - CG944**
10.85 miles away - Show on map
85 W. BURNSIDE AVE., Bronx, New York 10453
[Calendar icon] [Sun] [Mon] [Tue] [Select]

Map showing location options near 475 Market St, Elmwood Park, NJ 07407, USA.

BioReference®

Supplies


IMPORTANT REGULATORY NOTICE

The Federal Government, through the Centers for Medicare and Medicaid Services (CMS), has advised that laboratories may provide items, devices, or supplies if they are used solely to collect, transport, process, or store specimens. If any of these items, devices, or supplies may be used for purposes other than collection, transportation, processing, or storage of specimens, the laboratory may not provide them free of charge. An example of an item which a laboratory cannot provide free of charge is latex gloves. With regard to those items that are permitted, such as specimen collection devices (tubes, vials, needles, tourniquets etc.), the laboratory may provide only an amount that is reasonably related to the number of specimens referred to the laboratory by the healthcare provider.

NOTE: New York State law further restricts the provision of certain equipment and supplies that laboratories may provide to healthcare providers. New York permits laboratories to provide, at no charge, only those items, devices, or supplies that do not have any generally accepted use in healthcare practices other than to collect, transport, process, or store specimens. Further, these regulations specifically prohibit the distribution of additional items, such as adhesive bandages, alcohol prep pads, gauze pads, etc. Laboratories must provide permitted supplies, items, and devices of a size, type, and quantity reasonably related to the type and number of specimens being referred by the healthcare provider to the clinical laboratory.

Supplies and containers used for our laboratory services are provided at no additional charge. These include all blood collection tubes, needles, needle holders, slides, preservative solutions, and cytology fixative.

Supply request forms are available from the lab and should be used to request needed items. Please allow two to three days for delivery of your supplies.



BioReference
LABORATORIES
an **OPKO** Health Company

Fax Your Order To:
201-345-7048

Email Your Order To:
CustomerService@bioreference.com

Account Number:	Account Name:
Account Phone Number:	Account Email:

Date of Request: ____ - ____ - ____

Signature of Client: _____

Please ensure you keep enough supply on hand for delivery turn around:
For courier delivery, please allow 48-72 hours to process your order.
For FedEx delivery, please allow 3-5 days to process your order.

Important Notice: Federal and State Law mandates that Clinical Laboratories may give clients only those supplies directly related to the collection of specimens and in amounts proportionate to specimens received from the client.

CYTOLOGY			
Speedy#	Description	Units	Qty. Req.
325	SUREPATH BROOM	Pack/50	
368	SUREPATH BRUSHES/SCRAPER	Pack/25	
313	SUREPATH COLLECTION VIAL	Pack/25	
316	THIN PREP BROOM	Pack/50	
369	THIN PREP BRUSH/SPATULA	Pack/25	
305	THIN PREP COLLECTION VIAL	Pack/25	
KITS			
Speedy#	Description	Units	Qty. Req.
248	4KSCORE KIT	Each	
907	BREATH KITS - H. PYLORI KIT	Each	
404	CHAIN OF CUSTODY (SINGLE)	Each	
286	CLARITEST COLLECTION KIT	Each	
328	O&P (OVA & PARASITE) KIT	Pack/25	
412	PEDIATRIC URINE COLLECTION KIT	Each	
127	QUANTIFERON TB KIT	Each	
MICROBIOLOGY			

PRINTER SUPPLIES			
Speedy#	Description	Units	Qty. Req.
804	RESULT PAPER	Ream/500	
801	SMART PRINTER LABELS	Pack/2	
257	TONER - HP M452DN CF410X (BLACK)	Each	
256	TONER - HP M452DN CF411X (CYAN)	Each	
255	TONER - HP M452DN CF412X (YELLOW)	Each	
254	TONER - HP M452DN CF413X (MAGENTA)	Each	
852	TONER - HP PRO 400 CE410X (BLACK)	Each	
853	TONER - HP PRO 400 CE411A (CYAN)	Each	
854	TONER - HP PRO 400 CE412A (YELLOW)	Each	
855	TONER - HP PRO 400 CE413A (MAGENTA)	Each	
REQUISITIONS			
Speedy#	Description	Units	Qty. Req.
813	ADVANCE BENEFICIARY NOTICE (ABN) FORM	Each	
612	ALLERGEN REQUISITION	Each	
609	CHAIN OF CUSTODY	Each	
641	CYTOLOGY REQUISITION	Each	
601	GENERAL REQUISITION	Each	

COMMONLY USED COLLECTION SUPPLIES

A full list of supplies can be found on the supply order form



SST with Clot Activator Order #101:

This tube is the standard Serum Separator Tube. Please follow these instructions when using the SST Tube with Clot Activator in order to obtain the most accurate test results:

1. Collect blood specimen using standard venipuncture technique. Fill tube completely.
2. Gently invert tube five times to mix clot activator with blood.
3. Allow blood to clot for 30 minutes, do not exceed 60 minutes.
4. Centrifuge at high speed for 15 minutes.
5. Remove from centrifuge. Barrier will have formed, separating cells from serum. All of the separation gel should have moved from the bottom of the tube to form a barrier layer.
6. The sample is now ready to be transported to the laboratory. Do not remove stopper.




Lavender Top Order #102:

EDTA as anticoagulant. Avoid exposure to extreme temperatures. Available in microtainer sizes for capillary collections.



Tan Top Order #105:

K2EDTA (plasma), used for lead.




Yellow Top Order #106:

ACD A additive (plasma tube)



Royal Blue Top ADTA Order #107:

COMMONLY USED COLLECTION SUPPLIES



Royal Blue Top

Order #108:

No additive, trace metals free.



Grey Top

Order #112:

Potassium Oxalate as an anticoagulant. Sodium Fluoride as preservative.



Green Top

Order #116:

Sodium or Lithium Heparin as anticoagulant.



Red Top, Plain

Order #121:


Contains a clot activator. Used for tests where gel separator is not desirable, such as therapeutic drug monitoring.



Light Blue Top

Order #125:

Sodium citrate as anticoagulant. Available in 2.7mL size. This tube must be filled to capacity.



White Top

Order #126:

PPT gel tube containing EDTA (plasma tests).



ThinPrep Preservative

Fluid Collection Vial

Order #305:

Collection container for cervical and vaginal specimens. Used for GenPap testing.



Micro Test M4 Transport

Order #306:


Used for Chlamydia, Mycoplasma, and Ureaplasma cultures.

COMMONLY USED COLLECTION SUPPLIES



**SurePath Preservative Fluid
Collection Vial
Order #313:**

Collection container for cervical and vaginal specimens. Used for GenPap testing.




**ThinPrep Broom
Order #316:**

Used to collect a cytology sample.




**SurePath Broom
Order #325:**

Used to collect a cytology sample.



**O&P (Ova & Parasite)
Kit Order #328:**

For Faecal Concentration of Helminth Ova And Larvae / Protozoa Cysts and Oocysts.



**JEMBEC Plate
Order #354:**

Plastic compact container with chocolate agar and CO₂-generating tablet. Use for Gonorrhea cultures.




**SurePath Brush
Order #368:**

Used to collect a cytology sample.



**ThinPrep Spatula
Order #369:**

Used to collect a cytology sample.



**Sterile Container
Order #401:**

Plastic, wide-mouthed 5-ounce sterile container with graduated measure indicator on side. Supplied with cap.

COMMONLY USED COLLECTION SUPPLIES



Urine Container, for 24-hour Urine Collection
Order #406:

64-ounce plastic container with handle. Must specify: No Preservative, 6N HCL, or Boric Acid.



Urinalysis Tube
Order #409:

Plastic vial with graduated measurement indicator and yellow cap. Stabilizing preservative added. (Note: may be a coated tube, or contain tablet.)

***Do Not Use For Urine Cultures.**



BD Probetec Tube
Order #413:

Used for Male STI Urine Profile testing.




BBL Culture Swab
Order #501:

Used for Group B Strep by PCR, Group A Strep, DNA, Myobacterial (Acid Fast Bacillus) Culture/Stain, and Fungal Culture + Stain.



Culturette, White E-Swab
Order #502:

Used for aerobic/anaerobic cultures.



Aptima Swab
Order #503:

Used for Chlamydia and Gonorrhea testing.



Dacron Swab
Order #504:


Used for collection of anogenital samples.



Digene Swab
Order #505:

Used for HPV testing.

COMMONLY USED COLLECTION SUPPLIES



**GenProbe Urine Collection
Order #506:**

Plastic vial with white screw cap. Bacteriostatic preservative present. Use for Chlamydia and Gonorrhea testing.

***Do Not Use For Routine Urinalysis.**



**Meridian Viral Transport Vial
Order #507:**

Used for Herpes Select I/II cultures.




Swab Viral Culturette Order #509:

Used for Chlamydia, Mycoplasma, and Ureaplasma cultures.



**BD Affirm Swab
Order #511:**

Used for Candida, Trichomonas, and Gardnerella testing.



**Culturette, Blue E-Swab
Order #513:**

Used for aerobic/anaerobic cultures from nasopharynx and pediatric sample collection.



**Culturette, Green E-Swab
Order #514:**

Used for routine aerobic/anaerobic cultures from eye, ear, nasopharynx, and urogenital tracts.



**Charcoal Swab
Order #516:**

Used for Gonorrhea testing.

BioReference®

Specimen Preparation

SPECIMEN PREPARATION AND LABELING

All specimens may be potentially infectious material and should be handled, labeled, and transported accordingly.


It is essential that the following instructions are followed exactly to ensure delivery of a specimen that is adequate for testing. All specimens must be properly identified by using a minimum of two unique patient identifiers, per CLIA, that must include the patient’s full name and at least one additional identifier, such as: date of birth, chart number, or a patient ID. Also, the collector must be identified on all specimen containers submitted. The test requisition form requires completion and must include the time and date of the specimen collection, as well as the referring physician’s name and contact information, or other legally authorized person requesting the patient’s test(s) and such ordering provider’s signature. Electronic signature is acceptable as the order can be placed through an electronic medical record (EMR); or documentation needs to be in the patient medical record of the treating healthcare provider site clearly indicating intent to order the test(s) and signed by the ordering healthcare provider.

PROPER PHLEBOTOMY TECHNIQUES



1. Fill tubes to capacity required level

Always make sure to collect enough specimen.



2. Mix gently

Immediately following collection, all plastic tubes require mixing. Invert all tubes eight (8) times, except Light Blue tubes (invert four (4) times).



3. Separate

Do not use gel tubes (red/black SST) for toxicology/drug testing.

SPECIMEN SUBMISSION CRITERIA

Blood: When whole blood is requested, obtain the full amount in a vacuum tube as shown in the specimen requirement section of this compendium.

Lavender, Grey, Green, Yellow, Tan, and Royal Blue Top tubes contain different anticoagulants that inhibit blood coagulation. When drawing these specimens, immediately invert the tube 8-10 times. Invert Light Blue Top tubes 4 times.

Do not shake the tube, as this can cause hemolysis.

Serum: Obtain sufficient blood to yield the required volume of serum. A plain Red Top tube or Red/Serum Separator Tube (SST) should be used. (Refer to tests section for details.) A 10mL Red or SST tube will yield approximately 5 mL of serum.

When drawing specimens, immediately invert the tube eight (8) times. Allow the blood to clot upright in a rack for 30 minutes, not to exceed 60 minutes, then centrifuge for 15 minutes at 3135-3465 RPMs to separate the serum (4500-5000 RPMs for coagulation testing). If an SST tube is used, no other manipulations are required. Make sure that the gel has formed a thick, solid, intact barrier between the serum and the clotted cells. [If the gel doesn't form a barrier, transfer serum to the appropriately labeled aliquot/transfer tube.](#)

Plasma: When processing plasma, follow the instructions for each test.

Urine/Urinalysis: To adequately test urine specimens, the sample should be collected in a tube with a stabilizing chemical present. The tube provided contains a **yellow** "pop off" cap and a "Stabilur®" tablet or coated tube that preserves the formed elements, such as red cells, white cells, casts, and epithelial cells. For urinalysis, use a paper cup and transfer 8-10 mL of urine to the tube, replace cap and ensure proper labeling.

Urine Chemistry: Assays may require a 24-hour collection that may contain boric acid, hydrochloric acid, or sodium carbonate as a preservative or may be a non-additive collection. [Refer to the specific test in this compendium for test details.](#) Instruct the patient to discard the first urine voided upon arising in the morning thereafter, collect all urine specimens in a paper cup and transfer to the 24-hour container, including the first morning void of the following day. A normal intake of fluid is recommended. Measure the 24-hour volume and record it on the container and the test requisition form. Keep the specimen refrigerated until picked up by the laboratory. Urine chemistry tests that do not require a 24-hour collection container and can be submitted in the container specified in the compendium. Please keep refrigerated.

Urine, Drugs of Abuse (DAU): For routine DAU testing, submit a specimen in a sterile urine container with no preservative.

Urine, COC for Drugs of Abuse: The donor providing the urine specimen must be advised regarding the purpose of the specimen collection and the consequences it may have, legal or otherwise. The donor must sign the consent in chain of custody request form. Ask donor to remove unnecessary outer garments and personal belongings. The donor may retain his or her wallet. Check bathroom for any containers, waste, etc., before the donor enters and after the donor exits. Accompany the donor into the bathroom. Have the donor wash and dry hands. After washing hands, the donor shall remain in the presence of the collection site personnel and shall not have access to any water fountain, faucet, soap dispenser, cleaning agent, or any other material that could be used to adulterate the urine specimen. The donor is to be instructed not to flush the toilet until the specimen is delivered to the collector. The donor may be given a larger container from which the urine specimen is transferred to the sample container(s) for submission to the laboratory. The transfer and sealing of the sample must be done in the presence of the donor until such time as the sample bottle is sealed with a tamperproof seal. The donor is required to initial the seal.

DRUGS OF ABUSE TESTING

The following are intended to be used as general guidelines only. Many variables may affect duration of detection, such as drug metabolism and half-life, subject's physical condition, fluid balance and state of hydration, and route and frequency of ingestion.

DRUG NAME	SCREEN CUTOFF	CONFIRMATION CUTOFF
Amphetamines Also known as: Speed Pharmaceutical Names: Dexedrine, Benzedrine	1000 ng/mL	500 ng/mL
Benzodiazepines Pharmaceutical Names: Diazepam (Valium), Oxazepam (Serax), Chlordiazepoxide (Librium), Alprazolam (Xanax), Chlorazepate (Tranxene), Temazepam (Restoril)	100, 200 ng/mL	200 ng/mL
Clonazepam/7-aminoclonazepam		25 ng/mL
Cannabinoids Also known as: Dope, Weed, Hemp, Hash, Colombian, Sinsemilla Pharmaceutical Names: Marinol	100, 50, 20 ng/mL	15 ng/mL
Cocaine Also known as: Coke, Crack, Rock	100, 300 ng/mL	150 ng/mL
Methadone Also known as: Fizzies Pharmaceutical Names: Amidone, Dolophine	300 ng/mL	200 ng/mL
Phencyclidine Also known as: PCP, Angel Dust	25 ng/mL	25 ng/mL
Opiates Also known as: Cotton, Oxycotton, Percs Pharmaceutical Names: Oxycodone (OxyContin®), Hydrocodone (Vicodin®), Codeine, Morphine, etc.	200, 300 ng/mL	300 ng/mL
Propoxyphene Also known as: Biscuits, Chill Pills, Trammies Pharmaceutical Names: Darvon, Darvocet, Ultram	300 ng/mL	300 ng/mL
Methaqualone Also known as: Quaaludes, Sopors, Ludes, Mandrax	300 ng/mL	
Alcohol Biomarkers (ETG) Pharmaceutical Names: Diazepam (Valium)	500 ng/mL	500 ng/mL for ETG 250 ng/mL for ETS

DRUG NAME	SCREEN CUTOFF	CONFIRMATION CUTOFF
Alcohol, Urine	50, 20 mg/dL	20 mg/dL
Alcohol, Blood	10 mg/dL	10 mg/dL
Oxycodone Also known as: Oxycontin/Percocet	300 ng/mL	300 ng/mL
Buprenorphine Also known as: Subutex, Suboxone	5 ng/mL	5 ng/mL
Methadone Metabolite	300 ng/mL	
6-Monoacetyl Morphine (6-MAM)*	10 ng/mL	10 ng/mL
3, 4-Methylenedioxymethamphetamine (MDMA) Also known as: Ecstasy, Molly	300 ng/mL	
Meperidine (Demerol)	200 ng/mL	
Tramadol	200 ng/mL	
K2 Group 1 Also known as: Spice	20 ng/mL	
K2 Group 2 Also known as: Spice	10 ng/mL	
Fentanyl	2 ng/mL	
Serum Drug Screen 10 drug panel: Amphetamine, Methamphetamine, Opiates, Cocaine Benzodiazepines, THC, Barbiturates, Methadone, PCP, Oxycodone		

*6-Monoacetyl Morphine (6-MAM) is a specific and very short-lived metabolite of heroin. Confirmation of (6-MAM) by GC/MS in urine specimens serves as a positive indicator of recent heroin use.

Testing Methodologies: Specimen is initially screened using an enzyme immunoassay (EIA). Our initial screening assays are semi-quantitative and provide a numeric value as positive or negative result, indicating the presence or absence of a respective drug or drug class. Gas Chromatography/Mass Spectrometry (GC/MS) or Liquid chromatography/Mass Spectrometry (LC/MS) is recommended for any positive EIA results that will be used for legal or employment purposes, as this test can identify and quantify the amount of specific drug present in samples.

Chain of Custody Tracking Services: BioReference provides a custody and control form which is a six-part chain of custody document. A chain of custody document is a drug specimen collection tool that outlines the steps of the collection process to collectors. Also provided for chain of custody tests are TamperChecks, which are tests to confirm that specimens have not been adulterated. Validity testing also guards against adulteration, whereby you can guarantee the veracity of the urine drug test by accounting for pH, creatinine levels, and specific gravity, and testing for adulterants that may have been introduced to the urine specimen. The adulterants detected include bleach, nitrites, chromates, and other oxidizing agents.

INSTRUCTIONS FOR CLEAN CATCH URINE COLLECTIONS, URINE CULTURES, AND FROZEN SPECIMENS

Clean Catch Midstream Urine Collection:

1. Provide patient with a routine urine collection cup. To avoid contamination, instruct patient not to touch the inside of the cup, cap, or straw.
2. Patients should clean themselves with towelettes as follows:
 - a. **Males** – Wipe the head of the penis in a single motion with the first towelette, then repeat with the second. If uncircumcised, hold the foreskin back both before cleansing and while collecting the urine sample.
 - b. **Females** – Separate the labia, and wipe the inner folds of skin from front to back in a single motion with the first towelette, then down through the center of the labial folds with the second towelette. Make sure to keep the labia separated while collecting the urine sample.
3. Patients should urinate in a controlled stream into the toilet.
4. Patients should place the collection cup under the stream of urine and continue to urinate into the cup. Once the collection cup is full, they may finish urinating into the toilet.
5. Patient should then notify the healthcare practitioner of completed collection for next steps.

Urine Culture:

Specimen for Urine Culture should be transferred to a BD Vacutainer® Grey-Top urine culture tube. Instruct the patient not to urinate directly into the vial. The patient should urinate directly into a paper collection cup, and then immediately transfer the specimen into the tube. Refrigerate the specimen as soon as possible.

Frozen Specimens:

Certain tests must be submitted frozen because of the viability of the analyte being tested. Keep all frozen specimens separate from the routine tests and submit a separate test requisition form. As soon as possible, separate the serum or plasma and transfer to a plastic transfer tube. Place the specimen in the freezer and keep until it is solid. [Notify the Transportation Department as soon as possible that you have a frozen specimen for pick-up.](#)

STORE YOUR SPECIMEN IN THE FREEZER UNTIL PICK-UP.

CYTOLOGY INSTRUCTIONS

Use the Cytology requisition form for all cytology specimens. Relevant clinical information such as last menstrual period (LMP), prior diagnosis, etc., should be noted in the spaces provided.

Fine Needle Aspiration Biopsy (FNA):

1. Collect sample directly into 30 mL of CytoLyt® solution. If specimen must be collected in an intravenous solution, use a saline solution.
2. Label appropriately and send immediately in securely closed container.

Note: If possible, flush the needle and syringe with a sterile anticoagulant solution prior to aspirating the sample. Some anticoagulants may interfere with other cell processing techniques, so use caution if using the specimen for other testing.

Body Fluids (including urine, pleural fluid, ascites, CSF, and pericardial fluid):

1. Collect body fluids directly into 30 mL of Cytolyt solution. **Note:** Urine may be collected in PreservCyt® solution. A 2:1 urine-to-PreservCyt ratio is required. The UroCyto urine collection kit is recommended.
2. Label appropriately and send immediately in securely closed container.

Direct Smears:

1. Write patient's name and secondary identifier on frosted end of clean slide.
2. Spread material evenly over slide.
3. Fix immediately with cytology spray fixative from a distance of 10 to 12 inches until liquid droplets form.
4. Allow slide to dry before sending it out in designated slide holder.

Respiratory Specimens (including sputum, bronchial brush, or bronchial wash):

1. Collect sample directly into CytoLyt solution, or add 30 mL of CytoLyt solution to the fresh specimen as soon as possible.

Note: Large specimens (greater than 20 mL) should be concentrated before addition of CytoLyt solution to the sample.

2. Label appropriately and send immediately in securely closed container.

Ancillary Testing (including cell blocks for potential immunohistochemistry and special stains):

1. Collect specimen in 30 mL of CytoLyt solution.
2. Label appropriately and send immediately in securely closed container.

Halo Breast Pap:

1. If Nipple Aspirate Fluid is produced during procedure, specimens are to be collected with a custom, non-cell binding swab, which is placed in non-gyn Liquid Based Cytology (LBC) fixative, such as CytoLyt or CytoRich® Red.

Note: Use of gyn fixatives, such as PreservCyt, will not adversely affect the specimen, but their use is restricted (by FDA) to gyn applications.

2. Label appropriately and send immediately in securely closed container.

ThinPrep® Broom Collection Protocol:

1. Obtain an adequate sampling from the cervix using the broom-like device. Insert the central bristles of the broom into the endocervical canal deeply enough to allow the shorter bristles to fully contact the ectocervix. Push gently, and rotate the broom in a clockwise direction five (5) times.
2. Rinse the broom in the PreservCyt Solution vial by pushing the broom onto the bottom of the vial ten (10) times, forcing the bristles apart. As a final step, swirl the broom vigorously to further release material.

DISCARD THE BROOM.

3. Tighten the cap so that the torque line on the cap passes the torque line on the vial.
4. Record the patient's name and ID number on the vial. Record the patient information and medical history on the Cytology requisition form.
5. Place the vial and requisition in a specimen bag for transport to the laboratory.

ThinPrep Brush/Spatula Collection Protocol (for GenPap STI testing and Pap tests):

1. Obtain an adequate sampling from the ectocervix using a plastic spatula.
2. Rinse the spatula as quickly as possible into the PreservCyt Solution vial by swirling the spatula vigorously in the vial ten (10) times. **DISCARD THE SPATULA.**
3. Obtain an adequate sampling from the endocervix using an endocervical brush device. Insert the brush into the cervix until only the bottom-most fibers are exposed. Slowly rotate 1/4 or 1/2 turn in one direction. **DO NOT OVER-ROTATE.**
4. Rinse the brush as quickly as possible in the PreservCyt Solution by rotating the device in the solution ten (10) times while pushing against the vial wall. Swirl the brush vigorously to further release material. **DISCARD THE BRUSH.**
5. **After obtaining Pap sample, please use a sterile polyester tipped swab to isolate organisms from the vagina.**
6. Obtain an adequate sampling from vagina using a sterile polyester tipped swab.
7. Rinse the swab as quickly as possible into the PreservCyt Solution vial. **DISCARD THE SWAB.**
8. Tighten the cap so that the torque line on the cap passes the torque line on the vial.
9. Record the patient's name, date of birth, and/or ID number on the vial, and the patient information and medical history on the Cytology requisition form.
10. Place the vial and requisition in a specimen bag for transport to the laboratory.

SurePath Brush/Spatula Collection Protocol

- 1A. Insert the contoured end of the Pap Perfect® plastic spatula and rotate 360o around the entire exocervix.
- 2A. **SNAP THE DEVICE HANDLE AND DROP THE DETACHABLE HEAD OF THE DEVICE INTO THE SUREPATH™ VIAL.**
- 3A. Place cap on vial; do not tighten, go to Step 1B.
- 1B. Insert Cytobrush® Plus GT into the endocervix until only the bottom-most bristles are exposed at the os. Slowly rotate 1/4 to 1/2 turn in one direction.

To reduce unnecessary bleeding, do not over-rotate brush.

- 2B. **SNAP THE DEVICE HANDLE AND DROP THE DETACHABLE HEAD OF THE DEVICE INTO THE SUREPATH™ VIAL.**
- 3B. Place the cap on the vial and tighten. Send the SurePath™ vial to the lab for processing.

SurePath Broom Collection Protocol:

1. Position tip of longer bristles in cervical os (opening of the cervix). Begin rotating in clockwise direction. Bristles will begin to stiffen.
2. Continue rotating in a clockwise direction and gently push towards the cervix until the shorter bristles begin to bend extending over the ectocervix.
3. Complete five (5) 360 degree rotations. **SNAP THE DEVICE HANDLE AND DROP THE DETACHABLE HEAD OF THE DEVICE INTO THE SUREPATH™ VIAL.**
4. Cap the vial tightly.
5. Send the specimen containing the head(s) of the sampling device(s), with appropriate paperwork, to the laboratory for processing.

GenPap STI Lesion Profile:

Note: If an STI Lesion Profile is ordered with a Pap and/or vaginal STI sample, a separate SurePath or ThinPrep vial in which to collect the STI lesion sample is needed. **DO NOT** put the STI lesion sample in the vial with the Pap or with the STI vaginal sample. When taking sample from a lesion, please use the following instructions:

1. Vesicles/pustules should be unroofed. Vigorously swab the base of the lesion with a sterile polyester tipped swab.
2. If an ulcer is being tested, squeeze the lesion 3 or 4 times then vigorously rub the base of the ulcer with a sterile polyester tipped swab.
3. If the lesion is crusted, clean it first with a sterile saline solution and vigorously rub with a sterile polyester tipped swab.
4. Swirl the swab vigorously in the vial solution. **DISCARD THE SWAB.**
5. Tighten the cap so that the torque line on the cap passes the torque line on the vial.
6. Record the patient's name and ID number on the vial, and the patient information and medical history on the cytology requisition form.
7. Place the vial and requisition in a specimen bag for transport to the laboratory.

ANATOMIC PATHOLOGY (BIOPSIES AND SURGICAL SPECIMENS)

General Instructions:

1. Use the Surgical Pathology requisition form for all biopsies and surgical specimens. Relevant clinical information should be written down in the spaces provided.
2. Write patient's name and date of birth on specimen container.
3. Immediately place all tissues in 10% buffered formalin at ten times the volume of the specimen. Specimen containers with 10% formalin are provided by the laboratory.
4. Send immediately to the laboratory in a securely closed container. For additional information, please contact Client Services.

COLLECTION AND TRANSPORT OF SPECIMENS

Correct specimen collection and transport of clinical specimens to the laboratory are extremely important for rapid and accurate identification of significant microorganisms from patient samples.

General Considerations for Collection and Transport:

- Use sterile technique and transport to the laboratory as soon as possible.
- Close containers securely to prevent leaking of samples during transport. These specimens are biohazards.
- Whenever possible, obtain specimens prior to the administration of antibiotics.
- Do not use expired tubes or media for specimen collection.
- Write the patient's name and date of birth on each specimen container.
- If submitting two or more cultures, a separate test requisition is necessary for each culture.
- Send specimens in one of the following transport systems:
 1. E-Swabs with transport media: anaerobic and aerobic cultures.
 2. BBL Culture Swab: Mycobacterial (Acid Fast Bacillus); Fungal Culture & Stain.
 3. Sterile containers: sputum, body fluids (except blood and urine, see #4)
 4. Special transport systems:
 - **Blood:** Two 20 mL Vacutainer tubes with Supplemented Peptone Broth (SPB).
Note: E-Swabs cannot be used for these assays.
 - **Gonorrhea Specimen:** Urethral discharge or any source - ThinPrep vial, SurePath vial, Roche, Aptima collection tubes, BD, or E-Swab.
 - **Chlamydia Specimen:** Urethral discharge or any source - ThinPrep vial, SurePath vial, Roche, BD, Aptima collection tubes or Viral culturette.
 - **Urine:** BD Vacutainer Grey-Top tube for Urine Culture and Roche collection devices.
 - **Parasitology:** Special collection kits with formalin and PVA.

Body Fluid Culture (for PD fluid, See Peritoneal Fluid):

Pleural, pericardial, and synovial fluids must be aspirated aseptically. The skin at the body site should be disinfected with an iodophor prior to aspiration. Use sterile technique. Inoculate into sterile tube, container, or blood culture media.

Blood Culture:

Disinfect skin at the body site with iodophor or Chlorhexidine swabs which are 2% chlorhexidine Gluconate and 70% Isopropyl Alcohol, prior to venipuncture. Use sterile technique. Disinfect top of two (2) tubes with alcohol prior to inoculation. Two sets from separate venipuncture sites are recommended.

For adults, inoculate a set of blood culture bottles (1 Aerobic and 1 Anaerobic bottle) with 8-10 mL of blood in each bottle.

For pediatric collections, inoculate 1 pediatric blood culture bottle with 1-3 mL of blood.

Eye and Ear:

Use E-Swab (no Vial). Label with "LEFT" or "RIGHT," as applicable.

Fluid:

See Body Fluid or Peritoneal Fluid.

Fungi:

See Mycology.

Environmental (Water and Dialysate) Cultures:

See separate instructions for use of Millipore Sampler.

Genital Culture:

Collect these specimens using an E-Swab. Swabs must be stored at room temperature until transported to the laboratory.

Glucose Tolerance Testing (GTT):

For 2-hour and 3-hour panels: Collect blood sample in Grey Top tube upon patient arrival for the baseline glucose value. Dispense glucose beverage (generally 50g, 75g or 100g depending on test ordered). Glucose beverage should be consumed within 5 minutes. Collect subsequent blood samples at hourly intervals (1, 2, and/or 3 hours, depending on test ordered). All tubes must be labeled appropriately (Initial, 1-hour, 2-hour, and 3-hour, etc.). Some 1-hour or 2-hour tests do not require an initial baseline collection. (NOTE: 1, 2, 3 hr times start from when the patient finishes drinking the glucose beverage).

Gonorrhea/Chlamydia Endocervical NAAT Analysis:

1. Specimen collection swabs and transport media are supplied by the laboratory. Remove excess mucus from the cervical os and surrounding mucosa using one of the swabs provided and **DISCARD THE SWAB.**
2. Insert the second swab from the collection kit 1.0-1.5 cm into the endocervical canal.

3. Rotate the swab 30 seconds in the endocervical canal to ensure adequate sampling.
4. Withdraw the swab, avoiding any contact with the vaginal mucosa.
5. Insert this swab into the transport tube, or vial (ie. ThinPrep, Aptima) cap the tube, and store at 2 - 25°C until tested.

Group B Strep:

1. Collect specimen at 35 to 37 weeks gestation.
2. Using a single swab or two separate swabs, swab the distal vagina (vagina introitus), followed by the rectum (insert swab through the anal sphincter).
3. Insert swab in transport tube, cap securely, and label appropriately.

Fetal Fibronectin:

1. Submit one vaginal swab.

Maternal Serum Screening:

- Combined First, Combined First Plus, Integrated, Modified Sequential, Quad Screen, Serum Integrated, Triple Screen, and AFP:
 1. Draw one SST tube.
 2. Allow to clot for about 20 minutes.
 3. Centrifuge at 2500 RPM for 10 minutes.
 4. Refrigerate until shipment.

Molecular Genetics:

- **InheriGen® Plus:** Draw three lavender tubes. Invert tubes eight (8) times to prevent clotting.
- **Ashkenazi Jewish Tests:** Draw two lavender tubes. Invert tubes eight (8) times to prevent clotting.
- **Cystic Fibrosis, Fragile X, and Spinal Muscular Atrophy (SMA):** Draw one Lavender Top tube per test. Invert tube eight (8) times to prevent clotting.
- **Tay-Sachs Enzyme (Hexosaminidase A Enzyme):** Draw one Red Top and one Yellow Top tube. If patient is pregnant, only Yellow Top tube is required.
- **Chromosome Analysis/Karyotype, Blood:** Draw one green top sodium heparin tube. Invert gently 5-6 times.
- **Chromosomal Microarray Analysis, Blood:** Draw one EDTA lavender top tube. Invert 8-10 times. Do NOT shake tube.
- **Y-Chromosome Microdeletion:** Draw one EDTA lavender top tube. Invert 8-10 times. Do NOT shake tube.
- **ClariTest™ Core Non-Invasive Prenatal Screening:** Collect 2 full Cell-Free DNA White Top Plastic Tubes from mother. See kit for detailed collection instructions. It is recommended to use 21-gauge straight needle; avoid butterfly if possible. *****It is important to fill both tubes to full capacity.*****
*****Invert both tubes 8-10 times each.**

Mycology (Yeast):

Cultures for yeast can be submitted on a culturette. For fluids and sputum, best results are obtained by submitting the entire specimen. If dermatophytes are suspected, the specimen should be submitted in a dry, sterile tube.

Nasal Culture:

An E-Swab (**Blue Container Item # 24803**) is gently inserted through the nose to the posterior nasopharynx, where it is gently rotated. It should remain in this position for several seconds. The withdrawal should be slow to minimize irritation. Place the inoculated swab into the sterile plastic tube.

Parasitology Specimens:

Stools for Ova and Parasites should be shipped in Ova and Parasite kits. Transfer 5 gm minimum of stool in EACH of the paired vials. CLEAR tape preparation or pinworm paddle is appropriate for submission of specimens for pinworm examination. Submit intact parasites (insects or worms) in 70% alcohol.

Peritoneal Fluid or Dialysate Culture:

Disinfect bag's injection sampling port. Use sterile collection technique. Inoculate 3 mL each into four Vacutainer tubes with SPB (Supplemented Peptone Broth). Disinfect tops of tubes with alcohol prior to inoculation.

QuantiFERON®-TB GOLD Plus

- Collect one 6 mL Lithium Heparin (green top) tube. Allow tube to fill completely before removing.
- Gently mix by inverting the tube several times to dissolve the heparin. Label tube with two (2) patient identifiers and time of specimen collection.
- Place sample with completed requisition into a QTF-TB Gold Plus specimen bag, one (1) sample per bag.
- Seal bag and place in a refrigerator or cooler maintained at 2-8 degrees C. **Sample must be refrigerated within 3 hours of collection.**
- Ship to laboratory as soon as possible. **MUST REACH THE LAB** within 48 hours from collection.

Sputum Culture:

Instruct the patient to obtain material from a deep cough which is expectorated into a sterile container. Sputum containers are best suited for this collection. The volume of specimen need not be large (3 mL). Once collected, sputum should be refrigerated until transport. Be sure that the cap is tightly sealed on the container once the specimen is collected. A leaky container is a biohazard.

Stool Culture:

Submit stool in Cary-Blair transport media.

Throat Culture:

Use E-Swab to obtain all types of throat specimens. Rub the sterile swab firmly over the back of the throat (posterior pharynx), both tonsils or tonsillar fossa, and any area of inflammation. Once the specimen is collected, place the swab into the sterile plastic tube and break the ampule.

Throat, Group A Strep Screen:

Use BBL culturette swab or E-Swab to obtain specimen for Group A Strep only. Use Group A Strep, Rapid DNA test.

Wound Culture:

A superficial wound culture should be collected with an E-Swab. After collection, place swab back into plastic tube and refrigerate or leave at room temperature until transport. If the lesion is not open, a sterile needle and syringe should be employed to remove material. The swab may be inoculated with this sample.

NOTE:

- Please include site information on the requisition form and on swab with wound culture testing.
- CPT Codes shown with tests are to be used as guidelines only and may be subject to change. Their accuracy is neither expressed nor implied in this compendium. Please consult the AMA CPT Code Book for further information.

SPECIMEN QUALITY AND REJECTION

Hemolysis:

Some analytes may be reported erroneously if the serum is not promptly removed from the clot, or if the Barrier Gel tube is not centrifuged after the clot has formed. Possible discrepancies are low glucose, high potassium, and LDH. Hemolyzed hematological specimens are unsuitable for testing.

Inadequate Draw or Quantity Not Sufficient (QNS):

Hematology tests require a full tube of blood be obtained. This is to ensure the ratio of anticoagulant to the blood volume is exact for quality results. Particularly important are Sodium Citrate tubes (Light Blue Top) used for coagulation tests.

For CBCs, a "short draw" Lavender Top tube will result in red cell crenation, reduced MCV and Hematocrit, and possible changes in leukocyte morphology, platelets, and total leukocyte counts.

Clotted Specimens:

All hematological testing utilizes anticoagulated blood. For blood counts, a Lavender Top tube containing the anticoagulant EDTA is required. All specimens should be collected and the tube filled to the limit of the vacuum. Clotted samples, either macroscopic or microscopic in nature, cannot be processed for CBC testing, as results may produce false leukopenia, low red cell counts, and aberrant red cell indices. Similarly, small clots found in Blue Top tubes (for coagulation tests) will result in falsely-prolonged PT and PTT test results.

Alkaline Phosphatase:

If patient's age is not supplied, the reference range for adults will be used to analyze results. As the enzyme is increased in periods of bone growth, as well as in pathological bone disorders, the reference range for adolescents tends to be much higher than in adults. We will automatically provide adolescent reference ranges for all patients under 18 years of age if the age is clearly included on the test requisition form.

Decreased Bilirubin:

Bilirubin is photodegradable. Prolonged exposure of the specimen to bright light will produce decreased results.

Improperly Labeled Containers:

All specimen transport tubes and containers MUST be labeled properly to avoid specimen rejection and/or a delay in testing. All tubes and containers must be labeled clearly with patient name AND a second unique identifier, such as date of birth.

Decreased CO₂ Levels:

Carbon Dioxide (CO₂) levels are decreased if the specimen is not tested promptly. CO₂ escapes from serum in vitro, at a rate proportional to time. This can be minimized by keeping the stopper on the tube and by refrigeration of the specimen.

Icteric Specimens:

If the specimen is deeply icteric, falsely elevated cholesterol or bilirubin results may be obtained.

Lipemic Specimens:

Lipemia can falsely elevate AST and ALT. Additionally, it can indicate that the patient did not adequately fast for 12-18 hours before having the specimen collected. In this situation, triglycerides will be elevated.

Old Specimens:

Blood specimens older than 24 hours cannot be adequately tested for some analytes. Particularly sensitive are most hematology tests, including coagulation procedures.

Poor Cell Preservation:

Blood cells, particularly leukocytes, become fragile and can be distorted morphologically if the specimen is older than 24 hours. In such situations, a reliable differential white cell count cannot be carried out.

CYTOLOGY SPECIMEN QUALITY AND REJECTION

Unsatisfactory Specimen (UNSAT):

UNSATs are caused by a variety of factors including, but not limited to, patient biology, lubricants, collection technique, or laboratory interpretation. As such, one should investigate further to determine the root cause.

How to Reduce Frequency of UNSAT Results:

1. Review the rinsing technique. If using a brush to collect a cervical specimen, it is imperative to NOT over-rotate the device.
2. Rinsing technique is important to release the cells at the time of capture into the fixative solution. It is crucial to rinse the devices immediately after putting the head into the vial so that the cells do not stick to the device.
3. The use of lubricants is not recommended during Pap testing (Papanicolaou Technique Approved Guidelines [NCCLS Document GP15-A]). Lubricants increase the risk of contaminating or obscuring the cellular sample with both conventional Pap smears and all liquid-based methods. The Clinical and Laboratory Standards Institute (CLSI, formerly the NCCLS) recommends that lukewarm water be used to lubricate and warm the speculum.
 - If a lubricant must be used due to patient discomfort or other circumstances, it should be applied sparingly on the outer portion of the speculum with great care to avoid the tip, using a water-based lubricant.
 - The ONLY two lubricants recommended by Hologic are PAP Test Lubricating Jelly (Aseptic Control Products) and Surgilube Surgical Lubricant (HR Pharmaceuticals). It is critical to specimen quality that if a lubricant is used, only the recommended brands are used.

Lack of ECCs:

Lack of endocervical cells is common to both conventional and liquid-based Paps. Current image-guided technology is not designed to identify endocervical component.

Updated guidelines for lack of ECC results indicate yearly follow-up, provided the Pap results are normal and the patient is not in one of six high-risk categories (ASCCP 2003).

Collection is key to reducing this type of result: do not over-rotate the brush, and make sure to rinse immediately.

QNS for HPV Testing:

Depletion of the liquid is due to the amount of cells in the sample, not the amount of liquid. If the cell count is limited due to individual patient biology (e.g., menopausal) or a collection technique that does not capture enough cells, this may occur. **Again, proper collection is key to reducing this type of result.**

Expired Media:

Specimens collected/transported in expired vials may be rejected.

BioReference®

Connectivity Solutions

INSIGHTDX®: ONLINE TEST ORDERING AND RESULTS

BioReference's advanced digital solutions improve workflows for healthcare providers and help to improve diagnostic outcomes for patients. BioReference increases access to vital patient diagnostic information that can help healthcare providers reach a diagnosis quicker. InsightDx is a state-of-the-art, secure, HIPAA-compliant connectivity solution for easy test ordering and reporting. With the ability to interface with most electronic medical records (EMRs), InsightDx allows healthcare providers to easily access and manage patient and practice information from virtually any device.

Advantages of InsightDx:

- Online test directory and electronic ordering
- Custom test requisitions
- Prompt and easy retrieval of tests results
- Automated printing of test results
- Uni-directional or bi-directional interfacing
- Summary reporting on selected criteria
- Test utilization reports
- Cumulative patient reports
- Collection guidance for specimen requirements

InsightDx is easy to use, and your Account Representative can coordinate a demonstration for your office.

ORDER ENTRY

The test ordering feature has a step-by-step workflow designed with ease of use and test order information accuracy in mind.

Insight Dx

Announcements :: ENG | ESP :: Samuel Smith :: Logout

BioReference LABORATORIES an OPKO Health Company

Practice: Z1111 - Demo Practice Location: Z1112 - Demo Location

Results Orders Patients Admin Help

DOE, JANE Patient ID: DOB: 01/01/1990 Phone Number: 201-555-1111 Gender: F Primary Insurance: UHC

New Order 12345678 - Draw Now - Third Party Save Draft Cancel Order

1. Demographics 2. Tests 3. Diagnosis 4. Review Complete

Selected Tests Search Tests

Please select tests.

Favorite Test Groups

- Cholesterol / Lipids
- First Visit
- Heart Health Visit

Favorite Tests

- J148-8 4K Score
- J149-6 4K SCORE
- B655-2 BRCA-HEREDICV5-1
- B370-8 FAMILY MEMBER TESTING
- 3596-4 Fatty Acid Panel Non-Hy
- B415-1 URINE CYTO RFLX AUC-URO

Frequently Ordered

- 0286-5 25OH, VITAMIN D
- 0101-6 3hr.GTT (preg)
- J264-3 4K PSA rfx >3,0
- 0268-3 5-HIAA
- 0156-0 ABO/Rh BLOOD TYPE
- 0029-9 AFB SMEAR
- 0025-7 AFP, TUMOR MARKER
- 0037-2 Amylase, Urine Random
- 0038-0 ANA SCREEN
- 3747-3 ANTI-CCP
- 0364-0 ANTI-dsDNA
- 0041-4 ANTI-THYROGLOBULIN
- 0412-7 ANTI-TPO Ab
- 5030-2 AUTOMATED UROVYSION by FISH
- A519-1 B.BURGDORFERI DNA
- 2470-3 Babesia Microti IgG IgM
- 5000-5 BIOPSY POD @
- 4000-6 BIOPSY, GENERAL
- 0578-5 C-PEPTIDE
- 2130-3 CA 15-3
- 0823-5 CA 27.29
- 0536-3 CA-125
- 0051-3 CALCIUM, IONIZED
- 0034-9 CBC w/o DIFF
- 0053-9 CBC WITH DIFF
- 0055-4 CEA
- 1186-6 CELIAC DISEASE ABS.
- 0004-2 CHEM 12 PROFILE
- 0058-8 Cholesterol
- 0021-6 Citric Acid, 24 HR Ur
- 0068-7 CK
- J373-2 CLIENT ASSIST DRAW
- J374-0 CLIENT DRAW
- 3427-2 CMP, METABOLIC (CMP)

< Back Next >

Frequently ordered tests and panels can be combined into a single group to streamline ordering, based upon medical necessity

A real-time bridge with the EMR eliminates duplicate patient entries

User - defined "favorite" tests appear

Tests are alphabetically presented on the screen based on the provider's ordering history

COLLECTION GUIDANCE FOR SPECIMEN REQUIREMENTS




InsightDx® also features a collection guidance section to ensure the correct specimen types are drawn based on the tests you ordered.

COLLECTION GUIDANCE

(STANDARD)

03/14/2018 - J9999
10000010MN, 10000020MN

Total Draw Guidance

# TO DRAW	CONTAINER	TEST CODE(S)
2	 <p>SST SST Tube 8.5 ml</p>	0033; 0043; 0049; 0050; 0052; 0057; 0058; 0059; 0070; 0090; 0091; 0105; 0106; 0107; 0108; 0129; 0135; 0146; 0147; 0148; 0153; 0155; 0160; 0185; 0286; 0538; 0539; 0540; 0812; 1976; K135
1	 <p>LAV Lavender Top – EDTA 4 ml</p>	0053; 0086; 0102
1	 <p>USC Urine Cup</p>	0159

Identify the exact specimen requirements for orders placed

Processing Guidance

2 SST

- 1 = Refrigerate
- 1 = Pour off to 1 ALQS (Aliquot – Serum) - Refrigerate

1 LAV

- 1 = Refrigerate

1 USC


- 1 = Pour off to 1 UA (Urine Urinalysis Tube – Yellow) – Refrigerate

RESULTS REVIEW

Test results are easy to find in InsightDx®. Upon login, recent test results will be listed in the Results Inbox for quick review. Critical results will be flagged here so that you may triage your patients as needed. Results are also stored in the patients chart for retrieval any time. InsightDx provides two different views of test results, either a PDF format that is easily printed, or the discrete data (text only) format.

Name	Date	Accession	Flag	Alerts	Status	Provider	V	P
PATIENT_TEST	06/25/2019	123456789	N	N	F	Dr. John Smith	Y	N
PATIENT_TEST	06/25/2019	123456789	Y	N	F	Dr. John Smith	N	N
PATIENT_TEST	06/25/2019	123456789	N	N	F	Dr. John Smith	N	N
TEST_PATIENT	06/25/2019	123456789	N	N	F	Dr. John Smith	N	N
TEST_PATIENT	06/25/2019	123456789	N	N	F	Dr. John Smith	N	N
PATIENT_TEST	06/25/2019	123456789	N	Y	F	Dr. John Smith	N	N
PATIENT_TEST	06/25/2019	123456789	N	N	F	Dr. John Smith	N	N
PATIENT_TEST	06/25/2019	123456789	N	N	F	Dr. John Smith	N	N
PATIENT_TEST	06/25/2019	123456789	N	N	F	Dr. John Smith	N	N
PATIENT_TEST	06/25/2019	123456789	N	N	F	Dr. John Smith	N	N
PATIENT_TEST	06/25/2019	123456789	N	N	F	Dr. John Smith	N	N
PATIENT_TEST	06/25/2019	123456789	N	N	F	Dr. John Smith	N	N
PATIENT_TEST	06/25/2019	123456789	N	N	F	Dr. John Smith	N	N
PATIENT_TEST	06/25/2019	123456789	N	N	F	Dr. John Smith	N	N
PATIENT_TEST	06/25/2019	123456789	N	N	F	Dr. John Smith	N	N
PATIENT_TEST	06/25/2019	123456789	N	N	P	Dr. John Smith	N	N

InsightDx is compatible with any web browser, including Google Chrome, Mozilla Firefox, or Microsoft Internet Explorer. If you're using a smart phone or tablet, the touch-optimized software is easy to use and offers a seamless workflow similar to the desktop experience.



LABORATORY ANALYTICS

The value of laboratory analytics is immeasurable. We're transforming the laboratory experience by closing the gap between data and diagnosis, allowing healthcare providers to better anticipate patient needs, improve efficiencies, and provide more personalized care. Through Global CareIndex®, we aim to help in the treatment of patients with reliable, actionable, and understandable information.

Global CareIndex

Global CareIndex (GCI) is a proprietary data and analytics platform that provides practice management and patient action tools that:

- Allow healthcare providers to generate easy to understand, custom and cumulative reports. GCI provides graphic representations to compare and analyze laboratory test results in order to help ensure consistent patient care.
- Provide chronic disease management reports to track at-risk patients who have not been reevaluated for chronic disease based on abnormal test results.
- Help recognize an acute disease in patients to monitor treatment progress and identify need for an additional evaluation.
- Organize data benchmarks according to a healthcare provider's or practice's parameters to track and trend treatment data across multiple providers.

GCI Benefits:

- Help improve patient care and practice operations
- Facilitate patient follow-up
- Evaluate and track disease progression
- Benchmark, track and trend laboratory data
- Identify at-risk patients
- Help achieve practice performance goals



PATIENT PORTAL

With BioReference’s secure, HIPAA-compliant Patient Portal, patients have the laboratory at their fingertips. From a user-friendly website, patients can access their most recent laboratory test reports, review previous test results, pay their bills, and update health insurance information – all from their computer or mobile device.

Additionally, patients can view their laboratory test results in the Apple Health app on an iPhone (iOS 11.3 or later) or via the CommonHealth app on Android devices.

Parents or legal guardians are able to add dependents to their account in order to view a minor’s COVID-19 test results. Other test results for minors can be acquired by emailing patientportal@bioreference.com with their name, date of birth, and the date their specimen was collected. Minor children cannot register for the self-service portal.

For more information, visit: <https://www.bioreference.com/patient-portal/>
<https://www.personapay.com/bioreference/login>



BioReference®

Scarlet Health®



Scarlet Health® is BioReference's mobile phlebotomy service that provides specimen collection at a convenient location of your patients' choosing, such as their home or workplace.

WHY SCARLET®?

- Scarlet® can help improve testing compliance by:
 - Enabling more convenient specimen collection for your patients
 - Minimizing your patients' concerns over potential infectious disease exposure
 - Potentially reducing your patients' specimen collection anxiety through familiar and comfortable surroundings
- Phlebotomists are professionally trained BioReference employees
- The ever-increasing Scarlet-eligible test menu currently includes nearly 2,000 tests. Visit the BioReference Test Directory at www.bioreference.com/test-directory for the most up-to-date test listing
- Available for patients 10 years and older for all Scarlet-eligible tests
- Available for patients 4 years and older for COVID-19 PCR test nasal swab collections

HOW DOES SCARLET WORK?

Option 1: You can add Scarlet to your electronic BioReference lab test order:

- Add the Scarlet code TH31-8 to the test order in InsightDx or your EMR
- Scarlet will send your patient an email and a text to schedule their appointment
- A Scarlet phlebotomist goes directly to your patient to collect their specimen
- Collected specimens will be sent to BioReference for testing
- Lab test results will be sent to you, electronically
- Your patient may also access their test results via the secure, HIPAA-compliant BioReference Patient Portal

WOULD YOU LIKE TO LEARN MORE?

Contact your BioReference account representative or email hello@scarlethealth.com and our team will be in touch.

To watch a video on how Scarlet works, scan the code or go to scarlethealth.com/meetscarlet



BioReference®

Test Requisitions and Result Reports

TEST REQUISITION OVERVIEW

Requisitions are available through Customer Service by using the order number associated with each requisition. The order number can be found on supply forms or requested through your Account Representative. We offer a variety of unique test requisitions that are tailored to your medical specialty and testing needs.

The form is divided into three main vertical sections: **PATIENT** (top), **BILLING INFORMATION** (middle), and **PHYSICIANS ACCOUNT** (right). The **PATIENT** section includes fields for Patient ID/Room # (1), Comments (2), Fasting/Non-Fasting status (3), Name, Last (4), First (5), MI (6), Street (7), APT. # (8), City (9), State (10), ZIP (11), Date of Birth (12), Age (13), M/F (14), Patient Phone No. (15), Patient Cell No. (16), Patient Email/Fax (17), and Race/Ethnicity (18). The **BILLING INFORMATION** section includes Bill To (19), Insurance Carrier (20), Insurance ID # (21), Group # (22), Subscriber's Name (23), Date of Birth (24), Insurance Address (25), City (26), State (27), ZIP (28), Secondary Insurance Carrier (29), Insurance ID # (30), Group # (31), Diagnosis (32), DX Code (33), Referring Provider (34), Source of Referral (35), and Prior Authorization # (36). The **PHYSICIANS ACCOUNT** section includes the Ordering Provider (37), Send Results Copy (38), Fax Result (39), Name (40), Street (41), City (42), State (43), Zip (44), Fax (45), a signature line (46), and Date (47). At the bottom, there are sections for **COMMONLY ORDERED PANELS** and **CUSTOM DIAGNOSTIC PROFILES**.

1. Patient ID - Indicate the patient ID, if applicable. This will appear on results report.

2. Comments - Indicate any comments, if applicable. This will appear on results report.

3. Fasting - Check appropriate box to indicate if patient is fasting.

4. Account - Identifies the unique client number.

5. Physicians - Identifies the ordering physician(s) name/practice and contact information.

6. Ethnicity - Indicate patient's ethnicity by checking appropriate boxes.

NOTE: This is required for all reportable diseases.

7. Collected - Indicate date and time the specimen was collected.

8. Send Results - If a copy of results should be sent to another healthcare provider, include their contact information here.

9. Bill Information - All appropriate patient billing information is required.

10. Diagnosis Code Indicate the patient's diagnosis and appropriate ICD-10 diagnosis codes.

11. Patient Status - Required for testing on all patients.

12. Physician's Signature - Signature of physician or other authorized NPI provider is mandatory and required.

NOTE: If nothing is indicated, the admitting hospital may be billed for the charges.


TEST ADDITIONS AFTER SUBMISSION

The Customer Service Department can arrange for additional testing if the specimen is stable and the volume is adequate after your requested tests have been completed. We are required by federal mandate to ask for written approval for every test we perform. You will receive a request for written confirmation for verbal requests via hard-copy reporting or by facsimile. The physician or authorized healthcare provider must sign and return this written confirmation.

Change in Testing Authorization Form

Return Fax: 201-345-7048

Please review the information below and return with the required signature.



OPKO Health Companies

Patient Demographics →

Accession Number: _____	Date of Request: _____
Patient Name: _____	Form Initiated By: _____
Account: _____	Account Fax: _____
Date of Collection: _____	Requested By: _____
Patient DOB: _____	Attention To: _____
Patient ID: _____ Specimen Type: _____	Comments: _____

← **Origin of Request**

For Medicaid/Medicare patients, individual components must be ordered – **NO PROFILES/PANELS.**
 This form is not a guarantee testing will be performed, as testing is performed based on sample amount and availability.

ADD ON/DELETE REQUEST
 At your request, listed below are the tests to add or delete for this patient:

Test code(s) and test name(s) will be listed; Add/Delete must also be selected →

CODE	TEST NAME	ADD	DELETE
		<input type="checkbox"/>	<input type="checkbox"/>
		<input type="checkbox"/>	<input type="checkbox"/>
		<input type="checkbox"/>	<input type="checkbox"/>
		<input type="checkbox"/>	<input type="checkbox"/>
		<input type="checkbox"/>	<input type="checkbox"/>
		<input type="checkbox"/>	<input type="checkbox"/>
		<input type="checkbox"/>	<input type="checkbox"/>

****Additional Diagnosis (DX) Codes**

All additional DX codes should be emailed/scanned to BillingDXUpdate

Fax Attempts:

1st Request

2nd Request

← **Area to include additional diagnosis codes**

PLEASE SIGN BELOW:

I have authorized the above changes in testing.

Physician/authorized requestor's signature	Title	Date
--	-------	------

Signature is required within 30 days

← **Requesters signature, title, and date required**

Before sending CITA and creating a new add/delete incident for IHC and FISH, confirm material availability with lab.

PLEASE FAX COMPLETED FORM TO: 201-345-7048 OR EMAIL FORM TO: CustomerService@bioreference.com


481 Edward H. Ross Drive, Elmwood Park, NJ 07407 | T 1 800 229 5227 | F 1 201 791 1941 | www.BioReference.com
 GenPath is a division of BioReference © 2022 All Rights Reserved. BRD007286 V1 08/22

RESULTS REPORTS

Your reports through BioReference are designed for ease of use, including:

- Easy-to-interpret design
- Previous test results
- Detailed clinical abnormalities summary
- Abnormal test results clearly classified information

CLINICAL REPORT



an **OPKO** Health Company

FINAL REPORT

PHYSICIAN

SAMPLE PHYSICIAN
SAMPLE MEDICAL GROUP
123 Main Street
Anytown, USA 12345
Acct #: (12345-0) 20
P: (123) 456-7890

PATIENT

SAMPLE PATIENT
DOB: 01/01/1966 Age: 52 Sex: M
Address: 123 Main Street
Anytown, USA 12345
P: (123) 456-7890

SAMPLE

Specimen ID: 1000000
Date Of Report: 10/03/2017
Date Collected: 09/23/2017
Time Collected: 10:37 Date
Received: 09/23/2017 Time
Received: 18:31

North America Eastern Time

CLINICAL REPORT

CLINICAL ABNORMALITIES SUMMARY: (May not contain all abnormal results; narrative results may not have abnormal flags. Please review entire report.)

CO2	12 LO	LD	364 HI
Cholesterol	220 HI	LDL Cholesterol	122 HI
Milk	0.13 HI		
25OH, VITAMIN D	26.7 LO	ZnPP(EP),CHILD	37 HI

KIT, LAV, SST, U/A. NON FASTING

CHEMISTRY

Test	Result	Abnormal	Reference	Units	Rpt Date	Prior Result	Date
Total Protein	8.2		5.9-8.4	g/dL	09/24/17	7.4	09/08/16
Albumin	4.9		2.8-5.4	g/dL	09/24/17	4.5	09/08/16
Globulin	3.3		1.7-3.7	g/dL	09/24/17	2.9	09/08/16
A/G Ratio	1.5		1.1-2.9		09/24/17	1.6	09/08/16
Glucose	77		70-99	mg/dL	09/24/17	111 HI	09/08/16
Sodium	142		135-147	mmol/L	09/24/17	143	09/08/16
Potassium	4.6		3.5-5.5	mmol/L	09/24/17	4.5	09/08/16
Chloride	102		96-108	mmol/L	09/24/17	104	09/08/16
CO2		12 LO	22-29	mmol/L	09/24/17	22	09/08/16
BUN	14		5-18	mg/dL	09/24/17	12	09/08/16
Creatinine	0.62		0.30-0.70	mg/dL	09/24/17	0.50	09/08/16
BUN/Creat Ratio	22.6		10.0-28.0		09/24/17	24.0	09/08/16
Calcium	10.2		8.9-10.4	mg/dL	09/24/17	9.9	09/08/16
Uric Acid	4.6		3.4-8.5	mg/dL	09/24/17	4.1	09/08/16
Iron	108		59-158	ug/dL	09/24/17	104	09/08/16
Bilirubin, Total	<0.2		<1.2	mg/dL	09/24/17	<0.2	09/08/16
LD		364 HI	135-225	U/L	09/24/17	231 HI	09/08/16
Alk Phos	241		<330	U/L	09/24/17	200	09/08/16
AST	32		<40	U/L	09/24/17	28	09/08/16
Phosphorus	5.6		3.3-5.6	mg/dL	09/24/17	4.7	09/08/16
ALT	15		<41	U/L	09/24/17	15	09/08/16
GGTP	20		10-71	U/L	09/24/17	19	09/08/16

CARDIOVASCULAR/LIPIDS

Test	Result	Abnormal	Reference	Units	Rpt Date	Prior Result	Date
Cholesterol		220 HI	<170	mg/dL	09/24/17	172	09/08/16
Triglycerides	74		<90	mg/dL	09/24/17	194 HI	09/08/16
HDL CHOL., DIRECT	83		>or=60	mg/dL	09/24/17	68	09/08/16
HDL as % of Cholesterol	38		>14	%	09/24/17	40	09/08/16
Evaluation: BELOW AVERAGE RISK							
Chol/HDL Ratio	2.7		<7.4		09/24/17	2.5	09/08/16
Evaluation: BELOW AVERAGE RISK							
LDL/HDL Ratio	1.47		<3.56		09/24/17	0.96	09/08/16
LDL Cholesterol		122 HI	<110	mg/dL	09/24/17	65	09/08/16
VLDL, CALCULATED	15		7-32	mg/dL	09/24/17	39 HI	09/08/16

BioReference Laboratories, Inc. 481 Edward H. Ross Dr Elmwood Park, NJ 07407 (800) 229-5227	James Weisberger, M.D. Laboratory Director	Clinical Page 1 of 6
--	---	----------------------

CLARITEST CORE REPORT

WOMEN'S HEALTH
genpath
an **OPKO** Health Company

Final Report

PHYSICIAN

TEST, PROVIDER
SAMPLE MEDICAL GROUP
123 MAIN STREET
ANYTOWN, USA 12345
ACCT # (12345-0)
P: (123) 456-7890 F: (234) 567-8910

PATIENT


PATIENT, TEST
DOB: 01/01/1985 **AGE:** 34 **SEX:** F
PATIENT ID: 123
ADDRESS: 123 SAMPLE STREET,
ANYTOWN, NJ 12345

SAMPLE


Specimen ID: 123456789-0
Date of Report: 01/01/2021 00:00 PM
Date Collected: 01/01/2021 Time Unknown
Date Received: 01/01/2021 00:00 PM
Source: CFDNA White Top Tube

CLARITEST™ CORE NON-INVASIVE PRENATAL SCREEN
run utilizing the same core technology as **harmony**


RESULTS SUMMARY



Result:
Low Risk




Fetal Sex:
Consistent with Female



Fetal Fraction:
14.0%

CLINICAL INFORMATION

Pregnancy Type: Singleton
Maternal Weight: 150
EDD: 01/24/2022
Gestational Age: 14 weeks, 0 days
IVF Pregnancy: No **Donor:** N/A
Age at Egg Retrieval: N/A



RESULTS

CONDITION	RESULT	INTERPRETATION
Trisomy 21	Low Risk	Results consistent with two copies of chromosome 21
Trisomy 18	Low Risk	Results consistent with two copies of chromosome 18
Trisomy 13	Low Risk	Results consistent with two copies of chromosome 13
Sex Chromosome Aneuploidies	Low Risk	Results consistent with two sex chromosomes

Comments: This is a screening test. False positive and false negative results can occur. See Test Performance Data table below. The negative predictive value (NPV) is >99% for all low risk results. Results should be interpreted and genetic counseling considered based on clinical history. Invasive testing, such as chorionic villus sampling or amniocentesis, is necessary for diagnosis.

TESTING METHODOLOGY

ClariTest™ Core is a non-invasive prenatal screen that measures the relative proportion of chromosomes in maternal blood to assess the probability of fetal trisomies 21, 18, and 13, sex chromosome aneuploidies, and if ordered, 22q11.2 deletion syndrome. Evaluation of sex chromosome aneuploidies and 22q11.2 deletion is only available in singleton pregnancies. A directed analysis of cell-free DNA (cfDNA) is performed using Digital Analysis of Selected Regions (DANSR™) with microarray quantitation and Fetal Fraction Optimized Risk of Trisomy Evaluation (FORTE™) analysis. The result incorporates the fetal fraction of cfDNA, maternal age and gestational age.

Ayala Aviram, Ph.D.
Cytogenetic - Molecular Director, ex. 7510
This report was electronically signed.

BioReference Laboratories, Inc.
481 Edward H. Ross Drive
Elmwood Park, NJ 07407
(800) 633-4522

James Weisberger, M.D.
Laboratory Director

Page 1 of 2
Created: 01/01/2021 6:52 PM

PAP, HPV HIGH RISK & CT/GC REPORT

WOMEN'S HEALTH
genpath
an OPKO Health Company

Final Report

PHYSICIAN

SAMPLE PROVIDER
SAMPLE PROVIDER
123 MAIN STREET
ANYTOWN, USA 12345
ACCT #: **12345-0** **00**
P: (123) 456-7890

PATIENT

PATIENT, TEST
DOB: Age: 45 Y Sex: F
Surgical #: Patient ID:

SAMPLE

Specimen ID: **123456789**
Date Collected: 01/01/2021 00:00 AM
Date Received: 01/01/2021
Date of Report: 01/01/2021 00:00 PM
Source: VAGINAL
Specimen: THINPREP

GYN CYTOLOGY REPORT

Diagnosis
ATYPICAL SQUAMOUS CELLS OF UNDETERMINED SIGNIFICANCE (ASC-US)

Adequacy
Satisfactory for evaluation / Endocervical/transformation zone component present.


Comment
This Pap smear was screened with the assistance of the Hologic ThinPrep(TM) Imaging System and screened by a cytotechnologist.

Test Ordered
PAP, LIQUID-BASED

Clinical Information
LMP: 06/22/2021

Screened By: Rebecca Vergara, M.D. Case Electronically Signed 01/01/2017

Cervicovaginal cytology should be considered a screening procedure subject to false negatives and false positives. Results are more reliable when a satisfactory sample is obtained on a regular repetitive basis, and should be interpreted together with past and current clinical data.



STI MOLECULAR PROFILE

Test	Normal	Abnormal
HPV High Risk Non-16/18 ^{2,3,4,6}	Not Detected	
HPV 18 ^{2,4,6,7}		Detected
HPV 16 ^{2,3,4,6,7}	Not Detected	

Test	Normal	Abnormal
Chlamydia trachomatis ^{5,6,7}	Not Detected	
Neisseria gonorrhoea ^{1,5,6,7}		Detected

Notes

- Rare cross-reactivity may occur in this amplified DNA GC assay due to certain strains of *N. cinerea*, *N. subflava* and *N. lactamica*.
- The cobas® HPV test is FDA-cleared for ThinPrep® specimens and detects genomic HPV DNA in the polymorphic L1 region in 14 subtypes: Type 16, Type 18, and other high risk types (31,33,35,39,45,51,52,56,58,59,66,68). The test has been modified and validated for use in SurePath™ specimens.
- HPV types 16 and/or 18 that were Not Detected were undetectable or below the pre-set threshold.
- The non-repeat rate for HPV genotyping assays varies from 5 to 15%. In the NILM cytology category, there is a low positive predictive value (PPV = 15-20%) for CIN2+ with a positive high risk HPV result.
- Chlamydia trachomatis (CT) and Neisseria gonorrhoeae (NG) qualitative detection is performed on urogenital specimens or urine on one of the below platforms:
 - BD Probetec™ Qx Amplified DNA Assay tested with the BD Viper™ System using Strand Displacement Amplification technology. FDA cleared for both ThinPrep, SurePath and urine samples.
 - The Aptima combo 2® Assay detects ribosomal RNA (rRNA) using target capture and Transcription-Mediated Amplification (TMA) technology. FDA cleared for ThinPrep samples.
 - The cobas® CT/NG v2.0 using the cobas® 4800 or CT/NG on the cobas® 8800 detects DNA using Polymerase Chain Reaction (PCR) technology; the assay is FDA cleared for ThinPrep Samples.


BioReference Laboratories, Inc.
481 Edward H. Ross Dr | Elmwood Park, NJ 07407 | (800) 633-4522

James Weisberger, M.D.
Laboratory Director

Cytology Page 1 of 2
Created 01/01/21 00:00 PM EST

GenPath is a business unit of BioReference Laboratories, Inc.

INHERIGEN[®] CARRIER SCREENING REPORT



an **DPKO** Health Company

Final Report

DOCTOR

SMITH, JANE
 WOMAN'S HEALTH
 123 ANY ROAD
 ANYTOWN, NJ 012345
 ACCT # J0000
 P:(123) 456-7890 F:(555) 123-4567


PATIENT

DOE, JANE
 DOB: 01/01/1981 Age: 35 Y Sex: F
 Patient ID: 12345
 Address: 123 WINDING STREET
 ANYTOWN, NJ 012345

SAMPLE

Specimen ID: **123456789**
 Date Reported: 01/01/2017 00:00 PM
 Date Collected: 12/01/2016 Time Unknown
 Date Received: 12/01/2016 00:00 PM
 Source: Peripheral Blood
 Clinical Information: AMA,Other:ROUTINE SCREENING

InheriGen: Pan-Ethnic Carrier Screen



KEY RESULTS		
Disease	Results	Interpretation
Cerebrotendinous Xanthomatosis (CTX)	Positive for one copy of the c.1183C>T (p.Arg395Cys) mutation in CYP27A1	Carrier of Cerebrotendinous Xanthomatosis (CTX)
<ul style="list-style-type: none"> All other mutations tested were negative. Genetic counseling and partner testing are recommended. 		

This test was developed and its performance characteristics were determined by BioReference Laboratories. It has not been cleared by the U.S. Food and Drug Administration. The FDA has determined that such clearance or approval is not necessary. This test is used for clinical purposes. It should not be regarded as investigational or for research. This lab has been approved by CLIA 88 and designated as a high complexity laboratory and is qualified to perform this test.

M.L. Cremona Ph.D., FACMG
 Clinical Molecular Geneticist
 Electronically signed by M.L. Cremona Ph.D., FACMG

BioReference Laboratories, Inc.
 481 Edward H. Ross Drive
 Elmwood Park, NJ 07407
 (800) 633-4522

James Weisberger, M.D.
 Laboratory Director

Page 1 of 17
 Created 01/01/2017 00:00 PM

GenPath is a business unit of BioReference Laboratories, Inc.

AEROBIC VAGINITIS REPORT

WOMEN'S HEALTH
genpath
an *OPKO* Health Company

Final Report

PHYSICIAN

SAMPLE PROVIDER
SAMPLE PROVIDER
123 MAIN STREET
Anytown, USA 12345
ACCT #: **12345** **00**

P: (123) 456-7890

PATIENT

PATIENT, TEST M
DOB: **Age:** 45 Y **Sex:** F
Surgical #:
Patient ID:

SAMPLE

Specimen ID: 123456789
Date Collected: 01/01/2021 00:00AM
Date Received: 01/01/2021
Date of Report: 01/01/2021 00:00PM
Source: VAGINAL
Specimen: THINPREP

STI MOLECULAR PROFILE

Test	Normal	Abnormal	Test	Normal	Abnormal
L. reuteri/rhannosus ^{6,7,8}		Not Detected	Escherichia coli ^{2,7,8}	Not Detected	
Group A strep ^{5,7,8}	Not Detected		Staphylococcus aureus ^{2,7,8}		Detected
Group B strep ^{4,7,8}	Not Detected		Lactobacillus species ^{1,8,9}	Detected	
Enterococcus faecalis ^{7,8}		Detected	Lacto. Sp. Log Cells/mL ^{1,8,9}	>5.25	

Interpretation
RESULTS FAVOR AEROBIC VAGINITIS (AV).

Comment
No Lactobacillus DNA is detected. Aerobic bacterial DNA (from E. coli, Group B Strep, E. faecalis, S. aureus, Group A Strep) is detected.

Notes

1. This assay utilizes quantitative real time PCR (multiple detection temperatures technology, MuDT™) in liquid cytology specimens to quantify the levels of nucleic acid of Lactobacillus species, G. vaginalis and A. vaginae, and also detects nucleic acid from several obligate anaerobes associated with Bacterial Vaginosis (BV).
2. Escherichia coli (E.coli) is a gram-negative, rod-shaped coliform bacterium. It is a non-spore-forming, facultative anaerobic bacterium, making ATP by aerobic respiration if oxygen is present, but switching to fermentation or anaerobic respiration if oxygen is absent. Staphylococcus aureus is a facultative anaerobic, gram-positive coccus (round) bacterium that is non-motile and does not form spores. Aerobic vaginitis (AV) is caused by the presence of one or more aerobic pathogens, usually without the presence of Lactobacillus, causing a vaginal inflammatory immune response.
3. Enterococcus faecalis is a gram-positive, non-motile, commensal bacterium inhabiting the gastrointestinal tracts of humans. Aerobic vaginitis (AV) is caused by the presence of one or more aerobic pathogens, usually without the presence of Lactobacillus, causing a vaginal inflammatory immune response.
4. Group B streptococcus is a gram-positive, facultative anaerobic coccus bacterium. Aerobic vaginitis (AV) is caused by the presence of one or more aerobic pathogens, usually without the presence of Lactobacillus, causing a vaginal inflammatory immune response.
5. Group A streptococcus is a gram-positive, facultative anaerobic coccus bacterium. Aerobic vaginitis (AV) is caused by the presence of one or more aerobic pathogens, usually without the presence of Lactobacillus, causing a vaginal inflammatory immune response.
6. Lactobacillus rhamnosus is a short, gram-positive, facultative anaerobic, non-spore-forming rod that often appears in chains. Some strains of Lactobacillus rhamnosus bacteria are being used as probiotics. The Lactobacillus rhamnosus and Lactobacillus reuteri species are most commonly found in the healthy female genitourinary tract. Lactobacillus reuteri is a gram-positive, rod-shaped, lactic acid bacterium. These gram-positive bacteria are used as probiotics and as supplements for women who do not have the natural Lactobacillus species (Lactobacillus crispatus, Lactobacillus jensenii, and Lactobacillus gasseri)
7. The Allplex® Aerobic Vaginitis (AV) Panel is qPCR technology for the identification of Escherichia coli (EC), Staphylococcus aureus (SA), Enterococcus faecalis (EF), Group A Streptococcus (GAS), Group B Streptococcus (GBS), Lactobacillus spp. (Lacto), Lactobacillus reuteri/Lactobacillus rhamnosus (LR), and β-globin as an Internal Control (IC).
8. This test was evaluated and its performance characteristics determined by BioReference Laboratories. It has not been cleared or approved by the U.S. Food and Drug Administration. The FDA has determined that such clearance or approval is not necessary. BioReference Laboratories is certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA) as qualified to perform high complexity clinical testing. This test is used for clinical purposes. It should not be regarded as investigational or for research.
9. Results should be interpreted in light of current and past laboratory data and symptomatology.

BioReference Laboratories, Inc.
481 Edward H. Ross Dr | Elmwood Park, NJ 07407 | (800) 633-4522

James Weisberger, M.D.
Laboratory Director

Cytology Page 1 of 2
Created 01/01/2021 0:00:00 PM EST

GenPath is a business unit of BioReference Laboratories, Inc.

BioReference®

Billing and Payer Policies

BILLING POLICIES AND INSURANCE COVERAGE

Billing to Patient:

BioReference can bill your patients directly for our services. The patient's full name and street address (including apartment number, city, state, and zip code) must be clearly printed in the space provided on the test requisition. A complete address at the time the test is ordered is essential. Each test requisition will result in a separate bill to the patient. Payment of patient bills is due upon receipt and, if not paid, will be followed by subsequent reminders and normal collection activity.

BILLING INFORMATION	BILL TO: <input type="checkbox"/> INSURANCE <input type="checkbox"/> PATIENT <input type="checkbox"/> CLIENT <input type="checkbox"/> MEDICARE <input type="checkbox"/> MEDICAID <input type="checkbox"/> TRICARE			ALL INSURANCES RELATION TO SUBSCRIBER: <input type="checkbox"/> CHILD <input type="checkbox"/> SELF <input type="checkbox"/> SPOUSE <input type="checkbox"/> OTHER: _____		
	INSURANCE CARRIER		INSURANCE ID #	GROUP #		
	SUBSCRIBER'S NAME					DATE OF BIRTH
	INSURANCE ADDRESS		CITY	STATE	ZIP	
	SECONDARY INSURANCE CARRIER		INSURANCE ID #	GROUP #		
	DIAGNOSIS	DX CODE	DX CODE	DX CODE	DX CODE	DX CODE
	REFERRING PROVIDER		PRIORITY AUTHORIZATION #			
	PATIENT STATUS – ONE MUST BE CHECKED					
	<input type="checkbox"/> HOSPITAL INPATIENT <input type="checkbox"/> HOSPITAL OUTPATIENT <input type="checkbox"/> NOT A HOSPITAL PATIENT			HOSPITAL PATIENT DATE OF DISCHARGE ____/____/____		

↖ Bill To
↖ Insurance information

↖ Patient Status
↖ Diagnosis Code(s)

Billing to Third Party:

BioReference will bill insurance carriers directly, when eligible for coverage under the insurer's criteria, provided all the necessary information is supplied with the specimen. Please be certain to include all of the necessary information in the spaces provided on the test requisition at the time the test is ordered, including any and all diagnosis codes to the highest level of specificity. Complete information significantly reduces the need to interrupt and impose upon your staff or patients with requests for any missing information. BioReference will bill patients for all applicable copayments, coinsurance, and deductibles commensurate with the terms of their insurance policy and as indicated on the patient's Explanation of Benefits issued by the patient's health benefit plan.

Medicare / Medicaid:

For patients that are part of government healthcare programs, only AMA-approved panels should be ordered. All other tests must be ordered individually. A special test requisition is available upon request for all government sponsored health plans (**see page A77**).

Patient Status:

Required for testing on all patients. If nothing is indicated, the admitting hospital may be billed for the charges.

Reminder to Ordering Providers:

Diagnostic testing services should be ordered only when medically and reasonably necessary and when required for the diagnosis or treatment of a patient. Federal and state payers (Medicare and Medicaid) typically exclude most testing that is only ordered for screening purposes.

When ordering any of the above tests outside of their applicable utilization guidelines and without additional documentation, please have the patient sign the Advanced Beneficiary Notice (ABN), so they may be billed directly for the test.

PAYER POLICY

BioReference®

BioReference is committed to providing the highest quality laboratory services for your patients. An important part of that quality is helping to ensure our customers and patients understand all aspects of our services. The guiding philosophy of BioReference is to perform only medically appropriate testing as set out in the guidelines of leading medical associations and organizations, as well as leveraging our team of medical experts and key opinion leaders in the industry to deliver the best diagnostic tools available to healthcare providers. It is also essential that BioReference adhere to all applicable federal, state, and local laws and regulations regarding billing for laboratory services; this letter seeks to clarify the billing policies of the company.

In accordance with the Explanation of Benefits (EOB) provided by the patient's health insurance carrier, regardless of whether in-network or out-of-network, patients will be sent bills for the portion for which they are financially responsible in an effort to make a reasonable attempt to collect any and all outstanding balances. This will be followed by a pre-collection notice advising patients that if payment is not received in a timely manner, they may be outsourced to a collection agency. Financial Assistance Programs are available for qualified patients.

Uninsured patients or those who choose to self-pay are financially responsible for all services provided. All bills are expected to be paid in full. If payment is not received in a timely manner, they may be outsourced to a collection agency.

Patients may conveniently pay their bill online by visiting bioreference.com/patient-portal.

BioReference strives to provide unparalleled laboratory services for every one of your patients, and to enact the most patient-friendly billing policies in accordance with all laws and regulations. If you have any questions, please do not hesitate to contact the Billing Department at (833)-4MYLABS (833-469-5227).

For the most updated insurance information visit: bioreference.com/insurance

BioReference®

Allergy Evaluations

ALLERGY EVALUATION

All single allergens require 1 ml serum/spun barrier tube SST. If testing for any of the following profiles, please send 2 full SST tubes.

FOOD ALLERGY PROFILE TEST CODE: K426-7

ALLERGEN TYPE	INDIVIDUAL TEST CODE
Egg white	0671-8
Cow's milk	0672-6
Cod fish	0840-9
Wheat	0673-4
Sesame	0936-5
Peanut	0674-2
Soybean	0680-9
Hazelnut	1404-3
Almond	0815-1
Shrimp	0696-5
Tuna	0698-1
Salmon	0714-6
Cashew	1412-6
Walnut	1479-5
Scallop	2543-7
Total IgE	0996-9

CHILDHOOD ALLERGY PROFILE TEST CODE: K726-0

ALLERGEN TYPE	INDIVIDUAL TEST CODE
D. pteronyssinus	0715-3
D. farinae	1372-2
Cat dander	0665-0
Dog dander	0664-3
Mouse urine	1360-7
Egg white	0671-8
Milk	0672-6
Cod fish	0840-9
Wheat	0673-4
Peanut	0674-2
Soybean	0680-9
Shrimp	0696-5
Walnut	1479-5
Cockroach	0716-1
Cladosporium herbarum	0661-9
Alternaria alternate	0660-1
Total IgE	0996-9

**For Medicare/Medicaid patients, individual tests must be ordered.*

REGION 1 RESPIRATORY PROFILE (NORTHEAST)

TEST CODE: 6823-9

Includes: CT; MA; NJ; NY; PA; VT; ME; NH; RI

ALLERGEN TYPE	INDIVIDUAL TEST CODE
D. pteronyssinus	0715-3
D. farinae	1372-2
Cat dander	0665-0
Dog dander	0664-3
Mouse urine	1360-7
Bermuda grass (Cynodon dactylon)	1198-1
Timothy grass (Phleum pratense)	0683-3
Cockroach	0716-1
Penicillium notatum	2167-5
Cladosporium herbarum	0661-9
Aspergillus fumigatus	0694-0
Alternaria alternata	0660-1
Maple (box elder; Acer negundo)	0684-1
Birch (Betula verrucosa)	0685-8
Mountain Cedar (Juniperus sabinoides)	6334-7
Oak (Quercus alba)	0652-8
Elm (Ulmus americana)	0659-3
Walnut (Juglans californica)	0686-6
Maple leaf sycamore, London Plane	1305-2
Cottonwood (Populus deltoides)	0891-2
White ash (Fraxinus americana)	1306-0
Mulberry	3835-6
Common ragweed (short; Ambrosia elatior)	0717-9
Mugwort	0718-7
Rough pigweed (Amaranthus retroflexus)	1627-9
Sheep sorrel (Rumex acetosella)	0894-6
Total IgE	0996-9

REGION 2 RESPIRATORY PROFILE (MID-ATLANTIC)

TEST CODE: 6824-7

Includes: DE; MD; VA; DC; NC

ALLERGEN TYPE	INDIVIDUAL TEST CODE
D. pteronyssinus	0715-3
D. farinae	1372-2
Cat dander	0665-0
Dog dander	0664-3
Mouse urine	1360-7
Bermuda grass (Cynodon dactylon)	1198-1
Timothy grass (Phleum pratense)	0683-3
Johnson grass (Sorghum halepense)	1824-2
Cockroach	0716-1
Penicillium notatum	2167-5
Cladosporium herbarum	0661-9
Aspergillus fumigatus	0694-0
Alternaria alternata	0660-1
Maple (box elder; Acer negindo)	0684-1
Birch (Betula verrucosa)	0685-8
Mountain Cedar (Juniperus sabinoides)	6334-7
Oak (Quercus alba)	0652-8
Elm (Ulmus americana)	0659-3
Cottonwood (Populus deltoides)	0891-2
Pecan/Hickory (Carya soecue, pecan)	1938-0
Mulberry	3835-6
Common ragweed (short; Ambrosia elatior)	0717-9
Rough pigweed (Amaranthus retroflexus)	1627-9
Sheep sorrel (Rumex acetosella)	0894-6
Total IgE	0996-9

REGION 3 RESPIRATORY PROFILE (SOUTHEAST)

TEST CODE: 6825-4

Includes: GA; SC; North FL

ALLERGEN TYPE	INDIVIDUAL TEST CODE
D. pteronyssinus	0715-3
D. farinae	1372-2
Cat dander	0665-0
Dog dander	0664-3
Mouse urine	1360-7
Bermuda grass (Cynodon dactylon)	1198-1
Timothy grass (Phleum pratense)	0683-3
Bahia grass (Paspalum notatum)	1825-9
Cockroach	0716-1
Penicillium notatum	2167-5
Cladosporium herbarum	0661-9
Aspergillus fumigatus	0694-0
Alternaria alternata	0660-1
Maple (box elder; Acer negundo)	0684-1
Birch (Betula verrucosa)	0685-8
Mountain Cedar (Juniperus sabinooides)	6334-7
Oak (Quercus alba)	0652-8
Elm (Ulmus americana)	0659-3
Pecan/Hickory (Carya soecue, pecan)	1938-0
Common ragweed (short; Ambrosia elatior)	0717-9
Rough pigweed (Amaranthus retroflexus)	1627-9
Sheep sorrel (Rumex acetosella)	0894-6
Nettle (Urtica dioica)	6308-1
Total IgE	0996-9

REGION 4 RESPIRATORY PROFILE (SOUTH FLORIDA)

TEST CODE: 6826-2

Includes: FL, South of Orlando

ALLERGEN TYPE	INDIVIDUAL TEST CODE
D. pteronyssinus	0715-3
D. farinae	1372-2
Blomia tropicalis	B235-3
Cat dander	0665-0
Dog dander	0664-3
Mouse urine	1360-7
Bermuda grass (Cynodon dactylon)	1198-1
Timothy grass (Phleum pratense)	0683-3
Bahia grass (Paspalum notatum)	1825-9
Cockroach	0716-1
Penicillium notatum	2167-5
Cladosporium herbarum	0661-9
Aspergillus fumigatus	0694-0
Alternaria alternata	0660-1
Maple (box elder; Acer negundo)	0684-1
Mountain Cedar (Juniperus sabinooides)	6334-7
Oak (Quercus alba)	0652-8
Elm (Ulmus americana)	0659-3
Australian pine (Causarina equisetifolia)	3837-2
Common ragweed (short; Ambrosia elatior)	0717-9
Rough pigweed (Amaranthus retroflexus)	1627-9
Sheep sorrel (Rumex acetosella)	0894-6
Nettle (Urtica dioica)	6308-1
Total IgE	0996-9

REGION 10 RESPIRATORY PROFILE (WEST SOUTH CENTRAL) TEST CODE: 6827-0

Includes: TX; OK

ALLERGEN TYPE	INDIVIDUAL TEST CODE
D. pteronyssinus	0715-3
D. farinae	1372-2
Cat dander	0665-0
Dog dander	0664-3
Mouse urine	1360-7
Bermuda grass (Cynodon dactylon)	1198-1
Timothy grass (Phleum pratense)	0683-3
Cockroach	0716-1
Penicillium notatum	2167-5
Cladosporium herbarum	0661-9
Aspergillus fumigatus	0694-0
Alternaria alternata	0660-1
Maple (box elder; Acer negundo)	0684-1
Birch (Betula verrucosa)	0685-8
Mountain Cedar (Juniperous sabinoides)	6334-7
Oak (Quercus alba)	0652-8
Elm (Ulmus americana)	0659-3
Cottonwood (Populous deltoides)	0891-2
White Ash (Fraxinus americana)	1306-0
Pecan/Hickory (Carya soecue, pecan)	1938-0
Mulberry	3835-6
Common ragweed (short; Ambrosia elatior)	0717-9
Rough pigweed (Amaranthus retroflexus)	1627-9
Rough marsh elder (Iva)	1481-1
Sheep sorrel (Rumex acetosella)	0894-6
Nettle (Urtica dioica)	6308-1
Total IgE	0996-9

REGION 12 RESPIRATORY PROFILE (SOUTH CA, AZ DESERT) TEST CODE: 6828-8

Includes: South Ca Desert; Az Desert

ALLERGEN TYPE	INDIVIDUAL TEST CODE
D. pteronyssinus	0715-3
D. farinae	1372-2
Cat dander	0665-0
Dog dander	0664-3
Mouse urine	1360-7
Bermuda grass (Cynodon dactylon)	1198-1
Perennial rye grass (Lolium perenne)	0682-5
Johnson grass (Sorghum halepense)	1824-2
Cockroach	0716-1
Penicillium notatum	2167-5
Cladosporium herbarum	0661-9
Aspergillus fumigatus	0694-0
Alternaria alternata	0660-1
Olive (Olea europa)	1251-8
Mountain Cedar (Juniperous sabinoides)	6334-7
Oak (Quercus alba)	0652-8
Elm (Ulmus americana)	0659-3
Cottonwood (Populus deltoides)	0891-2
Acacia (Acacia species)	3832-3
Common ragweed (short; Ambrosia elatior)	0717-9
Mugwort (safebrush; Artemisia vulgaris)	0718-7
Russian Thistle (Saltwort, Salsola kali)	0732-8
Rough pigweed (Amaranthus retroflexus)	1627-9
Total IgE	0996-9

REGION 13 RESPIRATORY PROFILE (SOUTH CA COAST)

TEST CODE: 6820-5

Includes: South Ca Coast

ALLERGEN TYPE	INDIVIDUAL TEST CODE
D. pteronyssinus	0715-3
D. farinae	1372-2
Cat dander	0665-0
Dog dander	0664-3
Mouse urine	1360-7
Bermuda grass (Cynodon dactylon)	1198-1
Timothy grass (Phleum pratense)	0683-3
Johnson grass (Sorghum halepense)	1824-2
Cockroach	0716-1
Penicillium notatum	2167-5
Cladosporium herbarum	0661-9
Aspergillus fumigatus	0694-0
Alternaria alternata	0660-1
Alder (Alnus incana)	0561-1
Mountain Cedar (Juniperus sabinoides)	6334-7
Oak (Quercus alba)	0652-8
Elm (Ulmus americana)	0659-3
Olive (Olea europaea)	1251-8
Walnut (Juglans californica)	0686-6
Cottonwood (Populus deltoides)	0891-2
Mulberry	3835-6
Common ragweed (short; Ambrosia elatior)	0717-9
Mugwort (safebrush; Artemisia vulgaris)	0718-7
Russian Thistle (Saltwort, Salsola kali)	0732-8
Rough pigweed (Amaranthus retroflexus)	1627-9
Total IgE	0996-9

REGION 14 RESPIRATORY PROFILE (CENTRAL CA)

TEST CODE: 6821-3

Includes: Central Ca

ALLERGEN TYPE	INDIVIDUAL TEST CODE
D. pteronyssinus	0715-3
D. farinae	1372-2
Cat dander	0665-0
Dog dander	0664-3
Mouse urine	1360-7
Bermuda grass (Cynodon dactylon)	1198-1
Timothy grass (Phleum pratense)	0683-3
Cockroach	0716-1
Penicillium notatum	2167-5
Cladosporium herbarum	0661-9
Aspergillus fumigatus	0694-0
Alternaria alternata	0660-1
Alder (Alnus incana)	0561-1
Birch (Betula verrucosa)	0685-8
Mountain Cedar (Juniperus sabinoides)	6334-7
Oak (Quercus alba)	0652-8
Elm (Ulmus americana)	0659-3
Olive (Olea europaea)	1251-8
Sycamore (Plantanus acerfolia)	1305-2
Mulberry	3835-6
Common ragweed (short; Ambrosia elatior)	0717-9
Mugwort (safebrush; Artemisia vulgaris)	0718-7
Russian Thistle (Saltwort, Salsola kali)	0732-8
Rough pigweed (Amaranthus retroflexus)	1627-9
Total IgE	0996-9

REGION 17 RESPIRATORY PROFILE (NORTH CA AND PACIFIC NORTHWEST) TEST CODE: 6822-1

Includes: Northern Ca; Or; Wa

ALLERGEN TYPE	INDIVIDUAL TEST CODE
D. pteronyssinus	0715-3
D. farinae	1372-2
Cat dander	0665-0
Dog dander	0664-3
Mouse urine	1360-7
Timothy grass (Phleum pratense)	0683-3
Cockroach	0716-1
Penicillium notatum	2167-5
Cladosporium herbarum	0661-9
Aspergillus fumigatus	0694-0
Alternaria alternata	0660-1
Maple (box elder; Acer negundo)	0684-1
Alder (Alnus incana)	0561-1
Birch (Betula verrucosa)	0685-8
Mountain Cedar (Juniperus sabinoides)	6334-7
Oak (Quercus alba)	0652-8
Elm (Ulmus americana)	0659-3
Walnut (Juglans californica)	0686-6
Cottonwood (Populus deltoides)	0891-2
White ash (Fraxinus americana)	1306-0
Common ragweed (short; Ambrosia elatior)	0717-9
Rough pigweed (Amaranthus retroflexus)	1627-9
Sheep sorrel (Rumex acetosella)	0894-6
Nettle (Urtica dioica)	6308-1
Total IgE	0996-9

BioReference®

How to Use the Test Directory

HOW TO USE THE TEST DIRECTORY

This directory lists our testing services in alphabetical order using common nomenclature. An alphabetical index is provided at the end of the directory to assist you in finding the test you need.


Some tests may have alternate names, so please refer to the alternate name filed within the test listing. If you cannot find a particular test, please contact your Account Representative or our Client Services Department and they will be able to assist you.

Our Test Format:

↓ Test Name
Test Code ↓

Comprehensive Metabolic Panel

Alternate Name: COMP. METABOLIC (CMP)
Methodology:
Preferred Specimen: 1 mL SST Tube
Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes
Storage Instruction: Refrigerate
Turn Around Time: 1 days
Scarlet Eligible: Yes

 **3427-2**

CPT Code: 80053 (1)

↑ Test Information

The CPT codes provided are based on AMA guidelines and are for informational purposes only. CPT coding is the sole responsibility of the billing party. Please direct any questions regarding coding to the payer being billed.

Scarlet blood drop will appear next to test code if eligible for collection by a Scarlet phlebotomist.

BioReference®

AMA Disease-Specific Profiles

AMA DISEASE-SPECIFIC PROFILES

(APPROVED FOR GOVERNMENT HEALTH PLANS)

CODE	PROFILE
3283-9 Acute Hepatitis Panel	Hepatitis A Antibody (HAAb), Igm (86709)
	Hepatitis B Core Antibody (HBcAb), Igm Antibody (86705)
	Hepatitis B Surface Antigen (HBsAg) (87340)
	Hepatitis C Antibody (86803)
2555-1 Basic Metabolic Panel <i>(With regular Calcium; For BMP with Ionized Calcium, please use test code 5563-2)</i>	Calcium, Total (82310)
	Carbon Dioxide (Bicarbonate) (82374)
	Chloride (82435)
	Creatinine (82565)
	Glucose, Serum (82947)
	Potassium (84132)
	Sodium (84295)
	Blood Urea Nitrogen (BUN) (84520)
3427-2 Comprehensive Metabolic Panel	Albumin (82040)
	Bilirubin, Total (82247)
	Calcium, Total (82310)
	Carbon Dioxide (Bicarbonate) (82374)
	Chloride (82435)
	Creatinine (82565)
	Glucose (82947)
	Phosphatase, Alkaline (84075)
	Potassium, Serum (84132)
	Protein, Total (84155)
	Sodium, Serum (84295)
	Transferase, Alanine Amino (ALT) (SGPT) (84460)
	Transferase, Aspartate Amino (AST) (SGOT) (84450)
Blood Urea Nitrogen (BUN) (84520)	

AMA DISEASE-SPECIFIC PROFILES

(APPROVED FOR GOVERNMENT HEALTH PLANS)

CODE	PROFILE
0002-6 Electrolytes Panel	Carbon Dioxide (Bicarbonate) (82374)
	Chloride, Blood (82435)
	Potassium, Serum (84132)
	Sodium, Serum (84295)
3422-3 Hepatic Function Panel/Liver Panel	Albumin (82040)
	Bilirubin, Total (82247)
	Bilirubin, Direct (82248)
	Phosphatase, Alkaline (84075)
	Protein, Total (84155)
	Transferase, Alanine Amino (ALT) (SGPT) (84460)
	Transferase, Aspartate Amino (AST) (SGOT) (84450)
3284-7 Renal Function Panel	Albumin (82040)
	Calcium, Total (82310)
	Carbon Dioxide (Bicarbonate) (82374)
	Chloride (82435)
	Creatinine (82565)
	Glucose (82947)
	Phosphorus, Inorganic (Phosphate) (84100)
	Potassium, Serum (84132)
	Sodium, Serum (84295)
	Blood Urea Nitrogen (Bun) (84520)

BioReference®

481 EDWARD H. ROSS DRIVE
ELMWOOD PARK, NJ 07407
T 1-800-229-5227 • F 1-201-345-7048
WWW.BIOREFERENCE.COM

BioReference®

Laboratory Tests in Alphabetical Order

XCE 8 BY FISH

B514-1

Alternate Name: +8 BY FISH

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Place bone marrow into tube, label with patient's name

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 88373 (1)

Scarlet Eligible: No

-5/5Q- BY FISH-TC ONLY

A869-0

Alternate Name: -5/5Q- BY FISH-TC ONLY

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube. Also acceptable: Whole blood in green (sodium heparin) tube.

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88374 (1)

Scarlet Eligible: No

D5S2064 FOR -5/5q- BY FISH

B507-5

Alternate Name: -5/5q- BY FISH

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Place bone marrow into tube, label with patient's name

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 88374 (1)

Scarlet Eligible: No

-7/7Q- BY FISH-TC ONLY

A868-2

Alternate Name: -7/7Q- BY FISH-TC ONLY

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube. Also acceptable: Whole blood in green (sodium heparin) tube

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88374 (1)

Scarlet Eligible: No

D7S486/XCE 7 FOR -7/7q- BY FISH

B506-7

Alternate Name: -7/7q- BY FISH

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Place bone marrow into tube, label with patient's name

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 88374 (1)

Scarlet Eligible: No

GLYCOMARK

J114-0

Alternate Name: 1,5-ANHYDROGLUCITOL ,1,5-AG

Methodology: Enzymatic

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84378 (1)

Scarlet Eligible: Yes

10-Color Flow Cytometry Test Codes

Alternate Name: 10-Color Flow Cytometry Test Codes

Scarlet Eligible: No

10-Color Screening (Short) Panel

TK69

Alternate Name: 10-Color Screening (Short) Panel

Scarlet Eligible: No

10-Color Screening (Short) Panel TC Only

TK81

Alternate Name: 10-Color Screening (Short) Panel TC Only

Scarlet Eligible: No

ATM FOR 11q22.3 BY FISH

B512-5

Alternate Name: 11q22.3 BY FISH

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Place bone marrow into tube, label with patient's name

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 88373 (1)

Scarlet Eligible: No

13Q14 by FISH TC

A872

Alternate Name: 13Q14 by FISH TC

Scarlet Eligible: No

D13S319 (13q14.3) BY FISH

B510-9

Alternate Name: 13q14.3 BY FISH

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Place bone marrow into tube, label with patient's name

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 88373 (1)

Scarlet Eligible: No

Estradiol, Serum

0516-5

Alternate Name: 17 Beta-Estradiol, Free, Serum

Methodology: Electrochemiluminescence Immunoassay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82670 (1)

Scarlet Eligible: Yes

Hydroxyprogesterone(17), Alpha Serum/Plasma

0997-7

Alternate Name: 17-Hydroxyprogesterone

Methodology: Liquid Chromatography Mass Spectrometry

Preferred Specimen: 1 mL Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 83498 (1)

Scarlet Eligible: Yes

P53 by FISH

A314-7

Alternate Name: 17P53,TP53,17P

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube. Also acceptable: Whole blood in green (sodium heparin) tube.

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88367 (1)

Scarlet Eligible: No

1P/19Q,D(1;19)DELETION BY FISH

A311-3

Alternate Name: 1P/19Q,D(1;19)DELETION BY FISH

Methodology: FISH

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block from client with surgical number imprint.

Turn Around Time: 3 days

CPT Code: 88377 (2)

Scarlet Eligible: No

1q21/CKS1B by FISH

B509

Alternate Name: 1q21/CKS1B by FISH

Scarlet Eligible: No

20Q12 BY FISH-TC ONLY

A871-6

Alternate Name: 20Q12 BY FISH-TC ONLY

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube. Also acceptable: Whole blood in green (sodium heparin) tube.

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88367 (1)

Scarlet Eligible: No

D20C108 FOR DEL (20q) BY FISH

B508-3

Alternate Name: 20q12 BY FISH

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Place bone marrow into tube, label with patient's name

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 88367 (1)

Scarlet Eligible: No

Osmolality, Urine 24hr.

0988-6

Alternate Name: 24 Hours Urine Osmolality

Methodology: Freeze Point Depression

Preferred Specimen: 1 mL Urine Container - 24hr

Collection Instruction: Discard first morning sample, then collect all urines after for 24 hrs, on completion label with collection time and patient's name; record total volume on Requisition

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 83935 (1)

Scarlet Eligible: No

4K PSA rflx >3.0

J264-3

Alternate Name: 4K prostate screen, Prostate test, Kallikrein score, Four K

Methodology: Chemiluminescence

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 84153 (1)

Scarlet Eligible: No

4K PSA Total + Free Reflex >3.0 to 4Kscore

K135-4

Alternate Name: 4K prostate screen, Prostate test, Kallikrein score, Four K

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84153 (1), 84154 (1)

Scarlet Eligible: Yes

4K Prostate Score (Serum)

J148-8

Alternate Name: 4K prostate screen, Prostate test, Kallikrein score, Four K, 4Kscore

Methodology: Electrochemiluminescence Immunoassay

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 81539 (1)

Scarlet Eligible: Yes

5-Flucytosine, S

6180-4

Alternate Name: 5-FLUCYTOSINE

Methodology: Liquid Chromatography Tandem Mass Spectrometry

Preferred Specimen: Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 80299 (1)

Scarlet Eligible: Yes

A. Fumigatus Mold (m3), IgE

0694-0

Alternate Name: A. Fumigatus Mold

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

A. Vaginalis by RT-PCR

J212-2

Alternate Name: A. VAGINALIS BY PCR

Methodology: Polymerase Chain Reaction

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 5-7 days

CPT Code: 87799 (1)

Scarlet Eligible: No

Abetalipoproteinemia

B148-8

Alternate Name: ABETALIPOPROTEINEMIA

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81479 (1)

Scarlet Eligible: Yes

Onkosight NGS ABL1 Sequencing

J132-2

Alternate Name: ABL1 BY NGS

Preferred Specimen: Bone Marrow - Lavender Top

Collection Instruction: Place bone marrow into tube, label with patient's name

Storage Instruction: Room Temp

Turn Around Time: 7-14 days

CPT Code: 81170 (1)

Scarlet Eligible: No

ABO and Rh Blood Typing

0156-0

Alternate Name: ABO/Rh BLOOD TYPE

Methodology: ID-Micro Typing System gel test

Preferred Specimen: Pink Tube-EDTA

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86900 (1), 86901 (1)

Scarlet Eligible: Yes

Absolute Neutrophil Count

3617-8

Alternate Name: ABSOLUTE NEUTROPHIL, ANC

Methodology: Flow Cytometry

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 99999 (1)

Scarlet Eligible: Yes

Cardiolipin Antibodies, IgG, IgM, IgA

0431-7

Alternate Name: ACA, ANTI-CARDIOLIPIN ABS.

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86147 (3)

Scarlet Eligible: Yes

Adrenocorticotrophic Hormone (ACTH)

0507-4

Alternate Name: ACTH; Corticopin,Cosyntropin

Methodology: Chemiluminescence

Preferred Specimen: 1 mL Aliquot Plasma-EDTA-Lavender Top

Collection Instruction: Centrifuge lavender-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA and patient's name. Must submit frozen.

Storage Instruction: STRICT Frozen

Turn Around Time: 2 days

CPT Code: 82024 (1)

Scarlet Eligible: No

Acute Hepatitis Panel

3283-9

Alternate Name: ACUTE HEPATITIS PANEL

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80074 (1)

Scarlet Eligible: Yes

Acute Phase Reactant Profile

0584-3

Alternate Name: ACUTE PHASE REACTANT

Preferred Specimen: SST Tube,Lavender top- EDTA

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Turn Around Time: 11 days

CPT Code: 85651 (1), 85810 (1), 86140 (1)

Scarlet Eligible: Yes

Adenocarcinoma vs. Mesothelioma IHC Panel

5181-3

Alternate Name: ADENOCARCINOMA V.MESOTHEL.

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88341 (8), 88342 (1)

Scarlet Eligible: No

GP/AE1/AE3 (Stain Only)

5108-6

Alternate Name: AE1/AE3 (Stain Only)

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88342 (1)

Scarlet Eligible: No

Maternal Serum Alpha-Fetoprotein Screen (MSAFP) - ONTD Only

0274-1

Alternate Name: AFP

Methodology: Chemiluminescence

Preferred Specimen: 2 mL SST Tube

Collection Instruction: MUST BE 15-22WEEKS AND 6 DAYS. PROVIDE GESTATIONAL INFO.

Storage Instruction: Refrigerate

Turn Around Time: 4 days

CPT Code: 82105 (1)

Scarlet Eligible: Yes

Amniotic fluid AFP

2122-0

Alternate Name: AFP, Amniotic Fluid

Methodology: Chemiluminescence

Preferred Specimen: 20 mL Amniotic Fluid

Collection Instruction: 3mL Amniotic Fluid submitted in a sterile, leak-proof container

Storage Instruction: Room Temp

Turn Around Time: 4 days

CPT Code: 82106 (1)

Scarlet Eligible: No

Alpha Fetoprotein, Tumor Marker (AFP)

0025-7

Alternate Name: AFP, Tumor Marker

Methodology: Electrochemiluminescence Immunoassay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82107 (1)

Scarlet Eligible: Yes

AFP/HCG Tumor Marker

4514-6

Alternate Name: AFP/HCG Tumor Marker

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Turn Around Time: 1 days

CPT Code: 82105 (1), 84702 (1)

Scarlet Eligible: Yes

Ashkenazi Jewish Expanded Profile

J293-2

Alternate Name: AJ NGS 50

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 14 days

CPT Code: 81412 (1)

Scarlet Eligible: Yes

Ashkenazi Jewish DNA Panel (NON NY)

A282-6

Alternate Name: AJ8, Ashkenazi Jewish DNA

Preferred Specimen: 10 mL Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81200 (1), 81209 (1), 81242 (1), 81251 (1), 81255 (1), 81260 (1), 81290 (1), 81330 (1)

Scarlet Eligible: Yes

Post Dialysis Albumin

2256-6

Alternate Name: ALBUMIN POST DIALYSIS

Methodology: Colorimetry

Preferred Specimen: 1 mL Microtainer - Pediatric SST

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82040 (1)

Scarlet Eligible: No

Alcohol Biomarkers (ETG), Qual. (500)

6284-4

Alternate Name: ALCOHOL BIOMARKERS

Methodology: Enzyme Immunoassay

Preferred Specimen: 15 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80307 (1)

Scarlet Eligible: Yes

Aldosterone, Serum

0302-0

Alternate Name: ALDOSTERONE**Methodology:** Chemiluminescence**Preferred Specimen:** 2 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 4 days**CPT Code:** 82088 (1)**Scarlet Eligible:** Yes**ALK (EML4/ALK) BY IHC**

B565-3

Alternate Name: ALK (EML4/ALK) BY IHC**Methodology:** Immunohistochemistry**Preferred Specimen:** Formalin-fixed, Paraffin-embedded Tissue**Collection Instruction:** This comes in block form from client with surgical number imprint.**Storage Instruction:** Room Temp**Turn Around Time:** 1 days**CPT Code:** 88342 (1)**Scarlet Eligible:** No**ALK (EML4/ALK) by IHC TC**

B603

Alternate Name: ALK (EML4/ALK) by IHC TC**Scarlet Eligible:** No

ALK BY FISH (TC)

B236-1

Alternate Name: ALK BY FISH (TC)

Methodology: FISH

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block from client with surgical number imprint.

Storage Instruction: Room Temp

Turn Around Time: 3-5 days

CPT Code: 88374 (1)

Scarlet Eligible: No

ALK BY FISH, IF NEG REFLEX TO ROS1 BY FISH

J015-9

Alternate Name: ALK FISH,IF NEG RFX>ROS

Methodology: FISH

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Storage Instruction: Room Temp

Turn Around Time: 3-5 days

CPT Code: 88377 (1)

Scarlet Eligible: No

EML4-ALK by FISH (FDA): Lung

A241-2

Alternate Name: ALK Status

Methodology: FISH

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: Please write surgical number on block and on requisition. Provide fixation time on requisition as well

Storage Instruction: Room Temp

Turn Around Time: 3-5 days

CPT Code: 88377 (1)

Scarlet Eligible: No

ALK-2P23: Anaplastic Lymphoma

B684-2

Alternate Name: ALK-2P23-ALCL

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Place bone marrow into tube, label with patient's name

Storage Instruction: Room Temp

Turn Around Time: 3 days

CPT Code: 88377 (1)

Scarlet Eligible: No

ALL Prognosis Panel by FISH

P260-3

Alternate Name: ALL PROGNOSIS PNL/FISH

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube. Also acceptable: Whole blood in green (sodium heparin) tube.

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 88374 (2)

Scarlet Eligible: No

Allergen Cedar Red (t212), IgE

6305-7

Alternate Name: ALLERGEN CEDAR RED, IGE

Methodology: Immunocap

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Chestnut (t206), IgE

2420-8

Alternate Name: ALLERGEN CHESTNUT, IGE

Methodology: Immunocap

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Herring (f205), IgE

6267-9

Alternate Name: ALLERGEN HERRING, IGE

Methodology: Immunocap

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen White Mulberry (T70), IgE

3835-6

Alternate Name: ALLERGEN WHITE MULBERRY

Methodology: Immunocap

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergy Outdoor Panel

P267-8

Alternate Name: ALLERGY OUTDOOR PANEL

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Turn Around Time: 3 days

CPT Code: 86003 (20)

Scarlet Eligible: Yes

Alkaline Phosphatase, Serum

0185-9

Alternate Name: ALP; Phosphatase, Alkaline

Methodology: Colorimetry

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84075 (1)

Scarlet Eligible: Yes

Alpha-2-Macroglobulin

3882-8

Alternate Name: ALPHA-2-MACROGLOBULIN

Methodology: Nephelometry

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 83883 (1)

Scarlet Eligible: Yes

Alpha-Globin mut. Fetus

B669-3

Alternate Name: ALPHA-GLOBIN, MUT.FETUS

Methodology: Polymerase Chain Reaction

Preferred Specimen: 20 Amniotic Fluid

Collection Instruction: 20mL Amniotic Fluid submitted in a sterile, leak-proof container MATERNAL WHOLE BLOOD SHOULD BE SENT WITH ALL SAMPLE TYPES.

Storage Instruction: Room Temp

Turn Around Time: 23 days

CPT Code: 81257 (1)

Scarlet Eligible: No

Aluminum, Water, Dialysis

1817-6

Alternate Name: ALUMINUM,WATER,ANY

Preferred Specimen: 20 mL Cup Plain

Collection Instruction: Place sample in sterile cup label with patient's name and source

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 84999 (1)

Scarlet Eligible: No

Mitochondrial Antibody (AMA) with titer

0939-9

Alternate Name: AMA, Anti-Mitochondrial Ab

Methodology: Indirect Fluorescence Assay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 4 days

CPT Code: 86381 (1)

Scarlet Eligible: Yes

AML Diagnostic Panel by FISH

P261-1

Alternate Name: AML DIAGNOSTIC BY FISH

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube. Also acceptable: Whole blood in green (sodium heparin) tube.

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 88374 (5)

Scarlet Eligible: No

AML-M4Eo by FISH TC

J380

Alternate Name: AML-M4Eo by FISH TC

Scarlet Eligible: No

AML-M2 (AML/ETO) by Fish- TC ONLY

B891-3

Alternate Name: AML1/ETO,t(8;21) FISH

Methodology: FISH

Preferred Specimen: 5 mL Bone Marrow - Green Top

Collection Instruction: Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube. Also acceptable: Whole blood in green (sodium heparin) tube.

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 88374 (1)

Scarlet Eligible: No

Amylase, Pleural Fluid

6218-2

Alternate Name: AMYLASE, PLEURAL FLD.
Methodology: Spectrophotometry
Preferred Specimen: 1 mL Fluid-source required
Collection Instruction: Submit fluid in sterile cup
Storage Instruction: Refrigerate
Turn Around Time: 7 days
CPT Code: 82150 (1)
Scarlet Eligible: No

Amylase, Urine, Random

0037-2

Alternate Name: AMYLASE, URINE, RANDOM
Methodology: Enzyme Immunoassay
Preferred Specimen: 5 mL Urine Cup
Collection Instruction: Collect urine in cup, label with patients name
Storage Instruction: Refrigerate
Turn Around Time: 4 days
CPT Code: 82150 (1), 82570 (1)
Scarlet Eligible: Yes

Antinuclear Antibody Evaluation (ANA)

0594-2

Alternate Name: ANA
Preferred Specimen: 4 mL SST Tube
Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes
Storage Instruction: Refrigerate
Turn Around Time: 1 days
CPT Code: 86225 (1), 86235 (4)
Scarlet Eligible: Yes

ANA with Reflex to IFA

A233-9

Alternate Name: ANA w/Cascade

Methodology: Indirect Fluorescence Assay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86038 (1)

Scarlet Eligible: Yes

Antinuclear Ab. (ANA) W/Reflex

0038-0

Alternate Name: ANA; FANA

Methodology: Indirect Fluorescence Assay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86038 (1)

Scarlet Eligible: Yes

Anal GC/CT ThinPrep (Non-NY)

K987-8

Alternate Name: ANAL GC/CT THINPREP

Turn Around Time: 7 days

CPT Code: 87491 (1), 87591 (1)

Scarlet Eligible: No

Anemia I Profile

0998-5

Alternate Name: ANEMIA I

Preferred Specimen: SST Tube,Lavender top- EDTA

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82728 (1), 83540 (1), 83550 (1), 85025 (1), 85045 (1)

Scarlet Eligible: Yes

Anemia II Profile

0999-3

Alternate Name: ANEMIA II

Preferred Specimen: SST Tube,Lavender top- EDTA

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82607 (1), 82728 (1), 82746 (1), 83540 (1), 83550 (1), 83615 (1), 85025 (1), 85045 (1)

Scarlet Eligible: Yes

Anemia III Profile

2073-5

Alternate Name: ANEMIA III

Preferred Specimen: SST Tube,Lavender top- EDTA

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82607 (1), 82746 (1), 85025 (1)

Scarlet Eligible: Yes

Anemia Profile, Expanded

9893-9

Alternate Name: ANEMIA SCREEN-EXPANDED

Turn Around Time: 3 days

CPT Code: 82607 (1), 82668 (1), 82728 (1), 82746 (1), 83540 (1), 83550 (1), 85025 (1), 85045 (1)

Scarlet Eligible: Yes

Anthrax Culture Screen

3573-3

Alternate Name: ANTHRAX CULTURE SCRE

Methodology: Cell Culture

Preferred Specimen: Swab-E

Collection Instruction: Collect specimen with swab then place swab into transport carrier; Label with name and source

Turn Around Time: 2 days

CPT Code: 87081 (1)

Scarlet Eligible: No

Anti-Annexin V. Antibodies, IgG, IgM

5706-7

Alternate Name: ANTI-ANNEXIN V

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: Red Top

Collection Instruction: Fill tube, label with patient name

Turn Around Time: 5 days

CPT Code: 83516 (2)

Scarlet Eligible: Yes

Beta-2 Glycoprotein 1 Antibody, IgA

3905-7

Alternate Name: ANTI-BETA2 GLYCO,IgA

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86146 (1)

Scarlet Eligible: Yes

Beta-2 Glycoprotein 1 Antibody, IgG

3904-0

Alternate Name: ANTI-BETA2 GLYCO,IgG

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86146 (1)

Scarlet Eligible: Yes

Beta-2 Glycoprotein 1 Antibody, IgM

3906-5

Alternate Name: ANTI-BETA2 GLYCO,IgM

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86146 (1)

Scarlet Eligible: Yes

Anti CMV Avidity, IgG

J057-1

Alternate Name: ANTI-CMV IgG AVIDITY

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87332 (1)

Scarlet Eligible: Yes

Anti-Epidermal Antibody, S

0184-2

Alternate Name: ANTI-EPIDERMAL ABS.
Methodology: Immunofluorescence Assay
Preferred Specimen: 1 mL Red Top
Collection Instruction: Fill tube, label with patient name
Storage Instruction: Refrigerate
Turn Around Time: 14 days
CPT Code: 86255 (2)
Scarlet Eligible: Yes

Anti-Mullerian Hormone

3707-7

Alternate Name: ANTI-MULLERIAN (AMH)
Methodology: Chemiluminescence
Preferred Specimen: 1 mL SST Tube
Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes
Storage Instruction: Refrigerate
Turn Around Time: 1 days
CPT Code: 82397 (1)
Scarlet Eligible: Yes

Parietal Cell Antibody

0263-4

Alternate Name: ANTI-PARIETAL CELL Ab
Methodology: Indirect Fluorescence Assay
Preferred Specimen: 1 mL SST Tube
Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes
Storage Instruction: Refrigerate
Turn Around Time: 2 days
CPT Code: 86256 (1)
Scarlet Eligible: Yes

Antiphosphatidylserine Antibodies, IgG, IgM, Special Coagulation

5790-1

Alternate Name: ANTI-PHOS-SERINE ABS.

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: Aliquot Tube-Serum

Collection Instruction: Place 1-3 mL of serum in transport tube. Label as serum;refrigerate or freeze as required.

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 86148 (2)

Scarlet Eligible: Yes

Smith Antibody

0851-6

Alternate Name: ANTI-SM(Smith)

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86235 (1)

Scarlet Eligible: Yes

Antithrombin III Activity (AT III)

5714-1

Alternate Name: ANTI-THROMB.IIIfunct.

Methodology: Chromogenic

Preferred Specimen: 2 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 4 days

CPT Code: 85300 (1)

Scarlet Eligible: No

Antithrombin III Antigen (ATIII)

5715-8

Alternate Name: ANTI-THROMBIN III Ag

Methodology: Immunospectrometric

Preferred Specimen: 2 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 4 days

CPT Code: 85301 (1)

Scarlet Eligible: No

Arsenic, Blood

J019-1

Alternate Name: ARSENIC, BLOOD

Methodology: Inductively Coupled Plasma Mass Spectrometry

Preferred Specimen: Blue-Royal Blue Top K2 EDTA Plasma(*LAVENDER LINE)

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82175 (1)

Scarlet Eligible: Yes

Arthritis Profile

0325-1

Alternate Name: ARTHRITIC SCREEN

Preferred Specimen: SST Tube,Lavender top- EDTA

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 84550 (1), 85651 (1), 86038 (1), 86060 (1), 86431 (1)

Scarlet Eligible: Yes

Arthritis Profile - Expanded

9892-1

Alternate Name: ARTHRITIC SCREEN-EXPAND

Preferred Specimen: SST Tube,Lavender top- EDTA

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 84550 (1), 85651 (1), 86038 (1), 86140 (1), 86200 (1), 86225 (1), 86235 (4), 86431 (1)

Scarlet Eligible: Yes

Ashkenazi Jewish 11

P658-8

Alternate Name: ASHKENAZI JEWISH 11

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81412 (1)

Scarlet Eligible: Yes

Ashkenazi Jewish 18

P659-6

Alternate Name: ASHKENAZI JEWISH 18

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81412 (1)

Scarlet Eligible: Yes

Ashkenazi Jewish 25

Q161-2

Alternate Name: ASHKENAZI JEWISH 25

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81412 (1)

Scarlet Eligible: Yes

Ashkenazi Jewish DNA Carrier

3267-2

Alternate Name: ASHKENAZI JEWISH DNA

Methodology: Polymerase Chain Reaction

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 14 days

CPT Code: 81412 (1)

Scarlet Eligible: Yes

Ashkenazi Jewish DNA (NY patients only)

P223-1

Alternate Name: ASHKENAZI JEWISH DNA-NY

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81200 (1), 81209 (1), 81242 (1), 81251 (1), 81255 (1), 81260 (1), 81290 (1), 81330 (1)

Scarlet Eligible: Yes

Antistreptolysin O (ASO) Titer

0042-2

Alternate Name: ASO, Antistreptolysin, ASLO

Methodology: Immunoturbidimetric

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86060 (1)

Scarlet Eligible: Yes

Aspergillus Antibodies, Immunodiffusion, Serum

1741-8

Alternate Name: ASPERGILLUS Ab.PANEL

Methodology: Immunodiffusion

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 10 days

CPT Code: 86606 (3)

Scarlet Eligible: Yes

Salicylates, Serum

0526-4

Alternate Name: ASPIRIN

Methodology: Enzyme Immunoassay

Preferred Specimen: Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80179 (1)

Scarlet Eligible: Yes

Atenolol Level

2392-9

Alternate Name: ATENOLOL LEVEL

Methodology: Liquid Chromatography Tandem Mass Spectrometry

Preferred Specimen: 1 mL Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 7 days

CPT Code: 80375 (1)

Scarlet Eligible: Yes

Atherosclerotic Risk Expanded

6802-3

Alternate Name: ATHEROSCLEROTIC RISK EXPANDED

Preferred Specimen: Aliquot Plasma-Citrated-Blue Top,Red Top,Lavender top- EDTA,SST Tube

Collection Instruction: Submit 1 lav, 3 serum aliquots from 2 red tops & 4 frozen citrated plasma aliquots from 4 light blue tops.

Turn Around Time: 5 days

CPT Code: 80061 (1), 81240 (1), 81241 (1), 81291 (1), 83090 (1), 83695 (1), 83721 (1), 85240 (1), 85280 (1), 85379 (1), 85384 (1), 85598 (1), 85610 (1), 85613 (2), 85670 (1), 85730 (2), 85732 (8), 86141 (1), 86146 (3), 86147 (3)

Scarlet Eligible: Yes

Antithrombin III Anti-Xa Assay

5723-2

Alternate Name: ATIII ANTI-XA

Methodology: Chromogenic

Preferred Specimen: 2 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 5 days

CPT Code: 85300 (1)

Scarlet Eligible: No

11Q22.3 BY FISH-TC ONLY

A879-9

Alternate Name: ATM by FISH, TC Only

Methodology: FISH

Preferred Specimen: 5 mL Peripheral Blood - Green Top

Collection Instruction: Fill tube, invert 8-10 times (DO NOT SHAKE!)

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88373 (1)

Scarlet Eligible: No

Automated UroVysion By FISH

5030-2

Alternate Name: AUTOMATED UROVYSION by FISH

Methodology: FISH

Preferred Specimen: 50 mL ThinPrep Cytolyt Solution

Collection Instruction: Collect voided urine in Thinprep Cytolyte solution

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88121 (1)

Scarlet Eligible: No

AV Panel w/o Lacto Reuteri and Rhamnosis (Non NY)

L767-3

Alternate Name: AV Panel

Turn Around Time: 7 days

CPT Code: 87640 (1), 87653 (1), 87798 (1), 87799 (1), 87651 (1), 99999 (1)

Scarlet Eligible: No

Abacavir Genotype HLA B57

J096-9

Alternate Name: Abacavir Genotype HLA B57

Methodology: Polymerase Chain Reaction

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 13 days

CPT Code: 81381 (1)

Scarlet Eligible: No

Abnormal APTT Assessment

5851-1

Alternate Name: Abnormal Activated Partial Thromboplastin Time Assessment

Preferred Specimen: Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.

Turn Around Time: 4 days

CPT Code: 85240 (1), 85250 (1), 85270 (1), 85598 (1), 85610 (1), 85613 (2), 85730 (2), 85732 (8)

Scarlet Eligible: No

Abnormal PT and APTT Assessment

5853-7

Alternate Name: Abnormal Prothrombin Time/Activated Partial Thromboplastin Time Assessment

Preferred Specimen: Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 4 days

CPT Code: 85210 (1), 85220 (1), 85230 (1), 85240 (1), 85250 (1), 85260 (1), 85270 (1), 85379 (1), 85384 (1), 85598 (1), 85610 (1), 85611 (8), 85613 (2), 85670 (1), 85730 (2), 85732 (8)

Scarlet Eligible: No

Eosinophil Count, Absolute, Blood

0381-4

Alternate Name: Absolute Eosinophil Count; Eosinophil Count, EOS

Methodology: Flow Cytometry

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 85008 (1)

Scarlet Eligible: Yes

Activation Markers

5811-5

Alternate Name: Activation Markers

Preferred Specimen: 2.0 Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect 2 full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into 3 plastic transfer tubes. Label with date/time collected, PLASMA, patient's name and freeze.

Turn Around Time: 10 days

CPT Code: 83520 (2), 85379 (1)

Scarlet Eligible: No

Acute 10-Color Leukemia/Lymphoma Panel

TH61

Alternate Name: Acute 10-Color Leukemia/Lymphoma Panel

Scarlet Eligible: No

Acute 10-Color Leukemia/Lymphoma Panel TC Only

TH64

Alternate Name: Acute 10-Color Leukemia/Lymphoma Panel TC Only

Scarlet Eligible: No

Acute Leukemia + Myeloid Disorders Analysis (Technical Only) and StormPath BMM Report Template R735-2

Alternate Name: Acute Flow TC and BMM

Turn Around Time: 1 days

CPT Code: 99999 (1), 88184 (1), 88185 (27)

Scarlet Eligible: No

Acute Leukemia/MDS IHC Panel 5116-9

Alternate Name: Acute Leukemia/MDS IHC

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88341 (10), 88342 (1)

Scarlet Eligible: No

AML - M2 (AML/ETO) (RUNX1/RUNX1T1) by FISH 5024-5

Alternate Name: Acute Myeloid Leukemia with t(8;21)

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube. Also acceptable: Whole blood in green (sodium heparin) tube.

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 88374 (1)

Scarlet Eligible: No

AML - M2 (AML/ETO) by RT-PCR

5034-4

Alternate Name: Acute Myeloid Leukemia with t(8;21)

Methodology: Real-Time Polymerase Chain Reaction

Preferred Specimen: 2 mL Bone Marrow - Lavender Top

Collection Instruction: Place bone marrow into tube, label with patient's name. Send in bone marrow kit.

Turn Around Time: 5 days

CPT Code: 81401 (1)

Scarlet Eligible: No

AML - M4Eo by FISH

5025-2

Alternate Name: Acute Myelomonocytic Leukemia with Eosinophilia

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube. Also acceptable: Whole blood in green (sodium heparin) tube

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 88374 (1)

Scarlet Eligible: No

AML - M3 (PML/RARA) by RT-PCR

5261-3

Alternate Name: Acute Promyelocytic Leukemia with t(15;17)

Methodology: Real-Time Polymerase Chain Reaction

Preferred Specimen: 2 mL Bone Marrow - Lavender Top

Collection Instruction: Place bone marrow into tube, label with patient's name. Send in bone marrow kit.

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 81315 (1)

Scarlet Eligible: No

AML - M3 (PML/RARA) by FISH

5260-5

Alternate Name: Acute Promyelocytic Leukemia with t(15;17)

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube. Also acceptable: Whole blood in green (sodium heparin) tube.

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88374 (1)

Scarlet Eligible: No

Acylglycines, Quant. Urine

6310-7

Alternate Name: Acylglycines, Quant.

Methodology: Gas Chromat. Mass Spectrometry

Preferred Specimen: 13 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 16 days

CPT Code: 82542 (1)

Scarlet Eligible: Yes

OnkoSight Advanced CNS NGS Panel

TH57-3

Alternate Name: Advanced CNS NGS

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Storage Instruction: Room Temp

Turn Around Time: 10 days

CPT Code: 81445 (1)

Scarlet Eligible: No

OnkoSight Advanced Gastrointestinal Stromal Tumor NGS Panel

TH58-1

Alternate Name: Advanced GIST NGS Panel

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Storage Instruction: Room Temp

Turn Around Time: 10 days

CPT Code: 81445 (1)

Scarlet Eligible: No

OnkoSight Advanced GYN Tumor NGS

TH53-2

Alternate Name: Advanced GYN Tumor NGS

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Storage Instruction: Room Temp

Turn Around Time: 10 days

CPT Code: 81445 (1)

Scarlet Eligible: No

OnkoSight Advanced Lung Cancer NGS Panel

TH49-0

Alternate Name: Advanced Lung Cancer NGS

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Storage Instruction: Room Temp

Turn Around Time: 10 days

CPT Code: 81445 (1)

Scarlet Eligible: No

OnkoSight Advanced Melanoma NGS Panel

TH59-9

Alternate Name: Advanced Melanoma NGS
Methodology: Genotyping by Next Generation Sequencing
Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue
Storage Instruction: Room Temp
Turn Around Time: 10 days
CPT Code: 81445 (1)
Scarlet Eligible: No

OnkoSight Advanced Prostate NGS Panel

TH48-2

Alternate Name: Advanced Prostate NGS
Methodology: Genotyping by Next Generation Sequencing
Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue
Storage Instruction: Room Temp
Turn Around Time: 10 days
CPT Code: 81445 (1)
Scarlet Eligible: No

OnkoSight Advanced Solid Tumor NGS Panel

TJ16-5

Alternate Name: Advanced Solid TumorNGS
Methodology: Genotyping by Next Generation Sequencing
Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue
Storage Instruction: Room Temp
Turn Around Time: 10 days
CPT Code: 81445 (1)
Scarlet Eligible: No

Aerobic Vaginitis Profile (Non NY)

R732-9

Alternate Name: Aerobic Vaginitis Profile
Turn Around Time: 7 days
CPT Code: 87651 (1), 87653 (1), 87798 (1), 87640 (1), 87799 (1), 99999 (1)
Scarlet Eligible: No

Aluminum, Serum

0378-0

Alternate Name: Al, Serum

Methodology: Inductively Coupled Plasma Mass Spectrometry

Preferred Specimen: Blue-Royal Blue Trace Element Serum (**RED LINE**)

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 7 days

CPT Code: 82108 (1)

Scarlet Eligible: Yes

Albumin, Serum

0033-1

Alternate Name: Albumin

Methodology: Colorimetry

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82040 (1)

Scarlet Eligible: Yes

Alcohol, Methyl

A473-1

Alternate Name: Alcohol, Methyl

Methodology: Gas Chromatography/Flame Ionization

Preferred Specimen: Grey top- Sodium Fluoride

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 80320 (1)

Scarlet Eligible: Yes

Allergen Lime (RF306), IgE

3522-0

Alternate Name: Allergen Lime (Rf306)

Methodology: Immunocap

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Yeast Bakers, IgE

0813-6

Alternate Name: Allergen Yeast Bakers

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Digoxin, Serum

0083-6

Alternate Name: Allocar; Cardioreg; Digacin; Lanocor; Lanoxicaps; Lenoxin; Lanoxin; Purgoxin

Methodology: Enzyme Immunoassay

Preferred Specimen: Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80162 (1)

Scarlet Eligible: Yes

Allergen Almond Nut (f20), IgE

0815-1

Alternate Name: Almond Nut

Methodology: Immunocap

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Alpha Lactalbumin (f76), IgE

1406-8

Alternate Name: Alpha Lactalbumin

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86008 (1)

Scarlet Eligible: Yes

Vitamin E, Serum

0447-3

Alternate Name: Alpha Tocopherol;Tocopherol;a-Tocopherol

Methodology: High Pressure Liquid Chromatography

Preferred Specimen: 2 mL Amber Serum

Collection Instruction: Draw SST or RED, let clot for 8-10 minutes. Centrifuge and pour into either an aliquot tube (must be covered in foil) or poured into an amber tube and label with patient name and specimen type (serum).

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 84446 (1)

Scarlet Eligible: Yes

Triple Screen AFP3 (AFP, UE3, HCG)

0825-0

Alternate Name: Alpha-Fetoprotein 3 (AFP), including Estriol/Human Chorionic Gonadotropin

Methodology: Chemiluminescence

Preferred Specimen: 2 mL SST Tube

Collection Instruction: MUST BE 15-22 WEEKS AND 6 DAYS. PROVIDE GESTATIONAL INFO.

Storage Instruction: Refrigerate

Turn Around Time: 4 days

CPT Code: 81510 (1)

Scarlet Eligible: Yes

Alport Syndrome, Autosomal Recessive

J316-1

Alternate Name: Alport Syndrome, Autosomal Recessive

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81408 (1)

Scarlet Eligible: No

Allergen Alternaria Tenuis Mold (m6), IgE

0660-1

Alternate Name: Alternaria Tenuis Mold

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Theophylline, Serum

0380-6

Alternate Name: Aminophylline;Elixophyllin®;Ethylenediamine;Phyllcontin;Sio-Phyllin®;Sustaire®;Theo-Dur®;Theolair®,ç;Theospan®;Truphylline®

Methodology: Colorimetry

Preferred Specimen: Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 80198 (1)

Scarlet Eligible: Yes

Barbiturates, Urine, Screen (200)

0921-7

Alternate Name: Amobarb; Butalbital; Pentobarb; Secobarb;Phenobarbital

Methodology: Enzyme Immunoassay

Preferred Specimen: 5 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80307 (1)

Scarlet Eligible: Yes

Amylase, Serum

0036-4

Alternate Name: Amylase, Serum

Methodology: Enzyme Immunoassay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82150 (1)

Scarlet Eligible: Yes

Anal CT on ThinPrep (Non NY)

J792-3

Alternate Name: Anal CT

Methodology: Polymerase Chain Reaction

Preferred Specimen: 3 mL Anal ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 5-7 days

CPT Code: 87491 (1)

Scarlet Eligible: No

Anal GC on ThinPrep (Non NY)

J791-5

Alternate Name: Anal GC

Methodology: Polymerase Chain Reaction

Preferred Specimen: 3 mL Anal ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 5-7 days

CPT Code: 87591 (1)

Scarlet Eligible: No

Allergen Anchovy (rf313), IgE

1463-9

Alternate Name: Anchovy (Rf313)

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

PCP (6-Phencyclidine), Urine, Screen (25)

0926-6

Alternate Name: Angel Dust; Crystal Joint; Elephant Tranquilizers; Goon; Hog; Kay Jay; Killer weed; Mist; PCP; Peace Pills; Peace weed; Rocket Fuel; Sheets; Sherm; Snorts; Soma®; Supergrass; Wickistick; PCP

Methodology: Enzyme Immunoassay

Preferred Specimen: 5 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80307 (1)

Scarlet Eligible: Yes

Animal Mix (ex1), IgE

2180-8

Alternate Name: Animal Mix

Methodology: Immunocap

Preferred Specimen: 5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86005 (1)

Scarlet Eligible: Yes

Animals Screen(6)

0902-7

Alternate Name: Animals Screen(6)

Preferred Specimen: 5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (5), 82785 (1)

Scarlet Eligible: Yes

Endomysial Antibody (IgA) Screen with Titer

1206-2

Alternate Name: Anti-Endomysial Ab

Methodology: Indirect Fluorescence Assay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86231 (1)

Scarlet Eligible: Yes

Anti-Jo-1, IgG

T974-5

Alternate Name: Anti-Jo-1,IgG

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Room Temp

Turn Around Time: 1 days

CPT Code: 86235 (1)

Scarlet Eligible: Yes

Double-Stranded DNA (dsDNA) Antibody

0364-0

Alternate Name: Anti-dsDNA

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86225 (1)

Scarlet Eligible: Yes

Antibody Screen (Ind.Coombs)

0039-8

Alternate Name: Antiglobulin Test,Indirect

Methodology: Manual

Preferred Specimen: Pink Tube-EDTA

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86850 (1)

Scarlet Eligible: Yes

Apolipoprotein A-1 (Apo A-1)

0451-5

Alternate Name: Apo A-1

Methodology: Immunoturbidimetric

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82172 (1)

Scarlet Eligible: Yes

Apolipoprotein A-1 and B (Apo A1/B)

0895-3

Alternate Name: Apo A1/B

Methodology: Immunoturbidimetric

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82172 (2)

Scarlet Eligible: Yes

Apolipoprotein B (Apo B)

0457-2

Alternate Name: Apo B

Methodology: Immunoturbidimetric

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82172 (1)

Scarlet Eligible: Yes

Allergen Apple (f49), IgE

1134-6

Alternate Name: Apple

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Apricot (f237), IgE

2344-0

Alternate Name: Apricot

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Vaginitis/Vaginosis by RT-PCR, Aptima Tube

L050-4

Alternate Name: Aptima Swab Vaginitis

Preferred Specimen: Swab Aptima Orange label

Collection Instruction: Collect specimen-break swab shaft into transport tube, then label with name and source

Turn Around Time: 7 days

CPT Code: 87512 (1), 87799 (2), 87798 (2), 99999 (1), 87481 (2), 87661 (1)

Scarlet Eligible: No

Arthrogryposis, Mental Retardation, and Seizures

J297-3

Alternate Name: Arthrogryposis mental retardation and seizures

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81479 (1)

Scarlet Eligible: Yes

Arsenic, Urine, 24 Hours

0427-5

Alternate Name: As, Urine, 24 Hours

Methodology: Inductively Coupled Plasma Mass Spectrometry

Preferred Specimen: 10 mL Urine Container - 24hr

Collection Instruction: Discard first morning sample, then collect all urines after for 24 hrs. Collect urine in acid-washed or metal-free plastic container. On completion label with collection time and patient's name; record total volume on Requisition

Storage Instruction: Refrigerate

Turn Around Time: 8 days

CPT Code: 82175 (1)

Scarlet Eligible: No

Vitamin C

0291-5

Alternate Name: Ascorbic Acid

Methodology: High Pressure Liquid Chromatography

Preferred Specimen: 1 mL Amber Serum

Collection Instruction: Draw SST or RED, let clot for 8-10 minutes. Centrifuge and pour into either an aliquot tube (must be covered in foil) or poured into an amber tube and label with patient name and specimen type (serum). FREEZE.

Storage Instruction: STRICT Frozen

Turn Around Time: 3 days

CPT Code: 82180 (1)

Scarlet Eligible: No

Lorazepam (Ativan)), Urine

5924-6

Alternate Name: Ativan

Methodology: Liquid Chromatography Tandem Mass Spectrometry

Preferred Specimen: 1 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 8 days

CPT Code: 80346 (1)

Scarlet Eligible: No

Allergen Australian Pine, IgE (t73)

3837-2

Alternate Name: Australian Pine

Methodology: Immunocap

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

ISH Single Antibody Stain Only

5289-4

Alternate Name: Available probes: Kappa, Lambda, EBER, HPV high/low risk

Methodology: In Situ Hybridization

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88365 (1)

Scarlet Eligible: No

Allergen Avocado (f96), IgE

1945-5

Alternate Name: Avocado

Methodology: Immunocap

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

B. fragilis by RT PCR

J216-3

Alternate Name: B. FRAGILIS BY PCR

Methodology: Polymerase Chain Reaction

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 5-7 days

CPT Code: 87798 (1)

Scarlet Eligible: No

Beta-2-Microglobulin, Serum

0262-6

Alternate Name: B2-Microglobulin, Serum

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 82232 (1)

Scarlet Eligible: Yes

Beta-2-Microglobulin, Urine, Random

1754-1

Alternate Name: B2-Microglobulin, Urine, Random

Methodology: Chemiluminescence

Preferred Specimen: 4 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Turn Around Time: 2 days

CPT Code: 82232 (1)

Scarlet Eligible: Yes

Vaginitis/Vaginosis, DNA - BD Affirm

5093-0

Alternate Name: BACTERIAL VAGINOSIS SCR.N.,DNA

Methodology: Nucleic Acid Hybridization

Preferred Specimen: Affirm(BD) VPIII

Collection Instruction: Collect specimen-break swab shaft into BD Affirm, then label with name and source

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 87480 (1), 87510 (1), 87660 (1)

Scarlet Eligible: No

Basic Antiphospholipid/Lupus

5994-9

Alternate Name: BASIC ANTI-PHOS/LUPUS

Turn Around Time: 4 days

CPT Code: 85025 (1), 85597 (1), 85598 (1), 85610 (1), 85611 (8), 85613 (2), 85730 (2), 85732 (8), 86146 (3), 86147 (3)

Scarlet Eligible: No

Basic Food Allergy Panel

4633-4

Alternate Name: BASIC FOOD PANEL

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82785 (1), 86003 (10)

Scarlet Eligible: Yes

BCLC (3Q27) by FISH TC

B348

Alternate Name: BCLC (3Q27) by FISH TC

Scarlet Eligible: No

BCR/ABL by FISH TC t(9:22)

A880

Alternate Name: BCR/ABL by FISH TC t(9:22)

Scarlet Eligible: No

BCR/ABL Quantification by RT-PCR

5858-6

Alternate Name: BCR/ABL for CML Monitoring, Philadelphia Chromosome Analysis

Methodology: Polymerase Chain Reaction

Preferred Specimen: 5 mL Peripheral Blood - Lavender Top

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient's name

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 81206 (1), 81207 (1)

Scarlet Eligible: Yes

Vaginitis/Vaginosis by DNA probe, Thin Prep

P576-2

Alternate Name: BD AFFIRM THINPREP PNL

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 87480 (1), 87510 (1), 87660 (1)

Scarlet Eligible: No

Bean Profile

Q032-5

Alternate Name: BEAN PROFILE

Methodology: Immunocap

Preferred Specimen: 5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82785 (1), 86003 (4)

Scarlet Eligible: Yes

Benzodiazepines, Urine, Screen (100)

3233-4

Alternate Name: BENZODIAZEPINES,URN.SCRN.(100)

Methodology: Enzyme Immunoassay

Preferred Specimen: 5 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80307 (1)

Scarlet Eligible: Yes

Beta Carotene, Total

0512-4

Alternate Name: BETA CAROTENE,TOTAL

Methodology: High Pressure Liquid Chromatography

Preferred Specimen: 2 mL Amber Serum

Collection Instruction: Draw SST or RED, let clot for 8-10 minutes. Centrifuge and pour into either an aliquot tube (must be covered in foil) or poured into an amber tube and label with patient name and specimen type (serum).

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82380 (1)

Scarlet Eligible: Yes

Beta-2-Glycoprotein I (G,A,M)

2675-7

Alternate Name: BETA-2-GLYCOPROTEIN I (G,A,M)

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86146 (3)

Scarlet Eligible: Yes

Beta-Galactosidase, Leukocytes

3514-7

Alternate Name: BETA-GALACTOSIDASE

Methodology: Fluorometric

Preferred Specimen: Yellow top- ACD

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient's name

Storage Instruction: Refrigerate

Turn Around Time: 17 days

CPT Code: 82657 (1)

Scarlet Eligible: No

Biopsy Level II, PC ONLY

5233-2

Alternate Name: BIOPSY LEVEL II, PC ONLY

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Collect specimen place in formalin jar write patient name and site on jar

Turn Around Time: 1 days

CPT Code: 88302 (1)

Scarlet Eligible: No

Biopsy Level III, PC ONLY

5234-0

Alternate Name: BIOPSY LEVEL III, PC ONLY

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Collect specimen place in formalin jar write patient name and site on jar

Turn Around Time: 1 days

CPT Code: 88304 (1)

Scarlet Eligible: No

Biopsy Level IV, PC ONLY

5235-7

Alternate Name: BIOPSY LEVEL IV, PC ONLY

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Collect specimen place in formalin jar write patient name and site on jar

Turn Around Time: 1 days

CPT Code: 88305 (1)

Scarlet Eligible: No

Biopsy Level V, PC ONLY

5236-5

Alternate Name: BIOPSY LEVEL V, PC ONLY

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Collect specimen place in formalin jar write patient name and site on jar

Turn Around Time: 1 days

CPT Code: 88307 (1)

Scarlet Eligible: No

Biopsy Pod 2-@ _

5002-1

Alternate Name: BIOPSY POD 2-@(AREA)

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 1 days

CPT Code: 88305 (2), 88311 (2)

Scarlet Eligible: No

Biopsy Pod 3-@ _

5003-9

Alternate Name: BIOPSY POD 3-@(AREA)

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 1 days

CPT Code: 88305 (3), 88311 (3)

Scarlet Eligible: No

Biopsy Pod -@ _

5000-5

Alternate Name: BIOPSY POD @

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 2-3 days

CPT Code: 88305 (1)

Scarlet Eligible: No

Biopsy, general-@ (AREA)

4000-6

Alternate Name: BIOPSY, GENERAL

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Storage Instruction: Room Temp

Turn Around Time: 2 days

CPT Code: 88305 (1)

Scarlet Eligible: No

Biopsy Gen 2-@ (AREA)

4007-1

Alternate Name: BIOPSY, GENERAL X2**Methodology:** Histology**Preferred Specimen:** Biopsy In Formalin**Collection Instruction:** Submit in 10% Neutral Buffered Formalin, varied jar sizes.**Turn Around Time:** 1 days**CPT Code:** 88305 (2)**Scarlet Eligible:** No**Biopsy Gen 3@AREA**

4003-0

Alternate Name: BIOPSY, GENERAL X3**Methodology:** Histology**Preferred Specimen:** Biopsy In Formalin**Collection Instruction:** Submit in 10% Neutral Buffered Formalin, varied jar sizes.**Turn Around Time:** 1 days**CPT Code:** 88305 (3)**Scarlet Eligible:** No**Biopsy Gen 4-@ (AREA)**

4004-8

Alternate Name: BIOPSY, GENERAL X4**Methodology:** Histology**Preferred Specimen:** Biopsy In Formalin**Collection Instruction:** Submit in 10% Neutral Buffered Formalin, varied jar sizes.**Turn Around Time:** 1 days**CPT Code:** 88305 (4)**Scarlet Eligible:** No

Biopsy Gen 5-@ (AREA)

4005-5

Alternate Name: BIOPSY, GENERAL X5**Methodology:** Histology**Preferred Specimen:** Biopsy In Formalin**Collection Instruction:** Submit in 10% Neutral Buffered Formalin, varied jar sizes.**Turn Around Time:** 1 days**CPT Code:** 88305 (5)**Scarlet Eligible:** No**Biopsy Gen 6-@ (AREA)**

4006-3

Alternate Name: BIOPSY, GENERAL X6**Methodology:** Histology**Preferred Specimen:** Biopsy In Formalin**Collection Instruction:** Submit in 10% Neutral Buffered Formalin, varied jar sizes.**Turn Around Time:** 1 days**CPT Code:** 88305 (6)**Scarlet Eligible:** No**Biopsy, GI (1 Jar)**

5631-7

Alternate Name: BIOPSY, GI**Methodology:** Histology**Preferred Specimen:** Biopsy In Formalin**Collection Instruction:** Submit in 10% Neutral Buffered Formalin, varied jar sizes.**Turn Around Time:** 1 days**CPT Code:** 88305 (1)**Scarlet Eligible:** No

Biopsy, GI (10 JARS)

A296-6

Alternate Name: BIOPSY, GI (10 JARS)**Methodology:** Histology**Preferred Specimen:** Biopsy In Formalin**Collection Instruction:** Submit in 10% Neutral Buffered Formalin, varied jar sizes.**Turn Around Time:** 1 days**CPT Code:** 88305 (10)**Scarlet Eligible:** No**Biopsy, GI (2 Jars)**

5632-5

Alternate Name: BIOPSY, GI (2 JARS)**Methodology:** Histology**Preferred Specimen:** Biopsy In Formalin**Collection Instruction:** Submit in 10% Neutral Buffered Formalin, varied jar sizes.**Turn Around Time:** 1 days**CPT Code:** 88305 (2)**Scarlet Eligible:** No**Biopsy, GI (3 Jars)**

5633-3

Alternate Name: BIOPSY, GI (3 JARS)**Methodology:** Histology**Preferred Specimen:** Biopsy In Formalin**Collection Instruction:** Submit in 10% Neutral Buffered Formalin, varied jar sizes.**Turn Around Time:** 1 days**CPT Code:** 88305 (3)**Scarlet Eligible:** No

Biopsy, GI (4 Jars)

5634-1

Alternate Name: BIOPSY, GI (4 JARS)**Methodology:** Histology**Preferred Specimen:** Biopsy In Formalin**Collection Instruction:** Submit in 10% Neutral Buffered Formalin, varied jar sizes.**Turn Around Time:** 1 days**CPT Code:** 88305 (4)**Scarlet Eligible:** No**Biopsy, GI (5 Jars)**

5635-8

Alternate Name: BIOPSY, GI (5 JARS)**Methodology:** Histology**Preferred Specimen:** Biopsy In Formalin**Collection Instruction:** Submit in 10% Neutral Buffered Formalin, varied jar sizes.**Turn Around Time:** 1 days**CPT Code:** 88305 (5)**Scarlet Eligible:** No**Biopsy, GI (6 Jars)**

5636-6

Alternate Name: BIOPSY, GI (6 JARS)**Methodology:** Histology**Preferred Specimen:** Biopsy In Formalin**Collection Instruction:** Submit in 10% Neutral Buffered Formalin, varied jar sizes.**Turn Around Time:** 1 days**CPT Code:** 88305 (6)**Scarlet Eligible:** No

Biopsy, GI (7 JARS)

A293-3

Alternate Name: BIOPSY, GI (7 JARS)**Methodology:** Histology**Preferred Specimen:** Biopsy In Formalin**Collection Instruction:** Submit in 10% Neutral Buffered Formalin, varied jar sizes.**Turn Around Time:** 2 days**CPT Code:** 88305 (7)**Scarlet Eligible:** No**Biopsy, GI (8 JARS)**

A294-1

Alternate Name: BIOPSY, GI (8 JARS)**Methodology:** Histology**Preferred Specimen:** Biopsy In Formalin**Collection Instruction:** Submit in 10% Neutral Buffered Formalin, varied jar sizes.**Turn Around Time:** 1 days**CPT Code:** 88305 (8)**Scarlet Eligible:** No**Biopsy, GI (9 JARS)**

A295-8

Alternate Name: BIOPSY, GI (9 JARS)**Methodology:** Histology**Preferred Specimen:** Biopsy In Formalin**Collection Instruction:** Submit in 10% Neutral Buffered Formalin, varied jar sizes.**Turn Around Time:** 1 days**CPT Code:** 88305 (9)**Scarlet Eligible:** No

Biopsy, Level VI Global

5237-3

Alternate Name: BIOPSY, LEVEL VI

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Collect specimen place in formalin jar write patient name and site on jar

Turn Around Time: 1 days

CPT Code: 88309 (1)

Scarlet Eligible: No

Liver Biopsy

5009-6

Alternate Name: BIOPSY, LIVER

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 1 days

CPT Code: 88307 (1), 88313 (4)

Scarlet Eligible: No

Bleeding Diatheses Evaluation Profile I

5850-3

Alternate Name: BLEEDING DIATHESES EVAL1

Preferred Specimen: 2 mL Light blue top,Lavender top- EDTA

Turn Around Time: 1 days

CPT Code: 85025 (1), 85610 (1), 85730 (1)

Scarlet Eligible: Yes

Bleeding Diathesis Evaluation Profile II

5997-2

Alternate Name: BLEEDING DIATHESIS EVALUATION

Preferred Specimen: Aliquot Plasma-Citrated-Blue Top,Lavender top- EDTA

Collection Instruction: Collect one lavender top tube and maintain at room temperature. Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 9 days

CPT Code: 85025 (1), 85240 (1), 85245 (1), 85246 (1), 85250 (1), 85270 (1), 85290 (1), 85610 (1), 85730 (1)

Scarlet Eligible: No

Bloom Syndrome Carrier Test

A658-7

Alternate Name: BLOOM SYNDROME

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81209 (1)

Scarlet Eligible: Yes

Bone Marrow Morph.- Tech Only

B846-7

Alternate Name: BM MORP- TC ONLY

Methodology: Histology

Preferred Specimen: Bone Marrow - Clot in Formalin,Bone Marrow Core in Formalin,Unstained Slide

Collection Instruction: Submit in 10% Neutral Biffered Formalin, varied jar sixes. Submit core/clot in formalin. Send smear in slide holder.

Storage Instruction: Room Temp

Turn Around Time: 1-2 days

CPT Code: 88305 (2), 88311 (1), 88313 (4)

Scarlet Eligible: No

Bone Marrow Morph. Clot Only- Tech Only

B843-4

Alternate Name: BM MORP. CLOT ONLY-TC

Methodology: Histology

Preferred Specimen: Bone Marrow - Clot in Formalin

Collection Instruction: Collect clot section place in formalin, label with patient's name

Turn Around Time: 1-2 days

CPT Code: 88305 (1), 88313 (2)

Scarlet Eligible: No

Bone Marrow Morph. Core Only- Tech Only

B845-9

Alternate Name: BM MORP. CORE ONLY- TC

Methodology: Histology

Preferred Specimen: Bone Marrow Core in Formalin

Collection Instruction: Collect specimen place in formalin jar, label with patient's name and site

Storage Instruction: Room Temp

Turn Around Time: 1-2 days

CPT Code: 88305 (1), 88311 (1), 88313 (2)

Scarlet Eligible: No

Bone Marrow Morph. Smear Only- Tech Only

B844-2

Alternate Name: BM MORP. SMEAR ONLY-TC

Methodology: Histology

Preferred Specimen: Unstained Slide

Collection Instruction: Collect specimen and place on charged slide, label with pt name, place in slide holder for transport.

Storage Instruction: Room Temp

Turn Around Time: 1-2 days

CPT Code: 88313 (1)

Scarlet Eligible: No

Basic Metabolic Profile (BMP)

2555-1

Alternate Name: BMP, CHEM 7, SMA 7

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80048 (1)

Scarlet Eligible: Yes

B-Type Natriuretic Peptide (BNP)

3553-5

Alternate Name: BNP

Methodology: Chemiluminescence

Preferred Specimen: 1 mL White Top (PPT)

Collection Instruction: After collection centrifuge at high speed for 8-10 min. Remove plasma and submit in aliquot tube. Label it plasma

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 83880 (1)

Scarlet Eligible: No

Bone Marrow Morphology

5199-5

Alternate Name: BONE MARROW ANALYSIS

Methodology: Histology

Preferred Specimen: Bone Marrow - Clot in Formalin,Unstained Slide,Bone Marrow Core in Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 1-2 days

CPT Code: 85097 (1), 88305 (2), 88311 (1), 88313 (4)

Scarlet Eligible: No

Bone Marrow Analysis - Clot Only

5202-7

Alternate Name: BONE MARROW ANALYSIS (CLOT)

Methodology: Histology

Preferred Specimen: Bone Marrow - Clot in Formalin

Collection Instruction: Collect specimen place in formalin jar, label with patient's name and site

Turn Around Time: 1-2 days

CPT Code: 88305 (1), 88313 (2)

Scarlet Eligible: No

Bone Marrow Morphology - Core Only

5200-1

Alternate Name: BONE MARROW ANALYSIS (CORE)

Methodology: Histology

Preferred Specimen: Bone Marrow Core in Formalin

Collection Instruction: Collect specimen place in formalin jar, label with patient's name and site

Turn Around Time: 1-2 days

CPT Code: 88305 (1), 88311 (1), 88313 (2)

Scarlet Eligible: No

Bone Marrow Slide Consult

5207-6

Alternate Name: BONE MARROW SLIDE CONSULT.

Methodology: Histology

Preferred Specimen: Stained Slides

Turn Around Time: 1-2 days

CPT Code: 88321 (1)

Scarlet Eligible: No

X/Y Bone Marrow Transplant Monitoring

5066-6

Alternate Name: BONE MARROW X/Y CHROM(FISH)

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube. Also acceptable:Whole blood in green (sodium heparin) tube.

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88377 (1)

Scarlet Eligible: No

Bone Marrow Analysis - Smears

5211-8

Alternate Name: BONE MARROW/SMEARS ONLY

Methodology: Histology

Preferred Specimen: Unstained Slide

Collection Instruction: Client to collect specimen, prepare smears, label with patient's name, place in slide carrier.

Turn Around Time: 1-2 days

CPT Code: 85097 (1), 88313 (1)

Scarlet Eligible: No

Bone/Tissue Culture

P210-8

Alternate Name: BONE/TISSUE CULTURE

Methodology: Bacterial Culture

Preferred Specimen: Swab-E

Collection Instruction: Collect specimen with swab then place swab into transport carrier; Label with name and source

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87205 (1), 87070 (1), 87075 (1)

Scarlet Eligible: No

BRAF, if neg reflex to MLH1 Hypermethylation- Non NY

J107-4

Alternate Name: BRAF RFLX NEG TO MLH1

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 10-15 unstained slides at 5 microns (tumor content>10%) with H&E, Shavings (if tumor content >40%) with H&E, if sending DNA - Please call CS.

Storage Instruction: Room Temp

Turn Around Time: 5-10 days

CPT Code: 81210 (1)

Scarlet Eligible: No

BRAF, if neg reflex to MLH1 Hypermethylation NY Only

J463-1

Alternate Name: BRAF RFLX NEG TO MLH1

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 10-15 unstained slides at 5 microns (tumor content >10%) with H&E, Shavings (if tumor content >10%) with H&E, if sending DNA- Please call CS

Turn Around Time: 5-10 days

CPT Code: 81210 (1)

Scarlet Eligible: No

BRAF Sequencing- Non NY

J021-7

Alternate Name: BRAF Sequencing

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 10-15 unstained slides at 5 microns (tumor content>10%) with H&E, Shavings (if tumor content >40%) with H&E, if sending DNA - Please call CS.

Storage Instruction: Room Temp

Turn Around Time: 5-10 days

CPT Code: 81210 (1)

Scarlet Eligible: No

BRAF Sequencing NY Only

J456-5

Alternate Name: BRAF by NGS

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 10-15 unstained slides at 5 microns (tumor content >10%) with H&E, Shavings (if tumor content >10%) with H&E, if sending DNA- Please call CS

Turn Around Time: 5-10 days

CPT Code: 81210 (1)

Scarlet Eligible: No

OnkoSight BRAF Sequencing - BM/PB

T979-4

Alternate Name: BRAF by NGS - BM/PB

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: 2 mL Bone Marrow - Lavender Top

Collection Instruction: Place bone marrow into tube, label with patient's name

Storage Instruction: Room Temp

Turn Around Time: 5-10 days

CPT Code: 81210 (1)

Scarlet Eligible: No

BRAF/EGFR NGS if EGFR neg rfx to ALK FISH if neg rfx ROS1

TB34-4

Alternate Name: BRAF/EGFR rfx >ALK>ROS1

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 0-15 unstained slides at 5 microns (tumor content >10%) with H&E, Shavings (if tumor content >40%) with H&E, if sending DNA - Please call CS. Please give to Block Retrieval for processing. Alternative specimen: Extracted DNA is acceptable, providing that the isolation of nucleic acid occurs in a CLIA-certified laboratory or a laboratory meeting equivalent requirements as determined by CMS and/or the CAP.

Storage Instruction: Room Temp

Turn Around Time: 10 days

CPT Code: 81210 (1), 81235 (1)

Scarlet Eligible: No

BRCA1/2 Sequencing and Del/Dup Analysis

B362-5

Alternate Name: BRCA1/2 sequencing and del/dup

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: 2-5 mL Lavender top- EDTA

Collection Instruction: Lavender: 2 tubes (5-6ml each) Label with the patient name, date of birth and/or ID. Ideal handling is to ship immediately at ambient temperature for overnight delivery. Alternate sample type: DNA: Extract and submit DNA 50 ng/mL (total volume 400 mL) in sterile container. Ship DNA in sealed and protective packaging at ambient temperature, using a cool pack in hot weather. Labeling: Name of person or identifier, DNA concentration and volume MUST be provided on the specimen tube.

Storage Instruction: Room Temp

Turn Around Time: 14 days

CPT Code: 81162 (1)

Scarlet Eligible: Yes

Breast Carcinoma/Lymph Nodes/MM Metastasis

5170-6

Alternate Name: BREAST CA/LYMPH NODES/MM

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88342 (1)

Scarlet Eligible: No

Breast Carcinoma, R/O Microinvasion BY IHC-TC

A952-4

Alternate Name: BREAST MICROINVASION-TC

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88341 (2), 88342 (1)

Scarlet Eligible: No

BREAST SMEAR (slide)

1301-1

Alternate Name: BREAST SMEAR (slide)

Preferred Specimen: Slides

Collection Instruction: Label slides with patient name, place in slide holder for transport

Turn Around Time: 4 days

CPT Code: 88104 (1)

Scarlet Eligible: No

Blood Urea Nitrogen (BUN), Serum

0049-7

Alternate Name: BUN

Methodology: Colorimetry

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84520 (1)

Scarlet Eligible: Yes

BUN, Post Dialysis

1679-0

Alternate Name: BUN

Methodology: Urease

Preferred Specimen: 1 mL Microtainer - Pediatric SST

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84520 (1)

Scarlet Eligible: No

BUN + Creatinine with Ratio

4030-3

Alternate Name: BUN/Creat Ratio

Turn Around Time: 1 days

CPT Code: 82565 (1), 84520 (1)

Scarlet Eligible: Yes

Buprenorphine (Suboxone), Urine, Screen (5)

3956-0

Alternate Name: BUPRENORPHINE

Methodology: Enzyme Immunoassay

Preferred Specimen: 5 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80307 (1)

Scarlet Eligible: Yes

Bacterial Vaginosis and Aerobic Vaginitis (Non NY)

L571-9

Alternate Name: BV + AV

Turn Around Time: 7 days

CPT Code: 87512 (1), 87799 (1), 99999 (1), 87651 (1), 87653 (1), 87798 (1), 87640 (1)

Scarlet Eligible: No

BV PANEL BY RT-PCR

H870-9

Alternate Name: BV PANEL BY RT-PCR

Preferred Specimen: ThinPrep Vial

Turn Around Time: 7 days

CPT Code: 87512 (1), 87798 (2), 87799 (2)

Scarlet Eligible: No

BV Panel by RT-PCR

K835-9

Alternate Name: BV Panel by RT-PCR

Methodology: Real-Time Polymerase Chain Reaction

Preferred Specimen: Swab-Aptima Genprobe

Collection Instruction: Collect specimen-break swab shaft into transport tube, then label with name and source

Storage Instruction: Room Temp

Turn Around Time: 7 days

CPT Code: 87512 (1), 87799 (1), 87798 (1), 99999 (1)

Scarlet Eligible: No

Bacterial Vag. W/RFX TP-NY.

H575-4

Alternate Name: BV TP NY

Turn Around Time: 5 days

CPT Code: 87480 (1), 87510 (1), 87660 (1)

Scarlet Eligible: No

BV If Negative Reflex to Aerobic Vaginitis (Non NY)

L570-1

Alternate Name: BV rfx to Aerobic Vaginitis

Turn Around Time: 7 days

CPT Code: 87512 (1), 87798 (2), 87799 (2)

Scarlet Eligible: No

BVAB2 by RT-PCR

J215-5

Alternate Name: BVAB2 BY PCR

Methodology: Polymerase Chain Reaction

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 5-7 days

CPT Code: 87798 (1)

Scarlet Eligible: No

Allergen Bahia Grass (g17), IgE

1825-9

Alternate Name: Bahia Grass

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Bamboo Shoot (f51), IgE

3642-6

Alternate Name: Bamboo Shoot (F51)

Methodology: Immunocap

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Banana (f92), IgE

0863-1

Alternate Name: Banana

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Phenobarbital, Serum

0289-9

Alternate Name: Barbita®; Comizial®; Fenilcal®;Gardenal®;Luminal®; Phenemal; Phenemalum;Phenobarb;Phenobarbitone;PheNY patientslethylmalolurea;Solfoton®;Stental Extentabs®

Methodology: Enzyme Immunoassay

Preferred Specimen: Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80184 (1)

Scarlet Eligible: Yes

Bardet-Biedl Syndrome 2

J317-9

Alternate Name: Bardet-Biedl Syndrome BBS2-associated

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81406 (1)

Scarlet Eligible: Yes

Allergen Barley (f6), IgE

0816-9

Alternate Name: Barley

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Lipid Screen (Basic Lipid Profile)

0009-1

Alternate Name: Basic Lipid Panel; Lipid Screen

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80061 (1)

Scarlet Eligible: Yes

Allergen Basil (rf269), IgE

1461-3

Alternate Name: Basil (Rf269)

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Beech Tree, IgE (t5)

0948-0

Alternate Name: Beech Tree

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Beef (f27), IgE

0861-5

Alternate Name: Beef

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Beetroot (f319), IgE

3803-4

Alternate Name: Beetroot

Methodology: Immunocap

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Immunotyping, Urine

1644-4

Alternate Name: Bence Jones Protein, Immunofixation

Methodology: Electrophoresis

Preferred Specimen: 9 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86335 (1)

Scarlet Eligible: Yes

Amphetamines, Urine, Screen (1000)

0920-9

Alternate Name: Bennies; Crystal; Dexedrine®; Dexies; Ferndex®; Speed; Uppers

Methodology: Enzyme Immunoassay

Preferred Specimen: 5 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80307 (1)

Scarlet Eligible: Yes

Allergen Bermuda Grass (g2), IgE

1198-1

Alternate Name: Bermuda Grass

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Betalactoglobulin (f77), IgE

1268-2

Alternate Name: Beta Lactoglobulin

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86008 (1)

Scarlet Eligible: Yes

Pregnancy, Urine

0133-9

Alternate Name: Beta Subunit Human Chorionic Gonadotropin, Urine;hCG, Slide Test, Stat;hCG,Urine;Human Chorionic Gonadotropin, Urine

Methodology: Colorimetry

Preferred Specimen: 5 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 81025 (1)

Scarlet Eligible: Yes

Procainamide, Ser.(+napa #0525)

1846-5

Alternate Name: Bicoryl®;Novocainamidum;Novocamid®;Procaine Amide Hydrochloride;Procanbid®;Procan® SR;Pronestyl®;Pronestyl-SR®;Retard®;Rhythmin®

Methodology: Fluor. Polarization Immunoassay

Preferred Specimen: Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80192 (1)

Scarlet Eligible: Yes

Special Stain, Halls Bile - TC

T409-2

Alternate Name: Bile Halls SS- TC

Methodology: Histology

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint. 3-4 unstained positively charged slides(3 to 5 microns thick sections each)

Storage Instruction: Room Temp

Turn Around Time: 3-6 days

CPT Code: 88313 (1)

Scarlet Eligible: No

Bilirubin, Direct, Serum

0044-8

Alternate Name: Bilirubin, Conjugated; Direct Bilirubin

Methodology: Colorimetry

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Place 1-3 mL of serum in transport tube. Label as serum;refrigerate or freeze as required.

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82248 (1)

Scarlet Eligible: Yes

Bilirubin, Total + Direct, Serum

0045-5

Alternate Name: Bilirubin Fractionation

Methodology: Colorimetry

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Turn Around Time: 1 days

CPT Code: 82247 (1), 82248 (1)

Scarlet Eligible: Yes

Bilirubin, Total, Serum

0043-0

Alternate Name: Bilirubin, Total

Methodology: Colorimetry

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Place 1-3 mL of serum in transport tube. Label as serum;refrigerate or freeze as required.

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82247 (1)

Scarlet Eligible: Yes

Quinidine, Serum

0140-4

Alternate Name: Biquin®; Cardioquin®; Cin-Quin®; Kiditard®; Kiniden®; Quinaglute®; Dura-Tabs®; Quinalan®; Extentabs®; Quini®; Durules®; Quinora®; Systodin®

Methodology: Enzyme Immunoassay

Preferred Specimen: Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 80194 (1)

Scarlet Eligible: Yes

Allergen Black Pepper (f250), IgE

1422-5

Alternate Name: Black Pepper

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Blackberry (rf211), IgE

1486-0

Alternate Name: Blackberry

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

GP/Non-Automated Urovysion BY FISH

6274-5

Alternate Name: Bladder Cancer

Methodology: FISH

Preferred Specimen: 50 mL ThinPrep Cytolyt Solution

Collection Instruction: Collect voided urine in Thinprep Cytolyte solution

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88120 (1)

Scarlet Eligible: No

Allergen Blomia Tropicalis (d201), IgE

B235-3

Alternate Name: Blomia Tropicalis

Methodology: Immunocap

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Culture, Blood (spec.#2)

1350-8

Alternate Name: Blood Culture

Methodology: Bacterial Culture

Preferred Specimen: 10 mL Blood Culture Bottles-Aerobic/Anaerobic

Collection Instruction: Collect 5-10 mL of blood in each bottle, label with patient's name

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87040 (1)

Scarlet Eligible: No

Culture, Blood (spec.#3)

1351-6

Alternate Name: Blood Culture

Methodology: Bacterial Culture

Preferred Specimen: 10 mL Blood Culture Bottles-Aerobic/Anaerobic

Collection Instruction: Collect 5-10 mL of blood in each bottle, label with patient's name

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87040 (1)

Scarlet Eligible: No

Culture, Blood

0341-8

Alternate Name: Blood Culture

Methodology: Bacterial Culture

Preferred Specimen: Blood Culture Bottles-Aerobic/Anaerobic

Collection Instruction: Collect 5-10 mL of blood in each bottle, label with patient's name

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87040 (1)

Scarlet Eligible: No

Malaria/Babesia/Other Blood Parasites

J897-0

Alternate Name: Blood Parasites

Methodology: Microscopic Examination

Preferred Specimen: 2 mL Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Room Temp

Turn Around Time: 1 days

CPT Code: 87207 (1)

Scarlet Eligible: Yes

Occult Blood,Stool (1st Specimen)

0122-2

Alternate Name: Blood, Occult, Stool, Fecal Occult Blood Test; FOBT, Guaiac

Methodology: Colorimetry

Preferred Specimen: Stool-Hemoccult Card

Collection Instruction: Using applicator stick, place small amount of stool on card, close flap, label with patient name and date collected.

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82272 (1)

Scarlet Eligible: No

Occult Blood,Stool (3rd Spec.)

0323-6

Alternate Name: Blood, Occult, Stool; Fecal Occult Blood Test; FOBT; Guaiac

Methodology: Colorimetry

Preferred Specimen: Stool-Hemoccult Card

Collection Instruction: Using applicator stick, place small amount of stool on card, close flap, label with patient name and date collected.

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82272 (1)

Scarlet Eligible: No

Occult Blood,Stool (2nd Spec.)

0319-4

Alternate Name: Blood, Occult, Stool; Fecal Occult Blood Test; FOBT; Guaiac

Methodology: Colorimetry

Preferred Specimen: Stool-Hemoccult Card

Collection Instruction: Using applicator stick, place small amount of stool on card, close flap, label with patient name and date collected.

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82272 (1)

Scarlet Eligible: No

Allergen Blue Mussel (f37), IgE

0697-3

Alternate Name: Blue Mussel

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Blueberry (rf228), IgE

1487-8

Alternate Name: Blueberry

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Culture, Body Fluid W/Gram Stain

0556-1

Alternate Name: Body Fluid Culture with gram stain

Preferred Specimen: Fluid-source required

Collection Instruction: Submit fluid in sterile cup

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 87070 (1), 87205 (1), 87075 (1)

Scarlet Eligible: No

Ostase - Alkaline Phosphatase Bone Specific

T447-2

Alternate Name: Bone Specific Alkaline Phosphatase

Methodology: Chemiluminescence

Preferred Specimen: 2 SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 84080 (1)

Scarlet Eligible: Yes

Allergen Brazil Nut (f18), IgE

1405-0

Alternate Name: Brazil Nut

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Breast Basal Layer Triple Stain by IHC Global

TA18-9

Alternate Name: Breast Basal Layer Trip

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Storage Instruction: Room Temp

Turn Around Time: 1 days

CPT Code: 88344 (1)

Scarlet Eligible: No

ER / PR / HER2 (Tech Only)w/digital imaging

5127-6

Alternate Name: Breast Cancer IHC Tech Only: ER, PR, HER2

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88361 (3)

Scarlet Eligible: No

Breast Cancer Management Panel

J055-5

Alternate Name: Breast High/Moderate Risk Panel

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: 2-5 mL Lavender top- EDTA

Collection Instruction: 2 tubes (5-6mL each) Label with the patient name, date of birth, and/or ID.

Storage Instruction: Room Temp

Turn Around Time: 14 days

CPT Code: 81162 (1), 81321 (1), 81406 (1)

Scarlet Eligible: Yes

Allergen Broccoli (f260), IgE

1229-4

Alternate Name: Broccoli

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Brome Grass (G11), IgE

2394-5

Alternate Name: Brome Grass

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Breast Basal Layer Triple Stain by IHC Tech Only

TA17-1

Alternate Name: Brst Basal Layer x3-TC

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Storage Instruction: Room Temp

Turn Around Time: 1 days

CPT Code: 88344 (1)

Scarlet Eligible: No

Allergen Brussel Sprouts (rf217), IgE

1814-3

Alternate Name: Brussel Sprouts

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Buckwheat (f11), IgE

1629-5

Alternate Name: Buckwheat

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Clostridium difficile Toxin/GDH antigen with Reflex to PCR

0259-2

Alternate Name: C difficile

Methodology: Enzyme Immunoassay

Preferred Specimen: Stool Cont.-Cup

Collection Instruction: Using applicator stick, place 20-30 gms of stool into cup, label with patient's name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 87324 (1), 87449 (1)

Scarlet Eligible: No

C-Peptide, Serum

0578-5

Alternate Name: C-Peptide, Serum

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 84681 (1)

Scarlet Eligible: Yes

C-Reactive Protein (hs-CRP), High Sensitivity

3320-9

Alternate Name: C-RP, High Sensitivity, Cardiac CRP

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86141 (1)

Scarlet Eligible: Yes

C-Reactive Protein (CRP), Serum

0069-5

Alternate Name: C-RP, Serum

Methodology: Immunoturbidimetric

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86140 (1)

Scarlet Eligible: Yes

Candida glabrata by PCR

T783-0

Alternate Name: C. Glabrata by PCR

Methodology: Multiplex PCR

Preferred Specimen: Swab Aptima Orange label

Collection Instruction: Collect specimen-break swab shaft into transport tube, then label with name and source

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87481 (1)

Scarlet Eligible: No

Candida dubliniensis by PCR

T786-3

Alternate Name: C. dubliniensis by PCR

Methodology: Multiplex PCR

Preferred Specimen: Swab Aptima Orange label

Collection Instruction: Collect specimen-break swab shaft into transport tube, then label with name and source

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87481 (1)

Scarlet Eligible: No

Candida krusei

T449-8

Alternate Name: C. krusei

Methodology: Multiplex PCR

Preferred Specimen: Swab Aptima Orange label

Collection Instruction: Collect specimen-break swab shaft into transport tube, then label with name and source

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87481 (1)

Scarlet Eligible: No

Candida parapsilosis

T536-2

Alternate Name: C. parapsilosis

Methodology: Multiplex PCR

Preferred Specimen: Swab Aptima Orange label

Collection Instruction: Collect specimen-break swab shaft into transport tube, then label with name and source

Storage Instruction: Room Temp

Turn Around Time: 21 days

CPT Code: 87481 (1)

Scarlet Eligible: No

Candida tropicalis by PCR

T745-9

Alternate Name: C. tropicalis

Methodology: Multiplex PCR

Preferred Specimen: Swab Aptima Orange label

Collection Instruction: Collect specimen-break swab shaft into transport tube, then label with name and source

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87481 (1)

Scarlet Eligible: No

Chlamydia Trachomatis, Urine, rRNA

6368-5

Alternate Name: C.TRACHOMATIS, UR, rRNA

Methodology: Transcription Mediated Amplification

Preferred Specimen: Urine Container-Genprobe-Aptima

Collection Instruction: First morning sample, collect in sterile cup, transfer 2 ml into aptima holder. Correct level is between the black fill lines(window) on the holder.

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87491 (1)

Scarlet Eligible: Yes

Complement C3, Serum

0532-2

Alternate Name: C3, Serum

Methodology: Immunoturbidimetric

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86160 (1)

Scarlet Eligible: Yes

Complement C4, Serum

0533-0

Alternate Name: C4, Serum

Methodology: Immunoturbidimetric

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86160 (1)

Scarlet Eligible: Yes

CA 15-3, Serum

2130-3

Alternate Name: CA 15-3

Methodology: Chemiluminescence

Preferred Specimen: 1 mL Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86300 (1)

Scarlet Eligible: Yes

CA 19-9, Serum

0535-5

Alternate Name: CA 19-9

Methodology: Chemiluminescence

Preferred Specimen: 2 mL Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86301 (1)

Scarlet Eligible: Yes

CA 27.29, Serum

0823-5

Alternate Name: CA 27.29**Methodology:** Chemiluminescence**Preferred Specimen:** 1 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 1 days**CPT Code:** 86300 (1)**Scarlet Eligible:** Yes**CA 125, Serum**

0536-3

Alternate Name: CA-125**Methodology:** Electrochemiluminescence Immunoassay**Preferred Specimen:** 1 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 1 days**CPT Code:** 86304 (1)**Scarlet Eligible:** Yes**Cadmium, Blood**

J018-3

Alternate Name: CADMIUM, BLOOD**Methodology:** Inductively Coupled Plasma Mass Spectrometry**Preferred Specimen:** Blue-Royal Blue Top K2 EDTA Plasma(*LAVENDER LINE)**Collection Instruction:** Fill tube, invert gently 5-6 times, label with patient name**Storage Instruction:** Refrigerate**Turn Around Time:** 1 days**CPT Code:** 82300 (1)**Scarlet Eligible:** Yes

Vitamin D, 25-Hydroxy, Serum

0286-5

Alternate Name: CALCIFEDIOL

Methodology: Chemiluminescence

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 82306 (1)

Scarlet Eligible: Yes

CALL RESULT TO DOCTOR

9990-3

Alternate Name: CALL RESULT TO DOCTOR

CPT Code: 99999 (1)

Scarlet Eligible: Yes

Canavan Disease Test

A657-9

Alternate Name: CANAVAN DISEASE

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81200 (1)

Scarlet Eligible: Yes

Vulvovaginal Candidiasis Panel by PCR w/ reflex to species

B976-2

Alternate Name: CANDIDA GENUS PLUS ALBICANS

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 5-7 days

CPT Code: 87481 (1)

Scarlet Eligible: No

Breast Carcinoma, Rule Out Microinvasion, with Interp.

5176-3

Alternate Name: CARCINOMA PANEL

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block or slide form with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88341 (3), 88342 (1)

Scarlet Eligible: No

Cathepsin D by IHC

5151-6

Alternate Name: CATHEPSIN D by IHC

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 82387 (1)

Scarlet Eligible: No

CCP (CYCLIC CITRULLINATED PEPTIDE) IGG, SERUM

3747-3

Alternate Name: CCP, cyclic citrullinated peptide antibody

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86200 (1)

Scarlet Eligible: Yes

T-NK CD57 Flow Cytometry Analysis

5221-7

Alternate Name: CD57

Methodology: Flow Cytometry

Preferred Specimen: Peripheral Blood - Green Top

Collection Instruction: Place bone marrow into tube, label with patient's name. Send in bone marrow kit.

Storage Instruction: Room Temp

Turn Around Time: 1 days

CPT Code: 88184 (1), 88185 (2), 88187 (1)

Scarlet Eligible: Yes

Clostridium Difficile Toxins A+B Genes By DNA Amplification

B137-1

Alternate Name: CDIFF TOXINS A+B GENES

Methodology: LAMP DNA Amplification

Preferred Specimen: Stool Cont.-Cup

Collection Instruction: Using applicator stick, place 10 gms of stool in container, label with patient's name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 87493 (1)

Scarlet Eligible: No

Carcinoembryonic Antigen (CEA)

0055-4

Alternate Name: CEA

Methodology: Electrochemiluminescence Immunoassay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82378 (1)

Scarlet Eligible: Yes

CEA, Peritoneal Fluid

6169-7

Alternate Name: CEA PERITONEAL,FLD

Methodology: Immunoassay

Preferred Specimen: 3 mL Fluid-source required

Collection Instruction: Submit fluid in sterile cup

Storage Instruction: Refrigerate

Turn Around Time: 11 days

CPT Code: 82378 (1)

Scarlet Eligible: No

AML - CEBPA Mutational Analysis

8860-9

Alternate Name: CEBPA MUT. ANALYSIS

Methodology: Polymerase Chain Reaction

Preferred Specimen: 2 mL Bone Marrow - Lavender Top

Collection Instruction: Place bone marrow into tube, label with patient's name. Send in bone marrow kit.

Storage Instruction: Refrigerate

Turn Around Time: 6 days

CPT Code: 81218 (1)

Scarlet Eligible: No

Celiac Disease Antibody Evaluation1186-6

Alternate Name: CELIAC DISEASE ABS.**Preferred Specimen:** 2 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 3 days**CPT Code:** 83516 (2), 86255 (1)**Scarlet Eligible:** Yes

Celiac Disease Eval + HLA GeneticsK002-6

Alternate Name: CELIAC EVAL+HLA**Turn Around Time:** 12 days**CPT Code:** 81382 (2), 83516 (4), 86255 (1)**Scarlet Eligible:** No

Cell Count, Fluid0056-2

Alternate Name: CELL COUNT, FLUID**Preferred Specimen:** 2 mL Fluid-source required**Collection Instruction:** Submit fluid in sterile cup**Storage Instruction:** Refrigerate**CPT Code:** 89051 (1)**Scarlet Eligible:** No

Centromere Autoantibodies

0587-6

Alternate Name: CENTROMERE AUTOABS.

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86038 (1)

Scarlet Eligible: Yes

Cervicitis W/O PAP W/O HPV PLUS

Q265-1

Alternate Name: CERV W/O PAP W/O HPV+

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Turn Around Time: 5-7 days

CPT Code: 87491 (1), 87529 (2), 87591 (1), 87661 (1)

Scarlet Eligible: No

CF/Fragile X/SMA

2115-4

Alternate Name: CF/FRAGILE X/SMA

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10 days

CPT Code: 81220 (1), 81243 (1), 81400 (1)

Scarlet Eligible: Yes

Cystic Fibrosis (CFTR) Intron 9 Poly-T Analysis (5T/7T/9T)

J520-8

Alternate Name: CFTR Poly-T (5T/7T/9T)

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 7-10 days

CPT Code: 81224 (1)

Scarlet Eligible: Yes

Chem 12 Profile

0004-2

Alternate Name: CHEM 12 PROFILE

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80012 (1)

Scarlet Eligible: Yes

Chem 7 Panel

0768-2

Alternate Name: CHEM 7 PANEL

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80007 (1)

Scarlet Eligible: Yes

Chlamydia Trachomatis by Multiplex PCR

6337-0

Alternate Name: CHLAMYD.TRACH BY MULTIP

Methodology: Polymerase Chain Reaction

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 87491 (1)

Scarlet Eligible: No

Chloramphenicol

3792-9

Alternate Name: CHLORAMPHENICOL

Methodology: High Pressure Liquid Chromatography

Preferred Specimen: 1 mL Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 4 days

CPT Code: 82415 (1)

Scarlet Eligible: Yes

Cholinesterase, Serum

0545-4

Alternate Name: CHOLINESTERASE,S.

Methodology: Colorimetry

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82480 (1)

Scarlet Eligible: Yes

Chromosome Analysis- Breakage, Fanconi Anemia

B534-9

Alternate Name: CHROM. ANAL./FANCONI

Methodology: Cell Culture

Preferred Specimen: 4 mL Green top- Sodium Heparin

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient name

Storage Instruction: Room Temp, Refrigerate

Turn Around Time: 19 days

CPT Code: 88230 (1), 88249 (1)

Scarlet Eligible: No

Chromosome Analysis w/Reflex Microarray

5394-2

Alternate Name: CHROM.W/RFX MICROARRAY

Preferred Specimen: Microtainer - Pediatric Green Top,Lavender top- EDTA

Collection Instruction: Collect Green top and lavender- Room Temp.

Turn Around Time: 21 days

CPT Code: 81228 (1), 88230 (1), 88262 (1), 88291 (1)

Scarlet Eligible: Yes

CPK Isoenzymes

0558-7

Alternate Name: CK Isoenzymes; Ck-MB and Total CK; CPK Isoenzymes

Methodology: Electrophoresis

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 82550 (1), 82552 (1)

Scarlet Eligible: Yes

CK-903 for Prostate Cancer

5434-6

Alternate Name: CK-903, IHC

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88342 (1)

Scarlet Eligible: No

CPK (CK), Total

0068-7

Alternate Name: CK; CPK; Creatine Phosphokinase, Total, Serum

Methodology: Enzymatic Reference method with Hexokinase

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82550 (1)

Scarlet Eligible: Yes

Creatine Kinase-MB Fraction(CKMB)

2922-3

Alternate Name: CKMB

Methodology: Electrochemiluminescence Immunoassay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82553 (1)

Scarlet Eligible: Yes

1Q21 BY FISH-TC ONLY

A876-5

Alternate Name: CKS1B by FISH, TC Only

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube. Also acceptable: Whole blood in green (sodium heparin) tube.

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88374 (1)

Scarlet Eligible: No

OnkoSight Advanced Chronic Lymphoid Neoplasm NGS Panel

TH55-7

Alternate Name: CLL Cancer NGS Panel

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Bone Marrow - Lavender Top

Storage Instruction: Room Temp

Turn Around Time: 10 days

CPT Code: 81450 (1)

Scarlet Eligible: No

CLL FISH Panel - TC Only

P962-4

Alternate Name: CLL FISH PANEL- TC ONLY

Methodology: FISH

Preferred Specimen: 5 mL Peripheral Blood - Green Top

Collection Instruction: Fill tube, invert 8-10 times (DO NOT SHAKE!)

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88367 (1), 88373 (3)

Scarlet Eligible: No

Minimal Residual Disease for CLL/SLL

5155-7

Alternate Name: CLL MRD Panel

Methodology: Flow Cytometry

Preferred Specimen: 5 mL Peripheral Blood - Green Top

Collection Instruction: Fill tube, invert 8-10 times (DO NOT SHAKE!)

Storage Instruction: Room Temp

Turn Around Time: 24 hours

CPT Code: 88184 (1), 88185 (8), 88188 (1)

Scarlet Eligible: Yes

CMV Ab., IGG *POS. Reflex To IGM

2377-0

Alternate Name: CMV Ab. IgG(reflex)

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86644 (1)

Scarlet Eligible: Yes

Cytomegalovirus (CMV) by Multiplex PCR (NON-NY)

6353-7

Alternate Name: CMV BY MULTIPLEX PCR

Methodology: Polymerase Chain Reaction

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Turn Around Time: 5 days

CPT Code: 87496 (1)

Scarlet Eligible: No

Cytomegalovirus (CMV), DNA, Ultraquant

1161-9

Alternate Name: CMV,DNA PCR

Methodology: Polymerase Chain Reaction

Preferred Specimen: 2 mL Aliquot Plasma-EDTA-Lavender Top

Collection Instruction: Centrifuge lavender-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA and patient's name.

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87497 (1)

Scarlet Eligible: Yes

Carbon Dioxide, Serum (CO2)

0052-1

Alternate Name: CO2

Methodology: Automated Chemistry

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82374 (1)

Scarlet Eligible: Yes

COA/BMP

K935-7

Alternate Name: COA/BMP

Turn Around Time: 1 days

CPT Code: 80048 (1)

Scarlet Eligible: No

AtheroSclerotic Risk (Coag)

5874-3

Alternate Name: COAG/ATHEROSCLEROTIC RISK

Preferred Specimen: Red Top,Light blue top

Storage Instruction: Refrigerate

Turn Around Time: 4 days

CPT Code: 80061 (1), 83090 (1), 83695 (1), 85240 (1), 85384 (1), 85610 (1), 85670 (1), 85730 (1), 86141 (1)

Scarlet Eligible: No

Coagulation Profile

0994-4

Alternate Name: COAGULATION

Preferred Specimen: SST Tube,Lavender top- EDTA,Light blue top

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 83718 (1), 85025 (1), 85384 (1), 85610 (1), 85730 (1), G0059 (1)

Scarlet Eligible: Yes

Cocaine, Urine, Screen (100)

1478-7

Alternate Name: COCAINE(100), QUAL.

Methodology: Enzyme Immunoassay

Preferred Specimen: 5 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80307 (1)

Scarlet Eligible: Yes

Cocaine, Urine, Screen (300)

0923-3

Alternate Name: COCAINE, URINE SCREE**Methodology:** Enzyme Immunoassay**Preferred Specimen:** 25 mL Urine Cup**Collection Instruction:** Collect urine in cup, label with patients name**Storage Instruction:** Refrigerate**Turn Around Time:** 1 days**CPT Code:** 80307 (1)**Scarlet Eligible:** Yes

Collagen Binding Assay

5744-8

Alternate Name: COLLAGEN BINDING**Methodology:** Enzyme Linked Immunoabsorbance**Preferred Specimen:** 1 mL Aliquot Plasma-Citrated-Blue Top**Collection Instruction:** Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.**Storage Instruction:** STRICT Frozen**Turn Around Time:** 7 days**CPT Code:** 85246 (1)**Scarlet Eligible:** No

Colon Biopsy (2 Specimens)

5652-3

Alternate Name: COLON BIOPSY X2**Methodology:** Histology**Preferred Specimen:** Biopsy In Formalin**Collection Instruction:** Submit in 10% Neutral Buffered Formalin, varied jar sizes.**Turn Around Time:** 1 days**CPT Code:** 88305 (2)**Scarlet Eligible:** No

Comprehensive Heart Health Baseline

F306-6

Alternate Name: COMP HEARTHLTH BAS(ACT)

Preferred Specimen: 2 mL SST Tube,Lavender top- EDTA

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Turn Around Time: 2 days

CPT Code: 80053 (1), 82465 (1), 83036 (1), 83695 (1), 83704 (1), 83718 (1), 83721 (1), 84478 (1), 85025 (1), 86141 (1)

Scarlet Eligible: Yes

Comprehensive Heart Health Monitoring

F307-4

Alternate Name: COMP HRTHLTH MNITR(ACT)

Turn Around Time: 2 days

CPT Code: 80053 (1), 82465 (1), 83036 (1), 83704 (1), 83718 (1), 83721 (1), 84478 (1), 85025 (1), 86141 (1)

Scarlet Eligible: Yes

Comprehensive Metabolic Panel

3427-2

Alternate Name: COMP. METABOLIC (CMP)

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80053 (1)

Scarlet Eligible: Yes

Comp. Metabolic/GGTP

9150-4

Alternate Name: COMP. METABOLIC/GGTP

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80053 (1), 82977 (1)

Scarlet Eligible: Yes

Comprehensive Slide Consult

5256-3

Alternate Name: COMPREHENSIVE SLIDE CONSULT

Methodology: Morphology

Preferred Specimen: Stained Slides

Collection Instruction: Client to collect specimen, prepare smears, label with patient's name, place in slide carrier.

Turn Around Time: 1-2 days

CPT Code: 88325 (1)

Scarlet Eligible: No

Cone Biopsy/LEEP

5672-1

Alternate Name: CONE BIOPSY/LEEP

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 2 days

CPT Code: 88307 (1)

Scarlet Eligible: No

Congo Red Amyloid Stain

5011-2

Alternate Name: CONGO RED AMYLOID STAIN

Methodology: Histology

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Storage Instruction: Room Temp

Turn Around Time: 2 days

CPT Code: 88313 (1)

Scarlet Eligible: No

Coronary Risk Profile II

0993-6

Alternate Name: CORONARY RISK PROFIL

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80061 (1), 82172 (2), 83700 (1)

Scarlet Eligible: Yes

Cortisol, AM, baseline

0066-1

Alternate Name: CORTISOL, AM (baseline)

Methodology: Chemiluminescence

Preferred Specimen: 1 mL Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 82533 (1)

Scarlet Eligible: No

Novel Coronavirus COVID-19 Nasal/Nasopharynx

TH68-0

Alternate Name: COVID-19 Nasal/Nasopharynx**Methodology:** Real-Time Polymerase Chain Reaction**Preferred Specimen:** Swab-Viral Culturette**Storage Instruction:** Refrigerate**Turn Around Time:** 3 days**CPT Code:** U0003 (1)**Scarlet Eligible:** Yes

Novel Coronavirus COVID-19 Oropharynx

TH69-8

Alternate Name: COVID-19 Oropharynx**Methodology:** Real-Time Polymerase Chain Reaction**Preferred Specimen:** Swab-Viral Culturette**Storage Instruction:** Refrigerate**Turn Around Time:** 3 days**CPT Code:** U0003 (1)**Scarlet Eligible:** No

Novel Coronavirus COVID-19 Pooled N/NP/OP

TH71-4

Alternate Name: COVID-19 Pooled N/NP/OP**Methodology:** Real-Time Polymerase Chain Reaction**Preferred Specimen:** Swab-Viral Culturette**Storage Instruction:** Refrigerate**Turn Around Time:** 3 days**CPT Code:** U0003 (1)**Scarlet Eligible:** No

Carnitine Palmitoyltransferase Def. Type 2 Carrier Screening

J315-3

Alternate Name: CPT II carrier screen**Methodology:** Genotyping by Next Generation Sequencing**Preferred Specimen:** Lavender top- EDTA**Collection Instruction:** Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!**Storage Instruction:** Refrigerate**Turn Around Time:** 10-14 days**CPT Code:** 81404 (1)**Scarlet Eligible:** Yes

OnkoSight Advanced Colorectal Cancer NGS Panel

TH50-8

Alternate Name: CRC NGS Panel**Methodology:** Genotyping by Next Generation Sequencing**Preferred Specimen:** Formalin-fixed, Paraffin-embedded Tissue**Storage Instruction:** Room Temp**Turn Around Time:** 10 days**CPT Code:** 81445 (1)**Scarlet Eligible:** No

Post Dialysis Creatinine

1687-3

Alternate Name: CREATININE POST DIALYS.**Methodology:** Colorimetry**Preferred Specimen:** Microtainer - Pediatric SST**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 1 days**CPT Code:** 82565 (1)**Scarlet Eligible:** No

PDF Creatinine (Timed Test code)

2063-6

Alternate Name: CREATININE, 24hr.DIALYSATE FLD

Preferred Specimen: 10.0 mL Fluid-source required

Collection Instruction: Submit fluid in sterile cup

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82570 (1)

Scarlet Eligible: No

Creatinine, Arterial

4534-4

Alternate Name: CREATININE, ARTERIAL

Methodology: Colorimetry

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82565 (1)

Scarlet Eligible: No

Creatinine Clearance (24Hr. Urine)

0072-9

Alternate Name: CREATININE, CLEAR,(24

Preferred Specimen: 2 mL SST Tube,Urine Container - 24hr

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82575 (1)

Scarlet Eligible: No

Creatinine, Systemic

4535-1

Alternate Name: CREATININE, SYSTEMIC**Methodology:** Colorimetry**Preferred Specimen:** 1 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 1 days**CPT Code:** 82540 (1)**Scarlet Eligible:** Yes

Creatinine, Urine, 24 HR.

0071-1

Alternate Name: CREATININE, URINE 24 HR**Methodology:** Colorimetry**Preferred Specimen:** 10 mL Urine Container - 24hr**Collection Instruction:** Discard first morning sample, then collect all urines after for 24 hrs, on completion label with collection time and patient's name; record total volume on Requisition**Storage Instruction:** Refrigerate**Turn Around Time:** 1 days**CPT Code:** 81050 (1), 82570 (1)**Scarlet Eligible:** No

Cerebrospinal Fluid (CSF) Cytology

5112-8

Alternate Name: CSF Cytology**Methodology:** Immunohistochemistry**Preferred Specimen:** Fluid-source required**Collection Instruction:** Submit fluid in sterile cup**Storage Instruction:** Refrigerate**Turn Around Time:** 4 days**CPT Code:** 88108 (1)**Scarlet Eligible:** No

Chlamydia (CT) PCR UrineT407-6

Alternate Name: CT Urine PCR**Methodology:** Polymerase Chain Reaction**Preferred Specimen:** 10 mL Urine Cobas Vial**Collection Instruction:** Patients should not urinate 1 hour prior to specimen collection.

Approximately 10-15 mL of urine should be collected from the beginning of the stream. Once collected the urine should be transferred to the cobas PCR media tube within 24 hrs. Use transfer pipette to add urine to the cobas PCR tube. Final volume of the tube should be between fill lines. Replace and tighten the cap and invert the tube 5 times to mix.

Storage Instruction: Refrigerate**Turn Around Time:** 3 days**CPT Code:** 87491 (1)**Scarlet Eligible:** Yes

Chlamydia, Anal AptimaA114-1

Alternate Name: CT, ANAL APTIMA**Methodology:** Transcription Mediated Amplification**Preferred Specimen:** Swab-Aptima Genprobe**Collection Instruction:** Collect specimen-break swab shaft into transport tube, then label with name and source**Storage Instruction:** Refrigerate**Turn Around Time:** 3 days**CPT Code:** 87491 (1)**Scarlet Eligible:** No

Chlamydia, Oral Aptima

A113-3

Alternate Name: CT, ORAL APTIMA

Methodology: Transcription Mediated Amplification

Preferred Specimen: Swab-Aptima Genprobe

Collection Instruction: Collect specimen-break swab shaft into transport tube, then label with name and source

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87491 (1)

Scarlet Eligible: No

Chlamydia, swab, TMA, Aptima

1004-1

Alternate Name: CT, rRNA,TMA SWAB

Methodology: Transcription Mediated Amplification

Preferred Specimen: Swab-Aptima Genprobe

Collection Instruction: Collect specimen-break swab shaft into transport tube, then label with name and source

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87491 (1)

Scarlet Eligible: No

High Risk STI CT/GC

3851-3

Alternate Name: CT/GC

Methodology: Nucleic Acid Amplification Test

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin Prep or SurePath vial, label with patients name.

Turn Around Time: 2 days

CPT Code: 87491 (1), 87591 (1)

Scarlet Eligible: No

Chlamydia/N. Gonorrhoea CT/GC PCR UrineL344-1

Alternate Name: CT/GC PCR Urine**Methodology:** Real-Time Polymerase Chain Reaction**Preferred Specimen:** 10 Urine Cobas Vial**Storage Instruction:** Refrigerate**Turn Around Time:** 3 days**CPT Code:** 87491 (1), 87591 (1)**Scarlet Eligible:** Yes

Chlamydia + Gonorrhoea, Urine, Aptima0257-6

Alternate Name: CT/GC rRNA,APTIMA,URINE**Methodology:** Transcription Mediated Amplification**Preferred Specimen:** Urine Container-Genprobe-Aptima**Storage Instruction:** Refrigerate**Turn Around Time:** 3 days**CPT Code:** 87491 (1), 87591 (1)**Scarlet Eligible:** Yes

Chlamydia + Gonorrhoeae, Anal, Swab, AptimaP011-0

Alternate Name: CT/GC, ANAL, APTIMA, Chlamydia/gonococcus**Methodology:** Transcription Mediated Amplification**Preferred Specimen:** Swab-Aptima Genprobe**Collection Instruction:** Collect specimen-break swab shaft into transport tube, then label with name and source**Storage Instruction:** Refrigerate**Turn Around Time:** 3 days**CPT Code:** 87491 (1), 87591 (1)**Scarlet Eligible:** No

Chlamydia, Gonorrhea, and Trichomonas rRNA,TMA Aptima Urine

H001-1

Alternate Name: CT/GC/TV,rRNA,TMA Urine

Turn Around Time: 3 days

CPT Code: 87491 (1), 87591 (1), 87661 (1)

Scarlet Eligible: Yes

High Risk STI CT/GC/Trich/M.Gen.

P006-0

Alternate Name: CT/GC/Trich/MGenHR STI

Preferred Specimen: ThinPrep Vial

Turn Around Time: 5 days

CPT Code: 87491 (1), 87591 (1), 87661 (1), 87798 (1)

Scarlet Eligible: No

Culture, Catheter Tip

J010-0

Alternate Name: CULTURE, CATHETER TIP

Methodology: Bacterial Culture

Preferred Specimen: Cup Plain

Collection Instruction: Place sample in sterile cup label with patient's name and source

Storage Instruction: Room Temp

Turn Around Time: 1-4 days

CPT Code: 87070 (1)

Scarlet Eligible: No

Culture, Mouth

0960-5

Alternate Name: CULTURE, MOUTH

Methodology: Bacterial Culture

Preferred Specimen: Swab-E

Collection Instruction: Collect specimen with swab then place swab into transport carrier; Label with name and source

Turn Around Time: 3 days

CPT Code: 87070 (1)

Scarlet Eligible: No

Culture, Stool (2nd Specimen)

0226-1

Alternate Name: CULTURE, STOOL (2ND)

Methodology: Bacterial Culture

Preferred Specimen: Fecal Swab

Collection Instruction: Using applicator stick, place 10 gms of stool in container, label with patient's name

Storage Instruction: Refrigerate

Turn Around Time: 4 days

CPT Code: 87045 (1), 87046 (2)

Scarlet Eligible: No

Culture Stool (3rd Specimen)

0225-3

Alternate Name: CULTURE, STOOL (3RD SPECIMEN)

Methodology: Bacterial Culture

Preferred Specimen: Fecal Swab

Collection Instruction: Using applicator stick, place 10 gms of stool in container, label with patient's name

Storage Instruction: Refrigerate

Turn Around Time: 4 days

CPT Code: 87045 (1), 87046 (2)

Scarlet Eligible: No

Culture, ID Only (Any Source)

1043-9

Alternate Name: CULTURE, ID ONLY (ANY SOURCE)

Methodology: Bacterial Culture

Preferred Specimen: Cup Plain

Collection Instruction: Place sample in sterile cup label with patient's name and source

Turn Around Time: 4 days

CPT Code: 87086 (1)

Scarlet Eligible: No

Cyclospora ID, Modified Acid Fast

1468-8

Alternate Name: CYCLOSPORA ID

Methodology: Modified Acid-Fast Stain

Preferred Specimen: Ova & Parasite Mini

Collection Instruction: With applicator stick, place 10 gms of stool into each vial, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 87207 (1)

Scarlet Eligible: No

PAP Smear, 2 Slides, Autopap

1102-3

Alternate Name: CYTOLOGY, FEMALE, GENITAL

Preferred Specimen: Slides

Turn Around Time: 5 days

CPT Code: 88148 (1)

Scarlet Eligible: No

Cytology, Respiratory System

1402-7

Alternate Name: CYTOLOGY, RESPIRATOR

Preferred Specimen: Slides

Collection Instruction: Label slides with patient name, place in slide holder for transport

Turn Around Time: 4 days

CPT Code: 88161 (1)

Scarlet Eligible: No

Cytology, Thyroid By FNA

4300-0

Alternate Name: CYTOLOGY, THYROID

Preferred Specimen: Slides

Collection Instruction: Label slides with patient name, place in slide holder for transport.

Turn Around Time: 4 days

CPT Code: 88104 (1)

Scarlet Eligible: No

Cytology, Urine With Interpretation

2100-6

Alternate Name: CYTOLOGY, URINE W/INTERP

Preferred Specimen: 50% Alcohol Jar

Collection Instruction: Body fluid for cytology in 50% alcohol (ex...urine, sputum, nipple discharge). Add directly to an equal volume of 50% Alcohol Solution.

Storage Instruction: Room Temp

Turn Around Time: 4 days

CPT Code: 88108 (1)

Scarlet Eligible: No

Cytopatho. Urine-PC (Alcohol)

5314-0

Alternate Name: CYTOPATHO. URN-PC (ALCOHOL)

Methodology: Cytogenetics

Preferred Specimen: Slides

Collection Instruction: Label slides with patient name, place in slide holder for transport

Turn Around Time: 1 days

CPT Code: 88108 (1)

Scarlet Eligible: No

Cytopathology Urine-TC (Thinprep)

5218-3

Alternate Name: CYTOPATHOLOGY URINE-TC (THINPREP)

Methodology: Cytogenetics

Preferred Specimen: 50 mL ThinPrep Cytolyt Solution

Collection Instruction: Follow kit instructions

Storage Instruction: Refrigerate

Turn Around Time: 4 days

CPT Code: 88112 (1)

Scarlet Eligible: No

Anal Cytopathology, ThinPrep

6626-6

Alternate Name: CYTOPATHOLOGY, ANAL, THINPREP

Preferred Specimen: Anal ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Refrigerate

Turn Around Time: 4 days

CPT Code: 88112 (1)

Scarlet Eligible: No

Calcium, Ionized, Serum

0051-3

Alternate Name: Ca, Ionized, Serum

Methodology: Ion Selective Electrode

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82330 (1)

Scarlet Eligible: Yes

Calcium, Serum

0050-5

Alternate Name: Ca, Serum

Methodology: Colorimetry

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82310 (1)

Scarlet Eligible: Yes

Calcium, Urine, 24 Hr.

0359-0

Alternate Name: Ca, Urine, 24 Hours

Methodology: Colorimetry

Preferred Specimen: 10 mL Urine Container - 24hr

Collection Instruction: Discard first morning sample, then collect all urines after for 24 hrs, on completion label with collection time and patient's name; record total volume on requisition. 24 hr urine should be collected in containers with 20-30ml of 6N HCL.

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 81050 (1), 82340 (1)

Scarlet Eligible: No

Calcium, Urine, Timed Specimen

1660-0

Alternate Name: Ca, Urine, Time Specimen

Methodology: Colorimetry

Preferred Specimen: 5 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 82340 (1)

Scarlet Eligible: No

Vitamin D, 1, 25-Dihydroxy, Serum

0288-1

Alternate Name: Calcitriol

Methodology: Chemiluminescence

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 82652 (1)

Scarlet Eligible: Yes

Chromosome Analysis (Oncology)- Tissue

J913-5

Alternate Name: Cancer Cytogenetics Karyotype

Methodology: Chromosome Analysis

Preferred Specimen: Tissue in RPMI

Collection Instruction: Label specimen with 2 patient identification. Please include specimen collection date on requisition.

Storage Instruction: Refrigerate

Turn Around Time: 5-7 days

CPT Code: 88262 (1), 88239 (1)

Scarlet Eligible: No

Allergen C.Albicans (m5), IgE

0892-0

Alternate Name: Candida

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Candida Species by M. PCR - Swab

L416-7

Alternate Name: Candida Species by M. PCR - Swab

Collection Instruction: Collect specimen-break swab shaft into transport tube, then label with name and source

Turn Around Time: 21 days

CPT Code: 87481 (1)

Scarlet Eligible: No

Allergen Cardamon (rf267), IgE

2308-5

Alternate Name: Cardamon

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Homocysteine Level

1311-0

Alternate Name: Cardiac or Ultra sensitive Homocysteine

Methodology: Enzyme Immunoassay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 83090 (1)

Scarlet Eligible: Yes

Allergen Carrots (f31), IgE

0934-0

Alternate Name: Carrots

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Cat Dander (e1), IgE

0665-0

Alternate Name: Cat Dander

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Urine Culture, Catheterized

2603-9

Alternate Name: Catheter Urine Culture

Methodology: Bacterial Culture

Preferred Specimen: 5-10 mL Urine Tube - Grey Top

Collection Instruction: Perform a clean catch process, collect urine in a sterile container, transfer 10-15 ml into a Boricult, label with patient's name

Storage Instruction: Refrigerate

Turn Around Time: 4 days

CPT Code: 87086 (1)

Scarlet Eligible: No

Allergen Cauliflower (rf291), IgE

1815-0

Alternate Name: Cauliflower

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Celery (f85), IgE

1380-5

Alternate Name: Celery

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Cereal Mix (fx3), IgE

2185-7

Alternate Name: Cereal Mix

Methodology: Immunocap

Preferred Specimen: 5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86005 (1)

Scarlet Eligible: Yes

Culture, Cerebrospinal Fluid

0966-2

Alternate Name: Cerebrospinal Fluid Culture, CSF

Preferred Specimen: 1 mL Fluid-source required

Collection Instruction: Submit fluid in sterile cup

Turn Around Time: 1 days

CPT Code: 87070 (1), 87205 (1)

Scarlet Eligible: No

Protein Electrophoresis, CSF

3946-1

Alternate Name: Cerebrospinal Fluid Protein Electrophoresis

Methodology: Spectrophotometry

Preferred Specimen: 6 mL Fluid-source required

Collection Instruction: Submit fluid in sterile cup

Storage Instruction: Refrigerate

Turn Around Time: 9 days

CPT Code: 84157 (1), 84166 (1)

Scarlet Eligible: No

Cervical Biopsy

5660-6

Alternate Name: Cervical Biopsy

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 1 days

CPT Code: 88305 (1)

Scarlet Eligible: No

Allergen Cheddar Cheese (f81), IgE

0850-8

Alternate Name: Cheddar Cheese

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Cheesemolds (f82), IgE

1249-2

Alternate Name: Cheesemolds

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Chem 29 Profile

0165-1

Alternate Name: Chemprofile 24 and High-Density Lipoprotein

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 83718 (1), G0059 (1)

Scarlet Eligible: Yes

CH24/HDL,CBC/D/PLT

0855-7

Alternate Name: Chemprofile 24, High Density Lipoprotein,Complete Blood Count w/ Differential and Platelets

Preferred Specimen: SST Tube,Lavender top- EDTA

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 83718 (1), 85025 (1), G0059 (1)

Scarlet Eligible: Yes

Allergen Chick Pea (f309), IgE

1163-5

Alternate Name: Chick Pea

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Chicken Feathers (e85), IgE

2570-0

Alternate Name: Chicken Feathers

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Chicken Meat (f83), IgE

0817-7

Alternate Name: Chicken Meat

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Childhood Allergy Profile

K726-0

Alternate Name: Childhood Allergy Prof.

Preferred Specimen: SST Tube

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82785 (1), 86003 (16)

Scarlet Eligible: Yes

Allergen Chili Pepper (rf279), IgE

1423-3

Alternate Name: Chili Pepper

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Chlamydia NAAT

3852-1

Alternate Name: Chlamydia NAAT

Methodology: Nucleic Acid Amplification Test

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin Prep or SurePath vial, label with patients name.

Turn Around Time: 2 days

CPT Code: 87491 (1)

Scarlet Eligible: No

Chlamydia Trachomatis MultiPlex PCR - swab

J835-0

Alternate Name: Chlamydia PCR Swab

Methodology: Multiplex PCR

Preferred Specimen: Swab-Aptima Genprobe

Collection Instruction: Collect specimen-break swab shaft into transport tube, then label with name and source

Turn Around Time: 5 days

CPT Code: 87491 (1)

Scarlet Eligible: No

High Risk-STI CT/GC/Trich-TP/SP

6116-8

Alternate Name: Chlamydia, Gonorrhea, Trichomonas Vaginalis

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Turn Around Time: 3 days

CPT Code: 87491 (1), 87591 (1), 87661 (1)

Scarlet Eligible: No

Chlamydia + Gonorrhea ,Oral, Swab, Aptima

P010-2

Alternate Name: Chlamydia/gonococcus

Methodology: Transcription Mediated Amplification

Preferred Specimen: Swab-Aptima Genprobe

Collection Instruction: Collect specimen-break swab shaft into transport tube, then label with name and source

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87491 (1), 87591 (1)

Scarlet Eligible: No

Chlamydia and Gonorrhea, rRNA,TMA Aptima Swab

1005-8

Alternate Name: Chlamydia/gonococcus

Methodology: Transcription Mediated Amplification

Preferred Specimen: Swab-Aptima Genprobe

Collection Instruction: Collect specimen-break swab shaft into transport tube, then label with name and source

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87491 (1), 87591 (1)

Scarlet Eligible: No

Allergen Chocolate Cocoa (f93), IgE

1007-4

Alternate Name: Chocolate Cocoa

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Cholesterol

0058-8

Alternate Name: Cholesterol

Methodology: Enzyme Immunoassay

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82465 (1)

Scarlet Eligible: Yes

Chromium, Plasma/Serum

0277-4

Alternate Name: Chromium, Plasma/Serum**Methodology:** Inductively Coupled Plasma Mass Spectrometry**Preferred Specimen:** 2 mL Blue-Royal Blue Trace Element Serum (**RED LINE**)**Collection Instruction:** Fill tube, label with patient name**Storage Instruction:** Refrigerate**Turn Around Time:** 1-2 days**CPT Code:** 82495 (1)**Scarlet Eligible:** Yes

Chromogranin A

TA34-6

Alternate Name: Chromogranin A**Methodology:** Enzyme Linked Immunoabsorbance**Preferred Specimen:** Aliquot Tube-Serum**Collection Instruction:** Place 1-3 mL of serum in transport tube. Label as serum;refrigerate or freeze as required.**Storage Instruction:** STRICT Frozen**Turn Around Time:** 4 days**CPT Code:** 86316 (1)**Scarlet Eligible:** No

Karyotype CVS

A587-8

Alternate Name: Chromosome Analysis, CVS**Methodology:** Chromosome Analysis**Preferred Specimen:** 20 mg Chorionic Villi**Collection Instruction:** CVS specimen should be collected in a sterile, screw-top container with tissue culture transport medium (available upon Request). If transport medium is not avail., plain RPMI, Hanks solution or saline can be used**Storage Instruction:** Room Temp**Turn Around Time:** 7-14 days**CPT Code:** 88235 (1), 88261 (1), 88267 (1)**Scarlet Eligible:** No

Chromosome Analysis, POC, Tissue

1053-8

Alternate Name: Chromosome Analysis, POC, Tissue

Methodology: Chromosome Analysis

Preferred Specimen: 20-50 mL Products of Conception

Collection Instruction: Products of Conception (0.5-1.0 cm) specimen should be collected in a sterile, screw-top container filled with tissue culture transport medium (available upon Request from the laboratory). If transport medium is not available, plain RPMI, Hanks solution or saline can be used. Ship at room temperature or refrigerated.

Storage Instruction: Room Temp

Turn Around Time: 20 days

CPT Code: 88291 (1), 88230 (1), 88262 (1)

Scarlet Eligible: No

Chromosome Analysis,Amnio.FLD(IG)

2136-0

Alternate Name: Chromosome analysis, AF

Methodology: Chromosome Analysis

Preferred Specimen: 20 mL Amniotic Fluid

Collection Instruction: 20-25 ml Amniotic Fluid should be collected in a sterile, screw top leak-proof conical tube(s).

Storage Instruction: Room Temp

Turn Around Time: 7-14 days

CPT Code: 88269 (1), 88235 (1), 88280 (1)

Scarlet Eligible: No

Chromosome Analysis/Karyotype, Blood

0559-5

Alternate Name: Chromosome analysis, PB/Karotype

Methodology: Chromosome Analysis

Preferred Specimen: 2-5 mL Green top- Sodium Heparin

Collection Instruction: Fill tube, invert 8-10 times (DO NOT SHAKE!)

Storage Instruction: Room Temp

Turn Around Time: 7-14 days

CPT Code: 88230 (1), 88262 (1)

Scarlet Eligible: Yes

Chronic 10-Color Lymphoproliferative Panel

TH66

Alternate Name: Chronic 10-Color Lymphoproliferative Panel

Scarlet Eligible: No

Chronic 10-Color Lymphoproliferative Panel TC Only

TH65

Alternate Name: Chronic 10-Color Lymphoproliferative Panel TC Only

Scarlet Eligible: No

Chronic Fatigue Syndrome

F100-3

Alternate Name: Chronic Fatigue Syndrom

Preferred Specimen: SST Tube,Lavender top- EDTA

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 84436 (1), 84439 (1), 84443 (1), 85025 (1), 86038 (1), 86141 (1), 86618 (1), 86644 (1), 86645 (1), 86663 (1), 86664 (1), 86665 (2)

Scarlet Eligible: Yes

Lymphoproliferative Disorders Analysis (TC Only) and StormPath BMM ReportTemplate

R734-5

Alternate Name: Chronic Flow TC and BMM

Turn Around Time: 1 days

CPT Code: 99999 (1), 88184 (1), 88185 (17)

Scarlet Eligible: No

Hepatitis Comprehensive II

2302-8

Alternate Name: Chronic Hepatitis Profile

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86704 (1), 86706 (1), 86707 (1), 86708 (1), 86803 (1), 87340 (1), 87350 (1)

Scarlet Eligible: Yes

Lithium, Serum

0119-8

Alternate Name: Cibalith S; Eskalith; Lithane; Lithobid; Lithonate; Lithotabs; Pfl-Lith; Phasal

Methodology: Colorimetry

Preferred Specimen: Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80178 (1)

Scarlet Eligible: Yes

Cyclosporine, Whole Blood

0022-4

Alternate Name: Ciclosporin; CSA; Cya; Neoral®; Sandimmune®

Methodology: Enzyme Immunoassay

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Submit 2 lavenders when Cyclosporine is ordered in conjunction with any other test requiring lavender.

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80158 (1)

Scarlet Eligible: Yes

Allergen Coriander (rf317), IgE

2361-4

Alternate Name: Cilantro

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Cinnamon (rf220), IgE

2345-7

Alternate Name: Cinnamon

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Chloride, Serum

0057-0

Alternate Name: Cl, Serum

Methodology: Ion Selective Electrode

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82435 (1)

Scarlet Eligible: Yes

Chloride, Urine, Random

1794-7

Alternate Name: Cl, Urine Random

Methodology: Ion Selective Electrode

Preferred Specimen: 3 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82436 (1)

Scarlet Eligible: Yes

Chloride, Urine 24HR.

0439-0

Alternate Name: Cl, Urine, 24 Hours

Methodology: Ion Selective Electrode

Preferred Specimen: 3 mL Urine Container - 24hr

Collection Instruction: Discard first morning sample, then collect all urines after for 24 hrs, on completion label with collection time and patient's name; record total volume on Requisition

Storage Instruction: Refrigerate

CPT Code: 81050 (1), 82436 (1)

Scarlet Eligible: No

Allergen C. Herbarum (m2), IgE

0661-9

Alternate Name: Cladosporium herbarum (m2) IgE; HormoDendrum

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Clam (f207), IgE

1169-2

Alternate Name: Clam

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Drug Confirmation, Clonazepam, Urine (LC/MS)

L204-7

Alternate Name: Clonazepam & Metab Conf.

Methodology: Liquid Chromatography Tandem Mass Spectrometry

Preferred Specimen: 5 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 80346 (1)

Scarlet Eligible: Yes

Allergen Clove (rf268), IgE

2350-7

Alternate Name: Clove

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

CoEnzyme Q10 Deficiency, Primary, 7

J441-7

Alternate Name: CoEnzyme Q10 Def.

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81479 (1)

Scarlet Eligible: Yes

Allergen Cocklebur (w13), IgE

1939-8

Alternate Name: Cocklebur

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Cockroach (i6), IgE

0716-1

Alternate Name: Cockroach

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Coconut (f36), IgE

1318-5

Alternate Name: Coconut

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Codfish (f3), IgE

0840-9

Alternate Name: Codfish**Methodology:** Immunocap**Preferred Specimen:** 1 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 3 days**CPT Code:** 86003 (1)**Scarlet Eligible:** Yes**Coenzyme Q10**

J322-9

Alternate Name: Coenzyme Q10**Methodology:** Liquid Chromatography Tandem Mass Spectrometry**Preferred Specimen:** 1 mL Amber Serum**Collection Instruction:** Place 1-3 mL of serum in transport vial should be either an amber tube if not please wrap in foil to protect from light, label as serum and FREEZE.**Storage Instruction:** STRICT Frozen**Turn Around Time:** 3 days**CPT Code:** 82542 (1)**Scarlet Eligible:** No**Allergen Coffee (rf221), IgE**

1133-8

Alternate Name: Coffee**Methodology:** Immunocap**Preferred Specimen:** 3 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 3 days**CPT Code:** 86003 (1)**Scarlet Eligible:** Yes

Cold Agglutinins Titer

0063-8

Alternate Name: Cold Agglutinins

Methodology: Agglutination

Preferred Specimen: Aliquot Tube-Serum

Collection Instruction: Place 1-3 mL of serum in transport tube. Label as serum;refrigerate or freeze as required.

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86157 (1)

Scarlet Eligible: No

B-CrossLaps/Serum

T421-7

Alternate Name: Collagen Telopeptide

Methodology: Electrochemiluminescence Immunoassay

Preferred Specimen: Aliquot Tube-Serum

Collection Instruction: Place 1-3 mL of serum in transport tube. Label as serum;refrigerate or freeze as required.

Storage Instruction: STRICT Frozen

Turn Around Time: 7 days

CPT Code: 82523 (1)

Scarlet Eligible: No

Colorectal Cancer Panel

B274

Alternate Name: Colorectal Cancer Panel

Scarlet Eligible: No

Allergen Common Birch (t3), IgE

0685-8

Alternate Name: Common Birch**Methodology:** Immunocap**Preferred Specimen:** 3 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 2 days**CPT Code:** 86003 (1)**Scarlet Eligible:** Yes**Allergen Common Ragweed (w1) ,IgE**

0717-9

Alternate Name: Common Ragweed**Methodology:** Immunocap**Preferred Specimen:** 3 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 3 days**CPT Code:** 86003 (1)**Scarlet Eligible:** Yes**Allergen Common Reed (G7) IgE+**

2395-2

Alternate Name: Common Reed**Methodology:** Immunocap**Preferred Specimen:** 3 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 3 days**CPT Code:** 86003 (1)**Scarlet Eligible:** Yes

Comprehensive Thyroid-Expanded

9891-3

Alternate Name: Comp. Thyroid - Expand

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 84436 (1), 84439 (1), 84443 (1), 84479 (1), 84480 (1), 84481 (1), 86376 (1)

Scarlet Eligible: Yes

CBC w/Diff, Platelet Ct.

0053-9

Alternate Name: Complete Blood Count (CBC), Includes Differential and Platelets

Methodology: Flow Cytometry

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 85025 (1)

Scarlet Eligible: Yes

CBC w/o DIFF (Hemogram)/PLAT. CT

0034-9

Alternate Name: Complete Blood Count w/out Differential, Platelet Count

Methodology: Flow Cytometry

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 85027 (1)

Scarlet Eligible: Yes

Cortisol, Serum AM and PM

0374-9

Alternate Name: Compound F, Serum AM and PM; Hydrocortisone, Serum AM and PM

Methodology: Chemiluminescence

Preferred Specimen: 2 mL Red Top

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 82533 (2)

Scarlet Eligible: No

Cortisol, Random

0900-1

Alternate Name: Compound F, Serum, Random; Hydrocortisone, Serum, Random

Methodology: Chemiluminescence

Preferred Specimen: 1 mL Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 82533 (1)

Scarlet Eligible: Yes

Cortisol, PM

0067-9

Alternate Name: Compound F, Serum; Hydrocortisone, Serum

Methodology: Chemiluminescence

Preferred Specimen: 1 mL Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 82533 (1)

Scarlet Eligible: No

Comprehensive Common Cancer Panel

B275-9

Alternate Name: Comprehensive Cancer Panel

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: 2-5 mL Lavender top- EDTA

Collection Instruction: Lavender: 2 tubes (5-6ml each) Label with the patient name, date of birth and/or ID. Ideal handling is to ship immediately at ambient temperature for overnight delivery. Alternate sample type: DNA: Extract and submit DNA 50 ng/mL (total volume 400 mL) in sterile container. Ship DNA in sealed and protective packaging at ambient temperature, using a cool pack in hot weather. Labeling: Name of person or identifier, DNA concentration and volume MUST be provided on the specimen tube.

Storage Instruction: Room Temp

Turn Around Time: 14 days

CPT Code: 81162 (1), 81201 (1), 81203 (1), 81292 (1), 81294 (1), 81295 (1)

Scarlet Eligible: Yes

Congenital Amegakaryocytic Thrombocytopenia (CAMT)

J306-2

Alternate Name: Congenital Amegakaryocytic Thrombocytopenia (CAMT)

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81402 (1)

Scarlet Eligible: Yes

Conventional Pap Slide(s)

1100-7

Alternate Name: Conventional Pap
Methodology: Microscopic Examination
Preferred Specimen: Slide GYN
Collection Instruction: Label slides with patient name, place in slide holder for transport
Storage Instruction: Room Temp
Turn Around Time: 5 days
CPT Code: 88148 (1)
Scarlet Eligible: No

Allergen Corn (maize)(f8), IgE

0821-9

Alternate Name: Corn (Maize)
Methodology: Immunocap
Preferred Specimen: 3 mL SST Tube
Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes
Storage Instruction: Refrigerate
Turn Around Time: 3 days
CPT Code: 86003 (1)
Scarlet Eligible: Yes

Allergen Cotton-Fiber (o1), IgE

1436-5

Alternate Name: Cotton-Fiber
Methodology: Immunocap
Preferred Specimen: 3 mL SST Tube
Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes
Storage Instruction: Refrigerate
Turn Around Time: 3 days
CPT Code: 86003 (1)
Scarlet Eligible: Yes

Warfarin, Pharmacogenomic Assay

6261-2

Alternate Name: Coumadin Sensitivity, CYP2C9, VKORC1

Methodology: Polymerase Chain Reaction

Preferred Specimen: 2 mL Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 81355 (1)

Scarlet Eligible: No

Allergen Cow Dander (e4), IgE

0692-4

Alternate Name: Cow Dander

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Crab (f23), IgE

1168-4

Alternate Name: Crab

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Cranberry (rf341), IgE

3796-0

Alternate Name: Cranberry

Methodology: Immunocap

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Creatinine, Serum

0070-3

Alternate Name: Creatinine

Methodology: Colorimetry

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82565 (1)

Scarlet Eligible: Yes

Cryptosporidium EIA

0348-3

Alternate Name: Cryptosporidium EIA

Methodology: Enzyme Immunoassay

Preferred Specimen: Ova & Parasite Mini

Collection Instruction: With applicator stick, place 10 gms of stool into each vial, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 87328 (1)

Scarlet Eligible: No

Cryptosporidium ID, Modified Acid Fast

J745-1

Alternate Name: Cryptosporidium ID, (MAF)

Methodology: Modified Acid-Fast Stain

Preferred Specimen: Ova & Parasite Mini

Collection Instruction: With applicator stick, place 10 gms of stool into each vial, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 87015 (1), 87207 (1)

Scarlet Eligible: No

Crystal Identification, Body Fluid

0579-3

Alternate Name: Crystals, Fluid

Methodology: Microscopic Examination

Preferred Specimen: 1 mL Fluid-source required

Collection Instruction: Submit fluid in sterile cup

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 89060 (1)

Scarlet Eligible: No

Copper, Plasma/Serum

0479-6

Alternate Name: Cu, Plasma

Methodology: Inductively Coupled Plasma Mass Spectrometry

Preferred Specimen: Blue-Royal Blue Top K2 EDTA Plasma(*LAVENDER LINE)

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82525 (1)

Scarlet Eligible: Yes

Allergen Cucumber (f244), IgE

1664-2

Alternate Name: Cucumber

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Cultivated Rye (G12), IgE

2406-7

Alternate Name: Cultivated Rye

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Cultivated Wheat (g15), IgE

3221-9

Alternate Name: Cultivated Wheat (G15)

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Culture, Actinomyces

J011-8

Alternate Name: Culture Actinomyces, IUD Culture

Methodology: Bacterial Culture

Preferred Specimen: Swab-E

Collection Instruction: Place sample in sterile cup label with patient's name and source

Storage Instruction: Room Temp

Turn Around Time: 7 days

CPT Code: 87070 (1)

Scarlet Eligible: No

Medical Implants/Devices Culture

TG39-3

Alternate Name: Culture Implants/Device

Methodology: Bacterial Culture

Preferred Specimen: Cup Plain

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87070 (1)

Scarlet Eligible: No

Culture, Anal (Beta strep. Group B)

3773-9

Alternate Name: Culture, Anal for Beta Strep. Group B

Methodology: Bacterial Culture

Preferred Specimen: Swab-E

Collection Instruction: Collect specimen with swab then place swab into transport carrier; Label with name and source

Turn Around Time: 4 days

CPT Code: 87070 (1)

Scarlet Eligible: No

Culture, Environmental

0819-3

Alternate Name: Culture, Environmental

Methodology: Bacterial Culture

Preferred Specimen: Swab-E

Collection Instruction: Place sample in sterile cup, label with specimen source. DO NOT ADD fixative.

Turn Around Time: 7 days

CPT Code: 87070 (1)

Scarlet Eligible: No

Culture, Tracheal Airway

K909-2

Alternate Name: Culture, Tracheal Airway

Methodology: Bacterial Culture

Preferred Specimen: Cup-Sputum

Collection Instruction: Submit sputum in sterile container, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87205 (1), 87070 (1)

Scarlet Eligible: No

Allergen Curry (rf281), IgE

2356-4

Alternate Name: Curry

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Curvularia Lunata (m16), IgE

2476-0

Alternate Name: Curvularia Lunata

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Cystatin C

J013-4

Alternate Name: Cystatin C

Methodology: Immunoassay

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82610 (1)

Scarlet Eligible: Yes

Cystic Fibrosis 40- NGS

J435-9

Alternate Name: Cystic Fibrosis 40- NGS

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 7-10 days

CPT Code: 81220 (1)

Scarlet Eligible: Yes

Cytology, Bladder Washing (50% Alcohol)

J529-9

Alternate Name: Cyto. Bladder Washing

Methodology: Microscopic Examination

Preferred Specimen: 50% Alcohol Jar

Collection Instruction: Body fluid for cytology in 50% alcohol (ex...urine, sputum, nipple discharge). Add directly to an equal volume of 50% Alcohol Solution

Storage Instruction: Room Temp

Turn Around Time: 4 days

CPT Code: 88104 (1)

Scarlet Eligible: No

Cytology, Bladder Washing (ThinPrep)

J534-9

Alternate Name: Cyto. Bladder Washing-TP

Methodology: Microscopic Examination

Preferred Specimen: ThinPrep Cytolyt Solution

Collection Instruction: Collect specimen swish in Cytolyte vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 4 days

CPT Code: 88112 (1)

Scarlet Eligible: No

Oncology Chromosome Analysis- BM/PB

5250-6

Alternate Name: Cytogenetics,Karyotype, Cancer Cytogenetics Karyotype

Methodology: Chromosome Analysis

Preferred Specimen: 1-3 mL Bone Marrow - Green Top

Collection Instruction: Place bone marrow into tube, label with patient's name

Storage Instruction: Room Temp

Turn Around Time: 5-7 days

CPT Code: 88237 (1), 88262 (1), 88291 (1)

Scarlet Eligible: No

CMV IGM

0461-4

Alternate Name: Cytomegalovirus Antibody IgM
Methodology: Enzyme Linked Immunoabsorbance
Preferred Specimen: 1 mL SST Tube
Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes
Storage Instruction: Refrigerate
Turn Around Time: 2 days
CPT Code: 86645 (1)
Scarlet Eligible: Yes

CMV Ab. (IGG)

0400-2

Alternate Name: Cytomegalovirus Antibody, IgG
Methodology: Chemiluminescence
Preferred Specimen: 1 mL SST Tube
Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes
Storage Instruction: Refrigerate
Turn Around Time: 2 days
CPT Code: 86644 (1)
Scarlet Eligible: Yes

Cytopath Conc Tech - TC only

5210

Alternate Name: Cytopath Conc Tech - TC only
Scarlet Eligible: No

Cytopathology Non-Gyn TC only

B003

Alternate Name: Cytopathology Non-Gyn TC only
Scarlet Eligible: No

D-Dimer

5743-0

Alternate Name: D-DIMER

Methodology: Immunospectrometric

Preferred Specimen: 1 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA and patient's name.

Storage Instruction: STRICT Frozen

Turn Around Time: 2 days

CPT Code: 85379 (1)

Scarlet Eligible: No

Allergen D. Pteronyssinus (d1), IgE

0715-3

Alternate Name: D. Pteronyssinus

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

13Q14.3 BY FISH-TC ONLY

A878-1

Alternate Name: D13S319 by FISH, TC Only

Methodology: FISH

Preferred Specimen: 5 mL Peripheral Blood - Green Top

Collection Instruction: Fill tube, invert 8-10 times (DO NOT SHAKE!)

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88373 (1)

Scarlet Eligible: No

Decalcification (Nails/Bone)

3188-0

Alternate Name: DECALCIFICATION (NAILS/BONE)

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 2 days

CPT Code: 88311 (1)

Scarlet Eligible: No

Dehydroepiandrosterone (DHEA)

0178-4

Alternate Name: DHEA, Serum

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82626 (1)

Scarlet Eligible: Yes

DHEA Sulfate

0406-9

Alternate Name: DHEA-S, DHEA-SO4

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 82627 (1)

Scarlet Eligible: Yes

Dihydrotestosterone

J892-1

Alternate Name: DHT

Methodology: Liquid Chromatography Tandem Mass Spectrometry

Preferred Specimen: 2 mL Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 83789 (1)

Scarlet Eligible: Yes

Diabetes Profile

0391-3

Alternate Name: DIABETES PROFILE

Preferred Specimen: Lavender top- EDTA,SST Tube,Urine Urinalysis Tube - Yellow,Grey top- Sodium Fluoride

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 81001 (1), 82947 (1), 83036 (1), 85025 (1), G0059 (1)

Scarlet Eligible: Yes

Diabetic Screen

0222-0

Alternate Name: DIABETIC SCREEN

Preferred Specimen: Lavender top- EDTA,Grey top- Sodium Fluoride

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82947 (1), 83036 (1)

Scarlet Eligible: Yes

DIC Level 2

5872-7

Alternate Name: DIC Assesment Panel Level-1; Disseminated Intravascular Coagulation Panel Level-1

Preferred Specimen: 3 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA and patient's name.

Storage Instruction: STRICT Frozen

Turn Around Time: 10 days

CPT Code: 83520 (2)

Scarlet Eligible: No

DIC Assessment Panel Level-1

5839-6

Alternate Name: DIC Assesment Panel Level-2; Disseminated Intravascular Coagulation Panel Level-2

Preferred Specimen: Aliquot Plasma-Citrated-Blue Top,Lavender top- EDTA

Storage Instruction: Refrigerate

Turn Around Time: 4 days

CPT Code: 85025 (1), 85060 (1), 85379 (1), 85384 (1), 85610 (1), 85670 (1), 85730 (1)

Scarlet Eligible: No

DIC Assessment

5992-3

Alternate Name: DIC Assessment; DIC Screen

Preferred Specimen: Aliquot Plasma-Citrated-Blue Top,Lavender top- EDTA

Storage Instruction: STRICT Frozen, Refrigerate

Turn Around Time: 10 days

CPT Code: 83520 (2), 85025 (1), 85060 (1), 85379 (1), 85384 (1), 85610 (1), 85670 (1), 85730 (1)

Scarlet Eligible: No

Dihydropomide Dehydrogenase

A644-7

Alternate Name: DIHYDROLIPOMIDE DEHYDRO

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81400 (1)

Scarlet Eligible: Yes

DLBCL FISH PANEL: MYC BA, BCL2, BCL6

L218-7

Alternate Name: DLBCL FISH PANEL

Methodology: FISH

Preferred Specimen: Bone Marrow - Green Top

Collection Instruction: Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube. Also acceptable, whole blood in green (sodium heparin) tube. Label specimen with patient first and last names and date of birth.

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88374 (1), 88377 (2)

Scarlet Eligible: No

Lymphoma - Diffuse large B-Cell, BCL6 by FISH

5028-6

Alternate Name: DLBCL; BCL6, 3q27

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube. Also acceptable: Whole blood in green (sodium heparin) tube.

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88377 (1)

Scarlet Eligible: No

DNA Ploidy for Molar Pregnancy

5555

Alternate Name: DNA Ploidy for Molar Pregnancy

Scarlet Eligible: No

DNA Ploidy for Molar Pregnancy + P57

K735

Alternate Name: DNA Ploidy for Molar Pregnancy + P57

Scarlet Eligible: No

Dihydropyrimidine Dehydrogenase (DPD) Gene Mutation Analysis

6285-1

Alternate Name: DPD GENE MUTATION

Methodology: Polymerase Chain Reaction

Preferred Specimen: 5 mL Peripheral Blood - Lavender Top

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient's name

Storage Instruction: Refrigerate

Turn Around Time: 10-13 days

CPT Code: 81232 (1)

Scarlet Eligible: Yes

Drug Abuse Screen, Urine. (8 drug)

0628-8

Alternate Name: DRUG ABUSE SCREEN,URINE

Methodology: Enzyme Immunoassay

Preferred Specimen: 5 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80307 (1)

Scarlet Eligible: Yes

Allergen Casein (f78), IgE

1267-4

Alternate Name: Dairy

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86008 (1)

Scarlet Eligible: Yes

Allergen Dandelion (w8), IgE

1363-1

Alternate Name: Dandelion

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Deafness, Autosomal Recessive 1A (GJB2)

J308-8

Alternate Name: Deafness autosomal recessive 1A (GJB2)

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81252 (1)

Scarlet Eligible: Yes

Gliadin Antibodies, IgG, IgA

0608-0

Alternate Name: Deamidated Gliadin

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 83516 (2)

Scarlet Eligible: Yes

Meperidine, Urine, Screen, Qual. (200)

8545-6

Alternate Name: Demerol

Methodology: Enzyme Immunoassay

Preferred Specimen: 5 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80307 (1)

Scarlet Eligible: Yes

Valproic Acid

0542-1

Alternate Name: Depakote, Depacon, Depakene, Depakote, Depamidol Dipropylacet

Methodology: Enzyme Immunoassay

Preferred Specimen: Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80164 (1)

Scarlet Eligible: Yes

Allergen D. Farinae (Dust Mite) (d2), IgE1372-2

Alternate Name: Dermatophagoides farina**Methodology:** Immunocap**Preferred Specimen:** 3 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 3 days**CPT Code:** 86003 (1)**Scarlet Eligible:** Yes

Culture, Dialysate0857-3

Alternate Name: Dialysate Culture**Methodology:** Bacterial Culture**Preferred Specimen:** 2 mL Cup Plain**Collection Instruction:** Place sample in sterile cup label with patient's name and source**Storage Instruction:** Refrigerate**Turn Around Time:** 5 days**CPT Code:** 87070 (1)**Scarlet Eligible:** No

Culture, Dialysis Water0856-5

Alternate Name: Dialysis Water Culture**Methodology:** Bacterial Culture**Preferred Specimen:** 2 mL Cup Plain**Collection Instruction:** Place sample in sterile cup label with patient's name and source**Storage Instruction:** Refrigerate**Turn Around Time:** 5 days**CPT Code:** 87070 (1)**Scarlet Eligible:** No

Allergen Dill (f277), IgE

3605-3

Alternate Name: Dill

Methodology: Immunocap

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

DRVVT Complete

5730-7

Alternate Name: Dilute Russell's viper venom time

Methodology: Clot Detection

Preferred Specimen: 2 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and submit frozen.

Storage Instruction: STRICT Frozen

Turn Around Time: 3 days

CPT Code: 85613 (2)

Scarlet Eligible: No

Coombs Test, Direct

0064-6

Alternate Name: Direct Coombs

Methodology: Hemagglutination Assay

Preferred Specimen: Pink Tube-EDTA

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86880 (1)

Scarlet Eligible: Yes

Sickle Cell Screen

0366-5

Alternate Name: Dithionite Test; Itano Solubility Test; Murayama test; Sickle Cell Preparation, Metabisulfite Test; Sickle Cell Solubility Test; Sickledex; Sickle Screen; Sickle-Sol; Sickle-Stat

Methodology: Solubility

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 85660 (1)

Scarlet Eligible: Yes

Allergen Dog Dander (e5), IgE

0664-3

Alternate Name: Dog Dander

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Duck Feathers (e86), IgE

0846-6

Alternate Name: Duck Feathers

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Dust Mix (hx2), IgE

2182-4

Alternate Name: Dust Mix

Methodology: Immunocap

Preferred Specimen: 5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86005 (1)

Scarlet Eligible: Yes

Dyskeratosis Congenita, Autosomal Recessive 5

J299-9

Alternate Name: Dyskeratosis Congenita autosomal Recessive 5

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81479 (1)

Scarlet Eligible: Yes

Escherichia coli (Non NY)

J968-9

Alternate Name: E. COLI

Methodology: Real-Time Polymerase Chain Reaction

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87798 (1)

Scarlet Eligible: No

Enterococcus faecalis (Non NY)

J967-1

Alternate Name: E. faecalis

Methodology: Real-Time Polymerase Chain Reaction

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87798 (1)

Scarlet Eligible: No

EDDP (Methadone Metabolite), Urine, Screen, Qual. (300)

6286-9

Alternate Name: EDDP

Methodology: Enzyme Immunoassay

Preferred Specimen: 5 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80307 (1)

Scarlet Eligible: Yes

EGFR, if neg reflex to ALK- Non Ny

J105-8

Alternate Name: EGFR NGS RFX NEG>

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 10-15 unstained slides at 5 microns (tumor content>10%) with H&E, Shavings (if tumor content >40%) with H&E, if sending DNA - Please call CS.

Storage Instruction: Room Temp

Turn Around Time: 5-10 days

CPT Code: 81235 (1)

Scarlet Eligible: No

EGFR, if neg reflex to ALK NY Only

J461-5

Alternate Name: EGFR NGS RFX NEG>

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 10-15 unstained slides at 5 microns (tumor content >10%) with H&E, Shavings (if tumor content >10%) with H&E, if sending DNA- Please call CS

Turn Around Time: 5-10 days

CPT Code: 81235 (1)

Scarlet Eligible: No

EGFR, if neg reflex to ALK - TC only Non NY

J202-3

Alternate Name: EGFR NGS RFX NEG>ALK

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 10-15 unstained slides at 5 microns (tumor content>10%) with H&E, Shavings (if tumor content >40%) with H&E, if sending DNA - Please call CS.

Storage Instruction: Room Temp

Turn Around Time: 5-10 days

CPT Code: 81235 (1)

Scarlet Eligible: No

EGFR, if neg reflex to ALK-TC only- NY Only

J466-4

Alternate Name: EGFR NGS RFX NEG>ALK

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 10-15 unstained slides at 5 microns (tumor content >10%) with H&E, Shavings (if tumor content >10%) with H&E, if sending DNA- Please call CS

Turn Around Time: 5-10 days

CPT Code: 81235 (1)

Scarlet Eligible: No

EGFR, if neg reflex to ALK, if neg reflex to ROS1- Non NY

J106-6

Alternate Name: EGFR RFX >ALK>ROS1

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 10-15 unstained slides at 5 microns (tumor content>10%) with H&E, Shavings (if tumor content >40%) with H&E, if sending DNA - Please call CS.

Storage Instruction: Room Temp

Turn Around Time: 5-10 days

CPT Code: 81235 (1)

Scarlet Eligible: No

EGFR, if neg reflex to ALK, if neg reflex to ROS1 NY Only

J460-7

Alternate Name: EGFR RFX>ALK>ROS1

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 10-15 unstained slides at 5 microns (tumor content >10%) with H&E, Shavings (if tumor content >10%) with H&E, if sending DNA- Please call CS

Turn Around Time: 5-10 days

CPT Code: 81235 (1)

Scarlet Eligible: No

EGFR by FISH

A240-4

Alternate Name: EGFR by FISH (NON-NY)

Methodology: FISH

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: Please write surgical number on block and on requisition. Provide fixation time on requisition as well

Storage Instruction: Room Temp

Turn Around Time: 11 days

CPT Code: 88365 (1)

Scarlet Eligible: No

EGFR Sequencing NY Only

J454-0

Alternate Name: EGFR by NGS

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 10-15 unstained slides at 5 microns (tumor content >10%) with H&E, Shavings (if tumor content >10%) with H&E, if sending DNA- Please call CS

Turn Around Time: 5-10 days

CPT Code: 81235 (1)

Scarlet Eligible: No

EGFR and KRAS, if both neg, reflex ALK, if neg reflex ROS1 Non NY

H579-6

Alternate Name: EGFR+ KRAS>ALK>ROS1

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Turn Around Time: 5-10 days

CPT Code: 81235 (1), 81275 (1), 81276 (1)

Scarlet Eligible: No

EGFR by NGS, ALK, ROS1, and MET by FISH- Non NY

H581-2

Alternate Name: EGFR, ALK, ROS1, MET

Methodology: Genotyping by Next Generation Sequencing

Turn Around Time: 5-10 days

CPT Code: 81235 (1), 88377 (3)

Scarlet Eligible: Yes

EGFR by NGS, ALK, ROS1, and MET by FISH (Ny Only)

K335-0

Alternate Name: EGFR-NGS,ALK,ROS1,MET

Turn Around Time: 10 days

CPT Code: 81235 (1), 88377 (3)

Scarlet Eligible: No

EGFR and KRAS, if both neg, reflex ALK, if neg reflex ROS 1 NY Only

J465-6

Alternate Name: EGFR/KRAS>ALK>ROS1

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 10-15 unstained slides at 5 microns (tumor content >10%) with H&E, Shavings (if tumor content >10%) with H&E, if sending DNA- Please call CS

Turn Around Time: 5-10 days

CPT Code: 81235 (1), 81275 (1), 81276 (1)

Scarlet Eligible: No

Electrolytes (24 hr urine)

8017-6

Alternate Name: ELECTROLYTES (24 hr U)

Turn Around Time: 1 days

CPT Code: 81050 (1), 82436 (1), 84133 (1), 84300 (1)

Scarlet Eligible: No

Electrolytes, Urine, Random

0989-4

Alternate Name: ELECTROLYTES, URINE, RANDOM

Preferred Specimen: 10 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82436 (1), 84133 (1), 84300 (1)

Scarlet Eligible: Yes

Encainide + Metabolite Level

0410-1

Alternate Name: ENCAINIDE + METABOLITE LEVEL

Methodology: Liquid Chromatography Tandem Mass Spectrometry

Preferred Specimen: 2 mL Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 80375 (1)

Scarlet Eligible: Yes

Endo-Cervical Biopsy

5662-2

Alternate Name: ENDO-CERVICAL BIOPSY

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 1 days

CPT Code: 88305 (1)

Scarlet Eligible: No

Endometrial Biopsy

5664-8

Alternate Name: ENDOMETRIAL BIOPSY

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Storage Instruction: Room Temp

Turn Around Time: 48 hours

CPT Code: 88305 (1)

Scarlet Eligible: No

Enterococcus/VRE Screen-@

6108-5

Alternate Name: ENTEROCOCCUS/VRE SCREEN

Methodology: Bacterial Culture

Preferred Specimen: Stool Cont.-Cup

Collection Instruction: Using applicator stick, place 20-30 gms of stool into cup, label with patient's name

Storage Instruction: Refrigerate

Turn Around Time: 4 days

CPT Code: 87081 (1)

Scarlet Eligible: No

Enterovirus RNA, Qualitative

3788-7

Alternate Name: ENTEROVIRUS RNA, Qual

Methodology: Polymerase Chain Reaction

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 87498 (1)

Scarlet Eligible: Yes

Erythropoietin, Serum

0183-4

Alternate Name: EPO

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 2-3 times, let stand for 20 minutes, spin for 10-15 minutes label with patient name. This specimen must be sent to lab within 24 hours of draw time. If not being submitted same day as collection, specimen must be spun down, poured into an aliquot and frozen.

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82668 (1)

Scarlet Eligible: Yes

Epstein-Barr Virus By PCR

A315-4

Alternate Name: EPSTEIN-BARR VIRUS BY PCR**Methodology:** Polymerase Chain Reaction**Preferred Specimen:** Formalin-fixed, Paraffin-embedded Tissue**Collection Instruction:** This comes in block form from client with surgical number imprint.**Turn Around Time:** 7 days**CPT Code:** 87798 (1)**Scarlet Eligible:** No**ER,PR,KI-67,P53,HER2 FISH, TC**

A278-4

Alternate Name: ER,PR,KI-67,P53,HER2 FISH, TC**Methodology:** Immunohistochemistry**Preferred Specimen:** Formalin-fixed, Paraffin-embedded Tissue**Collection Instruction:** This comes in block form from client with surgical number imprint.**Turn Around Time:** 3 days**CPT Code:** 88361 (4), 88374 (1)**Scarlet Eligible:** No**ER / PR / Ki-67/ HER2 (Tech Only)**

5433-8

Alternate Name: ER/ PR/Ki-67/HER2 (Tech Only)**Methodology:** Immunohistochemistry**Preferred Specimen:** Formalin-fixed, Paraffin-embedded Tissue**Collection Instruction:** This comes in block form from client with surgical number imprint.**Turn Around Time:** 1 days**CPT Code:** 88361 (4)**Scarlet Eligible:** No

ER/PR (Tech Only)

5425-4

Alternate Name: ER/PR (TECH COMPONENT ONLY)

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88361 (2)

Scarlet Eligible: No

ER/PR - w/Interpretation w/digital imaging

5198-7

Alternate Name: ER/PR RECEP.PROFILE

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88361 (2)

Scarlet Eligible: No

ER/ PR/ DNA/ KI-67/ HER2 - Reflex to HER2 by FISH

5408-0

Alternate Name: ER/PR/DNA/KI-67/HER2 W/RFLX

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 3-5 days

CPT Code: 88182 (1), 88361 (4)

Scarlet Eligible: No

ER/PR /DNA/ Ki-67 - w/Interp w/digital imaging

5160-7

Alternate Name: ER/PR/DNA/Ki-67 W/INTERP.

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88182 (1), 88361 (3)

Scarlet Eligible: No

ER/PR/HER2 - w/Interpretation w/digital imaging

5161-5

Alternate Name: ER/PR/HER2-NEU W/INTERP.

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88361 (3)

Scarlet Eligible: No

ER/PR/HER2w/digital imaging- Rflex to HER2 by FISH

5406-4

Alternate Name: ER/PR/HER2-RFX BY FISH

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 3-5 days

CPT Code: 88361 (3)

Scarlet Eligible: No

ER/ PR/ KI-67/ HER2 - Reflex to HER2 by FISH

5407-2

Alternate Name: ER/PR/KI-67/HER2 W/REFLEX

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 3-5 days

CPT Code: 88361 (4)

Scarlet Eligible: No

ER/PR/ Ki-67/ HER2 - w/Interp w/digital imaging

5163-1

Alternate Name: ER/PR/KI-67/HER2/IHC/INTER.

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88361 (4)

Scarlet Eligible: No

ER/PRKI-67/HER2/P53 BY IHC with Interp. w/digital imaging

A954-0

Alternate Name: ER/PRKI-67/HER2/P53

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 2 days

CPT Code: 88361 (5)

Scarlet Eligible: No

ER/PRKI-67/HER2/P53 BY IHC,TC

A951-6

Alternate Name: ER/PRKI-67/HER2/P53-TC

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88361 (5)

Scarlet Eligible: No

ERG by IHC (EP111)

T411-8

Alternate Name: ERG by IHC

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Storage Instruction: Room Temp

Turn Around Time: 1-2 days

CPT Code: 88342 (1)

Scarlet Eligible: No

ERG by IHC (EP111) - Tech Only

T412-6

Alternate Name: ERG by IHC -TC

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Storage Instruction: Room Temp

Turn Around Time: 1-2 days

CPT Code: 88342 (1)

Scarlet Eligible: No

Esophageal Biospy X1

5630-9

Alternate Name: ESOPHAGEAL BIOPSY**Methodology:** Histology**Preferred Specimen:** Biopsy In Formalin**Collection Instruction:** Submit in 10% Neutral Buffered Formalin, varied jar sizes.**Turn Around Time:** 1 days**CPT Code:** 88305 (1)**Scarlet Eligible:** No**Estradiol (E2), Ultrasensitive, LC/MS, Serum**

J258-5

Alternate Name: ESTRADIOL (E2), SERUM**Methodology:** Liquid Chromatography Tandem Mass Spectrometry**Preferred Specimen:** 1 mL Red Top**Collection Instruction:** Fill tube, invert gently 2-3 times, let stand for 20 minutes, spin for 10-15 minutes label with patient name.**Storage Instruction:** Refrigerate**Turn Around Time:** 3 days**CPT Code:** 82670 (1)**Scarlet Eligible:** Yes**Estrogen Panel**

2317-6

Alternate Name: ESTROGEN PANEL**Collection Instruction:** Fill tube, invert gently 2-3 times, let stand for 20 minutes, spin for 10-15 minutes label with patient name.**Turn Around Time:** 3 days**CPT Code:** 82670 (1), 82677 (1), 82679 (1)**Scarlet Eligible:** Yes

Estrone (E1), Serum

J259-3

Alternate Name: ESTRONE (E1), SERUM

Methodology: Liquid Chromatography Tandem Mass Spectrometry

Preferred Specimen: 1 mL Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82679 (1)

Scarlet Eligible: Yes

Drug Confirmation, Ethanol, Urine (GC/HS)

6239-8

Alternate Name: ETHANOL URN.,CONFIRMATION

Methodology: Headspace Gas Chromatography

Preferred Specimen: 1 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80320 (1)

Scarlet Eligible: Yes

Expanded CF, Fragile X, SMA

P150-6

Alternate Name: EXPAN. CF, FRAG X, SMA

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Turn Around Time: 4-10 days

CPT Code: 81220 (1), 81243 (1), 81400 (1)

Scarlet Eligible: Yes

Expanded Anti-Phos/Lupus

5990-7

Alternate Name: EXPANDED ANTI-PHOS/LUPUS

Preferred Specimen: Aliquot Plasma-Citrated-Blue Top, Aliquot Serum- Red Top, Lavender top-EDTA

Collection Instruction: Draw 1 lavender keep room temp. Draw 2 plain red spin, transfer serum into 4 aliquots and freeze. Draw 3 lt blue spin, transfer plasma into 4 aliquots and freeze.

Storage Instruction: STRICT Frozen, Freeze, Refrigerate

Turn Around Time: 8 days

CPT Code: 83516 (2), 83520 (2), 85025 (1), 85597 (1), 85598 (1), 85610 (1), 85611 (8), 85613 (2), 85730 (2), 85732 (8), 86146 (3), 86147 (3), 86148 (2), 86849 (1)

Scarlet Eligible: No

ExtraTube(s) Drawn

B978-8

Alternate Name: EXTRA TUBE(S) DRAWN

CPT Code: 99999 (1)

Scarlet Eligible: No

Culture, Ear

0087-7

Alternate Name: Ear Culture

Methodology: Bacterial Culture

Preferred Specimen: Swab-E

Collection Instruction: Collect specimen with swab then place swab into transport carrier; Label with name and source

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87070 (1)

Scarlet Eligible: No

Culture, Ear (Left)

0462-2

Alternate Name: Ear Culture, Left

Methodology: Bacterial Culture

Preferred Specimen: Swab-E

Collection Instruction: Collect specimen with swab then place swab into transport carrier; Label with name and source

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87070 (1)

Scarlet Eligible: No

Culture, Ear (Right)

0634-6

Alternate Name: Ear Culture, Right

Methodology: Bacterial Culture

Preferred Specimen: Swab-E

Collection Instruction: Collect specimen with swab then place swab into transport carrier; Label with name and source

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87070 (1)

Scarlet Eligible: No

Egg Components IgE

H992-1

Alternate Name: Egg Components IgE

Methodology: Immunoassay

Preferred Specimen: 1 mL SST Tube

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82785 (1), 86003 (2)

Scarlet Eligible: Yes

Allergen Egg White (f1), IgE

0671-8

Alternate Name: Egg White

Methodology: Immunocap

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Egg White w/Reflex to Component Panel

M066-9

Alternate Name: Egg White w/Reflex to Component Panel

Turn Around Time: 3 days

Scarlet Eligible: Yes

Allergen Egg Yolk (f75), IgE

0849-0

Alternate Name: Egg Yolk

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86008 (1)

Scarlet Eligible: Yes

Allergen Eggplant (rf262), IgE

2265-7

Alternate Name: Eggplant

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Ehlers-Danlos Syndrome, Type VIIC

J294-0

Alternate Name: Ehlers-Danlos syndrome type VIIC

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient's name

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81479 (1)

Scarlet Eligible: Yes

Electrolytes

0002-6

Alternate Name: Electrolytes

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80051 (1)

Scarlet Eligible: Yes

Protein Electrophoresis

0085-1

Alternate Name: Electrophoresis, Protein Electrophoresis Serum, SPEP

Methodology: Electrophoresis

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 84155 (1), 84165 (1)

Scarlet Eligible: Yes

Allergen Elm Tree (t8), IgE

0659-3

Alternate Name: Elm Tree

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Eng. Plantain (w9), IgE

0650-2

Alternate Name: Eng. Plantain

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Nasal Smear (for Eosinophils)

0248-5

Alternate Name: Eosinophils, Smear**Methodology:** Microscopic Examination**Preferred Specimen:** Slides**Collection Instruction:** Label slides with patient name, place in slide holder for transport**Turn Around Time:** 2 days**CPT Code:** 89190 (1)**Scarlet Eligible:** No**Eosinophil Ct, Urine**

1183-3

Alternate Name: Eosinophils, Urine**Methodology:** Microscopic Examination**Preferred Specimen:** 5 mL Urine Cup**Collection Instruction:** Collect urine in cup, label with patients name**Storage Instruction:** Refrigerate**Turn Around Time:** 1 days**CPT Code:** 89050 (1)**Scarlet Eligible:** Yes**EBV Panel (Capsid/EBNA G+M)**

0572-8

Alternate Name: Epstein Barr virus**Preferred Specimen:** 2 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 1 days**CPT Code:** 86663 (1), 86664 (1), 86665 (2)**Scarlet Eligible:** Yes

EBV, Early Antigen Ab

0582-7

Alternate Name: Epstein-Barr Virus**Methodology:** Chemiluminescence**Preferred Specimen:** 1 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 1 days**CPT Code:** 86663 (1)**Scarlet Eligible:** Yes**EBV, Nuclear Antigen Ab, IgG**

0583-5

Alternate Name: Epstein-Barr Virus, IgG**Methodology:** Chemiluminescence**Preferred Specimen:** 1 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 1 days**CPT Code:** 86664 (1)**Scarlet Eligible:** Yes**EBV Capsid Ab, IgG**

0234-5

Alternate Name: Epstein-Barr Virus, IgG**Methodology:** Chemiluminescence**Preferred Specimen:** 1 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 1 days**CPT Code:** 86665 (1)**Scarlet Eligible:** Yes

EBV Capsid Ab, IgM

0580-1

Alternate Name: Epstein-Barr Virus, IgM

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86665 (1)

Scarlet Eligible: Yes

ESR (Sed-Rate)

0086-9

Alternate Name: Erythrocyte Sedimentation Rate (ESR); Sed-Rate

Methodology: Manual

Preferred Specimen: 1 mL Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 85651 (1)

Scarlet Eligible: Yes

Estriol, Unconjugated

3116-1

Alternate Name: Estriol, Free; 16-Hydroxyestradiol; unconjugated Estriol, Pregnancy; Unconjugated Estrogen, Serum

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 82677 (1)

Scarlet Eligible: Yes

Estrogens, Fractionated

H906-1

Alternate Name: Estrogens, Fractionated

Methodology: Liquid Chromatography Tandem Mass Spectrometry

Preferred Specimen: 1 mL Red Top

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82671 (1)

Scarlet Eligible: Yes

Drug Confirmation, Ethanol, Blood (GC/HS)

B784-0

Alternate Name: Ethanol (GC-HS), BLD

Methodology: Headspace Gas Chromatography

Preferred Specimen: Grey top- Sodium Fluoride

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80320 (1)

Scarlet Eligible: Yes

Alcohol (Ethanol), Blood

0985-2

Alternate Name: Ethanol, Blood; Ethyl Alcohol, Blood

Methodology: Enzyme Immunoassay

Preferred Specimen: Grey top- Sodium Fluoride

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82077 (1)

Scarlet Eligible: Yes

Alcohol (Ethanol), Urine, Screen (20)

0995-1

Alternate Name: Ethanol, Urine, Ethyl Alcohol, Urine

Methodology: Enzyme Immunoassay

Preferred Specimen: 5 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80307 (1)

Scarlet Eligible: Yes

Allergen Eucalyptus (t18), IgE

1480-3

Alternate Name: Eucalyptus

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

EWSR1(22Q12) Gene Rearrang. FISH

B556-2

Alternate Name: Ewing Sarcoma FISH

Methodology: FISH

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Storage Instruction: Room Temp

Turn Around Time: 3-5 days

CPT Code: 88377 (1)

Scarlet Eligible: No

Expanded Cystic Fibrosis Panel

J434-2

Alternate Name: Expanded Cystic Fibrosis Panel

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 7-10 days

CPT Code: 81220 (1)

Scarlet Eligible: Yes

Factor II Activity

5724-0

Alternate Name: FACT. II ACT.

Methodology: Clot Detection

Preferred Specimen: 2 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name, DOB and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 4 days

CPT Code: 85210 (1)

Scarlet Eligible: No

Factor IX Activity

5769-5

Alternate Name: FACT. IX ACT.

Methodology: Clot Detection

Preferred Specimen: 2 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 4 days

CPT Code: 85250 (1)

Scarlet Eligible: No

Factor V Activity5725-7

Alternate Name: FACT. V ACT.**Methodology:** Clot Detection**Preferred Specimen:** 2 mL Aliquot Plasma-Citrated-Blue Top**Collection Instruction:** Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name, DOB and freeze.**Storage Instruction:** STRICT Frozen**Turn Around Time:** 4 days**CPT Code:** 85220 (1)**Scarlet Eligible:** No**Factor VII Activity**5728-1

Alternate Name: FACT. VII ACT.**Methodology:** Clot Detection**Preferred Specimen:** 2 mL Aliquot Plasma-Citrated-Blue Top**Collection Instruction:** Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.**Storage Instruction:** STRICT Frozen**Turn Around Time:** 4 days**CPT Code:** 85230 (1)**Scarlet Eligible:** No

Factor VIII Activity5729-9

Alternate Name: FACT. VIII ACT.**Methodology:** Clot Detection**Preferred Specimen:** 2 mL Aliquot Plasma-Citrated-Blue Top**Collection Instruction:** Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name, DOB and freeze.**Storage Instruction:** STRICT Frozen**Turn Around Time:** 4 days**CPT Code:** 85240 (1)**Scarlet Eligible:** No

Factor X Activity5735-6

Alternate Name: FACT. X ACT.**Methodology:** Clot Detection**Preferred Specimen:** 2 mL Aliquot Plasma-Citrated-Blue Top**Collection Instruction:** Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.**Storage Instruction:** STRICT Frozen**Turn Around Time:** 4 days**CPT Code:** 85260 (1)**Scarlet Eligible:** No

Factor XI Activity

5737-2

Alternate Name: FACT. XI ACT.**Methodology:** Clot Detection**Preferred Specimen:** 2 mL Aliquot Plasma-Citrated-Blue Top**Collection Instruction:** Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.**Storage Instruction:** STRICT Frozen**Turn Around Time:** 4 days**CPT Code:** 85270 (1)**Scarlet Eligible:** No**Factor XII Activity**

5738-0

Alternate Name: FACT. XII ACT.**Methodology:** Clot Detection**Preferred Specimen:** 2 mL Aliquot Plasma-Citrated-Blue Top**Collection Instruction:** Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.**Storage Instruction:** STRICT Frozen**Turn Around Time:** 4 days**CPT Code:** 85280 (1)**Scarlet Eligible:** No**Factor V Mutation (Leiden)**

5726-5

Alternate Name: FACTOR V MUTATION**Methodology:** Polymerase Chain Reaction**Preferred Specimen:** 4 mL Lavender top- EDTA**Collection Instruction:** Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!**Storage Instruction:** Room Temp**Turn Around Time:** 5 days**CPT Code:** 81241 (1)**Scarlet Eligible:** Yes

Factor X Antigen

5736-4

Alternate Name: FACTOR X AG.

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 2 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.

Storage Instruction: STRICT Frozen

CPT Code: 85260 (1)

Scarlet Eligible: No

Factor X Chromogenic

5707-5

Alternate Name: FACTOR X CHROMOGENIC

Methodology: Chromogenic

Preferred Specimen: 1 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 6 days

CPT Code: 85260 (1)

Scarlet Eligible: No

Factor XIII Activity

5739-8

Alternate Name: FACTOR XIII ACT.

Methodology: Chromogenic

Preferred Specimen: Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 5 days

CPT Code: 85290 (1)

Scarlet Eligible: No

Familial Hyperinsulinism

A643-9

Alternate Name: FAM-HYPERINSULINISM-SENDOUT

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 7-10 days

CPT Code: 81400 (1), 81401 (1)

Scarlet Eligible: Yes

Familial Dysautonomia

A697-5

Alternate Name: FAMILIAL DYSAUTONOMIA

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81260 (1)

Scarlet Eligible: Yes

Fanconi Anemia, Type C

A656-1

Alternate Name: FANCC

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81242 (1)

Scarlet Eligible: Yes

Fax Result To Doctor

9993-7

Alternate Name: FAX RESULT TO DOCTOR

CPT Code: 99999 (1)

Scarlet Eligible: Yes

Glucose, Fasting

0095-0

Alternate Name: FBS, Plasma glucose

Methodology: Enzyme Immunoassay

Preferred Specimen: Grey top- Sodium Fluoride

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82947 (1)

Scarlet Eligible: Yes

Fecal Calprotectin

J920-0

Alternate Name: FECAL CALPROTECTIN

Methodology: Chemiluminescence

Preferred Specimen: Stool Cont.-Cup

Collection Instruction: Using applicator stick, place 20-30 gms of stool into cup, label with patient's name

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 83993 (1)

Scarlet Eligible: No

Fecal Lactoferrin

3374-6

Alternate Name: FECAL LACTOFERRIN

Methodology: Immunoassay

Preferred Specimen: Stool Cont.-Cup

Collection Instruction: Using applicator stick, place 20-30 gms of stool into cup, label with patient's name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 83630 (1)

Scarlet Eligible: No

Fentanyl, Urine, Screen, Qual. (2)

6194-5

Alternate Name: FENTANYL

Methodology: Enzyme Immunoassay

Preferred Specimen: 1 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80307 (1)

Scarlet Eligible: Yes

Ferritin

0088-5

Alternate Name: FERRITIN

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82728 (1)

Scarlet Eligible: Yes

Fetal Fibronectin

3576-6

Alternate Name: FFN
Methodology: Enzyme Immunoassay
Preferred Specimen: Fetal Fibronectin
Collection Instruction: Special swab-contact the lab
Storage Instruction: Refrigerate
Turn Around Time: 1 days
CPT Code: 82731 (1)
Scarlet Eligible: No

FGFR1 by FISH (Tech Only)

TG46-8

Alternate Name: FGFR1-TC
Methodology: FISH
Preferred Specimen: Bone Marrow - Green Top
Storage Instruction: Room Temp
Turn Around Time: 3 days
CPT Code: 88374 (1)
Scarlet Eligible: Yes

FGFR2/IGH t(4;14) by FISH

A291

Alternate Name: FGFR2/IGH t(4;14) by FISH
Scarlet Eligible: No

Fibrinogen, Activity

0089-3

Alternate Name: FIBRINOGEN

Methodology: Photometric

Preferred Specimen: 1 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA and patient's name.

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 85384 (1)

Scarlet Eligible: Yes

Fibrinolytic Evaluation - Standard

5860-2

Alternate Name: FIBRINOLYTIC EVALUATION STANDARD

Preferred Specimen: Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect 3 full blue-top (citrated) tubes. Centrifuge Blue-tops at high speed for 15 minutes. Transfer plasma into plastic 4 transfer tubes. Label with date/time collected, CITRATED PLASMA and patient's name. Must be submitted Frozen.

Turn Around Time: 4 days

CPT Code: 85379 (1), 85384 (1), 85420 (1), 85610 (1), 85670 (1), 85730 (1)

Scarlet Eligible: No

Chronic Eosinophilic Leukemia: FIP1L1-PDGFRα by FISH

5182-1

Alternate Name: FIP1L1-PDGFRα by FISH, CHIC2

Methodology: FISH

Preferred Specimen: 5 mL Peripheral Blood - Green Top

Collection Instruction: Fill tube, invert 8-10 times (DO NOT SHAKE!)

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88377 (1)

Scarlet Eligible: Yes

Lymphoma - Non-Hodgkin, BCL1, BCL2, BCL6 (FISH Panel)

5273-8

Alternate Name: FISH Panel: t(11;14), t(14;18), 3q27, t(8;14)

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube. Also acceptable: Whole blood in green (sodium heparin) tube.

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88374 (1), 88377 (2)

Scarlet Eligible: No

Multiple Myeloma FISH Panel

5282-9

Alternate Name: FISH Panel: 13q14, t(11;14) & +11, 17p13, t(4;14), 1q21, t(14;16)

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube.

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88367 (1), 88373 (1), 88374 (4)

Scarlet Eligible: No

FISH, FGFR1, 8p11-12

B557-0

Alternate Name: FISH, FGFR1, 8p11-12

Methodology: FISH

Preferred Specimen: 3 mL Peripheral Blood - Green Top

Collection Instruction: Fill tube, invert 8-10 times (do not shake). Clinical history/reason for referral is required with test order. Prior therapy and transplant history should be provided with test order.

Storage Instruction: Room Temp

Turn Around Time: 3 days

CPT Code: 88374 (1)

Scarlet Eligible: No

OnkoSight FLT3

A044-0

Alternate Name: FLT3 Mutation Analysis

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Peripheral Blood - Lavender Top

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient's name

Storage Instruction: Room Temp

Turn Around Time: 5-10 days

CPT Code: 81246 (1), 81245 (1)

Scarlet Eligible: Yes

Onkosight NGS FLT3 IF NEG. RFLX TO CEBPA+NPM1 NGS

J154-6

Alternate Name: FLT3 NGS RFX CEPBA+NPM1

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Peripheral Blood - Lavender Top

Collection Instruction: Place bone marrow into tube, label with patient's name

Storage Instruction: Room Temp

Turn Around Time: 7-14 days

CPT Code: 81246 (1), 81245 (1)

Scarlet Eligible: Yes

Onkosight NGS FLT3 AND NPM1 Sequencing

J155-3

Alternate Name: FLT3+NPM1 BY NGS

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Peripheral Blood - Lavender Top

Collection Instruction: Place bone marrow into tube, label with patient's name

Storage Instruction: Room Temp

Turn Around Time: 7-14 days

CPT Code: 81246 (1), 81245 (1), 81310 (1)

Scarlet Eligible: Yes

Influenza A/B, RSV BY PCR

1410-0

Alternate Name: FLU A/B, RSV BY PCR

Methodology: Polymerase Chain Reaction

Preferred Specimen: Swab-Viral Culturette

Collection Instruction: Collect sample, place in holder, label with patient name. NOT TO BE USED FOR BACTERIAL TRANSPORT

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87501 (2), 87807 (1)

Scarlet Eligible: No

Cytology, Alcohol-Based FNA

5243-1

Alternate Name: FNA Cytology, Alcohol

Preferred Specimen: 50% Alcohol Jar

Collection Instruction: Body fluid for cytology in 50% alcohol (ex...urine, sputum, nipple discharge). Add directly to an equal volume of 50% Alcohol Solution.

Turn Around Time: 4 days

CPT Code: 88108 (1)

Scarlet Eligible: No

Fine Needle Aspirate

5242-3

Alternate Name: FNA Cytology, ThinPrep Pap

Preferred Specimen: ThinPrep Cytolyt Solution

Collection Instruction: Collect specimen swish in Cytolyte vial, label with pt name.

Turn Around Time: 4 days

CPT Code: 88173 (1)

Scarlet Eligible: No

Folate, RBC

0873-0

Alternate Name: FOLATE, RBC

Preferred Specimen: Lavender top- EDTA,Amber Lavender

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate, STRICT Frozen

Turn Around Time: 1 days

CPT Code: 82747 (1), 85014 (1)

Scarlet Eligible: No

Food Allerpanel (4 mix)

2189-9

Alternate Name: FOOD ALLERPANEL (4 mix)

Methodology: Immunocap

Preferred Specimen: 5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82785 (1), 86005 (4)

Scarlet Eligible: Yes

Food Panel II

2797-9

Alternate Name: FOODS PANEL ALLERGEN

Preferred Specimen: 5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82784 (1), 82785 (1), 86003 (10)

Scarlet Eligible: Yes

FOXO1 (FKHR) (13q14) GENE BY FISH

B665-1

Alternate Name: FOXO1 (FKHR) BY FISH

Methodology: FISH

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Storage Instruction: Room Temp, Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88377 (1)

Scarlet Eligible: No

Fractional Excretion Of Sodium

P791-7

Alternate Name: FRAC. EXCRETION OF NA

Preferred Specimen: Urine Cup,SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82565 (1), 82570 (1), 84295 (1), 84300 (1)

Scarlet Eligible: Yes

Fragile X

2383-8

Alternate Name: FRAGILE X

Methodology: Polymerase Chain Reaction

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 81243 (1)

Scarlet Eligible: Yes

FSH (Follicle Stimulating Hormone)

0092-7

Alternate Name: FSH

Methodology: Electrochemiluminescence Immunoassay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 83001 (1)

Scarlet Eligible: Yes

Factor XI carrier screening (Hemophilia C)

J310-4

Alternate Name: Factor XI carrier screening

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81401 (1)

Scarlet Eligible: Yes

Allergen False Ragweed (w4), IgE

2407-5

Alternate Name: False Ragweed

Methodology: Immunocap

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Familial Hypercholesterolemia LDLR

B151-2

Alternate Name: Familial Hypercholester

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81479 (1)

Scarlet Eligible: Yes

Familial Mediterranean Fever

J307-0

Alternate Name: Familial Mediterranean Fever

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81404 (1)

Scarlet Eligible: Yes

Feather Mix (ex71), IgE

1300-3

Alternate Name: Feather Mix (Ex71)

Methodology: Immunocap

Preferred Specimen: 5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86005 (1)

Scarlet Eligible: Yes

Occult Blood, Stool, Immunochem (FIT)

6182-0

Alternate Name: Fecal Occult Blood Test; FOBT; Guaiac; FIT

Methodology: Immunochemical

Preferred Specimen: FOBT Sample Bottle

Collection Instruction: Add fresh stool to container, mix gently and submit

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82274 (1)

Scarlet Eligible: No

Female STI Urine Profile

P863-4

Alternate Name: Female STI Urine Profil

Preferred Specimen: Urine Tube BD

Collection Instruction: Collect Specimen- then label with name and source. Recommended first morning urine.

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 87491 (1), 87591 (1), 87661 (1), 87798 (1)

Scarlet Eligible: Yes

Feulgen Stain - Tech Only

T983-6

Alternate Name: Feulgen Stain - TC

Methodology: Special Stain

Preferred Specimen: ThinPrep Cytolyt Solution

Collection Instruction: Collect specimen swish in Cytolyte vial, label with pt name.

Storage Instruction: Refrigerate

Turn Around Time: 4 days

CPT Code: 88313 (1)

Scarlet Eligible: No

Fibrotest/Actitest

H510-1

Alternate Name: Fibrosure, Necroin inflammatory Activity Score

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Turn Around Time: 2 days

CPT Code: 82172 (1), 82247 (1), 82977 (1), 83010 (1), 83883 (1), 84460 (1)

Scarlet Eligible: Yes

Allergen Fig (f328), IgE

1937-2

Alternate Name: Fig

Methodology: Immunocap

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Fire Ant (i70), IgE

0904-3

Alternate Name: Fire Ant

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Combined First (NT, Total BHCG, PAPP-A)

6246-3

Alternate Name: First Trimester Screen w/ NT

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: MUST BE 11-13 WEEKS AND 6 DAYS. PROVIDE GESTATIONAL INFO.

Storage Instruction: Refrigerate

Turn Around Time: 4 days

CPT Code: 81508 (1)

Scarlet Eligible: Yes

Combined First Plus (NT, Total BHCG, Inhibin-A, PAPP-A)

8622-3

Alternate Name: First Trimester Screen w/ NT, inhibin A (DIA)

Methodology: Chemiluminescence

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 4 days

CPT Code: 81509 (1)

Scarlet Eligible: Yes

Folate (Folic Acid)

0090-1

Alternate Name: Folate Acid, Folic Acid

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82746 (1)

Scarlet Eligible: Yes

Food Allergy Profile w/ rfx to Milk and Peanut

K361-6

Alternate Name: Food Allergy Profile

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82785 (1), 86003 (16)

Scarlet Eligible: Yes

Food Allergy Profile

K426-7

Alternate Name: Food Allergy Profile

Preferred Specimen: 1 mL SST Tube

Turn Around Time: 3 days

CPT Code: 82785 (1), 86003 (15)

Scarlet Eligible: Yes

Foods I Screen (fish/6)

0743-5

Alternate Name: Foods I Screen (Fish/6)

Preferred Specimen: 5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82785 (1), 86003 (6)

Scarlet Eligible: Yes

Foods II Screen (10)

0745-0

Alternate Name: Foods II Screen (10)

Preferred Specimen: 5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82785 (1), 86003 (10)

Scarlet Eligible: Yes

Allergen Formaldehyde (k80), IgE

2106-3

Alternate Name: Formaldehyde

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Catecholamines, Fractionation, Plasma

0214-7

Alternate Name: Free Catecholamine Fractionation, Plasma**Methodology:** High Pressure Liquid Chromatography**Preferred Specimen:** 4 mL Aliquot Plasma Heparinized-green top**Collection Instruction:** Patient Preparation: Patients should be relaxed in either a supine or upright position before blood is drawn. States of anxiety and stress can cause fluctuations in the catecholamine levels. Draw specimen in a pre-chilled green-top vacutainer. Patient should avoid alcohol, coffee, tea, tobacco and strenuous exercise prior to collection. Overnight fasting is preferred. Collection Instructions: After drawing specimen, place green-top tube in ice bath. Centrifuge for 20 minutes 1000XG. The plasma ****must**** be free of RBCs. Freeze on dry ice to transport.**Storage Instruction:** STRICT Frozen**Turn Around Time:** 10 days**CPT Code:** 82384 (1)**Scarlet Eligible:** No**Thyroxine, Free (FT4)**

0091-9

Alternate Name: Free T4, FT4**Methodology:** Electrochemiluminescence Immunoassay**Preferred Specimen:** 1 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 1 days**CPT Code:** 84439 (1)**Scarlet Eligible:** Yes

GGTP (Gamma-GT)

0093-5

Alternate Name: G-GTP

Methodology: Enzyme Immunoassay

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82977 (1)

Scarlet Eligible: Yes

G. vaginalis by RT-PCR

J211-4

Alternate Name: G. VAGINALIS BY PCR

Methodology: Polymerase Chain Reaction

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 5-7 days

CPT Code: 87512 (1)

Scarlet Eligible: No

G.T.T. 2 hr (includ. 30 min)

Q605-8

Alternate Name: G.T.T. 2hr (w/30min)

Preferred Specimen: Grey top- Sodium Fluoride,Urine tube without preservative

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82945 (4), 82951 (1), 82952 (1)

Scarlet Eligible: No

Galactosemia Carrier Screen

B155-3

Alternate Name: GALACTOSEMIA SCREEN

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81401 (1)

Scarlet Eligible: Yes

Gastric Biopsy (1 Specimen)

5640-8

Alternate Name: GASTRIC BIOPSY X1

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 1 days

CPT Code: 88305 (1)

Scarlet Eligible: No

Gastric Washings/Brushings

5643-2

Alternate Name: GASTRIC WASHINGS/BRUSHINGS

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Collect specimen place in formalin jar write patient name and site on jar

Turn Around Time: 1 days

CPT Code: 88108 (1), 88305 (1)

Scarlet Eligible: No

GC MultiPlex PCR - swab

J834-3

Alternate Name: GC PCR on Swab

Methodology: Multiplex PCR

Preferred Specimen: Swab-Aptima Genprobe

Collection Instruction: Collect specimen-break swab shaft into transport tube, then label with name and source

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87591 (1)

Scarlet Eligible: No

GC, Anal Aptima

A112-5

Alternate Name: GC, ANAL APTIMA

Methodology: Transcription Mediated Amplification

Preferred Specimen: Swab-Aptima Genprobe

Collection Instruction: Collect specimen-break swab shaft into transport tube, then label with name and source

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87591 (1)

Scarlet Eligible: No

GC, Oral Aptima

A111-7

Alternate Name: GC, ORAL APTIMA

Methodology: Transcription Mediated Amplification

Preferred Specimen: Swab-Aptima Genprobe

Collection Instruction: Collect specimen-break swab shaft into transport tube, then label with name and source

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87591 (1)

Scarlet Eligible: No

GC,rRNA,TMA Aptima Swab

1003-3

Alternate Name: GC,rRNA,TMA SWAB

Methodology: Transcription Mediated Amplification

Preferred Specimen: Swab-Aptima Genprobe

Collection Instruction: Collect specimen-break swab shaft into transport tube, then label with name and source

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87591 (1)

Scarlet Eligible: No

GC/CT SWAB BY PCR

L583-4

Alternate Name: GC/CT SWAB BY PCR

Turn Around Time: 3 days

CPT Code: 87491 (1), 87591 (1)

Scarlet Eligible: No

General Inhibitor Screen

5800-8

Alternate Name: GENERAL INHIBITOR SCREEN

Preferred Specimen: Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Fill tube to line at top of tube, invert 8-10 times (DO NOT SHAKE!)

Turn Around Time: 4 days

CPT Code: 85610 (1), 85611 (8), 85670 (1), 85730 (2), 85732 (8)

Scarlet Eligible: No

General Thrombophilia Assessment

6800-7

Alternate Name: GENERAL THROMBOPHILIA ASSESSMENT

Preferred Specimen: Aliquot Plasma-Citrated-Blue Top, Aliquot Serum- Red Top, Lavender top-EDTA

Collection Instruction: Lavender must be sent separately at room temp.

Turn Around Time: 5 days

CPT Code: 81240 (1), 81241 (1), 83090 (1), 85025 (1), 85240 (1), 85250 (1), 85270 (1), 85280 (1), 85300 (1), 85303 (1), 85306 (2), 85307 (2), 85379 (1), 85384 (1), 85420 (1), 85598 (1), 85610 (1), 85613 (2), 85670 (1), 85730 (2), 86146 (3), 86147 (3)

Scarlet Eligible: No

Culture, Genital

0299-8

Alternate Name: GENITAL CULTURE

Methodology: Bacterial Culture

Preferred Specimen: Swab-E

Collection Instruction: Collect specimen with swab then place swab into transport carrier; Label with name and source

Turn Around Time: 4 days

CPT Code: 87070 (1)

Scarlet Eligible: No

Gentamicin, Peak

0730-2

Alternate Name: GENTAMICIN, PEAK

Methodology: Enzyme Immunoassay

Preferred Specimen: Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 80170 (1)

Scarlet Eligible: No

Gentamicin, Peak + Trough

9081-1

Alternate Name: GENTAMICIN, PEAK + TROUGH
Methodology: Fluor. Polarization Immunoassay
Preferred Specimen: 2 mL Red Top
Collection Instruction: Fill tubes, label with patient name.
Storage Instruction: Refrigerate
Turn Around Time: 3 days
CPT Code: 80170 (2)
Scarlet Eligible: No

Gentamicin (Random)

A616-5

Alternate Name: GENTAMICIN, RANDOM
Methodology: Fluor. Polarization Immunoassay
Preferred Specimen: Red Top
Collection Instruction: Fill tube, label with patient name
Storage Instruction: Refrigerate
Turn Around Time: 2 days
CPT Code: 80170 (1)
Scarlet Eligible: Yes

GENTAMICIN, TROUGH

1930-7

Alternate Name: GENTAMICIN, TROUGH
Methodology: Fluor. Polarization Immunoassay
Preferred Specimen: Red Top
Collection Instruction: Fill tube, label with patient name
Storage Instruction: Refrigerate
Turn Around Time: 2 days
CPT Code: 80170 (1)
Scarlet Eligible: Yes

Germ Cell Tumor IHC Panel

5194-6

Alternate Name: GERM CELL TUMOR, IHC

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88341 (9), 88342 (1)

Scarlet Eligible: No

Giardia Antigen EIA (2ND SPECIMEN)

3898-4

Alternate Name: GIARDIA ANTIGEN (2ND)

Methodology: Enzyme Immunoassay

Preferred Specimen: Ova & Parasite Mini

Collection Instruction: With applicator stick, place 10 gms of stool into each vial, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 87329 (1)

Scarlet Eligible: No

Giardia Antigen EIA (3RD SPECIMEN)

3899-2

Alternate Name: GIARDIA ANTIGEN (3RD)

Methodology: Enzyme Immunoassay

Preferred Specimen: Ova & Parasite Mini

Collection Instruction: With applicator stick, place 10 gms of stool into each vial, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 87329 (1)

Scarlet Eligible: No

GIST Profile (IHC)

5175-5

Alternate Name: GIST PROFILE

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block or slide form with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88341 (5), 88342 (1)

Scarlet Eligible: No

GIST Profile by IHC TC

J687

Alternate Name: GIST Profile by IHC TC

Scarlet Eligible: No

Gleason Score 3+4, Reflex To PTEN FISH

J143-9

Alternate Name: GLEASON SCORE REFLEX

Methodology: FISH

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Collect specimen place in formalin jar write patient name and site on jar

Turn Around Time: 3-5 days

CPT Code: 99999 (1)

Scarlet Eligible: No

Gliadin IGA, ABS.

5868-5

Alternate Name: GLIADIN(DP) IgA,Ab

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86258 (1)

Scarlet Eligible: Yes

Gliadin IgG, ABS

3822-4

Alternate Name: GLIADIN(DP) IgG,Ab

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86258 (1)

Scarlet Eligible: Yes

GLU.TOL.3HR.(Pregnancy)100GM

1015-7

Alternate Name: GLU.TOL.3HR.(PREGNANCY)

Preferred Specimen: Grey top- Sodium Fluoride

Collection Instruction: Label ALL specimens with time of collection (4 grey tops)

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82951 (1), 82952 (1)

Scarlet Eligible: No

Glucose, Post Dialysis

3110-4

Alternate Name: GLUCOSE POST DIALYSIS

Methodology: Enzyme Immunoassay

Preferred Specimen: Microtainer - Pediatric SST

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82947 (1)

Scarlet Eligible: No

Glucose Tolerance 2 Hour

8979-7

Alternate Name: GLUCOSE TOLERANCE 2 HOUR

Preferred Specimen: Grey top- Sodium Fluoride,Urine Urinalysis Tube - Yellow

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82945 (3), 82951 (1)

Scarlet Eligible: No

Glucose Tolerance 3 HRS.

0097-6

Alternate Name: GLUCOSE TOLERANCE 3 hrs

Preferred Specimen: Grey top- Sodium Fluoride,Urine Urinalysis Tube - Yellow

Collection Instruction: Label ALL specimens with time of collection (4 grey tops, 4 urine tubes)

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82945 (4), 82951 (1), 82952 (1)

Scarlet Eligible: No

Glucose Tolerance 4 HRS.

0098-4

Alternate Name: GLUCOSE TOLERANCE 4

Preferred Specimen: Grey top- Sodium Fluoride,Urine Urinalysis Tube - Yellow

Collection Instruction: Label ALL specimens with time of collection (6 grey tops, 6 urine tubes)

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82945 (6), 82951 (1), 82952 (3)

Scarlet Eligible: No

Glucose Tolerance 5 HRS.

0099-2

Alternate Name: GLUCOSE TOLERANCE 5

Preferred Specimen: Grey top- Sodium Fluoride,Urine tube without preservative

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82945 (7), 82951 (1), 82952 (4)

Scarlet Eligible: No

Glucose, Urine, Random

5042-7

Alternate Name: GLUCOSE, URINE, RANDOM

Methodology: Enzyme Immunoassay

Preferred Specimen: 5 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82945 (1)

Scarlet Eligible: Yes

Glucose, 2HR. PP W/Glucola

3157-5

Alternate Name: GLUCOSE,2HR PP (75g)

Methodology: Enzyme Immunoassay

Preferred Specimen: Grey top- Sodium Fluoride

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82950 (1)

Scarlet Eligible: No

Glucose, Urine 24HR Quantitative

0096-8

Alternate Name: GLUCOSE,URINE 24 HR

Methodology: Enzyme Immunoassay

Preferred Specimen: 5 mL Urine Container - 24hr

Collection Instruction: Discard first morning sample, then collect all urines after for 24 hrs, on completion label with collection time and patient's name; record total volume on Requisition

Storage Instruction: Refrigerate

CPT Code: 81050 (1), 82945 (1)

Scarlet Eligible: No

Glycogen Storage, Type IA (GSDIA)

A651-2

Alternate Name: GLYCOGEN STORAGE TYP 1A

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81250 (1)

Scarlet Eligible: Yes

Gonadotropin I

2318-4

Alternate Name: GONADOTROPIN I

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 83001 (1), 83002 (1), 84403 (1)

Scarlet Eligible: Yes

GP/Bile Duct Brushing

5119-3

Alternate Name: GP/BILE DUCT BRUSHING

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Collect specimen place in formalin jar write patient name and site on jar

Turn Around Time: 2 days

CPT Code: 88104 (1)

Scarlet Eligible: No

Bone Marrow Smear Interpretation

5285-2

Alternate Name: GP/BM SMEAR INTERPRETATION

Methodology: Histology

Preferred Specimen: Stained Slides

Collection Instruction: Collect specimen and place on charged slide, label with pt name, place in slide holder for transport.

Turn Around Time: 2 days

CPT Code: 85097 (1)

Scarlet Eligible: No

ISH Single Antibody, W/Interp.

5286-0

Alternate Name: GP/ISH SINGLE AB.,W/INTERP.

Methodology: In Situ Hybridization

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88365 (1)

Scarlet Eligible: No

GP/Level III Histo (2 Slides)

5159-9

Alternate Name: GP/LEVEL III HISTO TC

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 1 days

CPT Code: 88304 (2)

Scarlet Eligible: No

GP/Level III Histo. TC (1 Slide)

5158-1

Alternate Name: GP/LEVEL III HISTO,TC

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 1 days

CPT Code: 88304 (1)

Scarlet Eligible: No

Gram Stain

0103-2

Alternate Name: GRAM STAIN

Methodology: Microscopic Examination

Preferred Specimen: Slides

Collection Instruction: Label slides with patient name, place in slide holder for transport

Turn Around Time: 1 days

CPT Code: 87205 (1)

Scarlet Eligible: No

Group A Strep (Non NY)

J965-5

Alternate Name: GROUP A STREP

Methodology: Real-Time Polymerase Chain Reaction

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87651 (1)

Scarlet Eligible: No

Group B Strep (Non NY)

J966-3

Alternate Name: GROUP B STREP

Methodology: Real-Time Polymerase Chain Reaction

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87653 (1)

Scarlet Eligible: No

Growth Hormone

0510-8

Alternate Name: GROWTH HORMONE

Methodology: Chemiluminescence

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 83003 (1)

Scarlet Eligible: Yes

GTT,GEST.,2HR.(ADA) PANEL

P782-6

Alternate Name: GTT,GEST.,2HR.(ADA) PANEL

Preferred Specimen: Grey top- Sodium Fluoride

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82951 (1)

Scarlet Eligible: No

GTT,GEST,3HR.(ACOG,CARP/COUS)

6794-2

Alternate Name: GTT/ACOG,3hr

Preferred Specimen: Grey top- Sodium Fluoride

Collection Instruction: Submit 4 grey top tubes, each labeled with time collected.

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82951 (1), 82952 (1)

Scarlet Eligible: No

Allergen Garlic (f47), IgE

0965-4

Alternate Name: Garlic

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Gastric Distress

F101-1

Alternate Name: Gastric Distress

Preferred Specimen: SST Tube

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 83516 (4), 86003 (1), 86255 (1), 86677 (3)

Scarlet Eligible: Yes

Gastrin, Serum

0508-2

Alternate Name: Gastrin, Serum

Methodology: Immunoassay

Preferred Specimen: 1 mL Aliquot Tube-Serum

Collection Instruction: Place 1-3 mL of serum in transport tube. Label as serum and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 11 days

CPT Code: 82941 (1)

Scarlet Eligible: No

Gaucher Disease Carrier Screening

A655-3

Alternate Name: Gaucher Carrier Screen

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81251 (1)

Scarlet Eligible: Yes

Solid Tumor Gene Fusion Assay

J355-9

Alternate Name: Gene Fusion

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 10-15 unstained slides at 5 microns (tumor content>10%) with H&E, Shavings (if tumor content>40%) with H&E. If sending DNA- Please call Client Services.

Storage Instruction: Room Temp

Turn Around Time: 10 days

CPT Code: 81445 (1)

Scarlet Eligible: No

General Health Panel9149-6

Alternate Name: General Health Panel**Preferred Specimen:** SST Tube,Lavender top- EDTA**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 1 days**CPT Code:** 80053 (1), 84443 (1), 85025 (1)**Scarlet Eligible:** Yes

Chlamydia Swab by PCRT832-5

Alternate Name: Genital CT Swab PCR**Methodology:** Real-Time Polymerase Chain Reaction**Preferred Specimen:** Roche Collection Swab**Collection Instruction:** CLEAN: Using the woven swab (Swab A), remove excess mucus from the cervical os and surrounding mucosa. Discard the swab after cleaning. NOTE: Cleaning excess mucus from the cervical os is required to ensure an adequate sample is obtained for processing. COLLECT: To collect the specimen, hold flocked swab (Swab B) with the scoreline above your hand and insert into the endocervical canal. Gently rotate the swab 5 times in one direction in the endocervical canal. Do not over-rotate. Carefully withdraw the swab, avoiding any contact with the vaginal mucosa. ALIGN: Remove the cap from the cobasÂ® PCR Media Tube and lower the swab specimen into the tube until the visible scoreline on the swab shaft is aligned with the tube rim. The bud of the swab should not be submerged into the liquid prior to breaking the shaft.**Storage Instruction:** Room Temp**Turn Around Time:** 3 days**CPT Code:** 87491 (1)**Scarlet Eligible:** No

Genital Cult. R/O Yeast/Gard.

1295-5

Alternate Name: Genital Culture Rule Out Yeast/Gardnerella**Methodology:** Bacterial Culture**Preferred Specimen:** Swab-E**Collection Instruction:** Collect specimen with swab then place swab into transport carrier; Label with name and source**Turn Around Time:** 3 days**CPT Code:** 87070 (1)**Scarlet Eligible:** No**Gonorrhea Swab by PCR**

T833-3

Alternate Name: Genital GC Swab PCR**Methodology:** Real-Time Polymerase Chain Reaction**Preferred Specimen:** Roche Collection Swab

Collection Instruction: CLEAN: Using the woven swab (Swab A), remove excess mucus from the cervical os and surrounding mucosa. Discard the swab after cleaning. NOTE: Cleaning excess mucus from the cervical os is required to ensure an adequate sample is obtained for processing. COLLECT: To collect the specimen, hold flocked swab (Swab B) with the scoreline above your hand and insert into the endocervical canal. Gently rotate the swab 5 times in one direction in the endocervical canal. Do not over-rotate. Carefully withdraw the swab, avoiding any contact with the vaginal mucosa. ALIGN: Remove the cap from the cobas® PCR Media Tube and lower the swab specimen into the tube until the visible scoreline on the swab shaft is aligned with the tube rim. The bud of the swab should not be submerged into the liquid prior to breaking the shaft. BREAK: Carefully leverage the swab against the tube rim to break the swab shaft at the scoreline. CLOSE: Tightly re-cap the cobas® PCR Media Tube. The specimen is now ready for transport. Discard the top portion of the swab.

Storage Instruction: Room Temp**Turn Around Time:** 3 days**CPT Code:** 87591 (1)**Scarlet Eligible:** No

Culture, Genital (Beta strep. Group B)

3400-9

Alternate Name: Genital, Beta Strep, Group B Culture

Methodology: Bacterial Culture

Preferred Specimen: Swab-E

Collection Instruction: Collect specimen with swab then place swab into transport carrier; Label with name and source

Turn Around Time: 4 days

CPT Code: 87070 (1)

Scarlet Eligible: No

Allergen Giant Ragweed (w3), IgE

1197-3

Alternate Name: Giant Ragweed (W3)

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Giardia Antigen (EIA)

1396-1

Alternate Name: Giardia Antigen (EIA)

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: Ova & Parasite Mini

Collection Instruction: With applicator stick, place 10 gms of stool into each vial, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 87329 (1)

Scarlet Eligible: No

Allergen Ginger (rf270), IgE

2346-5

Alternate Name: Ginger

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Glioma: IDH1/IDH2 by NGS

J635-4

Alternate Name: Glioma: IDH1/IDH2 by NGS

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 10-15 unstained slides at 5 microns(tumor content>10%) with H&E, Shavings(if tumor content>40% with H&E, if sending DNA-Please call CS.

Storage Instruction: Room Temp

Turn Around Time: 5-10 days

CPT Code: 81120 (1), 81121 (1)

Scarlet Eligible: No

Glucose, 1hr, PP (Preg) w/Glucola

0368-1

Alternate Name: Glu.1hr(glucola)preg.

Methodology: Enzyme Immunoassay

Preferred Specimen: Grey top- Sodium Fluoride

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82950 (1)

Scarlet Eligible: Yes

Glucose

1976-0

Alternate Name: Glucose

Methodology: Enzymatic Reference method with Hexokinase

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82947 (1)

Scarlet Eligible: Yes

Glucose 2hr PP

0369-9

Alternate Name: Glucose 2hr PP

Methodology: Enzyme Immunoassay

Preferred Specimen: Grey top- Sodium Fluoride

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82947 (1)

Scarlet Eligible: No

Glucose, 1hr PP

0171-9

Alternate Name: Glucose, 1hr PP

Methodology: Enzyme Immunoassay

Preferred Specimen: Grey top- Sodium Fluoride

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82947 (1)

Scarlet Eligible: No

Allergen Gluten (f79), IgE

0935-7

Alternate Name: Gluten

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Fructosamine, Serum

0485-3

Alternate Name: Glycated Albumin;Glycated Proteins; Ketoamines, Plasma Protein

Methodology: Colorimetry

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82985 (1)

Scarlet Eligible: Yes

Glycemic Control Profile

H959-0

Alternate Name: Glycemic Control

Turn Around Time: 1 days

CPT Code: 82947 (1), 83036 (1), 84378 (1)

Scarlet Eligible: Yes

Hemoglobin A1c (Glycohgb)

0102-4

Alternate Name: Glycohgb, Glycated Hemoglobin

Methodology: Turbidimetric Inhibition Immunoassay

Preferred Specimen: 0.5 mL Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 83036 (1)

Scarlet Eligible: Yes

Allergen Glycyphagus Domesticus (d73), IgE

3813-3

Alternate Name: Glycyphagus Domest. IgE

Methodology: Enzyme Immunoassay

Preferred Specimen: 0.5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 6 days

CPT Code: 86003 (1)

Scarlet Eligible: No

Allergen Goat Milk (rf300), IgE

2263-2

Alternate Name: Goat Milk

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Golden Rod (w12), IgE

1364-9

Alternate Name: Golden Rod

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Culture, Neisseria gonorrhoeae (GC), Anal

1353-2

Alternate Name: Gonorrhea , Anal, Culture

Methodology: Bacterial Culture

Preferred Specimen: Swab-Amies with Charcoal-Black

Collection Instruction: Collect specimen with swab then place swab into transport carrier; Label with patient's name and source

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87081 (1)

Scarlet Eligible: No

Culture, Neisseria gonorrhoeae (GC), Throat

1352-4

Alternate Name: Gonorrhea , Throat, Culture

Methodology: Bacterial Culture

Preferred Specimen: Swab-E

Collection Instruction: Collect specimen with swab then place swab into transport carrier; Label with patient's name and source

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87081 (1)

Scarlet Eligible: No

Culture, Neisseria gonorrhoeae (GC)

0075-2

Alternate Name: Gonorrhea Culture

Methodology: Bacterial Culture

Preferred Specimen: Swab-E

Collection Instruction: Collect sample with bacterial swab, smear on plate, place plate in plastic bag with the pill provided, label with patient's name;

Turn Around Time: 4 days

CPT Code: 87081 (1)

Scarlet Eligible: No

Gonorrhea NAAT

3853-9

Alternate Name: Gonorrhea NAAT

Methodology: Nucleic Acid Amplification Test

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin Prep or SurePath vial, label with patients name.

Turn Around Time: 2 days

CPT Code: 87591 (1)

Scarlet Eligible: No

GC Smear

0094-3

Alternate Name: Gonorrhea smear

Methodology: Microscopic Examination

Preferred Specimen: Slides

Collection Instruction: Label slides with patient name, place in slide holder for transport

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 87205 (1)

Scarlet Eligible: No

Allergen Goose Feathers (e70), IgE

0847-4

Alternate Name: Goose Feathers

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Grape (f259), IgE

1225-2

Alternate Name: Grape (F259)

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Grapefruit (f209), IgE

1349-0

Alternate Name: Grapefruit

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Grass Mix (gx1), IgE

1841-6

Alternate Name: Grass Mix

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86005 (1)

Scarlet Eligible: Yes

Grasses Screen (6)

0741-9

Alternate Name: Grasses Screen (6)

Preferred Specimen: 5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82785 (1), 86003 (6)

Scarlet Eligible: Yes

Allergen Green (Bell) Pepper (rf263), IgE

1243-5

Alternate Name: Green (Bell) Pepper (RF263)

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen String Bean (rf315), IgE

1230-2

Alternate Name: Green Bean

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Green Pea (f12), IgE

1228-6

Alternate Name: Green Pea

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Grey Alder (t2), IgE

0561-1

Alternate Name: Grey Alder

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Culture, Genital - Beta strep. Group B (Penicillin Allergic)

A265-1

Alternate Name: Group B Strep Culture (Penicillin Allergic)

Methodology: Bacterial Culture

Preferred Specimen: Swab-E

Collection Instruction: Collect specimen with swab then place swab into transport carrier; Label with name and source

Turn Around Time: 4 days

CPT Code: 87070 (1)

Scarlet Eligible: No

Allergen Guinea Pig (e6), IgE

3566-7

Alternate Name: Guinea Pig

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

H. Pylori BY IHC Tech Only

J159-5

Alternate Name: H. PYLORI BY IHC - TC

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Storage Instruction: Room Temp

Turn Around Time: 1 days

CPT Code: 88342 (1)

Scarlet Eligible: No

H. Pylori Antibody Panel (G,A,M)

0953-0

Alternate Name: H. Pylori Antibody Profile, IgG, IgA, IgM

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86677 (3)

Scarlet Eligible: Yes

H. Pylori Antigen, Stool

2441-4

Alternate Name: H. Pylori Antigen,Stool

Methodology: Chemiluminescence

Preferred Specimen: 10 mL Stool Cont.-Cup

Collection Instruction: Using applicator stick, place 20-30 gms of stool into cup, label with patient's name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 87338 (1)

Scarlet Eligible: No

H. PYLORI IgM

7736-2

Alternate Name: H. pylori IgM

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86677 (1)

Scarlet Eligible: Yes

H. Pylori,IGA (For Panel #0953)

1766-5

Alternate Name: H. pylori, IgA

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86677 (1)

Scarlet Eligible: Yes

H. Pylori,IGG (For Panel #0953)

1765-7

Alternate Name: H. pylori, IgG

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86677 (1)

Scarlet Eligible: Yes

H.Pylori (Stain Ident.)

5645-7

Alternate Name: H.PYLORI (STAIN IDENT.)

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Collect specimen place in formalin jar write patient name and site on jar

CPT Code: 88312 (1)

Scarlet Eligible: No

H. Pylori BY IHC W/ Interp

J158-7

Alternate Name: H.PYLORI BY IHC W/INTER

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Storage Instruction: Room Temp

Turn Around Time: 1 days

CPT Code: 88342 (1)

Scarlet Eligible: No

H. Pylori by IHC. If+ rflx to H. pylori PCR with rflx to H. pylori Clarithromycin resistance detecti

L282-3

Alternate Name: H.pylori IHC+rflx Hp res

Turn Around Time: 1 days

CPT Code: 88342 (1)

Scarlet Eligible: No

Haptoglobin, Serum

0514-0

Alternate Name: HAPTOGLOBIN

Methodology: Immunoturbidimetric

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 83010 (1)

Scarlet Eligible: Yes

Hepatitis B Core Antibody, Total

0108-1

Alternate Name: HBcAB

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86704 (1)

Scarlet Eligible: Yes

Hepatitis B Core Antibody, IgM

0206-3

Alternate Name: HBcAB IgM

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86705 (1)

Scarlet Eligible: Yes

Hepatitis B Surface Antibody (Qual/Immunity)

0107-3

Alternate Name: HBsAB;Hepatitis B Surface Antibody

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86706 (1)

Scarlet Eligible: Yes

Hepatitis B Surface Antigen

0106-5

Alternate Name: HBsAG

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 87340 (1)

Scarlet Eligible: Yes

HCG Quant.,Progesterone

2379-6

Alternate Name: HCG QUANT.,PROGESTER

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84144 (1), 84702 (1)

Scarlet Eligible: Yes

Hematocrit

0019-0

Alternate Name: HCT

Methodology: Cumulative Pulse Height Detection

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 85014 (1)

Scarlet Eligible: Yes

HCV Viral Load reflex to NS3/NS5 (Non-NY)

K943-1

Alternate Name: HCV VL rfx NS3/NS5

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: 2 Aliquot Plasma-EDTA-Lavender Top

Collection Instruction: Centrifuge lavender-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA and patient's name.

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87522 (1)

Scarlet Eligible: Yes

Hepatitis C RNA, Quant., PCR

3376-1

Alternate Name: HCV, HEP C VIRAL LOAD

Methodology: Polymerase Chain Reaction

Preferred Specimen: 3 mL Aliquot Plasma-EDTA-Lavender Top

Collection Instruction: Seperate plasma within 4 hrs of collection. Centrifuge lavender-top at high speed for 15 minutes, transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA and patient's name

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87522 (1)

Scarlet Eligible: Yes

HDL Cholesterol

0059-6

Alternate Name: HDL Cholesterol

Methodology: Enzyme Immunoassay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 83718 (1)

Scarlet Eligible: Yes

Human Epididymal Protein 4

J032-4

Alternate Name: HE4 OVARIAN CANCER

Methodology: Electrochemiluminescence Immunoassay

Preferred Specimen: Aliquot Tube-Serum

Collection Instruction: Place 1-3 mL of serum in transport tube. Label as serum;refrigerate or freeze as required.

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86305 (1)

Scarlet Eligible: No

Heart Health Monitoring

Q472-3

Alternate Name: HEART HEALTH MNTR(ACT)

Turn Around Time: 2 days

CPT Code: 82465 (1), 83036 (1), 83704 (1), 83718 (1), 83721 (1), 84478 (1), 86141 (1)

Scarlet Eligible: Yes

Heart Health Baseline

Q471-5

Alternate Name: HEARTHLTH BASELINE(ACT)

Collection Instruction: Fill tube, allow to clot then centrifuge 8-10 min

Turn Around Time: 2 days

CPT Code: 82465 (1), 83036 (1), 83695 (1), 83704 (1), 83718 (1), 83721 (1), 84478 (1), 86141 (1)

Scarlet Eligible: Yes

Hemoglobin + Hematocrit

0104-0

Alternate Name: HEMOGLOBIN + HEMATOCRIT

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 85014 (1), 85018 (1)

Scarlet Eligible: Yes

Hepatitis A Antibody, IgM

0538-9

Alternate Name: HEP. A Ab., IgM

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86709 (1)

Scarlet Eligible: Yes

Hepatitis A Ab.(non-reflex)

2126-1

Alternate Name: HEP. A Ab., NO REFLEX

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86708 (1)

Scarlet Eligible: Yes

Hepatitis A Ab.,Total W/Reflex

0105-7

Alternate Name: HEP. A Ab., TOTAL

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86708 (1)

Scarlet Eligible: Yes

Hepatitis B Core W/Reflex To IgM

3138-5

Alternate Name: HEP. B CORE Ab w/REFLEX

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86704 (1)

Scarlet Eligible: Yes

Hepatitis B Surface Antigen, Quant.

1608-9

Alternate Name: HEP. B SURF Ag., QUANT.

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 87340 (1)

Scarlet Eligible: Yes

Hepatitis B Surface Antibody, Quant.

1606-3

Alternate Name: HEP. B SURF. Ab., QUANT

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86706 (1)

Scarlet Eligible: Yes

Hepatitis Be Antigen

0539-7

Alternate Name: HEP. Be Ag.

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 87350 (1)

Scarlet Eligible: Yes

Hepatitis C Genotype

2161-8

Alternate Name: HEP. C GENOTYPE

Methodology: Line Probe Assay utilizing reverse-hybridization

Preferred Specimen: 2 mL Aliquot Plasma-EDTA-Lavender Top

Collection Instruction: Centrifuge lavender-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA and patient's name.

Storage Instruction: Refrigerate

Turn Around Time: 5-7 days

CPT Code: 87902 (1)

Scarlet Eligible: Yes

Heparin Level (ANTI-XA) - LMWH

5761-2

Alternate Name: HEPARIN LEVEL

Methodology: Chromogenic

Preferred Specimen: 5 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 4 days

CPT Code: 85130 (1)

Scarlet Eligible: No

Hepatic Function Panel

3422-3

Alternate Name: HEPATIC FUNCTION

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80076 (1)

Scarlet Eligible: Yes

Hepatitis B Recovery

2075-0

Alternate Name: HEPATITIS B RECOVERY

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86706 (1), 87340 (1)

Scarlet Eligible: Yes

Hepatitis/ Comprehensive

0992-8

Alternate Name: HEPATITIS/COMPREHENS

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86704 (1), 86706 (1), 86707 (1), 86708 (1), 87340 (1), 87350 (1)

Scarlet Eligible: Yes

Hepatocellular vs. Cholangiocarcinoma vs. Metastasis IHC Pnl

5193-8

Alternate Name: HEPATOCELVCHOLANGVMETAS/IHC

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88341 (10), 88342 (1)

Scarlet Eligible: No

HER2 by IHC (Tech Only)w/digital imaging

5405-6

Alternate Name: HER-2 NEU, IHC TC Only

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88361 (1)

Scarlet Eligible: No

HER2 by IHC - reflex +2 to HER2 by FISH

5428-8

Alternate Name: HER-2 NEU, IHC w/reflex to FISH

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Turn Around Time: 3-5 days

CPT Code: 88361 (1)

Scarlet Eligible: No

HER-2, IHC

A301-4

Alternate Name: HER-2, IHC MANUAL

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88360 (1)

Scarlet Eligible: No

HER2 by IHC w/digital imaging

5171-4

Alternate Name: HER-2/neu, IHC

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88361 (1)

Scarlet Eligible: No

HER-2/neu,FISH (Manual)

5247-2

Alternate Name: HER-2/neu,FISH (Manual)

Methodology: FISH

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: Please write surgical number on block and on requisition. Provide fixation time on requisition as well.

Turn Around Time: 3-5 days

CPT Code: 88377 (1)

Scarlet Eligible: No

HER2 by FISH with Interpretation

5262-1

Alternate Name: HER-2/neu,FISH interp.**Methodology:** FISH**Preferred Specimen:** Formalin-fixed, Paraffin-embedded Tissue**Collection Instruction:** Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube. Also acceptable:Whole blood in green (sodium heparin) tube.**Storage Instruction:** Room Temp**Turn Around Time:** 3-5 days**CPT Code:** 88374 (1)**Scarlet Eligible:** No

HER2 FISH Global Gastric+GEJ

A427-7

Alternate Name: HER2 FISH-Gastric+GEJ**Methodology:** FISH**Preferred Specimen:** Formalin-fixed, Paraffin-embedded Tissue**Collection Instruction:** Please write surgical number on block and on requisition. Provide fixation time on requisition as well**Storage Instruction:** Room Temp**Turn Around Time:** 3-5 days**CPT Code:** 88374 (1)**Scarlet Eligible:** No

HER2 FISH-Reflex Equivocal to HER2 IHC

B812-9

Alternate Name: HER2 by FISH, if borderline reflex to HER2 by IHC**Methodology:** FISH**Preferred Specimen:** Formalin-fixed, Paraffin-embedded Tissue**Collection Instruction:** This comes in block form from client with surgical number imprint.**Storage Instruction:** Room Temp**Turn Around Time:** 3-5 days**CPT Code:** 88374 (1)**Scarlet Eligible:** No

Automated HER2 FISH TC Only

5259-7

Alternate Name: HER2 by FISH: TC

Methodology: FISH

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: Please write surgical number on block and requisition. Provide fixation time on requisition as well. Please include specimen collection date on requisition.

Turn Around Time: 3-5 days

CPT Code: 88374 (1)

Scarlet Eligible: No

HER2 BY IHC GI/ESO Carcinoma

A972-2

Alternate Name: HER2 by IHC in Gastric/Gastroesophageal Carcinoma

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88342 (1)

Scarlet Eligible: No

HER2 IHC Gastric/GEJ, TC Only w/digital imaging

B269-2

Alternate Name: HER2-IHC Gastric/GEJ,TC

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Storage Instruction: Room Temp

Turn Around Time: 1 days

CPT Code: 88361 (1)

Scarlet Eligible: No

Hereditary Hemochromatosis DNA Mutation Analysis (NY only)

2696-3

Alternate Name: HFE-C282Y-H63D-S65C

Methodology: Pyrosequencing

Preferred Specimen: 5 mL Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 3-4 days

CPT Code: 81256 (1)

Scarlet Eligible: No

Hereditary Hemochromatosis DNA Mutation Analysis (Non-NY)

3420-7

Alternate Name: HFE-C282Y-H63D-S65C

Methodology: Pyrosequencing

Preferred Specimen: 5 mL Peripheral Blood - Lavender Top

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient's name

Storage Instruction: Refrigerate

Turn Around Time: 3-4 days

CPT Code: 81256 (1)

Scarlet Eligible: Yes

Histo Chemical Micro Organism TCx2

B671-9

Alternate Name: HISTO CHEM MICRO ORGAN

Methodology: Histology

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88312 (2)

Scarlet Eligible: No

Histo Chemical Micro Organism x3 W/IN

B673-5

Alternate Name: HISTO CHEM MICRO ORGAN

Methodology: Histology

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88312 (3)

Scarlet Eligible: No

Histo Chemical Micro Organism TC

B672-7

Alternate Name: HISTO CHEM MICRO ORGAN

Methodology: Histology

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88312 (1)

Scarlet Eligible: No

Histo Chemical Micro Organism x2 W/IN

B674-3

Alternate Name: HISTO CHEM MICRO ORGAN

Methodology: Histology

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88312 (2)

Scarlet Eligible: No

Histo Chemical Micro Organism W/INTERP

B675-0

Alternate Name: HISTO CHEM MICRO ORGAN

Methodology: Histology

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88312 (1)

Scarlet Eligible: No

Histo Chemical Micro Organism TCx3

B670-1

Alternate Name: HISTO CHEM MICRO ORGAN

Methodology: Histology

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88312 (3)

Scarlet Eligible: No

Histology Processing

5416-3

Alternate Name: HISTOLOGY PROCESSING

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

CPT Code: 88305 (1)

Scarlet Eligible: No

HIT Profile-Level 1

5871-9

Alternate Name: HIT PROFILE-LEVEL 1

Preferred Specimen: Lavender top- EDTA,Red Top

Turn Around Time: 3 days

CPT Code: 85025 (1), 86022 (1)

Scarlet Eligible: Yes

HIV AG/AB 4th Generation

B688-3

Alternate Name: HIV P24

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 87389 (1)

Scarlet Eligible: Yes

HIV-1,RNA Viral Reflex Genotype

6364-4

Alternate Name: HIV-1,RNA VIRAL RFLX GENOTYPE

Preferred Specimen: 2 mL Aliquot Plasma-EDTA-Lavender Top

Collection Instruction: Centrifuge lavender-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA and patient's name.

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87536 (1)

Scarlet Eligible: Yes

HIV-1, RNA, Ultra/PCR, Viral Load

1010-8

Alternate Name: HIV-1,RNA,PCR,ULTRA

Methodology: Polymerase Chain Reaction

Preferred Specimen: Aliquot Plasma-EDTA-Lavender Top

Collection Instruction: Draw 2 lavender tubes, they must be centrifuged, poured off into transfer tube and label the pour off PLASMA EDTA with patient name. Can be submitted refrigerated or frozen.

Turn Around Time: 2-3 days

CPT Code: 87536 (1)

Scarlet Eligible: Yes

HLA Typing for Celiac Disease

J097-7

Alternate Name: HLA CELIAC GENETICS
Methodology: Polymerase Chain Reaction
Preferred Specimen: 3 mL Lavender top- EDTA
Storage Instruction: Refrigerate
Turn Around Time: 12 days
CPT Code: 81376 (2), 81382 (2)
Scarlet Eligible: No

HOLD BLOCK FOR ADDT'L TESTING

B700-6

Alternate Name: HOLD BLOCK
Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue
Collection Instruction: This comes in block form from client with surgical number imprint.
Storage Instruction: Room Temp
Turn Around Time: 8 days
CPT Code: 99999 (1)
Scarlet Eligible: No

Hold Cells For Further Testing

B232-0

Alternate Name: HOLD CELLS FURTHER TEST
Preferred Specimen: Green top- Sodium Heparin
CPT Code: 99999 (1)
Scarlet Eligible: No

HPV DNA 16, 18, NON 16/18

A531-6

Alternate Name: HPV DNA 16, 18, NON 16/18
Methodology: Real-Time Polymerase Chain Reaction
Preferred Specimen: ThinPrep Vial
Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.
Turn Around Time: 3-5 days
CPT Code: 87624 (1)
Scarlet Eligible: No

HPV High Risk by ISH

8714-8

Alternate Name: HPV HIGH by ISH

Methodology: In Situ Hybridization

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88365 (1)

Scarlet Eligible: No

HPV High by ISH Stain ONLY

8689-2

Alternate Name: HPV HIGH by ISH STAIN ONLY

Methodology: In Situ Hybridization

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88365 (1)

Scarlet Eligible: No

HPV High/Low Risk by ISH

3412-4

Alternate Name: HPV HIGH/LOW BY ISH

Methodology: Histology

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88364 (1), 88365 (1)

Scarlet Eligible: No

HPV High/Low Risk by ISH (Stain Only)

8696-7

Alternate Name: HPV HIGH/LOW by ISH STAIN ONLY

Methodology: In Situ Hybridization

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88364 (1), 88365 (1)

Scarlet Eligible: No

HPV PRIMARY SCREENING >=25 YRS, WHEN + RFX TO PAP

H852-7

Alternate Name: HPV PRIMARY SCREENING >=25 YRS, WHEN + RFX TO PAP

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87624 (1)

Scarlet Eligible: No

HPV mRNA screen, if pos, Reflex to 16, 18/45

B342-7

Alternate Name: HPV mRNA screen, if pos, Reflex to 16, 18/45

Methodology: Transcription Mediated Amplification

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 3-5 days

CPT Code: 87624 (1)

Scarlet Eligible: No

HPV DNA, High Risk Anal - TP

J025-8

Alternate Name: HPV, HIGH RISK ANAL

Methodology: Hybrid Capture

Preferred Specimen: Anal ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 3-5 days

CPT Code: 87624 (1)

Scarlet Eligible: No

HPV, High Risk

3910

Alternate Name: HPV, High Risk

Scarlet Eligible: No

HSV 1+2 Reflex to VZV

L009-0

Alternate Name: HSV

Methodology: Helicase Dependent Amplification

Preferred Specimen: Swab-Viral Culturette

Collection Instruction: Collect sample, place in holder, label with patient name. NOT TO BE USED FOR BACTERIAL TRANSPORT

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87529 (2)

Scarlet Eligible: No

Herpes Virus I/II DNA

H760-2

Alternate Name: HSV

Methodology: Helicase Dependent Amplification

Preferred Specimen: Swab-Viral Culturette

Collection Instruction: Collect sample, place in holder, label with patient name. NOT TO BE USED FOR BACTERIAL TRANSPORT

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87529 (2)

Scarlet Eligible: No

HSV 1+ 2/ VZV Assay

L008-2

Alternate Name: HSV 1 + 2/ VZV Assay

Preferred Specimen: Swab-Viral Culturette

Collection Instruction: Collect sample, place in holder, label with patient name. NOT TO BE USED FOR BACTERIAL TRANSPORT

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87529 (2), 87798 (1)

Scarlet Eligible: No

HERPES SIMPLEX VIRUS I DNA

J164-5

Alternate Name: HSV-1

Preferred Specimen: Swab-Viral Culturette

Collection Instruction: Collect sample, place in holder, label with patient name. NOT TO BE USED FOR BACTERIAL TRANSPORT

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87529 (1)

Scarlet Eligible: No

HERPES SIMPLEX VIRUS II DNA

J165-2

Alternate Name: HSV-2

Preferred Specimen: Swab-Viral Culturette

Collection Instruction: Collect sample, place in holder, label with patient name. NOT TO BE USED FOR BACTERIAL TRANSPORT

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87529 (1)

Scarlet Eligible: No

Herpes I/II Ab. Screen, IgG

0109-9

Alternate Name: HSV-I/II Ab. IgG

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86695 (1), 86696 (1)

Scarlet Eligible: Yes

HSV/LGV Lesion Profile (NON-NY)

P764-4

Alternate Name: HSV/LGV LESION PROFILE

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Turn Around Time: 7 days

CPT Code: 87529 (2), 87798 (1)

Scarlet Eligible: No

HSV 1+2 (NON NY)

L010-8

Alternate Name: HSV1+2
Methodology: Multiplex PCR
Preferred Specimen: Swab-Aptima Genprobe
Storage Instruction: Room Temp
Turn Around Time: 5-7 days
CPT Code: 87529 (1)
Scarlet Eligible: No

Hyperactive Coag. Panel

5993-1

Alternate Name: HYPERACTIVE COAG. PANEL
Preferred Specimen: Aliquot Plasma-Citrated-Blue Top
Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.
Storage Instruction: STRICT Frozen
Turn Around Time: 10 days
CPT Code: 83520 (2), 85379 (1)
Scarlet Eligible: No

Hypersensitivity Pneumonitis

1366-4

Alternate Name: HYPERSEN.PNEUMONITIS
Methodology: Immunodiffusion
Preferred Specimen: 1 mL SST Tube
Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes
Storage Instruction: Refrigerate
Turn Around Time: 8 days
CPT Code: 86331 (1), 86606 (1), 86609 (4)
Scarlet Eligible: Yes

Hyperthyroid Panel

0403-6

Alternate Name: HYPERTHYROID PANEL

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84436 (1), 84479 (1), 84480 (1)

Scarlet Eligible: Yes

Allergen Halibut (rf303), IgE

1485-2

Alternate Name: Halibut

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Hazelnut (f17), IgE

1404-3

Alternate Name: Hazelnut

Methodology: Immunocap

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Hemoglobin

1499-3

Alternate Name: Hb;Hgb

Methodology: Photometric

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 85018 (1)

Scarlet Eligible: Yes

Hemoglobin A1c with reflex to Glycomark

J218-9

Alternate Name: HbA1c Rfx To GlycoMark

Methodology: Turbidimetric Inhibition Immunoassay

Preferred Specimen: 0.5 mL Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 83036 (1)

Scarlet Eligible: Yes

Allergen Helminthosporium (m8), IgE

1894-5

Alternate Name: Helminthosporium

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Complement, Total (CH50), Serum

0534-8

Alternate Name: Hemolytic Complement

Methodology: Liposome

Preferred Specimen: 1 mL Aliquot Tube-Serum

Collection Instruction: Place 1-3 mL of serum in transport tube. Label as serum and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 6 days

CPT Code: 86162 (1)

Scarlet Eligible: No

Hepatitis B Panel

0197-4

Alternate Name: Hepatitis B Panel

Preferred Specimen: 5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86704 (1), 86706 (1), 87340 (1)

Scarlet Eligible: Yes

Hepatitis B Virus, DNA, Quantitative, RT-PCR

3389-4

Alternate Name: Hepatitis B Virus,Dna, Ultraquant

Methodology: Polymerase Chain Reaction

Preferred Specimen: 3 mL Aliquot Plasma-EDTA-Lavender Top

Collection Instruction: Centrifuge lavender-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA and patient's name.

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87517 (1)

Scarlet Eligible: Yes

Hepatitis C Antibody W/ Reflex RT PCR

B125-6

Alternate Name: Hepatitis C, HCV

Methodology: Chemiluminescence

Preferred Specimen: 1.0 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86803 (1)

Scarlet Eligible: Yes

Hepatitis Screen

0330-1

Alternate Name: Hepatitis Screen

Methodology: Chemiluminescence

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86704 (1), 86708 (1), 87340 (1)

Scarlet Eligible: Yes

Hepatitis Management Profile

0984-5

Alternate Name: Hepatitis/Management

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86704 (1), 86706 (1), 86707 (1), 87340 (1), 87350 (1)

Scarlet Eligible: Yes

HER2 FISH TC-Gastric+GEJ

A428-5

Alternate Name: Herceptin for Gastric

Methodology: FISH

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block from client with surgical number imprint.

Storage Instruction: Room Temp

Turn Around Time: 3-5 days

CPT Code: 88374 (1)

Scarlet Eligible: No

Hereditary MDS/Leukemia Panel

T830-9

Alternate Name: Hereditary Myeloid disorders

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: 2-5 mL Lavender top- EDTA

Collection Instruction: 2 tubes (5-6mL each) Label with the patient name, date of birth, and/or ID.

Storage Instruction: Room Temp

Turn Around Time: 21 days

CPT Code: 81405 (1), 81218 (1), 81479 (1)

Scarlet Eligible: Yes

Pancreatic Cancer Panel

B343-5

Alternate Name: Hereditary Pancreatic Cancer Panel

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: 2-5 mL Lavender top- EDTA

Collection Instruction: Lavender: 2 tubes (5-6ml each) Label with the patient name, date of birth and/or ID. Ideal handling is to ship immediately at ambient temperature for overnight delivery.

Alternate sample type: DNA: Extract and submit DNA 50 ng/mL (total volume 400 mL) in sterile container. Ship DNA in sealed and protective packaging at ambient temperature, using a cool pack in hot weather. Labeling: Name of person or identifier, DNA concentration and volume MUST be provided on the specimen tube.

Storage Instruction: Room Temp

Turn Around Time: 14 days

CPT Code: 81162 (1), 81201 (1), 81404 (1)

Scarlet Eligible: Yes

Hermansky-Pudlak Syndrome 3

B150-4

Alternate Name: Hermansky-Pudlak

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81479 (1)

Scarlet Eligible: Yes

Opiates, Urine, Screen (300)

0433-3

Alternate Name: Heroin; Narcotics; 6-O-Acetyl Morphine; Poppy Seeds

Methodology: Enzyme Immunoassay

Preferred Specimen: 5 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80307 (1)

Scarlet Eligible: Yes

Herpes I AB., IgG

0110-7

Alternate Name: Herpes I Ab.(IgG)

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86695 (1)

Scarlet Eligible: Yes

Herpes II Ab., IgG

0570-2

Alternate Name: Herpes II Ab.(IgG)

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86696 (1)

Scarlet Eligible: Yes

Herpes Simplex II by Multiplex (Non NY)

6352-9

Alternate Name: Herpes II by Multiplex

Methodology: Polymerase Chain Reaction

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Turn Around Time: 5-7 days

CPT Code: 87529 (1)

Scarlet Eligible: No

Herpes Simplex I by Multiplex (Non NY)

6351-1

Alternate Name: Herpes Simplex I Multi.

Methodology: Polymerase Chain Reaction

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Turn Around Time: 5-7 days

CPT Code: 87529 (1)

Scarlet Eligible: No

Varicella Zoster Ab., IgG0597-5

Alternate Name: Herpes Zoster Ab, IgG; Shingles**Methodology:** Chemiluminescence**Preferred Specimen:** 1 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 2 days**CPT Code:** 86787 (1)**Scarlet Eligible:** Yes

Histology Processing 8 H+EJ858-2

Alternate Name: Histo Processing 8 H+E**Methodology:** Histology**Preferred Specimen:** Biopsy In Formalin**Collection Instruction:** Collect specimen place in formalin jar write patient name and site on jar**Storage Instruction:** Room Temp**Turn Around Time:** 1 days**CPT Code:** 88305 (8)**Scarlet Eligible:** No

Histology Processing (GI)J862-4

Alternate Name: Histology Processing (GI)**Methodology:** Histology**Preferred Specimen:** Biopsy In Formalin**Collection Instruction:** Collect specimen place in formalin jar write patient name and site on jar**Storage Instruction:** Room Temp**Turn Around Time:** 1 days**CPT Code:** 88305 (11)**Scarlet Eligible:** No

Histology Processing 10 H+EJ861-6

Alternate Name: Histology Processing 10 H+E**Methodology:** Histology**Preferred Specimen:** Biopsy In Formalin**Collection Instruction:** Collect specimen place in formalin jar write patient name and site on jar**Storage Instruction:** Room Temp**Turn Around Time:** 1 days**CPT Code:** 88305 (10)**Scarlet Eligible:** No

Histology Processing 6 H+EJ857-4

Alternate Name: Histology Processing 6 H+E**Methodology:** Histology**Preferred Specimen:** Biopsy In Formalin**Collection Instruction:** Collect specimen place in formalin jar write patient name and site on jar**Storage Instruction:** Room Temp**Turn Around Time:** 1 days**CPT Code:** 88305 (6)**Scarlet Eligible:** No

Histology Processing 7 H+EJ859-0

Alternate Name: Histology Processing 7 H+E**Methodology:** Histology**Preferred Specimen:** Biopsy In Formalin**Collection Instruction:** Collect specimen place in formalin jar write patient name and site on jar**Storage Instruction:** Room Temp**Turn Around Time:** 1 days**CPT Code:** 88305 (7)**Scarlet Eligible:** No

Histology Processing 9 H+E

J860-8

Alternate Name: Histology Processing 9 H+E

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Collect specimen place in formalin jar write patient name and site on jar

Storage Instruction: Room Temp

Turn Around Time: 1 days

CPT Code: 88305 (9)

Scarlet Eligible: No

Histone Ab, IgG

T841-6

Alternate Name: Histone Ab., IgG

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 83516 (1)

Scarlet Eligible: Yes

Allergen Honey Bee (i1), IgE

0883-9

Alternate Name: Honey Bee

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Hops (RF324), IgE

A333-7

Alternate Name: Hops

Methodology: Immunocap

Preferred Specimen: 2.0 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Horse Dander (e3), IgE

0691-6

Alternate Name: Horse Dander

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Horsefly (i204), IgE

3556-8

Alternate Name: Horsefly (i204) IgE

Methodology: Enzyme Immunoassay

Preferred Specimen: 0.5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 6 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen House Dust (Greer) (h1), IgE

1631-1

Alternate Name: House Dust (Greer)

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen House Dust H/S (h2), IgE

0662-7

Alternate Name: House Dust H/S

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

House Dust Screen (5)

0910-0

Alternate Name: House Dust Screen

Preferred Specimen: 5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82785 (1), 86003 (5)

Scarlet Eligible: Yes

HLA-B27

0375-6

Alternate Name: Human Leukocyte Antigen B27
Methodology: Flow Cytometry
Preferred Specimen: 5 mL Green top- Sodium Heparin
Collection Instruction: Fill tube, invert gently 5-6 times, label with patient name
Storage Instruction: Refrigerate
Turn Around Time: 3 days
CPT Code: 86812 (1)
Scarlet Eligible: Yes

Prothrombin Fragment 1+2

5797-6

Alternate Name: Hypercoagulation Panel
Methodology: Enzyme Linked Immunoabsorbance
Preferred Specimen: 2 mL Aliquot Plasma-Citrated-Blue Top
Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.
Storage Instruction: STRICT Frozen
Turn Around Time: 10 days
CPT Code: 83520 (1)
Scarlet Eligible: No

IBD Differentiation Panel

H583-8

Alternate Name: IBD PANEL
Turn Around Time: 3 days
CPT Code: 86021 (2), 86255 (1), 86671 (2)
Scarlet Eligible: Yes

IGF Binding Protein-3

1316-9

Alternate Name: IGF Binding Protein-3; IGFB-3; Somatomedins

Methodology: Chemiluminescence

Preferred Specimen: 1 mL Aliquot Tube-Serum

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: STRICT Frozen

Turn Around Time: 2 days

CPT Code: 82397 (1)

Scarlet Eligible: Yes

Insulin-Like Growth Factor I

1072-8

Alternate Name: IGF-1;IGF 1;IGF-I; Sm-C; Somatomedin-C; Sulfation Factor

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 84305 (1)

Scarlet Eligible: Yes

Lymphoma - B-Cell, IGH by FISH

5032-8

Alternate Name: IGH, B-cell

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube. Also acceptable:Whole blood in green (sodium heparin) tube.

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88377 (1)

Scarlet Eligible: No

IGH/BCL2 t(14;18) by FISH TC

b350

Alternate Name: IGH/BCL2 t(14;18) by FISH TC

Scarlet Eligible: No

IGH/MAF t(14;16) by FISH

B564

Alternate Name: IGH/MAF t(14;16) by FISH

Scarlet Eligible: No

FISH, IGH/MAF, t(14;16)- TC

B760-0

Alternate Name: IGH/MAF- TC

Methodology: FISH

Preferred Specimen: Bone Marrow - Green Top

Collection Instruction: Place bone marrow into tube, label with patient's name

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88374 (1)

Scarlet Eligible: No

FISH, IGH/MAFB, t(14;20)

B545-5

Alternate Name: IGH/MAFB

Methodology: FISH

Preferred Specimen: 3 mL Bone Marrow - Green Top

Collection Instruction: Place bone marrow into tube, label with patient's name

Storage Instruction: Refrigerate

Turn Around Time: 9 days

CPT Code: 88374 (1)

Scarlet Eligible: No

IGH/MAFB T(14;20) - TC ONLY

TB61-7

Alternate Name: IGH/MAFB

Methodology: FISH

Preferred Specimen: 1 mL Bone Marrow - Green Top

Collection Instruction: Place bone marrow into tube, label with patient's name

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 88374 (1)

Scarlet Eligible: No

QFT - TB GOLD PLUS (QUANTIFERON)

T814-3

Alternate Name: IGRA

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 6 mL Green Top-Lithium Heparin

Collection Instruction: Please refer to special instruction sheet for instructions. Must be kept at Room temperature if not incubated. DO NOT SPIN. Must be received within 16 hrs of collection. If incubated stable for 3 days before centrifugation. After centrifugation stable for 7 days. Contact sales representative for account set-up. DO NOT REFRIGERATE OR FREEZE THE SAMPLE AT ANY TIME.

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86480 (1), 99999 (1)

Scarlet Eligible: Yes

IHC 1 Antibody Interp.

B024-1

Alternate Name: IHC 1 ANTIBODY INTERP.

Methodology: Immunohistochemistry

Preferred Specimen: Stained Slides

Storage Instruction: Room Temp

Turn Around Time: 1 days

CPT Code: 88342 (1)

Scarlet Eligible: No

IHC 1 Antibody/Stain Only

5144-1

Alternate Name: IHC 1 ANTIBODY/STAIN ONLY

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88342 (1)

Scarlet Eligible: No

IHC 2 Antibody/Stain Only

5145-8

Alternate Name: IHC 2 ANTIBODIES/STAIN ONLY

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88341 (1), 88342 (1)

Scarlet Eligible: No

IHC 3 Antibody/Stain Only

5146-6

Alternate Name: IHC 3 ANTIBODIES/STAIN ONLY

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88341 (2), 88342 (1)

Scarlet Eligible: No

IHC 4 Antibody/Stain Only

5147-4

Alternate Name: IHC 4 ANTIBODIES/STAIN ONLY

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88341 (3), 88342 (1)

Scarlet Eligible: No

IHC 5 Antibody/Stain Only

5148-2

Alternate Name: IHC 5 ANTIBODIES/STAIN ONLY

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88341 (4), 88342 (1)

Scarlet Eligible: No

IHC Antibody (2) With Interpretation

5137-5

Alternate Name: IHC ANTIBODY 2 W/INTERP.

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88341 (1), 88342 (1)

Scarlet Eligible: No

IHC Antibody (3) With Interpretation

5138-3

Alternate Name: IHC ANTIBODY 3 W/INTERP**Methodology:** Immunohistochemistry**Preferred Specimen:** Formalin-fixed, Paraffin-embedded Tissue**Collection Instruction:** This comes in a block or slide form from client with surgical number imprint.**Turn Around Time:** 1 days**CPT Code:** 88341 (2), 88342 (1)**Scarlet Eligible:** No**IHC Antibody (4) With Interpretation**

5139-1

Alternate Name: IHC ANTIBODY 4 W/INTERP.**Methodology:** Immunohistochemistry**Preferred Specimen:** Formalin-fixed, Paraffin-embedded Tissue**Collection Instruction:** This comes in block or slide from client with surgical number imprint.**Turn Around Time:** 1 days**CPT Code:** 88341 (3), 88342 (1)**Scarlet Eligible:** No**IHC Antibody (5) With Interpretation**

5141-7

Alternate Name: IHC ANTIBODY 5 W/INTERP**Methodology:** Immunohistochemistry**Preferred Specimen:** Formalin-fixed, Paraffin-embedded Tissue**Collection Instruction:** This comes in a block or slide form from client with surgical number imprint.**Turn Around Time:** 1 days**CPT Code:** 88341 (4), 88342 (1)**Scarlet Eligible:** No

IHC Single Antibody With Interpretation

5136-7

Alternate Name: IHC SINGLE AB W/INTER.

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Storage Instruction: Room Temp

Turn Around Time: 1 days

CPT Code: 88342 (1)

Scarlet Eligible: No

Interleukin 6

J785-7

Alternate Name: IL-6

Methodology: Chemiluminescence

Preferred Specimen: 1 mL Aliquot Tube-Serum

Collection Instruction: Place 1-3 mL of serum in transport tube. Label as serum;refrigerate or freeze as required.

Storage Instruction: STRICT Frozen

Turn Around Time: 48 hours

CPT Code: 83529 (1)

Scarlet Eligible: No

Imaged IHC (1) W/Interp.

5308-2

Alternate Name: IMAGED IHC (1) W/INTERP.

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88361 (1)

Scarlet Eligible: No

Imaged IHC, 1 Slide TC Only

5842-0

Alternate Name: IMAGED IHC, 1 SLIDE TC ONLY

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This specimen comes in block or slide form from client with surgical number imprint.

Turn Around Time: 2 days

CPT Code: 88361 (1)

Scarlet Eligible: No

Imaged IHC, 2 Slides TC Only

5838-8

Alternate Name: IMAGED IHC, 2 SLIDES TC ONLY

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This specimen comes in block or slide form from client with surgical number imprint.

Turn Around Time: 2 days

CPT Code: 88361 (2)

Scarlet Eligible: No

Imaged IHC, 2 Slides W/Interp.

5845-3

Alternate Name: IMAGED IHC, 2 SLIDES W/INTERP.

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This specimen comes in block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88361 (2)

Scarlet Eligible: No

Imaged IHC, 3 Slides TC Only

5837-0

Alternate Name: IMAGED IHC, 3 SLIDES TC ONLY

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This specimen comes in block or slide form from client with surgical number imprint.

Turn Around Time: 2 days

CPT Code: 88361 (3)

Scarlet Eligible: No

Imaged IHC, 3 Slides W/Interp.

5844-6

Alternate Name: IMAGED IHC, 3 SLIDES W/INTERP.

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This specimen comes in block or slide form from client with surgical number imprint.

Turn Around Time: 2 days

CPT Code: 88361 (3)

Scarlet Eligible: No

Imaged IHC, 4 Slides TC Only

5834-7

Alternate Name: IMAGED IHC, 4 SLIDES TC ONLY

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This specimen comes in block or slide form from client with surgical number imprint.

Turn Around Time: 2 days

CPT Code: 88361 (4)

Scarlet Eligible: No

Imaged IHC, 4 Slides W/Interp.

5843-8

Alternate Name: IMAGED IHC, 4 SLIDES W/INTERP.

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This specimen comes in block or slide form from client with surgical number imprint.

Turn Around Time: 2 days

CPT Code: 88361 (4)

Scarlet Eligible: No

Immunoglobulins, Serum

0520-7

Alternate Name: IMMUNOGLOBULINS, SER

Methodology: Immunoturbidimetric

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82784 (3)

Scarlet Eligible: Yes

Indoor Allergy Panel

P268-6

Alternate Name: INDOOR ALLERGY PANEL

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Turn Around Time: 3 days

CPT Code: 86003 (19)

Scarlet Eligible: Yes

Inherigen

A662-9

Alternate Name: INHERIGEN

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81443 (1)

Scarlet Eligible: Yes

Inherigen Plus

A659-5

Alternate Name: INHERIGEN PLUS

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81443 (1), 81329 (1), 81243 (1)

Scarlet Eligible: Yes

Inhibitor Titer

5757-0

Alternate Name: INHIBITOR TITER

Methodology: Clot Detection

Preferred Specimen: 2 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 3 days

CPT Code: 85335 (1)

Scarlet Eligible: No

Prenatal Rapid Aneuploidy FISH

3582-4

Alternate Name: INSIGHT BY FISH ANALYSIS; FLASH FISH

Methodology: Aneuvyvision by FISH

Preferred Specimen: 5 mL Amniotic Fluid

Collection Instruction: 5-6 ml Amniotic Fluid for OR 5-6 mg of Villi for CVS; please submit in a sterile, leak-proof tube(s)

Storage Instruction: Room Temp

Turn Around Time: 1-2 days

CPT Code: 88271 (5), 88274 (1), 88275 (1), 88291 (1)

Scarlet Eligible: No

Insulin Tolerance (3 Hour)

0024-0

Alternate Name: INSULIN TOLERANCE (3 hr)

Methodology: Electrochemiluminescence Immunoassay

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 83525 (4)

Scarlet Eligible: No

Insulin, 1st HR. Specimen

1649-3

Alternate Name: INSULIN, 1 HR.

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 83525 (1)

Scarlet Eligible: Yes

Insulin, Fasting

1648-5

Alternate Name: INSULIN, FASTING

Methodology: Indirect Fluorescence Assay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 83525 (1)

Scarlet Eligible: Yes

Interferon-ALPHA, EIA

3997-4

Alternate Name: INTERFERON-ALPHA

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 1 mL Aliquot Tube-Serum

Collection Instruction: Fill tube, invert gently 2-3 times, spin for 10-15 minutes label with patient name. Keep samples frozen. Freeze within 30 minutes of draw.

Storage Instruction: STRICT Frozen

Turn Around Time: 16 days

CPT Code: 83520 (1)

Scarlet Eligible: No

Interleukin 28B POLYMORPHISM (NON-NY)

A304-8

Alternate Name: INTERLEUKIN 28B

Methodology: Polymerase Chain Reaction

Preferred Specimen: 5 mL Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 8 days

CPT Code: 81283 (1)

Scarlet Eligible: Yes

Interleukin 5

J223-9

Alternate Name: INTERLEUKIN 5

Methodology: Electrochemiluminescence Immunoassay

Preferred Specimen: 1 mL Aliquot Tube-Serum

Collection Instruction: Place 1-3 mL of serum in transport tube. Label as serum;refrigerate or freeze as required.

Storage Instruction: STRICT Frozen

Turn Around Time: 9 days

CPT Code: 83520 (1)

Scarlet Eligible: No

Iodine, Serum/Plasma

J076-1

Alternate Name: IODINE S/P

Methodology: Inductively Coupled Plasma Mass Spectrometry

Preferred Specimen: Blue-Royal Blue Trace Element Serum (**RED LINE**)

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 83018 (1)

Scarlet Eligible: Yes

Iodine, Urine Panel

H474-0

Alternate Name: IODINE, URINE PANEL

Preferred Specimen: Urine Cup

Turn Around Time: 2 days

CPT Code: 82570 (1), 83789 (1)

Scarlet Eligible: Yes

Iron + TIBC

0250-1

Alternate Name: IRON + TIBC

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 83540 (1), 83550 (1)

Scarlet Eligible: Yes

Iron, Liver Studies, Quant.

5099-7

Alternate Name: IRON,LIVER QUANT.

Methodology: Inductively Coupled Plasma Mass Spectrometry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: Please write surgical number on block and on requisition. Provide fixation time on requisition as well

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 83540 (1)

Scarlet Eligible: No

ISH Single Antibody -PC Only

B960-6

Alternate Name: ISH SINGLE ANTIBODY-PC

Methodology: In Situ Hybridization

Preferred Specimen: Stained Slides

Collection Instruction: Client collect specimen, prepare smears, label with patient's name, and place in slide carrier.

Storage Instruction: Room Temp

Turn Around Time: 1-2 days

CPT Code: 88365 (1)

Scarlet Eligible: No

Isospora Belli ID, Modified Acid Fast

2670-8

Alternate Name: ISOSPORA BELLI ID, STOOL

Methodology: Modified Acid-Fast Stain

Preferred Specimen: Ova & Parasite Mini

Collection Instruction: With applicator stick, place 10 gms of stool into each vial, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 87015 (1), 87207 (1)

Scarlet Eligible: No

Accutype Ribavirin

B832-7

Alternate Name: ITPA

Methodology: Polymerase Chain Reaction

Preferred Specimen: 4.0 Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 87902 (1)

Scarlet Eligible: Yes

Immunoglobulin A, (IgA) Serum

0211-3

Alternate Name: IgA, Quantitative

Methodology: Immunoturbidimetric

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82784 (1)

Scarlet Eligible: Yes

Immunoglobulin G (IgG), Serum0541-3

Alternate Name: IgG, Quantitative**Methodology:** Immunoturbidimetric**Preferred Specimen:** 1 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 1 days**CPT Code:** 82784 (1)**Scarlet Eligible:** Yes

Lymphoma - B-Cell, IGH by PCR5278-7

Alternate Name: IgH Gene Rearrangement**Methodology:** Fluorescence**Preferred Specimen:** 2 mL Bone Marrow - Lavender Top**Collection Instruction:** Place bone marrow into tube, label with patient's name. Send in bone marrow kit.**Storage Instruction:** Refrigerate**Turn Around Time:** 6 days**CPT Code:** 81261 (1)**Scarlet Eligible:** No

Immunoglobulin M, (IgM) Serum0549-6

Alternate Name: IgM, Quantitative**Methodology:** Immunoturbidimetric**Preferred Specimen:** 1 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 1 days**CPT Code:** 82784 (1)**Scarlet Eligible:** Yes

CLL Prognosis: IgVH Mutation Analysis by PCR

5223-3

Alternate Name: IgVH

Methodology: Fluorescence

Preferred Specimen: 5 mL Peripheral Blood - Lavender Top

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient's name

Storage Instruction: Refrigerate

Turn Around Time: 7 days

CPT Code: 81263 (1)

Scarlet Eligible: Yes

Immunotyping, Serum

0413-5

Alternate Name: Immuno-electrophoresis, Immunofixation, IPEP

Methodology: Electrophoresis

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 4 days

CPT Code: 86334 (1)

Scarlet Eligible: Yes

Insulin, Total/Random

0113-1

Alternate Name: Immunoreactive Insulin

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 83525 (1)

Scarlet Eligible: Yes

InheriGen Plus w/o Fragile X

F021-1

Alternate Name: InheriGen Plus w/o Frag

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81443 (1), 81329 (1)

Scarlet Eligible: No

InheriGenTx

B735-2

Alternate Name: InheriGenTx

Methodology: Targeted Genotyping

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81205 (1), 81242 (1), 81250 (1), 81251 (1), 81400 (1), 81401 (1), 81402 (1)

Scarlet Eligible: Yes

Inhibin-A

TF79-1

Alternate Name: Inhibin-A

Methodology: Chemiluminescence

Preferred Specimen: 2 mL SST Tube

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86336 (1)

Scarlet Eligible: Yes

Phosphorus, Inorganic

0127-1

Alternate Name: Inorganic Phosphate, Blood

Methodology: Enzymatic Reference method with Hexokinase

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84100 (1)

Scarlet Eligible: Yes

Integrated Part 1 (NT, PAPP-A)

6375-0

Alternate Name: Integrated Screen First Stage

Methodology: Chemiluminescence

Preferred Specimen: 2 mL SST Tube

Collection Instruction: MUST BE 11-13 WEEKS AND 6 DAYS. PROVIDE GESTATIONAL INFO.

Storage Instruction: Refrigerate

Turn Around Time: 4 days

CPT Code: 99999 (1)

Scarlet Eligible: Yes

Integrated Part 2 (AFP, UE3, HCG, INHIBIN-A)

6376-8

Alternate Name: Integrated Screen Second Stage

Methodology: Chemiluminescence

Preferred Specimen: 3.0 mL SST Tube

Collection Instruction: MUST BE 15-22 WEEKS AND 6 DAYS. PROVIDE GESTATIONAL INFO.

Storage Instruction: Refrigerate

Turn Around Time: 4 days

CPT Code: 81511 (1), 84163 (1)

Scarlet Eligible: Yes

Intrinsic Factor AB

T835-8

Alternate Name: Intrinsic Factor AB

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 86340 (1)

Scarlet Eligible: Yes

Iron, Serum

0114-9

Alternate Name: Iron

Methodology: Colorimetry

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 83540 (1)

Scarlet Eligible: Yes

Iron Reflex to TIBC

M064-4

Alternate Name: Iron Reflex to TIBC

Turn Around Time: 1 days

CPT Code: 83540 (1)

Scarlet Eligible: Yes

Transferrin

0192-5

Alternate Name: Iron and Total Iron Total Binding Capacity/Transferrin, Serum

Methodology: Immunoturbidimetric

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84466 (1)

Scarlet Eligible: Yes

Iron/Ferritin Panel

1021-5

Alternate Name: Iron/Ferritin Panel

Turn Around Time: 1 days

CPT Code: 82728 (1), 83540 (1), 83550 (1)

Scarlet Eligible: Yes

Allergen Isocyanate MD1 (k76), IgE

2104-8

Alternate Name: Isocyanate MD1

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Isocyanate TD1 (k75), IgE

2103-0

Alternate Name: Isocyanate TD1

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Ispaghula (Psyllium) (k72), IgE

3968-5

Alternate Name: Ispaghula (K72)

Methodology: Enzyme Immunoassay

Preferred Specimen: 0.5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 6 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

JAK2 (including V617F and exon 12), MPL and CALR if negative reflex to MPN Panel

J632

Alternate Name: JAK2 (including V617F and exon 12), MPL and CALR if negative reflex to MPN Panel

Scarlet Eligible: No

JAK2 (including V617F and exon 12), MPL, and CALR if negative, reflex to MPN Panel J632-1

Alternate Name: JAK2/MPL/CALR by NGS rflx

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Bone Marrow - Lavender Top

Collection Instruction: Place bone marrow into tube, label with patient's name

Storage Instruction: Room Temp

Turn Around Time: 5-10 days

CPT Code: 81270 (1), 81219 (1), 81279 (1), 81338 (1)

Scarlet Eligible: No

JCV Ab w/ Index w/ RFX

B542-2

Alternate Name: JCV Ab. W/ Index w/ RFX

Methodology: Immunoassay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 16 days

CPT Code: 86711 (1)

Scarlet Eligible: Yes

Joubert Syndrome

A645-4

Alternate Name: JOUBERT SYNDROME

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81479 (1)

Scarlet Eligible: Yes

Allergen Johnson Grass (g10), IgE

1824-2

Alternate Name: Johnson Grass

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen June/Ken.Blue Grass (g8), IgE

0648-6

Alternate Name: June/Ken.Blue Grass

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Potassium, Serum

0129-7

Alternate Name: K+, Serum

Methodology: Ion Selective Electrode

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84132 (1)

Scarlet Eligible: Yes

Potassium, Urine 24 Hr.

0301-2

Alternate Name: K+, Urine, 24 Hours

Preferred Specimen: 10 mL Urine Container - 24hr

Collection Instruction: Discard first morning sample, then collect all urines after for 24 hrs, on completion label with collection time and patient's name; record total volume on Requisition

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 81050 (1), 84133 (1)

Scarlet Eligible: No

Cannabinoid, synthetic (K2,Spice), Urine, Screen, group 1, Qual.(20)

A292-5

Alternate Name: K2

Methodology: Enzyme Immunoassay

Preferred Specimen: 5 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 80307 (1)

Scarlet Eligible: Yes

Cannabinoid, synthetic (K2,Spice), Urine, Screen, group 2, Qual.(10)

J090-2

Alternate Name: K2 (SPICE)-GROUP2

Methodology: Enzyme Immunoassay

Preferred Specimen: 10 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 80307 (1)

Scarlet Eligible: Yes

Kappa Free Light Chain Serum

6241-4

Alternate Name: KAPPA FREE LIGHT CHAIN

Methodology: Immunoturbidimetric

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 83521 (1)

Scarlet Eligible: Yes

Free Kappa + Lambda Light Chain (serum)

3893-5

Alternate Name: KAPPA/LAMBDA LIGHT C

Methodology: Immunoturbidimetric

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 83520 (2)

Scarlet Eligible: Yes

Keppra Level (Levetiracetam)

3465-2

Alternate Name: KEPPRA LEVEL

Methodology: Immunoturbidimetric

Preferred Specimen: Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80177 (1)

Scarlet Eligible: Yes

KI-67 Imaged By IHC-Tech Only

A346-9

Alternate Name: KI-67 IMAGED

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block or slide form from client with surgical number imprint.

Storage Instruction: Room Temp

Turn Around Time: 1 days

CPT Code: 88361 (1)

Scarlet Eligible: No

KIT Sequencing (Solid Tumors) Non NY

J020-9

Alternate Name: KIT BY NGS

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 10-15 unstained slides at 5 microns (tumor content>10%) with H&E, Shavings (if tumor content >40%) with H&E, if sending DNA - Please call CS.

Storage Instruction: Room Temp

Turn Around Time: 5-10 days

CPT Code: 81272 (1)

Scarlet Eligible: No

Onkosight NGS KIT Sequencing (Non-NY)

J138-9

Alternate Name: KIT BY NGS

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Bone Marrow - Lavender Top

Collection Instruction: Place bone marrow into tube, label with patient's name

Storage Instruction: Room Temp

Turn Around Time: 5-10 days

CPT Code: 81273 (1)

Scarlet Eligible: No

KIT and PDGFRA Sequencing NY Only

J464-9

Alternate Name: KIT and PDGFRA Sequ.

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 10-15 unstained slides at 5 microns (tumor content >10%) with H&E, Shavings (if tumor content >10%) with H&E, if sending DNA- Please call CS

Turn Around Time: 5-10 days

CPT Code: 81272 (1), 81314 (1)

Scarlet Eligible: No

KIT and PDGFRA Sequencing- Non Ny

H578-8

Alternate Name: KIT and PDGFRA Sequ.

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Turn Around Time: 5-10 days

CPT Code: 81272 (1), 81314 (1)

Scarlet Eligible: No

KIT Sequencing (Solid Tumors) NY Only

J453-2

Alternate Name: KIT by NGS

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 10-15 unstained slides at 5 microns (tumor content >10%) with H&E, Shavings (if tumor content >10%) with H&E, if sending DNA- Please call CS

Turn Around Time: 5-10 days

CPT Code: 81272 (1)

Scarlet Eligible: No

KRAS Sequencing- Non NY

J024-1

Alternate Name: KRAS BY NGS

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 10-15 unstained slides at 5 microns (tumor content>10%) with H&E, Shavings (if tumor content >40%) with H&E, if sending DNA - Please call CS.

Storage Instruction: Room Temp

Turn Around Time: 5-10 days

CPT Code: 81275 (1), 81276 (1)

Scarlet Eligible: No

KRAS, if neg reflex to BRAF NY Only

J462-3

Alternate Name: KRAS NEG, RFLX TO BRAF

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 10-15 unstained slides at 5 microns (tumor content >10%) with H&E, Shavings (if tumor content >10%) with H&E, if sending DNA- Please call CS

Turn Around Time: 5-10 days

CPT Code: 81275 (1), 81276 (1)

Scarlet Eligible: No

KRAS, if neg reflex to BRAF- Non NY

J104-1

Alternate Name: KRAS NEG. RFLX TO BRAF

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 10-15 unstained slides at 5 microns (tumor content>10%) with H&E, Shavings (if tumor content >40%) with H&E, if sending DNA - Please call CS.

Storage Instruction: Room Temp

Turn Around Time: 5-10 days

CPT Code: 81275 (1), 81276 (1)

Scarlet Eligible: No

KRAS Sequencing NY Only

J458-1

Alternate Name: KRAS by NGS

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 10-15 unstained slides at 5 microns (tumor content >10%) with H&E, Shavings (if tumor content >10%) with H&E, if sending DNA- Please call CS

Turn Around Time: 5-10 days

CPT Code: 81275 (1), 81276 (1)

Scarlet Eligible: No

KRAS by NGS, if negative reflex to BRAF; + MSI by PCR, MSI-High, reflex to MLH1

H580-4

Alternate Name: KRAS>BRAF + MSI-(H)>MLH1

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Storage Instruction: Room Temp

Turn Around Time: 5-10 days

CPT Code: 81275 (1), 81301 (1), 81276 (1)

Scarlet Eligible: No

KRAS/NGS if neg rflx to BRAF/ MSI-PCR, if MSI-High, rflx MLH(NY)

J452-4

Alternate Name: KRAS>BRAF + MSI-(H)>MLH1

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 10-15 unstained slides at 5 microns (tumor content >10%) with H&E, Shavings (if tumor content >10%) with H&E, if sending DNA- Please call CS

Turn Around Time: 5-10 days

CPT Code: 81301 (1)

Scarlet Eligible: No

Peritoneal KT/V, Measured

6238-0

Alternate Name: KT/V, MEASURED

Preferred Specimen: 10 mL Fluid-source required,2 mL SST Tube

Collection Instruction: Record total volume on requisition. Include patient's weight and height. FLUID and SERUM are required.

Storage Instruction: Room Temp

Turn Around Time: 4 days

CPT Code: 84520 (2), 84540 (1)

Scarlet Eligible: No

Cabbage (f216), IgE

1244-3

Alternate Name: Kale

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Ketones (Acetoacetic), Serum

0028-1

Alternate Name: Ketones

Methodology: Manual

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 6 days

CPT Code: 82010 (1)

Scarlet Eligible: Yes

Ki-67 by IHC

5152-4

Alternate Name: Ki-67, IHC MANUAL

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88360 (1)

Scarlet Eligible: No

Allergen Kidney Bean (rf287), IgE

3122-9

Alternate Name: Kidney Bean

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Kidney Profile

L325-0

Alternate Name: Kidney Profile

Turn Around Time: 1 days

CPT Code: 82570 (1), 82043 (1), 99999 (1), 82565 (1)

Scarlet Eligible: Yes

Allergen Kiwi (f84), IgE

0842-5

Alternate Name: Kiwi

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Lupus-Hexagonal Phase

5749-7

Alternate Name: LA HEX DIFF

Methodology: Clot Detection

Preferred Specimen: 2 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 3 days

CPT Code: 85598 (1)

Scarlet Eligible: No

Lactobacillus (Species) by PCR

J213-0

Alternate Name: LACTOBACILLUS BY PCR

Methodology: Polymerase Chain Reaction

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 5-7 days

CPT Code: 87799 (1)

Scarlet Eligible: No

Lactose Tolerance Test1178-3

Alternate Name: LACTOSE TOLERANCE TEST**Methodology:** Enzyme Immunoassay**Preferred Specimen:** Grey top- Sodium Fluoride**Collection Instruction:** Submit one grey top (Sodium Fluoride) tube for EACH specimen, Label tubes with time of collection. Observe same rules as glucose tolerance. Collect all spec. (except fasting) after administration of 50gms of lactose in 8 ounce cup of water.**CPT Code:** 82951 (1), 82952 (3)**Scarlet Eligible:** No

LAL (Limulus Amebocyte Lysate)2475-2

Alternate Name: LAL**Methodology:** Chromogenic**Preferred Specimen:** 10 mL Other (forced comment) ?**Collection Instruction:** Collect water in special LAL tube. Order from Sendout Department.**Storage Instruction:** Refrigerate**Turn Around Time:** 5 days**CPT Code:** 87999 (1)**Scarlet Eligible:** No

Lambda Free Light Chain Serum6240-6

Alternate Name: LAMBDA FREE LIGHT CHAIN**Methodology:** Immunoturbidimetric**Preferred Specimen:** 1 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 2 days**CPT Code:** 83521 (1)**Scarlet Eligible:** Yes

Liver Cytosol (LC-1) Autoantibodies

B608-1

Alternate Name: LC-1 AUTOANTIBODIES

Methodology: Enzyme Immunoassay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 11 days

CPT Code: 86376 (1)

Scarlet Eligible: Yes

LDL Direct

2194-9

Alternate Name: LDL Direct

Methodology: Enzyme Immunoassay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 83721 (1)

Scarlet Eligible: Yes

Lead, Blood + E.P.

0588-4

Alternate Name: LEAD, BLOOD + E.P.

Preferred Specimen: Tan Top

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 83655 (1), 84202 (1)

Scarlet Eligible: Yes

Legionella Culture

3874-5

Alternate Name: LEGIONELLA CULTURE

Methodology: Bacterial Culture

Preferred Specimen: Cup-Sputum

Collection Instruction: Submit sputum in sterile container, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 12 days

CPT Code: 87081 (1)

Scarlet Eligible: No

Leptin (ELISA)

J220-5

Alternate Name: LEPTIN

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 83520 (1)

Scarlet Eligible: Yes

Level 3 Surgical Pathology

5695-2

Alternate Name: LEVEL 3 SURGICAL PATHOLOGY

Methodology: Histology

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 2 days

CPT Code: 88304 (1)

Scarlet Eligible: No

Level 4 Surgical Pathology

5696-0

Alternate Name: LEVEL 4 SURGICAL PATHOLOGY

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 2 days

CPT Code: 88305 (1)

Scarlet Eligible: No

Level 5 Surgical Pathology

5697-8

Alternate Name: LEVEL 5 SURGICAL PATHOLOGY

Methodology: Histology

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 2 days

CPT Code: 88307 (1)

Scarlet Eligible: No

Level II Histo,TC

A229-7

Alternate Name: LEVEL II HISTO,TC

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Collect specimen place in formalin jar write patient name and site on jar

Turn Around Time: 1 days

CPT Code: 88302 (1)

Scarlet Eligible: No

LEVEL II-PC SURGICAL PATHOLOGY

B012-6

Alternate Name: LEVEL II-PC SURGICAL PATHOLOGY

Methodology: Histology

Preferred Specimen: Slides

Collection Instruction: Label slides with patient name, place in slide holder for transport

Turn Around Time: 1 days

CPT Code: 88302 (1)

Scarlet Eligible: No

Level III Histo. (3 Slides)

5164-9

Alternate Name: LEVEL III HISTO TC

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 1 days

CPT Code: 88304 (3)

Scarlet Eligible: No

Level III Histo. (4 Slides)

5172-2

Alternate Name: LEVEL III HISTO. TC

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 1 days

CPT Code: 88304 (4)

Scarlet Eligible: No

Level III Histo.(5 Slides)

5173-0

Alternate Name: LEVEL III HISTO.(5 SLIDES)

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 1 days

CPT Code: 88304 (5)

Scarlet Eligible: No

LEVEL III-PC SURGICAL PATH.

B013-4

Alternate Name: LEVEL III-PC SURGICAL PATH.

Methodology: Histology

Preferred Specimen: Slides

Collection Instruction: Label slides with patient name, place in slide holder for transport

Turn Around Time: 1 days

CPT Code: 88304 (1)

Scarlet Eligible: No

Level IV Histo, TC

A231-3

Alternate Name: LEVEL IV HISTO, TC

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Collect specimen place in formalin jar write patient name and site on jar

Turn Around Time: 1 days

CPT Code: 88305 (1)

Scarlet Eligible: No

LEVEL IV-PC SURGICAL PATHOLOGY

B014-2

Alternate Name: LEVEL IV-PC SURGICAL PATHOLOGY

Methodology: Histology

Preferred Specimen: Slides

Collection Instruction: Label slides with patient name, place in slide holder for transport

Turn Around Time: 1 days

CPT Code: 88305 (1)

Scarlet Eligible: No

Level V Histo,TC

A228-9

Alternate Name: LEVEL V HISTO,TC

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Collect specimen place in formalin jar write patient name and site on jar

Turn Around Time: 1 days

CPT Code: 88307 (1)

Scarlet Eligible: No

LEVEL V-PC SURGICAL PATHOLOGY

B015-9

Alternate Name: LEVEL V-PC SURGICAL PATHOLOGY

Methodology: Histology

Preferred Specimen: Slides

Collection Instruction: Label slides with patient name, place in slide holder for transport

Turn Around Time: 1 days

CPT Code: 88307 (1)

Scarlet Eligible: No

Level VI Histo,TC

A230-5

Alternate Name: LEVEL VI HISTO,TC

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Collect specimen place in formalin jar write patient name and site on jar

Turn Around Time: 1 days

CPT Code: 88309 (1)

Scarlet Eligible: No

LEVEL VI-PC SURGICAL PATHOLOGY

B016-7

Alternate Name: LEVEL VI-PC SURGICAL PATHOLOGY

Methodology: Histology

Preferred Specimen: Slides

Collection Instruction: Label slides with patient name, place in slide holder for transport

Turn Around Time: 1 days

CPT Code: 88309 (1)

Scarlet Eligible: No

Lymphogranuloma Venereum (LGV) (Non-NY)

6356-0

Alternate Name: LGV, Non-NY patients, Chlamydia lesion

Methodology: Polymerase Chain Reaction

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Turn Around Time: 5-7 days

CPT Code: 87798 (1)

Scarlet Eligible: No

LH, FSH, Estradiol

0982-9

Alternate Name: LH, FSH, Estradiol

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82670 (1), 83001 (1), 83002 (1)

Scarlet Eligible: Yes

Luteinizing Hormone (LH)

0342-6

Alternate Name: LH; Follitropin; ICSH; Interstitial Cell Stimulating Hormone

Methodology: Electrochemiluminescence Immunoassay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 83002 (1)

Scarlet Eligible: Yes

Lipase, Peritoneal Fluid

6167-1

Alternate Name: LIPASE PERITONEAL,FLD

Methodology: Colorimetry

Preferred Specimen: 1 mL Fluid-source required

Collection Instruction: Submit fluid in sterile cup

Storage Instruction: Refrigerate

Turn Around Time: 4 days

CPT Code: 83690 (1)

Scarlet Eligible: No

Lipoprotein Electrophoresis w/Lipids

0551-2

Alternate Name: LIPOPROTEIN ELECTRO

Methodology: Electrophoresis

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80061 (1), 83700 (1)

Scarlet Eligible: Yes

Liver Panel

1074-4

Alternate Name: LIVER PANEL

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82040 (1), 82247 (1), 82248 (1), 82465 (1), 82977 (1), 83615 (1), 84075 (1), 84155 (1), 84450 (1), 84460 (1)

Scarlet Eligible: Yes

Lupus Antibody Panel (SLE)

4291-1

Alternate Name: LUPUS PANEL, SLE

Preferred Specimen: 1 mL SST Tube

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 83516 (1), 86038 (2), 86147 (3), 86160 (2), 86225 (1), 86235 (2)

Scarlet Eligible: Yes

Lyme AB West-Blot IgG/IgM

4033-7

Alternate Name: LYME AB WEST-BLOT IG

Methodology: Western Blot

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86617 (2)

Scarlet Eligible: Yes

Lyme Disease AB RFLX IGM TO WB G/M

2297-0

Alternate Name: LYME AB. RFLX WB G/M

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86618 (1)

Scarlet Eligible: Yes

Lyme Ab.,Total (No Reflex)

2382-0

Alternate Name: LYME Ab.(no reflex)

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86618 (1)

Scarlet Eligible: Yes

Lyme Disease Antibody, Serum

0568-6

Alternate Name: LYME DISEASE Ab.

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86618 (1)

Scarlet Eligible: Yes

Lyme Blot, IgG W/Bands

6272-9

Alternate Name: LYME WB, IGG W/BANDS

Methodology: Line Blot

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86617 (1)

Scarlet Eligible: Yes

Lyme Blot, IgM W/Bands

6273-7

Alternate Name: LYME WB, IGM W/BANDS

Methodology: Line Blot

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86617 (1)

Scarlet Eligible: Yes

Sjogren's Abs (SS-A, SS-B)

0668-4

Alternate Name: La Antibodies; Ro Antibodies; SS-A Antibodies;SS-B Antibodies

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86235 (2)

Scarlet Eligible: Yes

LabFinder Report

T420-9

Alternate Name: LabFinder Report

Collection Instruction: This code should not be reproted. It should auto complete in reporting.

CPT Code: 99999 (1)

Scarlet Eligible: No

Lactate Dehydrogenase (LDH)

0117-2

Alternate Name: Lactic Acid Dehydrogenase; LD; LDH

Methodology: Enzymatic Reference method with Hexokinase

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 83615 (1)

Scarlet Eligible: Yes

Lactobacillus reuteri/rhamnosus (Non NY)

J964-8

Alternate Name: Lactobacillus reuteri/rhamnosus

Methodology: Real-Time Polymerase Chain Reaction

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87798 (1)

Scarlet Eligible: No

Lamb/Mutton (f88) , IgE

0862-3

Alternate Name: Lamb/Mutton

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Latex (k82), IgE

1371-4

Alternate Name: Latex

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Lemon (f208), IgE

1224-5

Alternate Name: Lemon (F208)

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Lentil (f235), IgE

3521-2

Alternate Name: Lentil (F235)

Methodology: Immunocap

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Lettuce (f215), IgE

1409-2

Alternate Name: Lettuce

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Syphilis Reverse Algorithm

J275-9

Alternate Name: Liaison CIA (Antibody to T. pallidum)

Methodology: Chemiluminescence

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86780 (1)

Scarlet Eligible: Yes

Lipid Screen W/Reflex to Direct LDL

M063-6

Alternate Name: Lipid Screen W/Reflex to Direct LDL

Turn Around Time: 1 days

Scarlet Eligible: Yes

Liver Profile

0012-5

Alternate Name: Liver Profile

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80007 (1)

Scarlet Eligible: Yes

Allergen Lobster (f80), IgE

0841-7

Alternate Name: Lobster

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Lipoprotein-a (Lp-a)

1639-4

Alternate Name: Lp-A

Methodology: Immunoturbidimetric

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 83695 (1)

Scarlet Eligible: Yes

Lp-Pla2 (PLAC)

6155-6

Alternate Name: Lp-Pla2 , PLAC

Methodology: Enzyme Immunoassay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 83698 (1)

Scarlet Eligible: Yes

Lupus Anticoagulant Screen

5817-2

Alternate Name: Lupus Inhibitor Screen

Methodology: Clot Detection

Preferred Specimen: 3 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect 3 full blue-top (citrated) tubes. Centrifuge Blue-tops at high speed for 15 minutes. Transfer plasma into 5 plastic transfer tubes. Label with date/time collected, PLASMA, patient's name and FREEZE.

Storage Instruction: STRICT Frozen

Turn Around Time: 4 days

CPT Code: 85597 (1), 85598 (1), 85610 (1), 85613 (2), 85730 (2), 85732 (8)

Scarlet Eligible: No

Lyme Ab., IgM

1615-4

Alternate Name: Lyme Antibody, IgM; Lyme Disease Antibody Detection

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86618 (1)

Scarlet Eligible: Yes

Lynch/Colorectal High Risk Panel

B522

Alternate Name: Lynch/Colorectal High Risk Panel

Scarlet Eligible: No

Mycoplasma Genitalium by MultiPlex PCR- swab

J837-6

Alternate Name: M. Genitalium PCR Swab

Methodology: Multiplex PCR

Preferred Specimen: Swab-Aptima Genprobe

Collection Instruction: Collect specimen-break swab shaft into transport tube, then label with name and source

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87563 (1)

Scarlet Eligible: No

Multiple Myeloma FISH Panel - TC Only

P961-6

Alternate Name: M. MYELOMA FISH PANEL-TC ONLY

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube.

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88367 (1), 88373 (1), 88374 (4)

Scarlet Eligible: No

Magnesium, Urine, Random

3135-1

Alternate Name: MAGNESIUM, URINE, RANDOM

Methodology: Colorimetry

Preferred Specimen: 10 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 83735 (1)

Scarlet Eligible: Yes

Male STI Urine Profile (GP)

P310-6

Alternate Name: MALE STI URINE PROF(GP)

Preferred Specimen: 2 Urine Tube BD

Collection Instruction: Collect Specimen- then label with name and source. Recommended first morning urine.

Turn Around Time: 5 days

CPT Code: 87491 (1), 87591 (1), 87661 (1), 87798 (2)

Scarlet Eligible: Yes

Maple Syrup Urine Disease(Type 1A & 1B)

A650-4

Alternate Name: MAPLE SYRUP URN

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81205 (1)

Scarlet Eligible: Yes

MDM2 Gene Amplification by FISH

B544-8

Alternate Name: MDM2 Fish

Methodology: FISH

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Storage Instruction: Room Temp

Turn Around Time: 3 days

CPT Code: 88377 (1)

Scarlet Eligible: No

Ecstasy (MDMA), Urine, Screen, Qual. (300)

3174-0

Alternate Name: MDMA

Methodology: Enzyme Immunoassay

Preferred Specimen: 1 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80307 (1)

Scarlet Eligible: Yes

Megasphaera Type 1 by RT PCR

J214-8

Alternate Name: MEGASPHAERA TYPE 1 PCR

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 5-7 days

CPT Code: 87798 (1)

Scarlet Eligible: No

Metastatic Melanoma - Lymph Node

5126-8

Alternate Name: MELANOMA/LYMPH NODES

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88341 (2), 88342 (1)

Scarlet Eligible: No

Menopause/Perimenopause Risk Management Profile

P244-7

Alternate Name: MENOPAUSE/PERIMENOPAUSE

Preferred Specimen: SST Tube,Lavender top- EDTA

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 83718 (1), 83721 (1), 84439 (1), 84443 (1), 85025 (1), 86141 (1), G0059 (1)

Scarlet Eligible: Yes

Male High Risk Sexual Health Panel

F239-9

Alternate Name: MENS H/R SEXUAL PANEL

Collection Instruction: Submit 1 SST, 1 Viral Media, 1 Urine Aptima

Turn Around Time: 10 days

CPT Code: 86706 (1), 86803 (1), 87109 (2), 87389 (1), 87491 (1), 87591 (1), 87661 (1)

Scarlet Eligible: No

MET BY FISH-Global

B335-1

Alternate Name: MET BY FISH GLOBAL

Methodology: FISH

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block from client with surgical number imprint.

Storage Instruction: Room Temp

Turn Around Time: 3-5 days

CPT Code: 88377 (1)

Scarlet Eligible: No

MET BY FISH-Tech Only

B336-9

Alternate Name: MET BY FISH-TECH ONLY

Methodology: FISH

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Storage Instruction: Room Temp

Turn Around Time: 3-5 days

CPT Code: 88374 (1)

Scarlet Eligible: No

MET by IHC

A612-4

Alternate Name: MET amplification

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88342 (1)

Scarlet Eligible: No

Methadone, Urine, Screen (300)

0925-8

Alternate Name: METHADONE, URINE SCREEN

Methodology: Enzyme Immunoassay

Preferred Specimen: 5 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80307 (1)

Scarlet Eligible: Yes

MGMT Promoter Region Methylation

J876-4

Alternate Name: MGMT Promoter Methyla.

Methodology: Pyrosequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Storage Instruction: Room Temp

Turn Around Time: 7 days

CPT Code: 81287 (1)

Scarlet Eligible: No

Microalbumin, Urine, 24Hr.

3171-6

Alternate Name: MICROALBUMIN,URN. 24 HR

Methodology: Immunoturbidimetric

Preferred Specimen: 10 mL Urine Container - 24hr

Collection Instruction: Discard first morning sample, then collect all urines after for 24 hrs, on completion label with collection time and patient's name; record total volume on Requisition

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 81050 (1), 82043 (1)

Scarlet Eligible: No

Microsporidium ID, Stool

2673-2

Alternate Name: MICROSPORIDIUM ID, STOOL

Methodology: Microscopic Examination

Preferred Specimen: Ova & Parasite Mini

Collection Instruction: With applicator stick, place 10 gms of stool into each vial, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 87207 (1)

Scarlet Eligible: No

Milk Component Panel IgE

H191-0

Alternate Name: MILK COMPONENT PNL IGE

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82785 (1), 86003 (3)

Scarlet Eligible: Yes

Mixing Studies-APTT

5830-5

Alternate Name: MIXING STUDIES-APTT

Methodology: Clot Detection

Preferred Specimen: 3 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 3 days

CPT Code: 85732 (8)

Scarlet Eligible: No

Mixing Studies-PT

5840-4

Alternate Name: MIXING STUDIES-PT

Methodology: Clot Detection

Preferred Specimen: 3 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 3 days

CPT Code: 85611 (8)

Scarlet Eligible: No

MLH-1 Promo. Methylation

A313-9

Alternate Name: MLH-1 PROMO.METHYLATION

Methodology: Pyrosequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Turn Around Time: 7 days

CPT Code: 81288 (1)

Scarlet Eligible: No

AML-M4 and M5 by Fish- TC ONLY

B889-7

Alternate Name: MLL, 11q23

Methodology: FISH

Preferred Specimen: 5 mL Bone Marrow - Green Top

Collection Instruction: Place bone marrow into tube, label with patient's name

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 88374 (1)

Scarlet Eligible: No

AML - M4 + M5 by FISH

5427-0

Alternate Name: MLL, 11q23

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube. Also acceptable: Whole blood in green (sodium heparin) tube.

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 88374 (1)

Scarlet Eligible: No

Measles,Mumps,Rubella

2565-0

Alternate Name: MMR

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86735 (1), 86762 (1), 86765 (1)

Scarlet Eligible: Yes

Mobiluncus (Species) by RT PCR

J217-1

Alternate Name: MOBILUNCUS BY PCR

Methodology: Polymerase Chain Reaction

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 5-7 days

CPT Code: 87798 (1)

Scarlet Eligible: No

6-Monoacetyl Morphine (6-MAM), Urine, Screen (10)

3911-5

Alternate Name: MONOACETYLMORPHINE,6-QUANT

Methodology: Enzyme Immunoassay

Preferred Specimen: 3 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80307 (1)

Scarlet Eligible: Yes

Rule Out Mosaicism Chromosome Analysis, Blood

T982-8

Alternate Name: MOSAIC-CHROM-ANALYSIS-BLOOD-SEND

Methodology: Chromosome Analysis

Preferred Specimen: 5 mL Green top- Sodium Heparin

Collection Instruction: Fill tube, invert 8-10 times (DO NOT SHAKE!)

Storage Instruction: Room Temp

Turn Around Time: 7-14 days

CPT Code: 88263 (1), 88230 (1)

Scarlet Eligible: Yes

MPN: JAK2/MPL/CALR by NGS

J633-9

Alternate Name: MPN: JAK2/MPL/CALR by NGS

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Bone Marrow - Lavender Top

Collection Instruction: Place bone marrow into tube, label with patient's name

Storage Instruction: Room Temp

Turn Around Time: 5-10 days

CPT Code: 81270 (1), 81219 (1), 81279 (1), 81338 (1)

Scarlet Eligible: No

Myeloperoxidase Ab.

2192-3

Alternate Name: MPO

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86021 (1)

Scarlet Eligible: Yes

Minimal Residual Disease for CLL/SLL- Tech Only

J496-1

Alternate Name: MRD for CLL/SLL-TC Only

Methodology: Flow Cytometry

Preferred Specimen: 2 mL Peripheral Blood - Green Top

Collection Instruction: Fill tube, invert 8-10 times (DO NOT SHAKE!)

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 88184 (1), 88185 (8)

Scarlet Eligible: No

MRSA Screen

6119-2

Alternate Name: MRSA SCREEN

Methodology: Bacterial Culture

Preferred Specimen: Swab-E

Collection Instruction: Collect specimen with swab then place swab into transport carrier; Label with name and source

Turn Around Time: 3 days

CPT Code: 87070 (1)

Scarlet Eligible: No

Mismatch Repair(MMR) by IHC w/Interp

A943-3

Alternate Name: MSI BY IHC

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88341 (3), 88342 (1)

Scarlet Eligible: No

Mismatch Repair(MMR) by IHC, Tech Only

A944-1

Alternate Name: MSI BY IHC, TC**Methodology:** Immunohistochemistry**Preferred Specimen:** Formalin-fixed, Paraffin-embedded Tissue**Collection Instruction:** This comes in block or slide form from client with surgical number imprint.**Turn Around Time:** 1 days**CPT Code:** 88341 (3), 88342 (1)**Scarlet Eligible:** No**MSI: if MSI-H reflex to BRAF, if BRAF neg. reflex to MLH1**

J761-8

Alternate Name: MSI(H)>BRAF(neg)>MLH1**Methodology:** Polymerase Chain Reaction**Preferred Specimen:** Formalin-fixed, Paraffin-embedded Tissue**Collection Instruction:** This comes in block form from client with surgical number imprint.**Storage Instruction:** Room Temp**Turn Around Time:** 10 days**CPT Code:** 81301 (1)**Scarlet Eligible:** No**Microsatellite Instability (MSI) by PCR**

3371-2

Alternate Name: MSI, Mismatch Repair Protein Deficiency**Methodology:** Polymerase Chain Reaction**Preferred Specimen:** Formalin-fixed, Paraffin-embedded Tissue**Collection Instruction:** This comes in block from client with surgical number imprint.**Storage Instruction:** Refrigerate**Turn Around Time:** 7 days**CPT Code:** 81301 (1)**Scarlet Eligible:** No

MUC4 by IHC

TF93-2

Alternate Name: MUC4 by IHC

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Storage Instruction: Room Temp

Turn Around Time: 24-48 hours

CPT Code: 88342 (1)

Scarlet Eligible: No

MUC4 by IHC Tech Only

TF92-4

Alternate Name: MUC4 by IHC TC

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Storage Instruction: Room Temp

Turn Around Time: 24-48 hours

CPT Code: 88342 (1)

Scarlet Eligible: No

Mucopolidosis Type IV

A654-6

Alternate Name: MUCOLIPIDOSIS TYPE IV

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81290 (1)

Scarlet Eligible: Yes

Allergen Mucor Racemosus IgG

B770-9

Alternate Name: MUCOR RACEMOSUS IGG

Methodology: Immunoassay

Preferred Specimen: 0.5 mL Aliquot Tube-Serum

Collection Instruction: Collect 1-2 mL whole blood in red top tube, centrifuge and transfer 0.5 mL serum into a transfer tube. (All serum acceptable)

Storage Instruction: Refrigerate

Turn Around Time: 7 days

CPT Code: 86001 (1)

Scarlet Eligible: Yes

Mumps Virus Antibody, IgG

0316-0

Alternate Name: MUMPS VIRUS Ab.(IgG)

Methodology: Chemiluminescence

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86735 (1)

Scarlet Eligible: Yes

MYC-BA, 8q24 Rearrangement FISH- TECH ONLY

J144-7

Alternate Name: MYC-BA 8q24 BY FISH-TC

Methodology: FISH

Preferred Specimen: Bone Marrow - Green Top

Collection Instruction: Place bone marrow into tube, label with patient's name

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88374 (1)

Scarlet Eligible: No

MYC-BA, 8q24 Rearrangement FISH

J141-3

Alternate Name: MYC-MA 8q24 BY FISH

Methodology: FISH

Preferred Specimen: Bone Marrow - Green Top

Collection Instruction: Place bone marrow into tube, label with patient's name

Turn Around Time: 3-5 days

CPT Code: 88374 (1)

Scarlet Eligible: No

MYC/IGH t(8;14) by FISH TC

B439

Alternate Name: MYC/IGH t(8;14) by FISH TC

Scarlet Eligible: No

Lymphoma - Burkitt, MYC/IGH by FISH

5027-8

Alternate Name: MYC/IGH, t(8;14) for Burkitt Lymphoma, c-MYC

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube. Also acceptable: Whole blood in green (sodium heparin) tube.

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 88377 (1)

Scarlet Eligible: No

Mycoplasma Culture, Blood

3247-4

Alternate Name: MYCOPLASMA CULTURE,

Methodology: Bacterial Culture

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10 days

CPT Code: 87109 (1)

Scarlet Eligible: No

Mycoplasma Genitalium By Multiplex PCR

6340-4

Alternate Name: MYCOPLASMA GENITALIUM MULTIPLEX

Methodology: Polymerase Chain Reaction

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87563 (1)

Scarlet Eligible: No

Mycoplasma Hominis By Multiplex PCR (NON-NY)

6341-2

Alternate Name: MYCOPLASMA HOMINIS BY MULTIPLEX

Methodology: Polymerase Chain Reaction

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 87798 (1)

Scarlet Eligible: No

Mycoplasma pneumoniae IgG/IgM Antibody, IFA (CSF) (NON-NY)

3442-1

Alternate Name: MYCOPLASMA PNEUMONIA
Methodology: Indirect Fluorescence Assay
Preferred Specimen: 1 mL Fluid-source required
Collection Instruction: Submit fluid in sterile cup
Storage Instruction: Refrigerate
Turn Around Time: 7 days
CPT Code: 86738 (2)
Scarlet Eligible: No

Mycoplasma/Ureaplasma Culture

2523-9

Alternate Name: MYCOPLASMA/UREAPLASMA
Preferred Specimen: Swab-Viral Culturette
Collection Instruction: Collect sample, place in holder, label with patient name. NOT TO BE USED FOR BACTERIAL TRANSPORT
Storage Instruction: Refrigerate
Turn Around Time: 10 days
CPT Code: 87109 (1)
Scarlet Eligible: No

Myoglobin, Serum

0789-8

Alternate Name: MYOGLOBIN, SERUM
Methodology: Chemiluminescence
Preferred Specimen: 1 mL SST Tube
Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes
Storage Instruction: Refrigerate
Turn Around Time: 1 days
CPT Code: 83874 (1)
Scarlet Eligible: Yes

Allergen Macadamia Nuts (RF345), IgE

3644-2

Alternate Name: Macadamia Nut

Methodology: Immunocap

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Malt (f90), IgE

1630-3

Alternate Name: Malt

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Mandarin (rf302), IgE

3520-4

Alternate Name: Mandarin, Ige (Rf302)

Methodology: Immunocap

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Mango (f91), IgE

3456-1

Alternate Name: Mango (F91)**Methodology:** Immunocap**Preferred Specimen:** 3 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 3 days**CPT Code:** 86003 (1)**Scarlet Eligible:** Yes**Manual Diff/Platelet only**

1340-9

Alternate Name: Manual Diff/PLT only**Methodology:** Microscopic Examination**Preferred Specimen:** Lavender top- EDTA**Collection Instruction:** Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!**Storage Instruction:** Refrigerate**Turn Around Time:** 1 days**CPT Code:** 85007 (1), 85049 (1)**Scarlet Eligible:** Yes**Citrated Platelet Count**

J661-0

Alternate Name: Manual Platelet Count**Methodology:** Automated Hematology**Preferred Specimen:** 2.7 mL Light blue top**Collection Instruction:** Fill tube to line at top of tube, invert 8-10 times (DO NOT SHAKE!)**Storage Instruction:** Refrigerate**Turn Around Time:** 1 days**CPT Code:** 85049 (1)**Scarlet Eligible:** Yes

Allergen Maple (Box Elder) (t1), IgE

0684-1

Alternate Name: Maple (Box Elder)

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Marguerite/Ox-eye Daisy (w7), IgE

3879-4

Alternate Name: Marguerite/Ox-eye Daisy

Methodology: Enzyme Immunoassay

Preferred Specimen: 0.5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 4 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Cannabinoids (Marijuana), Urine, Screen (100)

0922-5

Alternate Name: Marijuana(THC)

Methodology: Enzyme Immunoassay

Preferred Specimen: 5 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80307 (1)

Scarlet Eligible: Yes

Allergen Marjoram (f274), IgE

6206-7

Alternate Name: Marjoram (f274) IgE

Methodology: Enzyme Immunoassay

Preferred Specimen: 0.5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 6 days

CPT Code: 86003 (1)

Scarlet Eligible: No

Allergen Marsh Elder (w16), IgE

1481-1

Alternate Name: Marsh Elder

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Meadow Fescue (g4), IgE

0681-7

Alternate Name: Meadow Fescue

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Rubeola/Measles IgG W/Reflex

T007-4

Alternate Name: Measles

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86765 (1)

Scarlet Eligible: Yes

Rubeola/Measles IgG Antibody

0567-8

Alternate Name: Measles Serology

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86765 (1)

Scarlet Eligible: Yes

Allergen Melaleuca Tree (t212), IgE

2642-7

Alternate Name: Melaleuca

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Melon (Cantaloupe) (f87), IgE

1227-8

Alternate Name: Melon (Cantaloupe)**Methodology:** Immunocap**Preferred Specimen:** 3 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 3 days**CPT Code:** 86003 (1)**Scarlet Eligible:** Yes**Meprobamate**

1936-4

Alternate Name: Meprobamate**Methodology:** Gas Chromatography/Flame Ionization**Preferred Specimen:** 1 mL Red Top**Collection Instruction:** Fill tube, label with patient name**Storage Instruction:** Refrigerate**Turn Around Time:** 12 days**CPT Code:** 80369 (1)**Scarlet Eligible:** Yes**Mesoridazine**

3709-3

Alternate Name: Mesoridazine**Methodology:** Liquid Chromatography Tandem Mass Spectrometry**Preferred Specimen:** 1 mL Red Top**Collection Instruction:** Fill tube, label with patient name**Storage Instruction:** Refrigerate**Turn Around Time:** 7 days**CPT Code:** 80342 (1)**Scarlet Eligible:** Yes

Vitamin B12

0160-2

Alternate Name: Methylcobalamin

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82607 (1)

Scarlet Eligible: Yes

Methylmalonic Acid, Serum

0349-1

Alternate Name: Methylmalonic Acid, Serum

Methodology: Liquid Chromatography Mass Spectrometry

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 83921 (1)

Scarlet Eligible: Yes

Magnesium, Serum

0120-6

Alternate Name: Mg, Serum

Methodology: Colorimetry

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 83735 (1)

Scarlet Eligible: Yes

Magnesium, Urine, 24Hr.

1171-8

Alternate Name: Mg, Urine, 24 Hours

Preferred Specimen: 10 mL Urine Cup 24 hour

Storage Instruction: Refrigerate

CPT Code: 81050 (1), 83735 (1)

Scarlet Eligible: No

Magnesium, RBC

1383-9

Alternate Name: Mg, Whole Blood

Methodology: Inductively Coupled Plasma Mass Spectrometry

Preferred Specimen: Blue-Royal Blue Top K2 EDTA Plasma(*LAVENDER LINE)

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 83735 (1)

Scarlet Eligible: No

Microalbumin, Urine, Random

0228-7

Alternate Name: Microalbumin/Creatinine Ratio

Methodology: Immunospectrometric

Preferred Specimen: 5 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82043 (1), 82570 (1)

Scarlet Eligible: Yes

MSI by PCR if MSI-H Reflex to MLH1

P264-5

Alternate Name: Microsatellite Instability, High Rflx to MLH1 Hypermethylation

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block or slide form from client with surgical number imprint.

Turn Around Time: 7 days

CPT Code: 81301 (1)

Scarlet Eligible: No

Mitochondrial Complex I Deficiency

J305-4

Alternate Name: Mitochondrial complex I deficiency

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81479 (1)

Scarlet Eligible: Yes

Mold Mix (mx1), IgE

2181-6

Alternate Name: Mold Mix

Methodology: Immunocap

Preferred Specimen: 5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86005 (1)

Scarlet Eligible: Yes

Molds Screen (6)

0901-9

Alternate Name: Molds Screen (6)

Preferred Specimen: 5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82785 (1), 86003 (6)

Scarlet Eligible: Yes

Skin Biopsy, (1 Jar)

5602-8

Alternate Name: Molluscum contagiosum

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes

Turn Around Time: 2-3 days

CPT Code: 88305 (1)

Scarlet Eligible: No

Mono Screen

0121-4

Alternate Name: Mono Screen

Methodology: Agglutination

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Freeze, Refrigerate

Turn Around Time: 1 days

CPT Code: 86308 (1)

Scarlet Eligible: Yes

Allergen Mosquito (i71), IgE

2179-0

Alternate Name: Mosquito

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Mountain Juniper (T6), IgE

6334-7

Alternate Name: Mountain Cedar

Methodology: Immunocap

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Mouse Epithelium (e71), IgE

1359-9

Alternate Name: Mouse Epithelium

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Mouse U. Protein (e72), IgE

1360-7

Alternate Name: Mouse U Protein

Methodology: Immunocap

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Mucor Racemosus (m4), IgE

1315-1

Alternate Name: Mucor Racemosus

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Mugwort (w6),IgE

0718-7

Alternate Name: Mugwort

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Multiple Double Stain by IHC Global

B962

Alternate Name: Multiple Double Stain by IHC Global

Scarlet Eligible: No

Multiple Double Stain by IHC TC Only

B963

Alternate Name: Multiple Double Stain by IHC TC Only

Scarlet Eligible: No

Multiple Sulfatase Deficiency

J296-5

Alternate Name: Multiple sulfatase deficiency

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81479 (1)

Scarlet Eligible: Yes

Mumps Virus Ab IgG W/Reflex

T009-0

Alternate Name: Mumps Virus Ab(IgG)Rfx

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86735 (1)

Scarlet Eligible: Yes

Allergen Mushroom (rf212), IgE

1462-1

Alternate Name: Mushroom (Rf212)

Methodology: Immunocap

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Mustard (f89), IgE

1823-4

Alternate Name: Mustard

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

MDS FISH Panel

5281-1

Alternate Name: Myelodysplastic Syndrome FISH Panel: +8, -5q, -7q, 20q-

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube. Also acceptable: Whole blood in green (sodium heparin) tube.

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 88367 (1), 88373 (1), 88374 (2)

Scarlet Eligible: No

MDS FISH Panel- TC Only

P960-8

Alternate Name: Myelodysplastic Syndrome FISH Panel: +8, -5q, -7q, 20q-

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube. Also acceptable: Whole blood in green (sodium heparin) tube.

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88367 (1), 88373 (1), 88374 (2)

Scarlet Eligible: No

Primidone (incl. Phenobarb)

0464-8

Alternate Name: Mysoline

Preferred Specimen: Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80184 (1), 80188 (1)

Scarlet Eligible: Yes

N. Gonorrhea, Urine, rRNA

6369-3

Alternate Name: N.GONORRHEA,URINE, rRNA

Methodology: Transcription Mediated Amplification

Preferred Specimen: Urine Container-Genprobe-Aptima

Collection Instruction: First morning sample, collect in sterile cup, transfer 2 ml into aptima holder. Correct level is between the black fill lines(window) on the holder.

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87591 (1)

Scarlet Eligible: Yes

Infant Tot.Bili (tot/dir=#0046)

1900-0

Alternate Name: N/A

Methodology: Colorimetry

Preferred Specimen: 1.0 mL Amber Serum

Collection Instruction: Draw SST or RED, let clot for 8-10 minutes. Centrifuge and pour into either an aliquot tube (must be covered in foil) or poured into an amber tube and label with patient name and specimen type (serum).

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82247 (1)

Scarlet Eligible: No

Chlamydia Direct Antigen

0782-3

Alternate Name: N/A

Methodology: Direct Immunofluorescence Staining

Preferred Specimen: Slides

Collection Instruction: Label slides with patient name, place in slide holder for transport

Storage Instruction: Refrigerate

Turn Around Time: 6 days

CPT Code: 87270 (1)

Scarlet Eligible: No

Infant Dir.Bili (tot/dir=#0046)

1901-8

Alternate Name: N/A

Methodology: Colorimetry

Preferred Specimen: 1.0 mL Amber Serum

Collection Instruction: Draw SST or RED, let clot for 8-10 minutes. Centrifuge and pour into either an aliquot tube (must be covered in foil) or poured into an amber tube and label with patient name and specimen type (serum).

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82248 (1)

Scarlet Eligible: No

Nefazodone Level

3147-6

Alternate Name: NEFAZODONE LEVEL

Methodology: High Pressure Liquid Chromatography

Preferred Specimen: 1 mL Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 12 days

CPT Code: 80338 (1)

Scarlet Eligible: Yes

Allergen Nettle (w20), IgE

6308-1

Alternate Name: NETTLE, IGE

Methodology: Immunocap

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Neuroendocrine Neoplasm IHC Panel

5195-3

Alternate Name: NEUROENDOCRINE NEOPLASM, IHC

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88341 (5), 88342 (1)

Scarlet Eligible: No

Neutrophil Function Panel

B124-9

Alternate Name: NEUTROPHIL FUNCTION PANEL

Methodology: Flow Cytometry

Preferred Specimen: 7 mL Green top- Sodium Heparin

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient name

Storage Instruction: Room Temp

Turn Around Time: 7 days

CPT Code: 82657 (1)

Scarlet Eligible: No

N. Gonorrhoea/NG (GC) PCR Urine

T408-4

Alternate Name: NG Urine PCR

Methodology: Polymerase Chain Reaction

Preferred Specimen: 10 mL Urine Cobas Vial

Collection Instruction: Patients should not urinate 1 hour prior to specimen collection.

Approximately 10-15 mL of urine should be collected from the beginning of the stream. Once collected the urine should be transferred to the cobas PCR media tube within 24 hrs. Use transfer pipette to add urine to the cobas PCR tube. Final volume of the tube should be between fill lines. Replace and tighten the cap and invert the tube 5 times to mix.

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87591 (1)

Scarlet Eligible: Yes

Ammonia, Plasma (As Nitrogen)

0035-6

Alternate Name: NH3, Plasma

Methodology: Enzyme Immunoassay

Preferred Specimen: 2 mL White Top (PPT)

Collection Instruction: After collection centrifuge at high speed for 8-10 min. Remove plasma and submit in aliquot tube. Label it plasma and freeze. If submitting a PPT, spin and freeze the entire tube immediately.

Storage Instruction: STRICT Frozen

Turn Around Time: 1 days

CPT Code: 82140 (1)

Scarlet Eligible: No

Niemann-Pick Type A+B Carrier Test

A653-8

Alternate Name: NIEMANN-PICK TYPE A+B

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81330 (1)

Scarlet Eligible: Yes

Nipple Discharge Smear

1303-7

Alternate Name: NIPPLE DISCHARGE SMR

Preferred Specimen: Slides

Collection Instruction: Label slides with patient name, place in slide holder for transport

Turn Around Time: 4 days

CPT Code: 88161 (1)

Scarlet Eligible: No

ClariTest Core Non-Invasive Prenatal Screen (NON-NY)

TH18-5

Alternate Name: NIPT, NIPS

Methodology: cfDNA is sequenced using microarray quantitation and DANSR and the FORTE algorithm.

Preferred Specimen: Cell-Free DNA White Top Plastic Tube

Storage Instruction: Room Temp

Turn Around Time: 5-7 days

CPT Code: 81507 (1)

Scarlet Eligible: Yes

ClariTest Core Non-Invasive Prenatal Screen with 22q11.2 Microdeletion (DiGeorge Syndrome)(NON-NY)

TH19-3

Alternate Name: NIPT, NIPS

Methodology: cfDNA is sequenced using microarray quantitation and DANSR and the FORTE algorithm.

Preferred Specimen: Cell-Free DNA White Top Plastic Tube

Storage Instruction: Room Temp

Turn Around Time: 5-7 days

CPT Code: 81507 (1)

Scarlet Eligible: Yes

Hodgkin Lymphoma IHC Panel

5101-1

Alternate Name: NON-HODGKINS V. HODGKINS

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 88341 (8), 88342 (1)

Scarlet Eligible: No

Northern Allerpanel (3 mix)

2187-3

Alternate Name: NORTHERN ALLERPANEL

Methodology: Immunocap

Preferred Specimen: 5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82785 (1), 86005 (3)

Scarlet Eligible: Yes

Onkosight NGS NPM1 Sequencing

J137-1

Alternate Name: NPM1 BY NGS

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Bone Marrow - Lavender Top

Collection Instruction: Place bone marrow into tube, label with patient's name

Turn Around Time: 7-14 days

CPT Code: 81310 (1)

Scarlet Eligible: No

NRAS Sequencing- Non NY

J023-3

Alternate Name: NRAS BY NGS

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 10-15 unstained slides at 5 microns (tumor content>10%) with H&E, Shavings (if tumor content >40%) with H&E, if sending DNA - Please call CS.

Storage Instruction: Room Temp

Turn Around Time: 5-10 days

CPT Code: 81311 (1)

Scarlet Eligible: No

NRAS Sequencing NY Only

J455-7

Alternate Name: NRAS by NGS

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 10-15 unstained slides at 5 microns (tumor content >10%) with H&E, Shavings (if tumor content >10%) with H&E, if sending DNA- Please call CS

Turn Around Time: 5-10 days

CPT Code: 81311 (1)

Scarlet Eligible: No

Sodium, Serum

0148-7

Alternate Name: Na, Serum

Methodology: Ion Selective Electrode

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84295 (1)

Scarlet Eligible: Yes

Sodium, Urine, 24 Hr.

0308-7

Alternate Name: Na, Urine, 24 Hours

Methodology: Ion Selective Electrode

Preferred Specimen: 10 mL Urine Container - 24hr

Collection Instruction: Discard first morning sample, then collect all urines after for 24 hrs, on completion label with collection time and patient's name; record total volume on Requisition

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 81050 (1), 84300 (1)

Scarlet Eligible: No

Sodium, Urine, Random

1793-9

Alternate Name: Na, Urine, Random
Methodology: Ion Selective Electrode
Preferred Specimen: 5 mL Urine Cup
Collection Instruction: Collect urine in cup, label with patients name
Storage Instruction: Refrigerate
Turn Around Time: 1 days
CPT Code: 84300 (1)
Scarlet Eligible: Yes

Naproxen (Aleve,Anaprox,Naprosyn)

6223-2

Alternate Name: Naproxen
Methodology: High Pressure Liquid Chromatography
Preferred Specimen: 1 mL Red Top
Collection Instruction: Fill tube, label with patient name
Storage Instruction: Refrigerate
Turn Around Time: 12 days
CPT Code: 80329 (1)
Scarlet Eligible: Yes

Culture, Nasopharyngeal/Nasal

0074-5

Alternate Name: Nasal Culture
Methodology: Bacterial Culture
Preferred Specimen: Swab-E
Collection Instruction: Collect specimen with swab then place swab into transport carrier; Label with name and source
Storage Instruction: Refrigerate
Turn Around Time: 3 days
CPT Code: 87070 (1)
Scarlet Eligible: No

Nasal Smear: Complete WBC w/ Diff

4009-7

Alternate Name: Nasal Smear

Methodology: Microscopic Examination

Preferred Specimen: Slides

Collection Instruction: Label slides with patient name, place in slide holder for transport

Turn Around Time: 2 days

CPT Code: 89190 (1)

Scarlet Eligible: No

GC By Multiplex PCR

6338-8

Alternate Name: Neisseria Gonorrhoeae

Methodology: Polymerase Chain Reaction

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 87591 (1)

Scarlet Eligible: No

Nemaline Myopathy 2

A646-2

Alternate Name: Nemaline Myopathy 2

Methodology: Polymerase Chain Reaction

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81479 (1)

Scarlet Eligible: Yes

New PD-L1

Alternate Name: New PD-L1

Scarlet Eligible: No

Cytology, Thinprep Fluids

5244-9

Alternate Name: Nipple Discharge, Bronchial Wash,Paracentesis,Pleural Fl,Sputum

Preferred Specimen: ThinPrep Cytolyt Solution

Collection Instruction: Collect specimen swish in Cytolyte vial, label with pt name.

Turn Around Time: 4 days

CPT Code: 88112 (1)

Scarlet Eligible: No

Cytology, Alcohol-Based Fluids

5245-6

Alternate Name: Nipple Discharge, Bronchial Washings, Paracentesis, Pleural Fluid

Preferred Specimen: 50% Alcohol Jar

Collection Instruction: Body fluid for cytology in 50% alcohol (ex...urine, sputum, nipple discharge).
Add directly to an equal volume of 50% Alcohol Solution.

Turn Around Time: 4 days

CPT Code: 88108 (1)

Scarlet Eligible: No

Non-Gyn Other (Slides) - Interp Only

J949-9

Alternate Name: Non-Gyn (Slides) -PC

Methodology: Microscopic Examination

Preferred Specimen: Slides

Collection Instruction: Label slides with patient name, place in slide holder for transport

Storage Instruction: Room Temp

Turn Around Time: 3 days

CPT Code: 88160 (1)

Scarlet Eligible: No

Allergen Profile, Northeast

0740-1

Alternate Name: Northeast Allergen Profile

Preferred Specimen: 5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82785 (1), 86003 (10)

Scarlet Eligible: Yes

Allergen Nut Mix (fx1), IgE

2183-2

Alternate Name: Nut Mix

Methodology: Immunocap

Preferred Specimen: 5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86005 (1)

Scarlet Eligible: Yes

Allergen Nutmeg (rf282), IgE

2352-3

Alternate Name: Nutmeg

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Ova and Parasite (O&P) Panel, Stool W/ Trichrome(3 specimens)

2326-7

Alternate Name: O + P Panel (3 Specimens), Ova and Parasites,

Methodology: Microscopic Examination

Preferred Specimen: Ova & Parasite Mini,O&P Kit (Total Fix)

Collection Instruction: With applicator stick, place 10 gms of stool into each vial, label with patient name

Storage Instruction: Room Temp

Turn Around Time: 2 days

CPT Code: 87177 (3), 87209 (3)

Scarlet Eligible: No

Ova and Parasite (O&P), Sputum

J002-7

Alternate Name: O&P Sputum, Ova and Parasites

Methodology: Microscopic Examination

Preferred Specimen: Cup-Sputum

Collection Instruction: Submit sputum in sterile container, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 48 hours

CPT Code: 87177 (1)

Scarlet Eligible: No

Ova and Parasite (O&P), Urine

2534-6

Alternate Name: O&P, Urine Ova and Parasites,

Methodology: Microscopic Examination

Preferred Specimen: Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Turn Around Time: 48 hours

CPT Code: 87177 (1)

Scarlet Eligible: Yes

Obstetric/Prenatal I

0008-3

Alternate Name: OBSTETRIC/PRENATAL I

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Turn Around Time: 1 days

CPT Code: 85025 (1), 86592 (1), 86762 (1), 86850 (1), 86900 (1), 86901 (1)

Scarlet Eligible: Yes

Obstetric/Prenatal II

0981-1

Alternate Name: OBSTETRIC/PRENATAL II

Preferred Specimen: SST Tube,Lavender top- EDTA,Pink Tube-EDTA

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate, Room Temp

Turn Around Time: 1 days

CPT Code: 85025 (1), 86592 (1), 86762 (1), 86850 (1), 86900 (1), 86901 (1), 87340 (1), G0059 (1)

Scarlet Eligible: Yes

OnkoSight Lung Cancer Panel, ALK and ROS1 by FISH

H125-8

Alternate Name: ONK LNG CA+ALK+ROS1

Methodology: Genotyping by Next Generation Sequencing

Collection Instruction: This comes in block form from client with surgical imprint. 10-15 unstained slides at 5 microns (tumor content>10%) with H&E, Shavings (if tumor content>40%) with H&E. If sending DNA- Please call Client Services.

Turn Around Time: 14 days

CPT Code: 81210 (1), 81235 (1), 81275 (1), 81276 (1), 81404 (1), 81405 (1), 81479 (1), 88377 (2)

Scarlet Eligible: No

Drug Confirmation, Opiates, Urine (GC/MS)

2596-5

Alternate Name: OPIATE CONFIRM. (GC/MS) 4 METAB.

Methodology: Gas Chromat. Mass Spectrometry

Preferred Specimen: 5 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 80363 (1)

Scarlet Eligible: Yes

Opiates, Urine, Screen (200)

6258-8

Alternate Name: OPIATES URINE SCREEN (200)

Methodology: Enzyme Immunoassay

Preferred Specimen: 4 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80307 (1)

Scarlet Eligible: Yes

Allergen Oak Tree (t7), IgE

0652-8

Alternate Name: Oak Tree

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Oat (f7), IgE

0814-4

Alternate Name: Oat

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Culture, Eye

0076-0

Alternate Name: Ocular/Eye Culture

Methodology: Bacterial Culture

Preferred Specimen: Swab-E

Collection Instruction: Collect specimen with swab then place swab into transport carrier; Label with name and source

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87070 (1)

Scarlet Eligible: No

Culture, Eye (Left)

0565-2

Alternate Name: Ocular/Eye Culture, Left

Methodology: Bacterial Culture

Preferred Specimen: Swab-E

Collection Instruction: Collect specimen with swab then place swab into transport carrier; Label with name and source

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87070 (1)

Scarlet Eligible: No

Culture, Eye (Right)

0633-8

Alternate Name: Ocular/Eye Culture, Right

Methodology: Bacterial Culture

Preferred Specimen: Swab-E

Collection Instruction: Collect specimen with swab then place swab into transport carrier; Label with name and source

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87070 (1)

Scarlet Eligible: No

Olive Tree (t9), IgE

1251-8

Alternate Name: Olive Tree

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

OneCheck Hematopathology

5500-4

Alternate Name: OneCheck

Preferred Specimen: Bone Marrow Core in Formalin

Storage Instruction: Refrigerate

Turn Around Time: 7 days

CPT Code: 99999 (1)

Scarlet Eligible: No

Allergen Onion (f48), IgE

0864-9

Alternate Name: Onion**Methodology:** Immunocap**Preferred Specimen:** 3 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 3 days**CPT Code:** 86003 (1)**Scarlet Eligible:** Yes**OnkoSight AML Panel Sequencing**

B894-7

Alternate Name: OnkoSight AML Panel**Methodology:** Genotyping by Next Generation Sequencing**Preferred Specimen:** 2 mL Bone Marrow - Lavender Top**Collection Instruction:** Place blood/bone marrow into tube, label with patient's name. Send in provided kit. If sending DNA - Please call C/S.**Storage Instruction:** Refrigerate**Turn Around Time:** 3-5 days**CPT Code:** 81450 (1)**Scarlet Eligible:** No**OnkoSight Colorectal Cancer Panel**

B822-8

Alternate Name: OnkoSight Colorectal Ca Pnl**Methodology:** Genotyping by Next Generation Sequencing**Preferred Specimen:** Formalin-fixed, Paraffin-embedded Tissue**Collection Instruction:** This comes in block form from client with surgical imprint. 10-15 unstained slides at 5 microns (tumor content>10%) with H&E, Shavings (if tumor content>40%) with H&E. If sending DNA- Please call Client Services.**Turn Around Time:** 5-10 days**CPT Code:** 81445 (1)**Scarlet Eligible:** No

OnkoSight Glioma Panel

J634-7

Alternate Name: OnkoSight Glioma Panel

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 10-15 unstained slides at 5 microns(tumor content>10%) with H&E, Shavings(if tumor content>40% with H&E, if sending DNA-Please call CS.

Storage Instruction: Room Temp

Turn Around Time: 5-10 days

CPT Code: 81445 (1)

Scarlet Eligible: No

OnkoSight Lung Cancer Panel

B819-4

Alternate Name: OnkoSight Lung Ca Panel

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical imprint. 10-15 unstained slides at 5 microns (tumor content>10%) with H&E, Shavings (if tumor content>40%) with H&E. If sending DNA- Please call Client Services.

Turn Around Time: 7-14 days

CPT Code: 81445 (1)

Scarlet Eligible: No

OnkoSight Lung Comprehensive

K924-1

Alternate Name: OnkoSight Lung Comprehensive

Turn Around Time: 14 days

CPT Code: 81210 (1), 81235 (1), 81275 (1), 81276 (1), 81404 (1), 81405 (1), 81479 (2), 88377 (1), 88342 (1)

Scarlet Eligible: No

OnkoSight Lung Panel, if EGFR- neg. rflx to (J015) ALK by FISH, if neg. rflx to ROS1 by FISH L113-0

Alternate Name: OnkoSight Lung Panel

Turn Around Time: 14 days

CPT Code: 81210 (1), 81235 (1), 81275 (1), 81276 (1), 81404 (1), 81405 (1), 81479 (2)

Scarlet Eligible: No

OnkoSight MDS Panel Sequencing B893-9

Alternate Name: OnkoSight MDS Panel

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: 2 mL Bone Marrow - Lavender Top

Collection Instruction: Place bone/blood marrow into tube, label with patient's name. Send in provided kit. If sending DNA - Please call CS.

Storage Instruction: Refrigerate

Turn Around Time: 5-10 days

CPT Code: 81450 (1)

Scarlet Eligible: No

OnkoSight MPN Panel Sequencing B818-6

Alternate Name: OnkoSight MPN Panel

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: 2 mL Peripheral Blood - Lavender Top

Collection Instruction: Place bone marrow into tube, label with patients name. Send in provided kit. If sending DNA-Please call Client Service.

Turn Around Time: 5-10 days

CPT Code: 81450 (1)

Scarlet Eligible: Yes

OnkoSight MYD88 Sequencing (Non-NY)

J130-6

Alternate Name: OnkoSight MYD88

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Bone Marrow - Lavender Top

Collection Instruction: Place bone marrow into tube, label with patient's name

Storage Instruction: Room Temp

Turn Around Time: 5-10 days

CPT Code: 81305 (1)

Scarlet Eligible: No

OnkoSight Melanoma Panel

B820-2

Alternate Name: OnkoSight Melanoma Pnl

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical imprint. 10-15 unstained slides at 5 microns (tumor content>10%) with H&E, Shavings (if tumor content>40%) with H&E. If sending DNA- Please call Client Services.

Turn Around Time: 5-10 days

CPT Code: 81445 (1)

Scarlet Eligible: No

OnkoSight Myeloid Disorder Panel

B823-6

Alternate Name: OnkoSight Myeloid Disorder Panel

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: 2 mL Bone Marrow - Lavender Top

Collection Instruction: Place bone marrow into tube, label with patients name. Send in provided kit. If sending DNA-Please call Client Service.

Turn Around Time: 5-10 days

CPT Code: 81450 (1)

Scarlet Eligible: No

OnkoSight Solid Tumor Panel

B821-0

Alternate Name: OnkoSight Solid Tumor Pnl

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical imprint. 10-15 unstained slides at 5 microns (tumor content>10%) with H&E, Shavings (if tumor content>40%) with H&E. If sending DNA- Please call Client Services.

Turn Around Time: 5-10 days

CPT Code: 81445 (1)

Scarlet Eligible: No

Opiates, Serum Screen

J865-7

Alternate Name: Opiates

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 6 mL Red Top

Collection Instruction: Fill tube, label with patient name

Turn Around Time: 3 days

CPT Code: 80307 (1)

Scarlet Eligible: Yes

Allergen Orange (f33), IgE

0835-9

Alternate Name: Orange

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Orchard Grass (g3), IgE

0649-4

Alternate Name: Orchard Grass

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Oregano (rf283), IgE

2351-5

Alternate Name: Oregano

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Ova and Parasite (O&P), Stool W/ Trichrome (2nd specimen)

0355-8

Alternate Name: Ova and Parasites

Methodology: Microscopic Examination

Preferred Specimen: Ova & Parasite Mini

Collection Instruction: With applicator stick, place 10 gms of stool into each vial, label with patient name

Storage Instruction: Room Temp

Turn Around Time: 2 days

CPT Code: 87177 (1), 87209 (1)

Scarlet Eligible: No

Ova and Parasite (O&P), Stool W/ Trichrome (3rd specimen)

0227-9

Alternate Name: Ova and Parasites,
Methodology: Microscopic Examination
Preferred Specimen: Ova & Parasite Mini
Storage Instruction: Room Temp
Turn Around Time: 2 days
CPT Code: 87177 (1), 87209 (1)
Scarlet Eligible: No

Allergen, Ovalbumin (F232), IgE

J281-7

Alternate Name: Ovalbumin IgE
Methodology: Immunocap
Preferred Specimen: 1.0 mL SST Tube
Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes
Storage Instruction: Refrigerate
Turn Around Time: 3 days
CPT Code: 86008 (1)
Scarlet Eligible: Yes

Allergen Ovomuroid (f233), IgE

J280-9

Alternate Name: Ovomuroid IgE
Methodology: Immunocap
Preferred Specimen: 1.0 mL SST Tube
Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes
Storage Instruction: Refrigerate
Turn Around Time: 3 days
CPT Code: 86008 (1)
Scarlet Eligible: Yes

Oxycodone, Urine, Screen (300)

A595-1

Alternate Name: Oxycontin

Methodology: Enzyme Immunoassay

Preferred Specimen: 5 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Turn Around Time: 1 days

CPT Code: 80307 (1)

Scarlet Eligible: Yes

Allergen Oyster (f290), IgE

3121-1

Alternate Name: Oyster

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Vitamin B6

0296-4

Alternate Name: P-5'-Pyridoxal-5-Phosphate;Pyridoxine

Methodology: High Pressure Liquid Chromatography

Preferred Specimen: 2 mL Amber Serum

Collection Instruction: Draw SST or RED, let clot for 8-10 minutes. Centrifuge and pour into either an aliquot tube (must be covered in foil) or poured into an amber tube and label with patient name and specimen type (serum).

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 84207 (1)

Scarlet Eligible: Yes

Neutrophil Cytoplasmic AB/ANCA

0636-1

Alternate Name: P-ANCA/C-ANCA

Methodology: Indirect Fluorescence Assay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86036 (1)

Scarlet Eligible: Yes

Prostate Specific Acid Phos (P.A.P.)

0030-7

Alternate Name: P.A.P.

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 84066 (1)

Scarlet Eligible: Yes

PTT (Partial Thromboplastin Time)

0139-6

Alternate Name: P.T.T.

Methodology: Photometric

Preferred Specimen: Light blue top

Collection Instruction: Fill tube to line at top of tube, invert 8-10 times (DO NOT SHAKE!)

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 85730 (1)

Scarlet Eligible: Yes

PAI 4G/5G Polymorphism

A269-3

Alternate Name: P4G/5G

Methodology: Polymerase Chain Reaction

Preferred Specimen: 2.0 mL Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 7 days

CPT Code: 81479 (1)

Scarlet Eligible: Yes

Chronic Lymphocytic Leuk.(CLL): FISH Pnl(+12, 13q, 17p, 11q)

5280-3

Alternate Name: P53/D13S319/ATM/+12

Methodology: FISH

Preferred Specimen: 5 mL Peripheral Blood - Green Top

Collection Instruction: Fill tube, invert 8-10 times (DO NOT SHAKE!)

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88367 (1), 88373 (3)

Scarlet Eligible: Yes

PANEL L202

L202-1

Alternate Name: PANEL L202

Turn Around Time: 4 days

CPT Code: 87086 (1), 83718 (1), 84443 (1), 85651 (1), 85025 (1)

Scarlet Eligible: Yes

PAP (RFLX HPV HIGH ASC-US or >)

6630-8

Alternate Name: PAP (RFLX HPV HIGH ASC-US or >)

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 88175 (1)

Scarlet Eligible: No

PAP (RFLX HPV HIGH ASC-US)

6629-0

Alternate Name: PAP (Rflx HPV Hi/Lo ASC-US)

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 88175 (1)

Scarlet Eligible: No

PAP + CT/GC

9241-1

Alternate Name: PAP + CT/GC

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Turn Around Time: 5 days

CPT Code: 87491 (1), 87591 (1), 88175 (1)

Scarlet Eligible: No

PAP + HPV HI

9102-5

Alternate Name: PAP + HPV HI

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87624 (1), 88175 (1)

Scarlet Eligible: No

PAP + HPV HI + CT/GC

9235-3

Alternate Name: PAP + HPV HI/LO + CT/GC

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 87491 (1), 87591 (1), 87624 (1), 88175 (1)

Scarlet Eligible: No

PAP RFLX HPV mRNA ASCUS,RFX 16/18/45 w/CT/GC

A253-7

Alternate Name: PAP CT/GC RFX HPVmRNA

Methodology: Transcription Mediated Amplification

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87491 (1), 87591 (1), 88175 (1)

Scarlet Eligible: No

PAP Reflex to HPV mRNA when ASCUS or greater,RFX 16,18/45

A254-5

Alternate Name: PAP RLX HPVmRNA, ASCUS>

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 88175 (1)

Scarlet Eligible: No

PAP W/O Image Reflex HR HPV

5240-7

Alternate Name: PAP W/O IMAGE RFLX HPV

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Turn Around Time: 5-7 days

CPT Code: 88142 (1)

Scarlet Eligible: No

PAP(RFX HPV Plus)

A489-7

Alternate Name: PAP(RFX HPV PLUS)

Preferred Specimen: ThinPrep Vial

Turn Around Time: 5 days

CPT Code: 88175 (1)

Scarlet Eligible: No

PAP+CT+GC(RF HPV+ LO ASC OR >)

P383-3

Alternate Name: PAP+CT+GC(RF HPV+ LO ASC OR >)

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Turn Around Time: 5 days

CPT Code: 87491 (1), 87591 (1), 88175 (1)

Scarlet Eligible: No

PAP+CT+GC(RFX HPV PLUS)

P374-2

Alternate Name: PAP+CT+GC(RFX HPV PLUS)

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Turn Around Time: 5 days

CPT Code: 87491 (1), 87591 (1), 88175 (1)

Scarlet Eligible: No

PAP + HPV Plus + CT + GC

P372-6

Alternate Name: PAP+HPV PLUS +CT +GC

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Turn Around Time: 5 days

CPT Code: 87491 (1), 87591 (1), 87624 (1), 88175 (1)

Scarlet Eligible: No

PAP, Liquid-Based

1962-0

Alternate Name: PAP, LIQUID-BASED

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Turn Around Time: 5 days

CPT Code: 88175 (1)

Scarlet Eligible: No

StormPath PAP, Liquid based

B352-6

Alternate Name: PAP-STORMPATH

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 4 days

CPT Code: 88175 (1)

Scarlet Eligible: No

PAP/CT/GC (Reflex HPV PLUS when ASC-US)

6431-1

Alternate Name: PAP/CT/GC (RFLX HPV)

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 87491 (1), 87591 (1), 88175 (1)

Scarlet Eligible: No

PAP/CT/GC (Reflex HPV PLUS when ASC-US or UNSAT)

L327-6

Alternate Name: PAP/CT/GC (Reflex HPV PLUS)

Turn Around Time: 5 days

CPT Code: 88175 (1), 87624 (1), 87491 (1), 87591 (1)

Scarlet Eligible: No

PAP RFX HPV Plus mRNA ASCUS>,RFX 16/18/45 w/CT/GC

A255-2

Alternate Name: PAP/CT/GC RFLXmRNA,ASC>

Preferred Specimen: ThinPrep Vial

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87491 (1), 87591 (1), 88175 (1)

Scarlet Eligible: No

PAP/CT/GC/ HPV mRNA,RFX 16/18/45 GENOTYPE

A252-9

Alternate Name: PAP/CT/GC/HPV mRNA

Methodology: Transcription Mediated Amplification

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87491 (1), 87591 (1), 87624 (1), 88175 (1)

Scarlet Eligible: No

PAP/HPVHI/GC RFX>GENCERV-P+HPV

P699-2

Alternate Name: PAP/HPVHI/GC RFX>GC-P+H

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Turn Around Time: 5 days

CPT Code: 87491 (1), 87591 (1), 87624 (1), 88175 (1), 88184 (1), 88185 (1), 88365 (1)

Scarlet Eligible: No

PAP/PSA

0956-3

Alternate Name: PAP/PSA

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 84066 (1), 84153 (1)

Scarlet Eligible: Yes

PAP/HPV+/GC RFX>GENCERV-PAP+HP (NON-NY)

P701-6

Alternate Name: PAPHPV+GC RFX>GC-PAP+HP

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Turn Around Time: 5 days

CPT Code: 87491 (1), 87591 (1), 87624 (1), 88175 (1), 88184 (1), 88185 (1), 88365 (1)

Scarlet Eligible: No

Parainfluenza Antibody

0946-4

Alternate Name: PARAINFLUENZA 1-3

Methodology: Complement Fixation

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 9 days

CPT Code: 86790 (3)

Scarlet Eligible: Yes

Path Cons On Ref.BLK/SLD GLOB.

5283-7

Alternate Name: PATH CONS ON REF.BLK/SLD GLOB.

Methodology: Morphology

CPT Code: 88323 (1)

Scarlet Eligible: No

Pathology Slide Consultation

5111-0

Alternate Name: PATHOLOGY CONSULT

Methodology: Histology

Preferred Specimen: Stained Slides

Collection Instruction: Client to collect specimen, prepare smears, label with patient's name, place in slide carrier.

Storage Instruction: Room Temp

Turn Around Time: 1-2 days

CPT Code: 88321 (1)

Scarlet Eligible: No

PBG Deaminase, WB

3245-8

Alternate Name: PBG Deaminase, WB

Methodology: Quantitative Spectrophotometry

Preferred Specimen: Green top- Sodium Heparin

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 10 days

CPT Code: 82657 (1)

Scarlet Eligible: Yes

Procalcitonin

T969-5

Alternate Name: PCT

Methodology: Chemiluminescence

Preferred Specimen: 1 mL Aliquot Tube-Serum

Collection Instruction: Place 1-3 mL of serum in transport tube, Label as serum and freeze as required.

Storage Instruction: STRICT Frozen

Turn Around Time: 4 days

CPT Code: 84145 (1)

Scarlet Eligible: No

PD-1 BY IHC

B994-5

Alternate Name: PD-1 BY IHC

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Storage Instruction: Room Temp

Turn Around Time: 1 days

CPT Code: 88342 (1)

Scarlet Eligible: No

PD-L1 (SP142) , TECENTRIQ.

T413-4

Alternate Name: PD-L1 (SP142)

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Storage Instruction: Room Temp

Turn Around Time: 1-2 days

CPT Code: 88360 (1)

Scarlet Eligible: No

PD-L1 22C3 for Keytruda Other Tumor Types

TK85

Alternate Name: PD-L1 22C3 for Keytruda Other Tumor Types

Collection Instruction: Not orderable as of 2/1 but launching within the month of Feb

Scarlet Eligible: No

PD-L1 IHC, Keytruda for Cervical Cancer W/INTERP.

T988-5

Alternate Name: PD-L1 Cervical Cancer

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Storage Instruction: Room Temp

Turn Around Time: 4 days

CPT Code: 88360 (1)

Scarlet Eligible: No

PD-L1 (Clone E1L3N) IHC w/Interp.

J257-7

Alternate Name: PD-L1 Generic

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Alternate specimen is 5 unstained slides cut @ 3 to 5 microns.

Turn Around Time: 4 days

CPT Code: 88360 (1)

Scarlet Eligible: No

Keytruda PD-L1 Head & Neck Squamous Carcinoma (HNSCC)

TH63-1

Alternate Name: PD-L1 HNSCC Keytruda

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Storage Instruction: Room Temp

Turn Around Time: 2 days

CPT Code: 88360 (1)

Scarlet Eligible: No

PD-L1 IHC 22C3 pharmDx, Keytruda W/INTERP.

B993-7

Alternate Name: PD-L1 Keytruda

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Storage Instruction: Room Temp

Turn Around Time: 4 days

CPT Code: 88360 (1)

Scarlet Eligible: No

PD-L1 IHC 28-8 pharmDx (OPDIVO) W/INTERP.

J221-3

Alternate Name: PD-L1 Opdivo

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Turn Around Time: 4 days

CPT Code: 88360 (1)

Scarlet Eligible: No

PD-L1 (Clone E1L3N) IHC Tech Only

TB35-1

Alternate Name: PD-L1(E1L3N)-Tech Only

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Storage Instruction: Room Temp

Turn Around Time: 2 days

CPT Code: 88360 (1)

Scarlet Eligible: No

PD-L1 IHC , Keytruda for Urothelial Carcinoma W/INTERP.

TA31-2

Alternate Name: PD-L1-Urothelial Carcin

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Storage Instruction: Room Temp

Turn Around Time: 4 days

CPT Code: 88360 (1)

Scarlet Eligible: No

PDGFRA Sequencing- Non NY

J045-6

Alternate Name: PDGFRA Sequencing

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 10-15 unstained slides at 5 microns (tumor content >10%) with H&E, Shavings (if tumor content >40%) with H&E, if sending DNA - Please call CS.

Storage Instruction: Room Temp

Turn Around Time: 5-10 days

CPT Code: 81314 (1)

Scarlet Eligible: No

PDGFRA Sequencing

J459-9

Alternate Name: PDGFRA Sequencing

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 10-15 unstained slides at 5 microns (tumor content >10%) with H&E, Shavings (if tumor content >10%) with H&E, if sending DNA- Please call CS

Turn Around Time: 5-10 days

CPT Code: 81314 (1)

Scarlet Eligible: No

PDGFRB by FISH (Global)- Non NY

T781-4

Alternate Name: PDGFRB by FISH

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Place bone marrow into tube, label with patient's name

Storage Instruction: Room Temp

Turn Around Time: 3 days

CPT Code: 88374 (1)

Scarlet Eligible: No

PDGFRB by FISH (Tech Only) Non NY

T782-2

Alternate Name: PDGFRB by FISH- TC

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Place bone marrow into tube, label with patient's name

Storage Instruction: Room Temp

Turn Around Time: 3 days

CPT Code: 88374 (1)

Scarlet Eligible: No

Keytruda PD-L1 Esophageal Squamous Cell Carcinoma by IHC

TH62-3

Alternate Name: PDL1 Esoph-SCC Keytruda

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Storage Instruction: Room Temp

Turn Around Time: 2 days

CPT Code: 88360 (1)

Scarlet Eligible: No

PDL1 (22C3 Keytruda) IHC + EGFR NGS + ALK FISH + ROS1 FISH

L393-8

Alternate Name: PDL1+EGFR- NGS+ALK+ROS1

Turn Around Time: 10 days

CPT Code: 88360 (1), 81235 (1), 88377 (1)

Scarlet Eligible: No

Pediatric Arthritis Profile

0596-7

Alternate Name: PEDIATRIC ARTHRITIS

Preferred Specimen: SST Tube,Lavender top- EDTA

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 82784 (3), 85651 (1), 86038 (1), 86431 (1), 86618 (1)

Scarlet Eligible: No

Pediatric Blood Culture #1

B331-0

Alternate Name: PEDIATRIC BL CULT #1

Methodology: Bacterial Culture

Preferred Specimen: Pediatric Blood Culture Bottle

Collection Instruction: Collect 1-3 mL of blood in bottle, label with patient's name.

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87040 (1)

Scarlet Eligible: No

Pediatric Blood Culture #2

B332-8

Alternate Name: PEDIATRIC BL CULT #2

Methodology: Bacterial Culture

Preferred Specimen: Pediatric Blood Culture Bottle

Collection Instruction: Collect 1-3 mL of blood in bottle, label with patient's name.

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87040 (1)

Scarlet Eligible: No

Pediatric Blood Culture #3

B333-6

Alternate Name: PEDIATRIC BL CULT #3

Methodology: Bacterial Culture

Preferred Specimen: Pediatric Blood Culture Bottle

Collection Instruction: Collect 1-3 mL of blood in bottle, label with patient's name.

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87040 (1)

Scarlet Eligible: No

Pediatric Profile I

0980-3

Alternate Name: PEDIATRIC PROFILE I

Preferred Specimen: Lavender top- EDTA,Tan Top

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 83655 (1), 84202 (1), 85025 (1), 85660 (1)

Scarlet Eligible: No

Pediatric Profile II

0979-5

Alternate Name: PEDIATRIC PROFILE II

Preferred Specimen: SST Tube,Lavender top- EDTA,Tan Top

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 83655 (1), 84202 (1), 85025 (1), 85660 (1), G0059 (1)

Scarlet Eligible: No

Perennial Allerpanel (3 mix)

2188-1

Alternate Name: PERENNIAL ALLERPANEL

Methodology: Immunocap

Preferred Specimen: 5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82785 (1), 86005 (3)

Scarlet Eligible: Yes

PET II (Dialysis Fluid)

4447-9

Alternate Name: PET II (DIALYSIS FLU

Preferred Specimen: Fluid-source required

Collection Instruction: Submit 1 SST (master tube) and 3 separate containers of peritoneal fluid (timed 0hr, 2hr, 4hr) Refrigerate

Turn Around Time: 4 days

CPT Code: 82565 (1), 82570 (3), 82945 (3), 82947 (1), 84520 (4)

Scarlet Eligible: No

Phenosense Entry (incl.Fuzeon)

6235-6

Alternate Name: PHENOSENSE ENTRY

Preferred Specimen: 10 mL Aliquot Plasma-EDTA-Lavender Top

Collection Instruction: Centrifuge lavender-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA and patient's name. Must submit frozen.

Storage Instruction: STRICT Frozen

Turn Around Time: 24 days

CPT Code: 87903 (1)

Scarlet Eligible: No

Phlebotomy Charge

0305-3

Alternate Name: PHLEBOTOMY CHARGE

CPT Code: 36415 (1)

Scarlet Eligible: Yes

Phosphatidyl Ethanolamine

5722-4

Alternate Name: PHOSPHATIDYL ETHANOL

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 8 days

CPT Code: 83520 (2)

Scarlet Eligible: Yes

Phosphorus, Urine, Random

3134-4

Alternate Name: PHOSPHORUS, URINE, RANDOM

Methodology: Enzymatic Reference method with Hexokinase

Preferred Specimen: 10 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 84105 (1)

Scarlet Eligible: Yes

Phosphates, Urine 24 Hr.

0411-9

Alternate Name: PHOSPHORUS, URINE 24 HR

Preferred Specimen: 10 mL Urine Container - 24hr

Collection Instruction: Discard first morning sample, then collect all urines after for 24 hrs, on completion label with collection time and patient's name; record total volume on Requisition

Storage Instruction: Refrigerate

CPT Code: 81050 (1), 84105 (1)

Scarlet Eligible: No

PID/Infertility W/O PAP Expanded

F368-6

Alternate Name: PID/INFRT NO PAP EXPAND

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Turn Around Time: 7 days

CPT Code: 87491 (1), 87512 (1), 87591 (1), 87661 (1), 87798 (4), 87799 (2)

Scarlet Eligible: No

Pinworm Tape Test

0204-8

Alternate Name: PINWORM TAPE TEST

Methodology: Microscopic Examination

Preferred Specimen: Pinworm slide w/tape, paddle

Collection Instruction: Press paddle or sticky side of tape without smearing on to the fold of the perianal area (press the sticky side of the paddle against the skin around the rectum 3-4 times).

Collect specimen first thing in the morning before bathing or bowel movement.

Turn Around Time: 1 days

CPT Code: 87172 (1)

Scarlet Eligible: No

Pituitary I

2319-2

Alternate Name: PITUITARY I

Preferred Specimen: 4 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 83001 (1), 83002 (1), 84146 (1), 84443 (1)

Scarlet Eligible: Yes

Plasma Cell Dyscrasia/Multiple Myeloma IHC Panel

5579-8

Alternate Name: PLASMA CELL DYSCRASIA/MM

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88341 (10), 88342 (1)

Scarlet Eligible: No

Plasma Cell Analysis by Flow Cytometry (Tech Only) and StormPath BMM Report Template

R736-0

Alternate Name: PLASMA CELL FlowTC and BMM

Turn Around Time: 1 days

CPT Code: 99999 (1), 88184 (1), 88185 (11)

Scarlet Eligible: No

Plasminogen Quant.

5775-2

Alternate Name: PLASMINOGEN, QUANT.

Methodology: Chromogenic

Preferred Specimen: 2 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 3 days

CPT Code: 85420 (1)

Scarlet Eligible: No

Platelet Antibody Identification5776-0

Alternate Name: PLATELET AB. IDENT.**Methodology:** Enzyme Linked Immunoabsorbance**Preferred Specimen:** 5 mL Aliquot Tube-Serum**Collection Instruction:** Fill tube, invert gently 2-3 times, let stand for 20 minutes, spin for 10-15 minutes label with patient name.**Storage Instruction:** STRICT Frozen**Turn Around Time:** 15 days**CPT Code:** 86022 (7)**Scarlet Eligible:** No

PNH by Flow Cytometry5564-0

Alternate Name: PNH Test**Methodology:** Flow Cytometry**Preferred Specimen:** 5 mL Peripheral Blood - Green Top**Collection Instruction:** Fill tube, invert gently 5-6 times, label with patient's name**Storage Instruction:** Refrigerate**Turn Around Time:** 1 days**CPT Code:** 88184 (1), 88185 (7), 88187 (1)**Scarlet Eligible:** Yes

Platelet Neutralizing Procedure5784-4

Alternate Name: PNP**Methodology:** Clot Detection**Preferred Specimen:** 2 mL Aliquot Plasma-Citrated-Blue Top**Collection Instruction:** Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.**Storage Instruction:** STRICT Frozen**Turn Around Time:** 4 days**CPT Code:** 85597 (1)**Scarlet Eligible:** No

Polycyst.Ovary Synd. Treat. Panel

P259-5

Alternate Name: POLYCYST.OVARY SYND.

Preferred Specimen: SST Tube,Lavender top- EDTA

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 83036 (1), 83718 (1), G0059 (1)

Scarlet Eligible: Yes

Polycyst.Ovary Synd. Conf.Panel

P258-7

Alternate Name: POLYCYST.OVARY SYND.CON

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 82040 (1), 82627 (1), 83001 (1), 83002 (1), 84146 (1), 84270 (1), 84305 (1), 84403 (1), 84439 (1), 84443 (1)

Scarlet Eligible: Yes

Post Dialysis Potassium

1688-1

Alternate Name: POTASSIUM POST DIALYSIS

Methodology: Ion Selective Electrode

Preferred Specimen: Microtainer - Pediatric SST

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84132 (1)

Scarlet Eligible: No

Potassium, Urine, Timed Specimen

1659-2

Alternate Name: POTASSIUM,URINE TIMED

Methodology: Ion Selective Electrode

Preferred Specimen: 3 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84133 (1)

Scarlet Eligible: Yes

Proteinase-3 AutoAbs

2663-3

Alternate Name: PR3

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86021 (1)

Scarlet Eligible: Yes

15q Del Characterization, FISH

B796-4

Alternate Name: PRADER-WILLI SYND CVS

Methodology: In Situ Hybridization

Preferred Specimen: 20-25 mL Chorionic Villi

Collection Instruction: CVS specimen should be collected in a sterile, screw-top container with tissue culture transport medium (available upon Request). If transport medium is not available, plain RPMI, Hanks solution or saline can be used

Storage Instruction: Room Temp

Turn Around Time: 12 days

CPT Code: 88271 (2), 88291 (1)

Scarlet Eligible: No

Preg.+Throm.Assessment

6801-5

Alternate Name: PREG.+THROM.ASSESSMENT

Preferred Specimen: Aliquot Plasma-Citrated-Blue Top,Aliquot Serum- Red Top,Lavender top-EDTA

Collection Instruction: Lavender must be sent separately at room temp. Aliquot samples must be frozen.

Turn Around Time: 5 days

CPT Code: 81240 (1), 81241 (1), 83090 (1), 85025 (1), 85305 (1), 85306 (2), 85598 (1), 85610 (1), 85613 (2), 85730 (2), 86146 (3), 86147 (3)

Scarlet Eligible: No

Pregnenolone, Serum

J150-4

Alternate Name: PREGNENOLONE, SERUM

Methodology: Liquid Chromatography Tandem Mass Spectrometry

Preferred Specimen: 1.0 Aliquot Tube-Serum

Collection Instruction: Place 1-3 mL of serum/plasma in transport tube. Label as serum or plasma and freeze as required.

Storage Instruction: STRICT Frozen

Turn Around Time: 3 days

CPT Code: 84140 (1)

Scarlet Eligible: No

Prenatal + HIV

2765-6

Alternate Name: PRENATAL + HIV

Preferred Specimen: SST Tube,Lavender top- EDTA,Pink Tube-EDTA

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 85025 (1), 86592 (1), 86762 (1), 86850 (1), 86900 (1), 86901 (1), 87340 (1), 87389 (1), G0059 (1)

Scarlet Eligible: Yes

Procainamide/NAPA

0525-6

Alternate Name: PROCAINAMIDE/NAPA

Preferred Specimen: Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80192 (1)

Scarlet Eligible: Yes

PROCESS AND HOLD FOR FISH

B664-4

Alternate Name: PROCESS/HOLD FOR FISH

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Test code is used when physician is waiting for DX on other testing. We will prep the sample and hold it until the client calls to proceed or cancel. If no request receives within 7 days of holding, we will call to confirm if FISH testing is needed.

Turn Around Time: 3 days

CPT Code: 88237 (1)

Scarlet Eligible: No

Progesterone, Serum

0335-0

Alternate Name: PROGESTERONE

Methodology: Electrochemiluminescence Immunoassay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84144 (1)

Scarlet Eligible: Yes

Prolactin, Serum

0134-7

Alternate Name: PROLACTIN, SERUM

Methodology: Electrochemiluminescence Immunoassay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84146 (1)

Scarlet Eligible: Yes

Propoxyphene, Urine, Screen (300)

0963-9

Alternate Name: PROPOXYPHENE, URINE, QUAL.

Methodology: Enzyme Immunoassay

Preferred Specimen: 5 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80307 (1)

Scarlet Eligible: Yes

Prostate Cancer Triple Stain By IHC (PC)

B964-8

Alternate Name: PROSTATE CANCER PC ONLY

Methodology: Immunohistochemistry

Preferred Specimen: Stained Slides

Collection Instruction: Client to collect specimen, prepare smears, label with patient's name, and place in slide carrier.

Storage Instruction: Room Temp

Turn Around Time: 1-2 days

CPT Code: 88344 (1)

Scarlet Eligible: No

Prostate Cancer/ Lymph Nodes

5122-7

Alternate Name: PROSTATE CANCER/LYMPH NODES

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88341 (2), 88342 (1)

Scarlet Eligible: No

Prostate Core Biopsy (1 Spec.)

5611-9

Alternate Name: PROSTATE CORE BIOPSY

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes

Turn Around Time: 1 days

CPT Code: 88305 (1)

Scarlet Eligible: No

Prostate Core Biopsy (8 Spec.)

5618-4

Alternate Name: PROSTATE CORE BIOPSY

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 1 days

CPT Code: 88305 (8)

Scarlet Eligible: No

Prostate Core Biopsy (4 Spec.)

5614-3

Alternate Name: PROSTATE CORE BIOPSY

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 1 days

CPT Code: 88305 (4)

Scarlet Eligible: No

Prostate Core Biopsy (3 Spec.)

5613-5

Alternate Name: PROSTATE CORE BIOPSY

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 1 days

CPT Code: 88305 (3)

Scarlet Eligible: No

Prostate Core Biopsy (20 Spec.)

5174-8

Alternate Name: PROSTATE CORE BIOPSY

Methodology: Immunohistochemistry

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 1 days

CPT Code: 88305 (20)

Scarlet Eligible: No

Prostate Core Biopsy (2 Spec.)

5612-7

Alternate Name: PROSTATE CORE BIOPSY

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 1 days

CPT Code: 88305 (2)

Scarlet Eligible: No

Prostate Core Biopsy (7 Spec.)

5617-6

Alternate Name: PROSTATE CORE BIOPSY

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 1 days

CPT Code: 88305 (7)

Scarlet Eligible: No

Prostate Core Biopsy (5 Spec.)

5615-0

Alternate Name: PROSTATE CORE BIOPSY

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 1 days

CPT Code: 88305 (5)

Scarlet Eligible: No

Prostate Core Biopsy (6 Spec.)

5616-8

Alternate Name: PROSTATE CORE BIOPSY

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 1 days

CPT Code: 88305 (6)

Scarlet Eligible: No

Prostate Core Biopsy (12 Spec.)

5622-6

Alternate Name: PROSTATE CORE BIOPSY 12

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

CPT Code: 88305 (12)

Scarlet Eligible: No

Activated Protein C Resistance Modified

5705-9

Alternate Name: PROT.C RES.MODIFIED

Methodology: Clot Detection

Preferred Specimen: 2 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 4 days

CPT Code: 85307 (1)

Scarlet Eligible: No

Protein C Antigen

5786-9

Alternate Name: PROTEIN C AG.

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 1 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 5 days

CPT Code: 85302 (1)

Scarlet Eligible: No

Activated Protein C Resistance

5704-2

Alternate Name: PROTEIN C RESIS.

Methodology: Clot Detection

Preferred Specimen: 2 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA and patient's name and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 4 days

CPT Code: 85307 (1)

Scarlet Eligible: No

Protein S Activity

5787-7

Alternate Name: PROTEIN S ACT.

Methodology: Colorimetry

Preferred Specimen: 2 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 4 days

CPT Code: 85306 (1)

Scarlet Eligible: No

Protein S Antigen Total

5789-3

Alternate Name: PROTEIN S AG.,TOTAL

Methodology: Colorimetry

Preferred Specimen: 1 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 5 days

CPT Code: 85305 (1)

Scarlet Eligible: No

Protein S Antigen, Free

5788-5

Alternate Name: PROTEIN S AG.FREE

Methodology: Immunoturbidimetric

Preferred Specimen: 1 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 5 days

CPT Code: 85306 (1)

Scarlet Eligible: No

Protein, Urine, 24 Hr.

0136-2

Alternate Name: PROTEIN, 24HR. URINE

Preferred Specimen: 10 mL Urine Container - 24hr

Collection Instruction: Discard first morning sample, then collect all urines after for 24 hrs, on completion label with collection time and patient's name; record total volume on Requisition

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 81050 (1), 84156 (1)

Scarlet Eligible: No

Protein/Creatinine w/Ratio

3329-0

Alternate Name: PROTEIN/CREAT. w/RATIO

Preferred Specimen: Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Turn Around Time: 2 days

CPT Code: 82570 (1), 84156 (1)

Scarlet Eligible: Yes

Prothrombin ABS IGG

5820-6

Alternate Name: PROTHROMBIN ABS.IGG
Methodology: Enzyme Linked Immunoabsorbance
Preferred Specimen: Aliquot Tube-Serum
Collection Instruction: Fill tube, label with patient name
Storage Instruction: Refrigerate
Turn Around Time: 5 days
CPT Code: 86849 (1)
Scarlet Eligible: Yes

COAG/Prothrombin Time

5798-4

Alternate Name: PROTHROMBIN TIME (PT)
Methodology: Clot Detection
Preferred Specimen: Aliquot Plasma-Citrated-Blue Top
Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.
Storage Instruction: STRICT Frozen
Turn Around Time: 4 days
CPT Code: 85610 (1)
Scarlet Eligible: No

PSA Free + Total

2088-3

Alternate Name: PSA FREE/TOTAL
Methodology: Chemiluminescence
Preferred Specimen: 1 mL SST Tube
Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes
Storage Instruction: Refrigerate
Turn Around Time: 1 days
CPT Code: 84153 (1), 84154 (1)
Scarlet Eligible: Yes

PSA Reflex to Free PSA

M067-7

Alternate Name: PSA Reflex to Free PSA

Turn Around Time: 1 days

Scarlet Eligible: Yes

PT + PTT

2237-6

Alternate Name: PT + PTT

Preferred Specimen: Light blue top

Collection Instruction: Fill tube to line at top of tube, invert 8-10 times (DO NOT SHAKE!)

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 85610 (1), 85730 (1)

Scarlet Eligible: Yes

Prothrombin G20210A Mutation

5795-0

Alternate Name: PT G20210A MUTATION

Methodology: Polymerase Chain Reaction

Preferred Specimen: 5 mL Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 81240 (1)

Scarlet Eligible: Yes

Prothrombin Time/INR (PT)

0137-0

Alternate Name: PT/INR

Methodology: Photometric

Preferred Specimen: Light blue top

Collection Instruction: Fill tube to line at top of tube, invert 8-10 times (DO NOT SHAKE!)

Storage Instruction: Refrigerate

Turn Around Time: 8 hours

CPT Code: 85610 (1)

Scarlet Eligible: Yes

PTEN BY FISH

B854-1

Alternate Name: PTEN BY FISH

Methodology: FISH

Preferred Specimen: Unstained Slide

Collection Instruction: Collect specimen and place on charged slide, label with pt name, place in slide holder for transport.

Storage Instruction: Room Temp

Turn Around Time: 3-4 days

CPT Code: 88377 (1)

Scarlet Eligible: No

PTEN BY IHC With Interpretation

A985-4

Alternate Name: PTEN BY IHC- W/INTERP

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88342 (1)

Scarlet Eligible: No

PTH, Intact w/ Calcium

1689-9

Alternate Name: PTH, INTACT

Preferred Specimen: Lavender top- EDTA,SST Tube

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82310 (1), 83970 (1)

Scarlet Eligible: Yes

PTH, Intact (use 1689 if with Calcium)

0598-3

Alternate Name: PTH, INTACT (use 1689), iPTH

Methodology: Chemiluminescence

Preferred Specimen: 1 mL Lavender top- EDTA

Collection Instruction: Centrifuge lavender-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA and patient's name.

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 83970 (1)

Scarlet Eligible: Yes

PTH, Intact, FNA

J047-2

Alternate Name: PTH, INTACT, FNA

Methodology: Chemiluminescence

Preferred Specimen: 1 mL Fluid-source required

Collection Instruction: Collect sample in sterile container.

Storage Instruction: STRICT Frozen

Turn Around Time: 9 days

CPT Code: 83970 (1)

Scarlet Eligible: No

Part.Thromb.Time(PTT)-LA

5771-1

Alternate Name: PTT-LA sens.

Methodology: Clot Detection

Preferred Specimen: 2 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 3 days

CPT Code: 85730 (1)

Scarlet Eligible: No

Panel A (Inhalants)

1632-9

Alternate Name: Panel A (Inhalants)

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82784 (1), 82785 (1), 86003 (15)

Scarlet Eligible: Yes

Pap + Pap Dependent HPV

B975-4

Alternate Name: Pap Dep HPV

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87624 (1), 88175 (1)

Scarlet Eligible: No

Pap + Pap Dependent HPV + CT/GC

P734-7

Alternate Name: Pap Dep HPV+CT/GC

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Turn Around Time: 5 days

CPT Code: 87491 (1), 87591 (1), 87624 (1), 88175 (1)

Scarlet Eligible: No

Pap (Reflex to HPV DNA Genotyping 16,18 when ASC-US or UNSAT)

L326-8

Alternate Name: Pap Reflex to HPV DNA Gen

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 88175 (1)

Scarlet Eligible: No

Pap (Reflex to HPV DNA Genotyping 16,18 when ASC-US)

6195-2

Alternate Name: Pap(rfl HPV DNA ASCUS)

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 88175 (1)

Scarlet Eligible: No

Pap (Reflex to HPV mRNA Genotyping 16, 18/45 when ASC-US)

A251-1

Alternate Name: Pap(rfl HPV mRNA ASCUS)

Methodology: Transcription Mediated Amplification

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 88175 (1)

Scarlet Eligible: No

Pap and HPV DNA Genotyping 16,18

P079-7

Alternate Name: Pap+HPV DNA

Preferred Specimen: ThinPrep Vial

Turn Around Time: 5 days

CPT Code: 87624 (1), 88175 (1)

Scarlet Eligible: No

Pap and HPV mRNA Genotyping 16, 18/45

A250-3

Alternate Name: Pap+HPV mRNA

Methodology: Transcription Mediated Amplification

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87624 (1), 88175 (1)

Scarlet Eligible: No

PAP/CT/GC(RFLX HPV PLUS ASC-US or >)

9985-3

Alternate Name: Pap/Ct/Gc(Rflx Hpv Hi Asc-Us Or >)

Methodology: Polymerase Chain Reaction

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Turn Around Time: 5 days

CPT Code: 87491 (1), 87591 (1), 88175 (1)

Scarlet Eligible: No

Allergen Papaya (f293), IgE

1942-2

Alternate Name: Papaya

Methodology: Immunocap

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Paper Wasp (i4), IgE

0886-2

Alternate Name: Paper Wasp

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Paprika (f218), IgE

1424-1

Alternate Name: Paprika

Methodology: Immunocap

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Ova and Parasite (O&P), Stool W/ Trichrome (1st specimen)

0377-2

Alternate Name: Parasites, Stool; Parasitology Examination, Stool; Stool Exam for Ova and Parasites

Methodology: Microscopic Examination

Preferred Specimen: Ova & Parasite Mini

Collection Instruction: With applicator stick, place 10 gms of stool into each vial, label with patient name

Storage Instruction: Room Temp

Turn Around Time: 2 days

CPT Code: 87177 (1), 87209 (1)

Scarlet Eligible: No

Allergen Parsley (f86), IgE

1355-7

Alternate Name: Parsley

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Parvovirus B19 Ab, IgG

T790-5

Alternate Name: Parvovirus B19 Ab, IgG

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86747 (1)

Scarlet Eligible: Yes

Parvovirus B19 Ab, IgM

T791-3

Alternate Name: Parvovirus B19 Ab, IgM

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86747 (1)

Scarlet Eligible: Yes

Parvovirus B19 Ab (IgG/IgM)

L399-5

Alternate Name: Parvovirus B19 G/M

Preferred Specimen: SST Tube

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86747 (1)

Scarlet Eligible: Yes

Allergen Passion Fruit (rf294), IgE

3480-1

Alternate Name: Passion Fruit (Rf294)**Methodology:** Immunocap**Preferred Specimen:** 2 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 3 days**CPT Code:** 86003 (1)**Scarlet Eligible:** Yes**Lead, Blood (Adult)**

1438-1

Alternate Name: Pb, Blood (Adult)**Methodology:** Inductively Coupled Plasma Mass Spectrometry**Preferred Specimen:** Tan Top**Collection Instruction:** Fill tube, invert gently 5-6 times, label with patient name**Storage Instruction:** Refrigerate**Turn Around Time:** 2 days**CPT Code:** 83655 (1)**Scarlet Eligible:** Yes**Lead, Blood (Child)**

0398-8

Alternate Name: Pb, Blood (Child)**Methodology:** Inductively Coupled Plasma Mass Spectrometry**Preferred Specimen:** Tan Top**Collection Instruction:** Fill tube, invert gently 5-6 times, label with patient name**Storage Instruction:** Refrigerate**Turn Around Time:** 2 days**CPT Code:** 83655 (1)**Scarlet Eligible:** Yes

Allergen Peach (f95), IgE

1226-0

Alternate Name: Peach (F95)

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Peanut (f13),IgE

0674-2

Alternate Name: Peanut

Methodology: Immunocap

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Peanut Component Ara (h1), IgE

J284-1

Alternate Name: Peanut Ara (h1)

Methodology: Immunocap

Preferred Specimen: 1.0 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86008 (1)

Scarlet Eligible: Yes

Peanut Component Ara (h2), IgE

J283-3

Alternate Name: Peanut Ara (h2)**Methodology:** Immunocap**Preferred Specimen:** 1.0 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 3 days**CPT Code:** 86008 (1)**Scarlet Eligible:** Yes**Peanut Component Ara (h3), IgE**

J285-8

Alternate Name: Peanut Ara (h3)**Methodology:** Immunocap**Preferred Specimen:** 1 SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 3 days**CPT Code:** 86008 (1)**Scarlet Eligible:** Yes**Peanut Component Ara (h6), IgE**

TG35-1

Alternate Name: Peanut Ara (h6)**Methodology:** Immunocap**Preferred Specimen:** 1 mL SST Tube**Storage Instruction:** Refrigerate**Turn Around Time:** 3 days**CPT Code:** 86008 (1)**Scarlet Eligible:** Yes

Peanut Component Ara (h8), IgE

J287-4

Alternate Name: Peanut Ara (h8)

Methodology: Immunocap

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86008 (1)

Scarlet Eligible: Yes

Peanut Component Ara (h9), IgE

J286-6

Alternate Name: Peanut Ara (h9)

Methodology: Immunocap

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86008 (1)

Scarlet Eligible: Yes

Peanut Components IgE

H993-9

Alternate Name: Peanut Components IgE

Preferred Specimen: 1 mL SST Tube

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82785 (1), 86003 (5)

Scarlet Eligible: Yes

Peanut w/Reflex to Component Panel

M065-1

Alternate Name: Peanut w/Reflex to Component Panel

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Pear (f94), IgE

0930-8

Alternate Name: Pear (F94)

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Pecan Nuts (f201), IgE

1592-5

Alternate Name: Pecan

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Pecan/Hickory Tree (t22), IgE

1938-0

Alternate Name: Pecan/Hickory

Methodology: Immunocap

Preferred Specimen: 3.0 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Pediatric Mix (fx5), IgE

2186-5

Alternate Name: Pediatric Mix

Methodology: Immunocap

Preferred Specimen: 5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86005 (1)

Scarlet Eligible: Yes

Allergen Penicillium Notatum Mold (m1), IgE

2167-5

Alternate Name: Penicillium Mold

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Perennial Rye Grass (g5),IgE0682-5

Alternate Name: Perennial Rye Grass**Methodology:** Immunocap**Preferred Specimen:** 3 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 3 days**CPT Code:** 86003 (1)**Scarlet Eligible:** Yes

Perimenopausal Hormonal Profile8431-9

Alternate Name: Perimenopausal Hormonal**Preferred Specimen:** SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 1 days**CPT Code:** 82670 (1), 83001 (1)**Scarlet Eligible:** Yes

Peroxisome Biogenesis Disorder 5A (Zellweger)J303-9

Alternate Name: Peroxisome Biogenesis Disorder 5A (Zellweger)**Methodology:** Genotyping by Next Generation Sequencing**Preferred Specimen:** Lavender top- EDTA**Collection Instruction:** Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!**Storage Instruction:** Refrigerate**Turn Around Time:** 10-14 days**CPT Code:** 81479 (1)**Scarlet Eligible:** No

Culture, Throat, Routine

0078-6

Alternate Name: Pharyngeal Culture

Methodology: Bacterial Culture

Preferred Specimen: Swab-E

Collection Instruction: Collect specimen with swab then place swab into transport carrier; Label with name and source

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87070 (1)

Scarlet Eligible: No

Phenylketonuria (PKU) Carrier Screening

J304-7

Alternate Name: Phenylketonuria (PKU) Carrier Screening

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81406 (1)

Scarlet Eligible: No

Dilantin

0084-4

Alternate Name: Phenytoin, Serum

Methodology: Enzyme Immunoassay

Preferred Specimen: Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80185 (1)

Scarlet Eligible: Yes

BCR/ABL: t(9;22) by FISH

5265-4

Alternate Name: Philadelphia Chromosome Analysis

Methodology: FISH

Preferred Specimen: 5 mL Bone Marrow - Green Top

Collection Instruction: Place bone marrow (2 ml req.) into tube, label with patient's name. Send in bone marrow kit

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88374 (1)

Scarlet Eligible: No

Allergen Phoma Betae (m13), IgE

2477-8

Alternate Name: Phoma Betae

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Phosphoglycerate dehydrogenase deficiency

J302-1

Alternate Name: Phosphoglycerate dehydrogenase deficiency

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81479 (1)

Scarlet Eligible: Yes

Allergen Pignolia (Pine) Nuts (t213), IgE

3638-4

Alternate Name: Pignolia Nuts

Methodology: Immunocap

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergens Pineapple (f210), IgE

1389-6

Alternate Name: Pineapple

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Pistachio (rf203), IgE

1195-7

Alternate Name: Pistachio (Rf203)

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Plasma 10-Color (Short) Panel

TK94

Alternate Name: Plasma 10-Color (Short) Panel

Collection Instruction: Not orderable as of 2/1 but launching within the month of Feb

Scarlet Eligible: No

Plasma 10-Color (Short) Panel TC Only

TK95

Alternate Name: Plasma 10-Color (Short) Panel TC Only

Collection Instruction: Not orderable as of 2/1 but launching within the month of Feb

Scarlet Eligible: No

Plasma 10-Color Myeloma/Monoclonal Gammopathy

TK97

Alternate Name: Plasma 10-Color Myeloma/Monoclonal Gammopathy

Collection Instruction: Not orderable as of 2/1 but launching within the month of Feb

Scarlet Eligible: No

Plasma 10-Color Myeloma/Monoclonal Gammopathy TC Only

TK96

Alternate Name: Plasma 10-Color Myeloma/Monoclonal Gammopathy TC Only

Collection Instruction: Not orderable as of 2/1 but launching within the month of Feb

Scarlet Eligible: No

Glucose, Non-Fasting

0223-8

Alternate Name: Plasma glucose

Methodology: Enzyme Immunoassay

Preferred Specimen: Grey top- Sodium Fluoride

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82947 (1)

Scarlet Eligible: Yes

Heparin Platelet Antibody

5746-3

Alternate Name: Platelet factor IV, PF4
Methodology: Enzyme Linked Immunoabsorbance
Preferred Specimen: Aliquot Serum- Red Top
Collection Instruction: Fill tube, label with patient name
Turn Around Time: 3 days
CPT Code: 86022 (1)
Scarlet Eligible: Yes

Allergen Plum (f255), IgE

1584-2

Alternate Name: Plum
Methodology: Immunocap
Preferred Specimen: 3 mL SST Tube
Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes
Storage Instruction: Refrigerate
Turn Around Time: 3 days
CPT Code: 86003 (1)
Scarlet Eligible: Yes

Pneumococcal Antibody Panel (14 Serotypes)

1250-0

Alternate Name: Pneumococcal Abs. IGG (14);Strep pneumoniae
Methodology: Immunoassay
Preferred Specimen: 2 mL SST Tube
Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes
Storage Instruction: Refrigerate
Turn Around Time: 7 days
CPT Code: 86317 (14)
Scarlet Eligible: Yes

Polycystic Kidney Disease, Autosomal Recessive (ARPKD)

J301-3

Alternate Name: Polycystic Kidney Disease Autosomal Recessive (ARPKD)

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81408 (1)

Scarlet Eligible: Yes

Allergen Poplar White IgE

3704-4

Alternate Name: Poplar White IgE

Methodology: Enzyme Immunoassay

Preferred Specimen: 0.5 mL Aliquot Tube-Serum

Collection Instruction: Collect 1-2 mL whole blood in red top tube, centrifuge and transfer 0.5 mL serum into a transfer tube. (All serum acceptable)

Storage Instruction: Refrigerate

Turn Around Time: 7 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Poppy Seed (rf224), IgE

1031-4

Alternate Name: Poppy Seed (Rf224)

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Pork (f26), IgE

0836-7

Alternate Name: Pork

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Potassium, RBC

J571-1

Alternate Name: Potassium, RBC

Methodology: Inductively Coupled Plasma Mass Spectrometry

Preferred Specimen: 2 mL Green top- Sodium Heparin

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 84132 (1)

Scarlet Eligible: Yes

Allergen Potato (f35), IgE

0844-1

Alternate Name: Potato

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Obstetric (AMA) Panel

0010-9

Alternate Name: Prenatal Panel

Preferred Specimen: SST Tube,Lavender top- EDTA,Pink Tube-EDTA

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 85025 (1), 86592 (1), 86762 (1), 86850 (1), 86900 (1), 86901 (1), 87340 (1)

Scarlet Eligible: Yes

Ciliary Dyskinesia, Primary 1

J311-2

Alternate Name: Primary Ciliary Dyskinesia DNAI1-associated

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81479 (1)

Scarlet Eligible: No

Allergen Privet Tree (T210), IgE

A321-2

Alternate Name: Privet Tree (T210)

Methodology: Immunocap

Preferred Specimen: 2.0 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Process and Hold for Cytogenetics

J655-2

Alternate Name: Process/Hold >Cytog

CPT Code: 99999 (1)

Scarlet Eligible: No

Tacrolimus (FK506)

1871-3

Alternate Name: Prograf

Methodology: Enzyme Immunoassay

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Submit 2 lavenders when Tacrolimus is ordered in conjunction with any other test requiring lavender.

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80197 (1)

Scarlet Eligible: Yes

Prostate Cancer (Stain Only)

5104-5

Alternate Name: Prostate Cancer (Stain Only)

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88344 (1)

Scarlet Eligible: No

Hereditary Prostate Cancer Panel

J665-1

Alternate Name: Prostate Cancer Panel**Methodology:** Genotyping by Next Generation Sequencing**Preferred Specimen:** 2-5 mL Lavender top- EDTA**Collection Instruction:** Lavender: 2 tubes (5-6 mL each) Label with the patient name, date of birth and/or ID. Ideal handling is to ship immediately at ambient temperature for overnight delivery. Alternate specimen type: DNA: Extract and submit DNA 50 ng/mL (total volume 400 mL) in sterile container. Ship DNA in sealed and protective packaging at ambient temperature, using a cool pack in hot weather. Labeling: Name of person or identifier, DNA concentration and volume MUST be provided on the specimen tube.**Storage Instruction:** Room Temp**Turn Around Time:** 14 days**CPT Code:** 81405 (1), 81162 (1)**Scarlet Eligible:** Yes

Prostate Carcinoma IHC Panel (CK903, p63, P504S)

5097-1

Alternate Name: Prostate Cancer Profile**Methodology:** Immunohistochemistry**Preferred Specimen:** Formalin-fixed, Paraffin-embedded Tissue**Collection Instruction:** This comes in block or slide form with surgical number imprint.**Storage Instruction:** Refrigerate**Turn Around Time:** 1 days**CPT Code:** 88344 (1)**Scarlet Eligible:** No

PSA Total

0190-9

Alternate Name: Prostate Specific Antigen (PSA)

Methodology: Electrochemiluminescence Immunoassay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84153 (1)

Scarlet Eligible: Yes

Protein C Activity

5785-1

Alternate Name: Protein C functional

Methodology: Chromogenic

Preferred Specimen: 3 mL Aliquot Plasma-Citrate-Blue Top

Collection Instruction: Collect full blue-top (citrate) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 4 days

CPT Code: 85303 (1)

Scarlet Eligible: No

Protein Electrophoresis, Urine, Random

0404-4

Alternate Name: Protein Electrophoresis, Urine

Methodology: Electrophoresis

Preferred Specimen: 15 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 84156 (1), 84166 (1)

Scarlet Eligible: Yes

Protein,Urine,Timed or Random

1662-6

Alternate Name: Protein,Urn.Timed/Rand

Methodology: Colorimetry

Preferred Specimen: 10 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 84156 (1)

Scarlet Eligible: Yes

Allergen Pumpkin (f225), IgE

1488-6

Alternate Name: Pumpkin

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Pumpkin Seed (rf226), IgE

1484-5

Alternate Name: Pumpkin Seed

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Pyruvic Acid (Pyruvate), Blood

3767-1

Alternate Name: Pyruvic Acid (Pyruvate)

Methodology: Enzymatic

Preferred Specimen: 5 mL PYR- Pyruvate tube

Collection Instruction: Collect 4 mL random whole blood in an EDTA (lavender-top), sodium or lithium heparin (green-top), or fluoride oxalate (gray-top) tube and immediately mix with 4 mL ice cold 7% or 8% Perchloric Acid. Let mixture stand for 10 minutes; then, centrifuge. Separate and submit 4.5 mL supernatant fluid for assay, refrigerated 2-8 degrees C. Transport refrigerated (cold packs). Please specify on the request form the volumes of blood and Perchloric Acid used.

Storage Instruction: Refrigerate

Turn Around Time: 10 days

CPT Code: 84210 (1)

Scarlet Eligible: No

Methaqualone, Urine, Screen (300)

3239-1

Alternate Name: Quaalude

Methodology: Enzyme Multiplied Immunoassay

Preferred Specimen: 5 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80307 (1)

Scarlet Eligible: Yes

Alpha-Fetoprotein (AFP 4) includes HCG, UE3, Inhibin-A

3158-3

Alternate Name: Quad Screen (AFP 4)

Methodology: Chemiluminescence

Preferred Specimen: 3 mL SST Tube

Collection Instruction: MUST BE 15-22 WEEKS AND 6 DAYS. PROVIDE GESTATIONAL INFO.

Storage Instruction: Refrigerate

Turn Around Time: 4 days

CPT Code: 81511 (1)

Scarlet Eligible: Yes

Pregnancy, Serum (Qualitative)

T824-2

Alternate Name: Qualitative HCG

Methodology: Immunoassay

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84703 (1)

Scarlet Eligible: Yes

Allergen R. Nigricans (m11), IgE

1373-0

Alternate Name: R. Nigricans

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Rf (Rheumatoid Arthritis) Titer

0796-3

Alternate Name: RA Latex; RA Latex, Turbidimetry; RF Titer; Rheumatoid Factor (RF)

Methodology: Immunosubimetric

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86431 (1)

Scarlet Eligible: Yes

Allergen Red Raspberry (rf343), IgE

3508-9

Alternate Name: RASPBERRY RED

Methodology: Immunocap

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

RB1 for 13q14 by FISH

B511

Alternate Name: RB1 for 13q14 by FISH

Scarlet Eligible: No

RBC

1498-5

Alternate Name: RBC

Methodology: Flow Cytometry

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 85041 (1)

Scarlet Eligible: Yes

Reactive Hyperplasia vs. Lymphoma (NHL) IHC Panel

5103-7

Alternate Name: REACT.HYPERPLASIAV.LYMPHOMA

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88341 (12), 88342 (1)

Scarlet Eligible: No

Renal Function (Medicare)

3284-7

Alternate Name: RENAL FUNCTION (MEDICARE)

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80069 (1)

Scarlet Eligible: Yes

Reptilase Time

5799-2

Alternate Name: REPTILASE TIME

Methodology: Clot Detection

Preferred Specimen: 1 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 4 days

CPT Code: 85635 (1)

Scarlet Eligible: No

Rheumatoid Factor IgA

TA08-0

Alternate Name: RF IgA

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 83516 (1)

Scarlet Eligible: Yes

Rheumatoid Factor IgM

TA09-8

Alternate Name: RF IgM

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 83516 (1)

Scarlet Eligible: Yes

ROMA

L012-4

Alternate Name: ROMA

Methodology: Electrochemiluminescence Immunoassay

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 81500 (1)

Scarlet Eligible: Yes

ROS-1 BY FISH-TECH Only

B337-7

Alternate Name: ROS-1 BY FISH-TECH ONLY

Methodology: FISH

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: Please write surgical number on block and on requisition. Provide fixation time on requisition as well

Storage Instruction: Room Temp

Turn Around Time: 3-5 days

CPT Code: 88374 (1)

Scarlet Eligible: No

ROS1 BY FISH Global

B334-4

Alternate Name: ROS1 BY FISH GLOBAL

Methodology: FISH

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block from client with surgical number imprint.

Storage Instruction: Room Temp

Turn Around Time: 3-5 days

CPT Code: 88377 (1)

Scarlet Eligible: No

POST TREATMENT SYPHILIS (RPR)

R733-7

Alternate Name: RPR

Turn Around Time: 1 days

CPT Code: 86592 (1)

Scarlet Eligible: Yes

RPR Serology

0142-0

Alternate Name: RPR Serology, Syphilis,Rapid Plasma Reagin Test for Syphilis

Methodology: Flocculation

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86592 (1)

Scarlet Eligible: Yes

Resp. Syncytial Virus

1397-9

Alternate Name: RSV

Methodology: Polymerase Chain Reaction

Preferred Specimen: Swab-Viral Culturette

Collection Instruction: Collect: Respiratory specimen: Nasopharyngeal washes, aspirates and swab specimens in M4 Transport Media or UTM. Alternately, an ESwab can be used. Specimen Preparation: Transfer the collected respiratory sample in 1 mL of viral transport media (Min: 0.5 mL)

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87634 (1)

Scarlet Eligible: No

Rubella Ab, IgG, Qual, Immune Status

3151-8

Alternate Name: RUBELLA IgG Ab(IMMUNITY)

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86762 (1)

Scarlet Eligible: Yes

Sirolimus, Blood

3481-9

Alternate Name: Rapamycin

Methodology: Enzyme Immunoassay

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Submit 2 lavenders when Sirolimus is ordered in conjunction with any other test requiring lavender.

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80195 (1)

Scarlet Eligible: Yes

Allergen Redtop (Bent) Grass (g9), IgE

1199-9

Alternate Name: Redtop (Bent) Grass

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Rest of common cancer panel after Lynch Syndrome Panel B985

B984-6

Alternate Name: Reflex To Remainder of Cancer Panel (After B985)

Preferred Specimen: Peripheral Blood - Lavender Top

Collection Instruction: 2 tubes (5-6mL) Label with the patient name, date of birth and/or ID.

Storage Instruction: Refrigerate

Turn Around Time: 14 days

CPT Code: 81405 (1), 81307 (1)

Scarlet Eligible: Yes

Rest of Comprehensive Cancer Panel

B363-3

Alternate Name: Reflex to Comprehensive Cancer Panel

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: 2-5 mL Lavender top- EDTA

Collection Instruction: Lavender: 2 tubes (5-6ml each) Label with the patient name, date of birth and/or ID. Ideal handling is to ship immediately at ambient temperature for overnight delivery. Alternate sample type: DNA: Extract and submit DNA 50 ng/mL (total volume 400 mL) in sterile container. Ship DNA in sealed and protective packaging at ambient temperature, using a cool pack in hot weather. Labeling: Name of person or identifier, DNA concentration and volume MUST be provided on the specimen tube.

Storage Instruction: Room Temp

Turn Around Time: 14 days

CPT Code: 81479 (1)

Scarlet Eligible: Yes

Rest of Comprehensive Cancer Panel (after Prostate Panel)

J664-4

Alternate Name: Reflex to Remainder of Cancer Panel (after J665))

Methodology: Full Gene Sequencing

Preferred Specimen: 5 mL Peripheral Blood - Lavender Top

Collection Instruction: Lavender: 2 tubes (5-6 mL each) Label with the patient name, date of birth and/or ID. Ideal handling is to ship immediately at ambient temperature for overnight delivery. Alternate specimen type: DNA: Extract and submit DNA 50 ng/mL (total volume 400 mL) in sterile container. Ship DNA in sealed and protective packaging at ambient temperature, using a cool pack in hot weather. Labeling: Name of person or indentifer, DNA concentration and volume MUST be provided on the specimen tube.

Storage Instruction: Room Temp

Turn Around Time: 14 days

CPT Code: 81408 (1)

Scarlet Eligible: Yes

Reflex to Remainder of Common Cancer Panel (after J055)

J054-8

Alternate Name: Reflex to Remainder of Cancer Panel(after J055),Rest of Cancer P

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: 2-5 mL Lavender top- EDTA

Collection Instruction: 2 tubes (5-6mL each) Label with the patient name, date of birth, and/or ID.

Storage Instruction: Room Temp

Turn Around Time: 8-10 days

CPT Code: 81404 (1), 81406 (2)

Scarlet Eligible: Yes

Region 1 Allergy Profile (North East)

6823-9

Alternate Name: Region 1 Allergy Prof.

Preferred Specimen: 2 mL SST Tube

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 82785 (1), 86003 (26)

Scarlet Eligible: Yes

Region 10 Allergy Profile (TX,OK)

6827-0

Alternate Name: Region 10 Allergy Prof.

Preferred Specimen: SST Tube

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 82785 (1), 86003 (26)

Scarlet Eligible: Yes

Region 11 Allergy Profile (AZ, ID, NM, WY, CO, UT)

K908-4

Alternate Name: Region 11 Allergy Prof.

Turn Around Time: 5 days

CPT Code: 86003 (25), 82785 (1)

Scarlet Eligible: Yes

Region 13 Allergy Profile (S. CA Coast)

6820-5

Alternate Name: Region 13 Allergy Prof.

Preferred Specimen: SST Tube

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 82785 (1), 86003 (25)

Scarlet Eligible: Yes

Region 14 Allergy Profile (Central CA)

6821-3

Alternate Name: Region 14 Allergy Prof.

Preferred Specimen: SST Tube

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 82785 (1), 86003 (24)

Scarlet Eligible: Yes

Region 17 Allergy Profile (North West)

6822-1

Alternate Name: Region 17 Allergy Prof.

Preferred Specimen: SST Tube

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82785 (1), 86003 (24)

Scarlet Eligible: Yes

Region 2 Allergy Profile (Mid Atlantic)

6824-7

Alternate Name: Region 2 Allergy Prof.

Preferred Specimen: SST Tube

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 82785 (1), 86003 (24)

Scarlet Eligible: Yes

Region 3 Allergy Profile (South East)

6825-4

Alternate Name: Region 3 Allergy Prof.

Preferred Specimen: SST Tube

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82785 (1), 86003 (23)

Scarlet Eligible: Yes

Region 4 Allergy Profile (S.Florida)

6826-2

Alternate Name: Region 4 Allergy Prof.

Preferred Specimen: SST Tube

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82785 (1), 86003 (23)

Scarlet Eligible: Yes

Region 5 Allergy Profile (IN, OH, TN, WV, KY)

6829-6

Alternate Name: Region 5 Allergy Profile

Turn Around Time: 5 days

CPT Code: 82785 (1), 86003 (27)

Scarlet Eligible: Yes

Renal Cancer Panel

B394

Alternate Name: Renal Cancer Panel

Scarlet Eligible: No

Renal Cancer Panel

B394

Alternate Name: Renal Cancer Panel

Scarlet Eligible: No

Renal Cancer Panel

B394

Alternate Name: Renal Cancer Panel

Scarlet Eligible: No

Rest of Demyelinating CMT Panel (after PMP22 Del/Dup)

J795-6

Alternate Name: Rest of Demyelinating CMT Panel

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 28 days

CPT Code: 81404 (1), 81405 (2), 81406 (1)

Scarlet Eligible: Yes

Reticulocyte Count

0141-2

Alternate Name: Retic Count

Methodology: Fluorescence

Preferred Specimen: 4 mL Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 85045 (1)

Scarlet Eligible: Yes

RET-HE

5848-7

Alternate Name: Reticulocyte Hemoglobin Equivalent

Methodology: Calculated

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 85046 (1)

Scarlet Eligible: Yes

Retinisis Pigmentosa 59

J313-8

Alternate Name: Retinitis pigmentosa 59

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81479 (1)

Scarlet Eligible: Yes

Vitamin A (Retinol), Serum

0337-6

Alternate Name: Retinol

Methodology: High Pressure Liquid Chromatography

Preferred Specimen: 2 mL Amber Serum

Collection Instruction: Draw SST or RED, let clot for 8-10 minutes. Centrifuge and pour into either an aliquot tube (must be covered in foil) or poured into an amber tube and label with patient name and specimen type (serum).

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 84590 (1)

Scarlet Eligible: Yes

Reflex to GenCerv

A780-9

Alternate Name: Rfx GenCerv

Methodology: Flow Cytometry

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Turn Around Time: 5 days

CPT Code: 88184 (1), 88365 (1)

Scarlet Eligible: No

Rheumatoid Factor IgG, IgA, IgM

L626-1

Alternate Name: Rheumatoid Factor

Turn Around Time: 3 days

CPT Code: 83516 (1), 86431 (1)

Scarlet Eligible: Yes

Vitamin B2

3325-8

Alternate Name: Riboflavin

Methodology: High Pressure Liquid Chromatography

Preferred Specimen: 2 mL Amber Serum

Collection Instruction: Preferred sample type plasma or serum , In either an aliquot tube (must be covered in foil) or poured into an amber tube and label with patient name and specimen type (plasma from Lavender), serum can be from a Red or SST.

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 84252 (1)

Scarlet Eligible: Yes

Allergen Rice (f9), IgE

0843-3

Alternate Name: Rice

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Rough Careless Pigweed (w14), IgE

1627-9

Alternate Name: Rough Careless Pigweed

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Rubella Ab IgG W/Reflex

T385-4

Alternate Name: Rubella Ab (IgG) Rfx

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86762 (1)

Scarlet Eligible: Yes

RUBELLA IGM

J900-2

Alternate Name: Rubella IgM

Methodology: Chemiluminescence

Preferred Specimen: 1.0 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 24 hours

CPT Code: 86762 (1)

Scarlet Eligible: Yes

Rubella Ab, IgG, Quant

0973-8

Alternate Name: Rubella Serology

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86762 (1)

Scarlet Eligible: Yes

Allergen Russian Thistle (w11), IgE

0732-8

Alternate Name: Russian Thistle

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Rye (food) (f5), IgE

0822-7

Alternate Name: Rye (Food) (F5)

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Hemoglobin Fractionation, Capillary Electrophoresis

0216-2

Alternate Name: S Solubility;Sickle Cell Conf.,Hemoglobin Electrophoresis

Methodology: Electrophoresis

Preferred Specimen: 2 mL Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 83021 (1)

Scarlet Eligible: Yes

Staphylococcus aureus (Non NY)

J970-5

Alternate Name: S. AUREUS

Methodology: Real-Time Polymerase Chain Reaction

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87640 (1)

Scarlet Eligible: No

S. Cerevisiae Antibodies (IGG/IGA)

H940-0

Alternate Name: SACCHAROMYCES, ASCA

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86671 (2)

Scarlet Eligible: Yes

Allergen Sardine, IgG (NON-NY)

B424-3

Alternate Name: SARDINE IGG

Methodology: Enzyme Immunoassay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 7 days

CPT Code: 86001 (1), 83520

Scarlet Eligible: Yes

Scleroderma Profile

0593-4

Alternate Name: SCLERODERMA PROFILE

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86038 (2), 86235 (1)

Scarlet Eligible: Yes

Selenium Serum/Plasma

0232-9

Alternate Name: SELENIUM,SER/PLASMA

Methodology: Inductively Coupled Plasma Mass Spectrometry

Preferred Specimen: 2 mL Blue-Royal Blue Trace Element Serum (**RED LINE**)

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1-2 days

CPT Code: 84255 (1)

Scarlet Eligible: Yes

Send Copy Of Results To:

9991-1

Alternate Name: SEND COPY OF RESULTS TO:

CPT Code: 99999 (1)

Scarlet Eligible: Yes

Sequential Screen Part 1 (NT, Total BHCG, PAPP-A)

8602-5

Alternate Name: SEQ. SCREEN PART 1

Methodology: Chemiluminescence

Preferred Specimen: 2 mL SST Tube

Collection Instruction: MUST BE 11-13 WEEKS AND 6 DAYS. PROVIDE GESTATIONAL INFO.

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 81508 (1)

Scarlet Eligible: Yes

Sequential Screen Part 2 (AFP, UE3, HCG, Inhibin-A)

8601-7

Alternate Name: SEQ. SCREEN PART 2

Methodology: Chemiluminescence

Preferred Specimen: 3 mL SST Tube

Collection Instruction: MUST BE 15-22 WEEKS AND 6 DAYS. PROVIDE GESTATIONAL INFO.

Storage Instruction: Refrigerate

Turn Around Time: 4 days

CPT Code: 81511 (1)

Scarlet Eligible: Yes

Serum Integrated Part 1 (PAPP-A)

A335-2

Alternate Name: SERUM INTEGRATED PART 1

Methodology: Chemiluminescence

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 4 days

CPT Code: 99999 (1)

Scarlet Eligible: Yes

Serum Integrated Part 2 (AFP, UE3, HCG, Inhibin-A)

A334-5

Alternate Name: SERUM INTEGRATED PART 2

Methodology: Chemiluminescence

Preferred Specimen: 3.0 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 4 days

CPT Code: 81511 (1), 84163 (1)

Scarlet Eligible: Yes

Aspartate Aminotransferase (AST/SGOT)

0146-1

Alternate Name: SGOT (AST)

Methodology: Colorimetry

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84450 (1)

Scarlet Eligible: Yes

SGPT (ALT)

0147-9

Alternate Name: SGPT (ALT)

Methodology: Colorimetry

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84460 (1)

Scarlet Eligible: Yes

Sex-Hormone Binding Globulin

0658-5

Alternate Name: SHBG

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 84270 (1)

Scarlet Eligible: Yes

Silicate, Water

6219-0

Alternate Name: SILICATE, WATER

Preferred Specimen: 25 mL Cup Plain

Collection Instruction: Place sample in sterile cup label with patient's name and source

Storage Instruction: Refrigerate

Turn Around Time: 10 days

CPT Code: 99999 (1)

Scarlet Eligible: No

Sinus Culture W/ Gram Stain

F245-6

Alternate Name: SINUS

Methodology: Bacterial Culture

Preferred Specimen: Swab-E

Collection Instruction: Submit sinus drainage in sterile cup.

Turn Around Time: 5 days

CPT Code: 87070 (1), 87075 (1), 87205 (1)

Scarlet Eligible: No

Skin Biopsy, (2 Jars)

5603-6

Alternate Name: SKIN BIOPSY, 2 JARS

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes

Turn Around Time: 2-3 days

CPT Code: 88305 (2)

Scarlet Eligible: No

Skin Biopsy, (3 Jars)

5604-4

Alternate Name: SKIN BIOPSY, 3 JARS

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes

Turn Around Time: 2-3 days

CPT Code: 88305 (3)

Scarlet Eligible: No

Skin Biopsy, (4 Jars)

5605-1

Alternate Name: SKIN BIOPSY, 4 JARS

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes

Turn Around Time: 2-3 days

CPT Code: 88305 (4)

Scarlet Eligible: No

Skin Biopsy, (5 Jars)

5606-9

Alternate Name: SKIN BIOPSY, 5 JARS

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes

Turn Around Time: 2-3 days

CPT Code: 88305 (5)

Scarlet Eligible: No

Skin Biopsy, (6 Jars)

5607-7

Alternate Name: SKIN BIOPSY, 6 JARS

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes

Turn Around Time: 2-3 days

CPT Code: 88305 (6)

Scarlet Eligible: No

Spinal Muscular Atrophy (SMA) Carrier (SMN1 Dosage)

6380-0

Alternate Name: SMA

Methodology: Reverse transcription polymerase chain reaction

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 5-7 days

CPT Code: 81329 (1)

Scarlet Eligible: Yes

Chem 24 (SMAC)

0005-9

Alternate Name: SMAC

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: G0059 (1)

Scarlet Eligible: Yes

Small Round Cell Tumor IHC Panel

5177-1

Alternate Name: SMALL BLUE ROUND TUMORS, IHC

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88341 (12), 88342 (1)

Scarlet Eligible: No

Smith-Lemli-Opitz Synd.Prenat (Non-NY)

B139-7

Alternate Name: SMITH-LEMLI-OPTIZ SYND.

Preferred Specimen: 5 mL Amniotic Fluid

Collection Instruction: 3mL Amniotic Fluid submitted in a sterile, leak-proof container

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 82542 (1)

Scarlet Eligible: No

SOX-10 BY IHC W/Interp

J147-0

Alternate Name: SOX-10 BY IHC W/INTERP

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block or slide form with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88342 (1)

Scarlet Eligible: No

SOX-10 BY IHC TECH Only

J146-2

Alternate Name: SOX-10 BY IHC-TC

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block or slide form with surgical number print.

Storage Instruction: Room Temp

Turn Around Time: 1 days

CPT Code: 88342 (1)

Scarlet Eligible: No

Special Stain Enzyme Constituents - PC

B756-8

Alternate Name: SPEC STAIN ENZY CON-PC

Methodology: Histology

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Storage Instruction: Room Temp

Turn Around Time: 1 days

CPT Code: 88319 (1)

Scarlet Eligible: No

Special Stain Enzyme Constituents -TC

B758-4

Alternate Name: SPEC STAIN ENZY CON-TC

Methodology: Histology

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Storage Instruction: Room Temp

Turn Around Time: 1 days

CPT Code: 88319 (1)

Scarlet Eligible: No

Special Stain Enzyme Constituents - Global

B757-6

Alternate Name: SPEC STN ENZ CON-GLOBAL

Methodology: Histology

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Storage Instruction: Room Temp

Turn Around Time: 1 days

CPT Code: 88319 (1)

Scarlet Eligible: No

Specific Factor Inhibitor

5818-0

Alternate Name: SPECIFIC FACTOR INHIBITOR

Preferred Specimen: Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 4 days

CPT Code: 85240 (1), 85335 (1), 85610 (1), 85611 (8), 85730 (2), 85732 (8)

Scarlet Eligible: No

Sperm Count (Post Vasectomy)

0145-3

Alternate Name: SPERM CT.POST VASECTOMY

Methodology: Microscopic Examination

Preferred Specimen: Semen Cup

Collection Instruction: For FERTILITY must be scheduled. For POST VASECTOMY, sample can be dropped off at any PSC for testing, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 89310 (1)

Scarlet Eligible: No

Spindle Cell Neoplasm IHC Panel

5191-2

Alternate Name: SPINDLE CELL NEOPLASM, IHC

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88341 (10), 88342 (1)

Scarlet Eligible: No

Spore Strips

0402-8

Alternate Name: SPORE STRIPS

Methodology: Bacterial Culture

Preferred Specimen: Duo Spore Test Kit

Collection Instruction: Submit control and sample spore at the same time, label each with source

Turn Around Time: 8 days

CPT Code: 87081 (1)

Scarlet Eligible: No

Lung - Adenocarcinoma vs. Squamous cell carcinoma

A130-7

Alternate Name: SQUAMOS V. ADENOCARCINOMA

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88341 (3), 88342 (1)

Scarlet Eligible: No

Synovial Sarcoma 18q11.2 FISH

J331-0

Alternate Name: SS18 or SYT

Methodology: FISH

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint. Or submit 5 consecutively cut unstained slides cut @ 5microns

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88377 (1)

Scarlet Eligible: No

Allergen StemphyliumBotryosum(m10),IgE

3828-1

Alternate Name: STEMPHYLIUM BOTRY. IgE

Methodology: Enzyme Immunoassay

Preferred Specimen: 0.5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 6 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

STI Lesion (NON-NY)

Q301-4

Alternate Name: STI LESION (NON-NY)

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Turn Around Time: 7 days

CPT Code: 87529 (2), 87798 (1)

Scarlet Eligible: No

STI Serology

P216-5

Alternate Name: STI SEROLOGY

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86592 (1), 86695 (1), 86696 (1), 86704 (1), 86803 (1), 87340 (1), 87389 (1)

Scarlet Eligible: Yes

Stormpath PAP + HPV HR W/Rfx to HPV 16 and 18/45

Q723-9

Alternate Name: STORMPATH PAP W/HPV PLUS

Preferred Specimen: ThinPrep Vial

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 87624 (1), 88175 (1)

Scarlet Eligible: No

Stormpath PAP+CT/GC+HPV HR w/RFX to HPV 16 and 18/45.

Q724-7

Alternate Name: STORMPATH PAP W/HPV+GC/CT

Preferred Specimen: ThinPrep Vial

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 87591 (1), 87491 (1), 87624 (1), 88175 (1)

Scarlet Eligible: No

Strep Group B (Amplified Probe)

0202-2

Alternate Name: STREP GRP.B (AMPLIFIED PROBE)

Methodology: Amplified Probe

Preferred Specimen: Swab-E

Collection Instruction: Vaginal/rectal specimens for this nucleic acid amplification assay should be collected at 35-37 weeks gest. using standard E-swab.

Turn Around Time: 2 days

CPT Code: 87653 (1)

Scarlet Eligible: No

Pneumococcal Antibody Panel (23 Serotype)

J048-0

Alternate Name: STREP PNEUMONIAE IGG ABS

Methodology: Quantitative Multiplex Bead Assay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 7 days

CPT Code: 86317 (23)

Scarlet Eligible: Yes

Surgical Pathology, Level IV

5619-2

Alternate Name: SURG.PATH.IV

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 2-3 days

CPT Code: 88305 (1)

Scarlet Eligible: No

SURGICAL PATHOLOGY LEVEL 2

5224-1

Alternate Name: SURGICAL PATHOLOGY LEVEL 2

Methodology: Histology

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Storage Instruction: Room Temp

Turn Around Time: 2 days

CPT Code: 88302 (1)

Scarlet Eligible: No

Allergen Salmon (f41), IgE

0714-6

Alternate Name: Salmon

Methodology: Immunocap

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Beta-Hexosaminidase B Enzyme Analysis (Sandhoff Disease)

B761-8

Alternate Name: Sandhoff Enzyme Testing

Preferred Specimen: Yellow top- ACD

Collection Instruction: Must fill 2 ACD tubes. Invert gently 5-6 times. Label with patients name.

Storage Instruction: Room Temp

Turn Around Time: 14 days

CPT Code: 82657 (1), 83080 (1), 84155 (1)

Scarlet Eligible: Yes

Allergen Sardine (Pilchard) (f308), IgE

6266-1

Alternate Name: Sardine (Pilchard)

Methodology: Immunocap

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Scallop (f338), IgE

2543-7

Alternate Name: Scallop

Methodology: Immunocap

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Hypercoagulation

2321-8

Alternate Name: Screen for Hypercoagulation; Thrombophilia Panel; Thrombotic Disease Screen

Preferred Specimen: 3 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Draw 4 light blue tops. Centrifuge blue-tops at high speed for 15 minutes. Transfer plasma into 4 plastic transfer tubes. Label with date/time collected, PLASMA, patient's name and DOB. SPECIMEN MUST BE SUBMITTED FROZEN.

Storage Instruction: STRICT Frozen

Turn Around Time: 4 days

CPT Code: 85300 (1), 81405 (1), 85303 (1), 85306 (1), 85384 (1), 85420 (1), 85613 (2)

Scarlet Eligible: No

Rast/Seafood Mix(fx2) IgE +

1843-2

Alternate Name: Seafood Mix

Methodology: Immunocap

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86005 (1)

Scarlet Eligible: Yes

Culture, Seminal Fluid

1827-5

Alternate Name: Seminal Fluid Culture, semen, sperm

Preferred Specimen: Semen Cup

Collection Instruction: For FERTILITY must be scheduled. For POST VASECTOMY, sample can be dropped off at any PSC for testing, label with patient name

Turn Around Time: 5 days

CPT Code: 87070 (1), 87205 (1)

Scarlet Eligible: No

Serum Drug Screen

L063-7

Alternate Name: Serum Drug Screen

Preferred Specimen: 6 mL Red Top

Collection Instruction: Fill tube, label with patient name and ID

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 80324 (1), 80361 (1), 80353 (1), 80365 (1), 80345 (1), 80358 (1), 83992 (1), 80346 (1)

Scarlet Eligible: Yes

Allergen Sesame Seed (f10), IgE

0936-5

Alternate Name: Sesame Seed

Methodology: Immunocap

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Sheep Epithelium (e81), IgE

1948-9

Alternate Name: Sheep Epithelium**Methodology:** Immunocap**Preferred Specimen:** 3 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 3 days**CPT Code:** 86003 (1)**Scarlet Eligible:** Yes**Allergen Sheep Sorrel (w18), IgE**

0894-6

Alternate Name: Sheep Sorrel (W18)**Methodology:** Immunocap**Preferred Specimen:** 3 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 2 days**CPT Code:** 86003 (1)**Scarlet Eligible:** Yes**Allergen Shrimp (f24), IgE**

0696-5

Alternate Name: Shrimp**Methodology:** Immunocap**Preferred Specimen:** 1 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 2 days**CPT Code:** 86003 (1)**Scarlet Eligible:** Yes

Allergen Silk (k74), IgE

1253-4

Alternate Name: Silk

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Extractable Nuclear Antigen (ENA) Antibodies (Sm/RnP)

0318-6

Alternate Name: Sm/RnP, Anti U1RNP, Anti Ribonuclear Protein, U1 RNP

Methodology: Multiplex Flow Immunoassay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86235 (2)

Scarlet Eligible: Yes

Smith-Lemli-Opitz Syndrome

J314-6

Alternate Name: Smith-Lemli-Opitz Syndrome

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81405 (1)

Scarlet Eligible: Yes

Smooth Muscle Antibody (ASM)

0419-2

Alternate Name: Smooth Muscle Antibody, Anti SMA**Methodology:** Indirect Fluorescence Assay**Preferred Specimen:** 1 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 4 days**CPT Code:** 86015 (1)**Scarlet Eligible:** Yes**Allergen Soybean (f14), IgE**

0680-9

Alternate Name: Soy Bean**Methodology:** Immunocap**Preferred Specimen:** 1 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 2 days**CPT Code:** 86003 (1)**Scarlet Eligible:** Yes**Spastic Tetraplegia, Thin Corpus Callosum and Progressive Microcephaly**

J298-1

Alternate Name: Spastic Tetraplegia, Thin Corpus Callosum and Progressive Mic**Methodology:** Genotyping by Next Generation Sequencing**Preferred Specimen:** Lavender top- EDTA**Collection Instruction:** Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!**Storage Instruction:** Refrigerate**Turn Around Time:** 10-14 days**CPT Code:** 81479 (1)**Scarlet Eligible:** Yes

Allergen Trichoderma Viride (m15), IgE

3829-9

Alternate Name: Specimen is accepted ambient, refrigerated or frozen.

Methodology: Immunoassay

Preferred Specimen: 0.5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 6 days

CPT Code: 86003 (1)

Scarlet Eligible: No

Culture, Sputum W/ Gram Stain

0382-2

Alternate Name: Sputum Culture

Methodology: Bacterial Culture

Preferred Specimen: Cup-Sputum

Collection Instruction: Submit sputum in sterile container, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 87070 (1), 87205 (1)

Scarlet Eligible: No

Culture, Stool (1st Specimen)

0077-8

Alternate Name: Stool Culture incl. Salmonella, Shigella, Campy, Yersinia, E Coli

Methodology: Bacterial Culture

Preferred Specimen: Fecal Swab

Collection Instruction: Using applicator stick, place 10 gms of stool in container, label with patient's name

Storage Instruction: Refrigerate

Turn Around Time: 4 days

CPT Code: 87045 (1), 87046 (2)

Scarlet Eligible: No

Allergen Storage Mite (Acarus) (d70), IgE

1605-5

Alternate Name: Storage Mite (Acarus)

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

StormPath PAP(RFLX to HPV mRNA GENOTYPING 16,18/45 w/ASC-US)

TB33-6

Alternate Name: StormPap>rflx HPV mRNA

Methodology: Microscopic Examination

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 3 days

CPT Code: 88175 (1)

Scarlet Eligible: No

StormPath PAP, liquid based + HPV DNA Genotyping 16/18

L079-3

Alternate Name: StormPath PAP+ PD HPV

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Turn Around Time: 3-5 days

CPT Code: 88175 (1), 87624 (1)

Scarlet Eligible: No

StormPath PAP, liquid based + HPV DNA + CT/GC

L080-1

Alternate Name: StormPath PAP, PD HPV + CT/GC

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 3-5 days

CPT Code: 88175 (1), 87491 (1), 87624 (1), 87591 (1)

Scarlet Eligible: No

Stormpath Pap/GC/CT reflex HPV when ASC-US

Q725

Alternate Name: Stormpath Pap/GC/CT reflex HPV when ASC-US

Scarlet Eligible: No

Allergen Strawberry (f44), IgE

0931-6

Alternate Name: Strawberry

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Strep. Group A Screen DNA

0079-4

Alternate Name: Strep Screen, Group A, Beta-Hemolytic Molecular

Methodology: LAMP DNA Amplification

Preferred Specimen: Swab-E

Collection Instruction: Collect specimen with swab then place swab into transport carrier; Label with name and source

Turn Around Time: 1 days

CPT Code: 87651 (1)

Scarlet Eligible: No

Allergen Sunflower Seed (k84), IgE

1414-2

Alternate Name: Sunflower Seed

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

High Risk STI- CT/GC- Aptima Swab

L062-9

Alternate Name: Swab CT/GC

Methodology: Multiplex PCR

Preferred Specimen: Swab-Aptima Genprobe

Collection Instruction: Collect specimen-break swab shaft into transport tube, then label with name and source

Storage Instruction: Room Temp

Turn Around Time: 2 days

CPT Code: 87491 (1), 87591 (1)

Scarlet Eligible: No

High Risk STI-CT/GC/Trich-Aptima Swab

L047-0

Alternate Name: Swab CT/GC/Trich

Preferred Specimen: Swab-Aptima Genprobe

Collection Instruction: Collect specimen-break swab shaft into transport tube, then label with name and source

Storage Instruction: Room Temp

Turn Around Time: 3 days

CPT Code: 87661 (1), 87491 (1), 87591 (1)

Scarlet Eligible: No

High Risk STI- CT/GC/Trich/M.Gen-Aptima Swab

L048-8

Alternate Name: Swab CT/GC/Trich/MGen

Preferred Specimen: Swab-Aptima Genprobe

Collection Instruction: Collect specimen-break swab shaft into transport tube, then label with name and source

Storage Instruction: Room Temp

Turn Around Time: 3 days

CPT Code: 87661 (1), 87491 (1), 87591 (1), 87798 (1)

Scarlet Eligible: No

Aptima Swab Cervicitis (NON-NY)

L049-6

Alternate Name: Swab Cervicitis

Preferred Specimen: Swab-Aptima Genprobe

Collection Instruction: Collect specimen-break swab shaft into transport tube, then label with name and source

Turn Around Time: 5-7 days

CPT Code: 87661 (1), 87491 (1), 87591 (1), 87529 (2), 87798 (1)

Scarlet Eligible: No

Allergen Sweet Potato (f54), IgE

1231-0

Alternate Name: Sweet Potato (F54)

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Sweet Vernal (g1), IgE

1196-5

Alternate Name: Sweet Vernal Grass (g1) IgE

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Sycamore (t11), IgE

1305-2

Alternate Name: Sycamore

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

C-KIT (D816V)

5179-7

Alternate Name: Systemic Mastocytosis

Methodology: Pyrosequencing

Preferred Specimen: 2 mL Bone Marrow - Lavender Top

Collection Instruction: Place bone marrow into tube, label with patient's name. Send in bone marrow kit.

Storage Instruction: Refrigerate

Turn Around Time: 3-4 days

CPT Code: 81273 (1)

Scarlet Eligible: No

Immune Deficiency Panel CD4/CD8 w/CBC

1763-2

Alternate Name: T & B Cell, T cell subset, T-cell subset,Lymphocyte enumeration

Methodology: Flow Cytometry

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Turn Around Time: 2 days

CPT Code: 86355 (1), 86357 (1), 86359 (1), 86360 (1)

Scarlet Eligible: Yes

T(11;14)+11 BY FISH-TC Only

A873-2

Alternate Name: T(11;14)+11 BY FISH-TC

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube. Also acceptable:Whole blood in green (sodium heparin) tube.

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88374 (1)

Scarlet Eligible: No

T(4;14) BY FISH-TC Only

A875-7

Alternate Name: T(4;14) BY FISH-TC ONLY

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube. Also acceptable:Whole blood in green (sodium heparin) tube.

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88374 (1)

Scarlet Eligible: No

T.Pallidum By Multiplex PCR (NON-NY)

6354-5

Alternate Name: T.PALLIDUM BY MULTIPLEX PCR

Methodology: Polymerase Chain Reaction

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Turn Around Time: 5-7 days

CPT Code: 87798 (1)

Scarlet Eligible: No

T3, Total

0150-3

Alternate Name: T3 (THYRONINE), TOTAL

Methodology: Electrochemiluminescence Immunoassay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84480 (1)

Scarlet Eligible: Yes

T3 Uptake

0152-9

Alternate Name: T3 Uptake (T3U)

Methodology: Electrochemiluminescence Immunoassay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84479 (1)

Scarlet Eligible: Yes

T3 Free, Serum

0271-7

Alternate Name: T3, FREE (FT3)

Methodology: Electrochemiluminescence Immunoassay

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84481 (1)

Scarlet Eligible: Yes

T4 (Thyroxine) Total

0151-1

Alternate Name: T4;Tetraiodothyronine;

Methodology: Electrochemiluminescence Immunoassay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84436 (1)

Scarlet Eligible: Yes

EGFR Sequencing- Non NY

J022-5

Alternate Name: T790M

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 10-15 unstained slides at 5 microns (tumor content>10%) with H&E, Shavings (if tumor content >40%) with H&E, if sending DNA - Please call CS.

Storage Instruction: Room Temp

Turn Around Time: 5-10 days

CPT Code: 81235 (1)

Scarlet Eligible: No

TAO Brush Endometrial Cytology

B685-9

Alternate Name: TAO BRUSH CYTOLOGY

Methodology: Microscopic Examination

Preferred Specimen: ThinPrep Cytolyt Solution

Collection Instruction: Collect specimen swish in Cytolyte vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 88305 (1)

Scarlet Eligible: No

Tay-Sachs Disease Carrier Test

A652-0

Alternate Name: TAY-SACHS DNA

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81255 (1)

Scarlet Eligible: Yes

Osteopetrosis, Autosomal Recessive 1

J295-7

Alternate Name: TCIRG1

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81479 (1)

Scarlet Eligible: No

Lymphoma - T-Cell, TCR Gamma by PCR

5031-0

Alternate Name: TCR, T-Cell Gene Rearrangement, CLONALITY IN T-CELLS

Methodology: Fluorescence

Preferred Specimen: 2 mL Bone Marrow - Lavender Top

Collection Instruction: Place bone marrow into tube, label with patient's name. Send in bone marrow kit.

Storage Instruction: Refrigerate

Turn Around Time: 6 days

CPT Code: 81342 (1)

Scarlet Eligible: No

TEL/AML1: t(12;21) by FISH

5038-5

Alternate Name: TEL/AML1 Gene Rearrangement by FISH, t(12;21)

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube. Also acceptable: Whole blood in green (sodium heparin) tube.

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 88377 (1)

Scarlet Eligible: No

Testosterone Free + Total

7301-5

Alternate Name: TESTOSTERONE FREE + TOTAL

Preferred Specimen: 2.5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 82040 (1), 84270 (1), 84403 (1)

Scarlet Eligible: Yes

Testosterone, Serum

0379-8

Alternate Name: TESTOSTERONE, TOT.,S.

Methodology: Electrochemiluminescence Immunoassay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84403 (1)

Scarlet Eligible: Yes

Cannabinoids (Marijuana), Urine, Screen (50)

3971-9

Alternate Name: THC

Methodology: Enzyme Immunoassay

Preferred Specimen: 5 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80307 (1)

Scarlet Eligible: Yes

Cannabinoids (Marijuana), Urine, Screen (20)

1576-8

Alternate Name: THC

Methodology: Enzyme Immunoassay

Preferred Specimen: 5 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80307 (1)

Scarlet Eligible: Yes

Thrombin Time

5807-3

Alternate Name: THROMBIN TIME(TT)

Methodology: Clot Detection

Preferred Specimen: 2 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 4 days

CPT Code: 85670 (1)

Scarlet Eligible: No

Thrombin/Antithrombin

5808-1

Alternate Name: THROMBIN/ANTITHROMBIN

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 2 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 10 days

CPT Code: 83520 (1)

Scarlet Eligible: No

Thyroglobulin, Serum

0577-7

Alternate Name: THYROGLOBULIN

Methodology: Chemiluminescence

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 84432 (1)

Scarlet Eligible: Yes

Thyroid Comprehensive Panel

0867-2

Alternate Name: THYROID COMPREHENSIVE

Preferred Specimen: 5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84436 (1), 84439 (1), 84443 (1), 84479 (1), 84480 (1)

Scarlet Eligible: Yes

Thyroid PNL-T3U,T4,TSH

0007-5

Alternate Name: THYROID PNL-T3U,T4,T

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84436 (1), 84443 (1), 84479 (1)

Scarlet Eligible: Yes

Thyroid Screen

0006-7

Alternate Name: THYROID SCREEN

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84436 (1), 84479 (1)

Scarlet Eligible: Yes

Thyroiditis

2322-6

Alternate Name: THYROIDITIS

Preferred Specimen: 5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 84439 (1), 84443 (1), 86376 (1), 86800 (1)

Scarlet Eligible: Yes

Tissue Path.(POC)

3001-5

Alternate Name: TISSUE PATH.(POC)

Methodology: Histology

Preferred Specimen: Products of Conception

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 1 days

CPT Code: 88304 (1)

Scarlet Eligible: No

Tissue Pathology, Gross Only

4002-2

Alternate Name: TISSUE PATHOLOGY, GROSS

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 1 days

CPT Code: 88300 (1)

Scarlet Eligible: No

Tobramycin, Peak

1235-1

Alternate Name: TOBRAMYCIN, PEAK
Methodology: Enzyme Immunoassay
Preferred Specimen: Red Top
Collection Instruction: Fill tube, label with patient name
Storage Instruction: Refrigerate
Turn Around Time: 3 days
CPT Code: 80200 (1)
Scarlet Eligible: No

Tobramycin (Random)

A615-7

Alternate Name: TOBRAMYCIN, RANDOM
Methodology: Enzyme Immunoassay
Preferred Specimen: Red Top
Collection Instruction: Fill tube, label with patient name
Storage Instruction: Refrigerate
Turn Around Time: 3 days
CPT Code: 80200 (1)
Scarlet Eligible: Yes

Tobramycin, Trough

2132-9

Alternate Name: TOBRAMYCIN, TROUGH
Methodology: Fluor. Polarization Immunoassay
Preferred Specimen: Red Top
Collection Instruction: Fill tube, label with patient name
Storage Instruction: Refrigerate
Turn Around Time: 1 days
CPT Code: 80200 (1)
Scarlet Eligible: No

Toxoplasma Antibodies, (IgG, IgM)

2721-9

Alternate Name: TOXO IGG/IGM**Preferred Specimen:** 3 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 1 days**CPT Code:** 86777 (1), 86778 (1)**Scarlet Eligible:** Yes**Toxoplasma IgG W/Reflex**

2376-2

Alternate Name: TOXO.Ab.,IgG(reflex)**Methodology:** Chemiluminescence**Preferred Specimen:** 1 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 1 days**CPT Code:** 86777 (1)**Scarlet Eligible:** Yes**Toxoplasma Ab. (IgG), EIA**

0552-0

Alternate Name: TOXOPLASMA Ab.(IgG)**Methodology:** Chemiluminescence**Preferred Specimen:** 1 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 1 days**CPT Code:** 86777 (1)**Scarlet Eligible:** Yes

Toxoplasma Ab. (IgM)

0877-1

Alternate Name: TOXOPLASMA Ab.(IgM)

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86778 (1)

Scarlet Eligible: Yes

TP53 Sequencing

J044-9

Alternate Name: TP53 Sequencing

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Bone Marrow - Lavender Top

Collection Instruction: Place bone marrow into tube, label with patient's name

Storage Instruction: Room Temp

Turn Around Time: 10 days

CPT Code: 81352 (1)

Scarlet Eligible: No

OnkoSight NGS TP53 Sequencing for Solid Tumors

T980-2

Alternate Name: TP53 Sequencing

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: 10-15 unstained slides at 5 microns (tumor content >10%) with H&E, Shavings (if tumor content >40%) with H&E, if sending DNA - Please call CS. Please give to Block Retrieval for processing. Alternative specimen: Extracted DNA is acceptable, providing that the isolation of nucleic acid occurs in a CLIA-certified laboratory or a laboratory meeting equivalent requirements as determined by CMS and/or the CAP.

Storage Instruction: Room Temp

Turn Around Time: 10 days

CPT Code: 81352 (1)

Scarlet Eligible: No

Tramadol, Urine, Screen, Qual. (200)

5917-0

Alternate Name: TRAMADOL, QUAL

Methodology: Enzyme Immunoassay

Preferred Specimen: 5 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80307 (1)

Scarlet Eligible: Yes

Transglutaminase, IgA Abs.

3472-8

Alternate Name: TRANSGLUTAMINASE,IgA,Ab

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 1 SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86364 (1)

Scarlet Eligible: Yes

Transglutaminase, IGG ABS.

5869-3

Alternate Name: TRANSGLUTAMINASE,IgG,Ab

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86364 (1)

Scarlet Eligible: Yes

Trichomonas Vaginalis, Probe rRNA

A861-7

Alternate Name: TRICHOMONAS VAG rRNA

Methodology: Transcription Mediated Amplification

Preferred Specimen: Urine Container-Genprobe-Aptima

Collection Instruction: First morning sample, collect in sterile cup, transfer 2 ml into aptima holder. Correct level is between the black fill lines(window) on the holder.

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87661 (1)

Scarlet Eligible: Yes

Trimipramine, Serum

0613-0

Alternate Name: TRIMIPRAMINE, SERUM

Methodology: Liquid Chromatography Tandem Mass Spectrometry

Preferred Specimen: Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 9 days

CPT Code: 80335 (1)

Scarlet Eligible: Yes

Trisomy 12 by FISH-TC Only

A877-3

Alternate Name: TRISOMY 12 BY FISH-TC

Methodology: FISH

Preferred Specimen: 5 mL Peripheral Blood - Green Top

Collection Instruction: Fill tube, invert 8-10 times (DO NOT SHAKE!)

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88373 (1)

Scarlet Eligible: No

Trisomy 8 By FISH -TC Only

A870-8

Alternate Name: TRISOMY 8 BY FISH-TC

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube. Also acceptable: Whole blood in green (sodium heparin) tube.

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88373 (1)

Scarlet Eligible: No

Trofile DNA

A205-7

Alternate Name: TROFILE DNA

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: STRICT Frozen

Turn Around Time: 14-21 days

CPT Code: 87999 (1)

Scarlet Eligible: No

Troponin I

2163-4

Alternate Name: TROPONIN-I

Methodology: Chemiluminescence

Preferred Specimen: Aliquot Tube-Serum

Collection Instruction: Place 1-3 mL of serum in transport tube. Label as serum and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 1 days

CPT Code: 84484 (1)

Scarlet Eligible: No

TSH W/ Reflex To Free T4

A518-3

Alternate Name: TSH W/RFX TO FREE T4

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84443 (1)

Scarlet Eligible: Yes

Thyroid Stimulating Hormone (TSH)

0153-7

Alternate Name: TSH,Thyrotropin;Ultrasensitive TSH

Methodology: Electrochemiluminescence Immunoassay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84443 (1)

Scarlet Eligible: Yes

TSH (Pregnancy)

J265-0

Alternate Name: TSH;Thyroidopin;Ultrasensitive TSH

Methodology: Electrochemiluminescence Immunoassay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84443 (1)

Scarlet Eligible: Yes

Thyroid Stimulating Immunoglobulin

J510-9

Alternate Name: TSI

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84445 (1)

Scarlet Eligible: Yes

Transglutaminase IGG/A,AB, Evaluation

3473-6

Alternate Name: TTG

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 83516 (2)

Scarlet Eligible: Yes

Tumor of Unknown Primary IHC Panel

5190-4

Alternate Name: TUMOR OF UNK.PRIM.ORIGIN

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88341 (13), 88342 (1)

Scarlet Eligible: No

Tyrosinemia, Type 1

B149-6

Alternate Name: TYROSINEMIA TYPE 1

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81479 (1)

Scarlet Eligible: Yes

Tysabri Antibodies

6293-5

Alternate Name: TYSABRI AB.

Methodology: Enzyme Immunoassay

Preferred Specimen: 2.0 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 16 days

CPT Code: 83516 (1)

Scarlet Eligible: Yes

Hexosaminidase Leukocytes (Tay-Sachs Enzyme)

0865-6

Alternate Name: Tay-Sachs Enzyme

Preferred Specimen: Yellow top- ACD

Collection Instruction: 2 full ACD tubes or 2 full Lavenders

Storage Instruction: Room Temp

Turn Around Time: 15 days

CPT Code: 82657 (1), 83080 (1), 84155 (1)

Scarlet Eligible: Yes

Allergen Tea (rf222), IgE

1390-4

Alternate Name: Tea

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Carbamazepine

0154-5

Alternate Name: Tegretol

Methodology: Enzyme Immunoassay

Preferred Specimen: 1 mL Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80156 (1)

Scarlet Eligible: Yes

Total Testosterone by LC/MS/MS

J451-6

Alternate Name: Testosterone,Tot.LC/MS/MS

Methodology: Liquid Chromatography Tandem Mass Spectrometry

Preferred Specimen: 1 mL Red Top

Collection Instruction: Place 1-3 mL of serum in transport tube. Label as serum;refrigerate or freeze as required.

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 84403 (1)

Scarlet Eligible: Yes

Throat Culture, R/O Group A Strep. ONLY

TJ06-6

Alternate Name: Throat Culture, R/O GAS

Methodology: Bacterial Culture

Preferred Specimen: Swab-E

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87070 (1)

Scarlet Eligible: No

Culture, Fungal (Yeast)

A502-7

Alternate Name: Thrush, Monilia Culture

Methodology: Fluorescence

Preferred Specimen: Swab-E

Collection Instruction: Collect specimen with swab then place swab into transport carrier; Label with name and source

Storage Instruction: Refrigerate

Turn Around Time: 12 days

CPT Code: 87102 (1), 87206 (1)

Scarlet Eligible: No

Calcitonin, Serum

0196-6

Alternate Name: Thyrocalcitonin

Methodology: Chemiluminescence

Preferred Specimen: 2 mL Aliquot Tube-Serum

Collection Instruction: Place 1-3 mL of serum in transport tube. Label as serum and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 2 days

CPT Code: 82308 (1)

Scarlet Eligible: No

Anti-Thyroglobulin Ab. (ATG)

0041-4

Alternate Name: Thyroglobulin Antibody, ATA,ATG

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86800 (1)

Scarlet Eligible: Yes

Anti-Microsomal Ab(Thyroid Peroxidase)TPO

0412-7

Alternate Name: Thyroid Peroxidase Abs; Thyroid Peroxidase Antibodies

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86376 (1)

Scarlet Eligible: Yes

Tick Identification

J291-6

Alternate Name: Tick ID

Methodology: Microscopic Examination

Preferred Specimen: Tick

Collection Instruction: Place tick in sterile container, label with patient name.

Storage Instruction: Room Temp

Turn Around Time: 3 days

CPT Code: 87168 (1)

Scarlet Eligible: Yes

Allergen Tilapia (RF414), IgE

8963-1

Alternate Name: Tilapia, Ige (Rf414)

Methodology: Immunocap

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Timothy Grass (g6), IgE

0683-3

Alternate Name: Timothy Grass

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Tissue Pathology One Check

T793

Alternate Name: Tissue Pathology One Check

Scarlet Eligible: No

Allergen Tomato (f25), IgE

0837-5

Alternate Name: Tomato

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

SCL70 (Scleroderma Ab.), Serum

0315-2

Alternate Name: Topoisomerase I Antibody,SCL70

Methodology: Chemiluminescence

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86235 (1)

Scarlet Eligible: Yes

TORCH Panel (IgG/IgM)

1297-1

Alternate Name: Torch Panel, IgG, IgM

Preferred Specimen: 4 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86644 (1), 86645 (1), 86695 (1), 86762 (1), 86777 (1), 86778 (1)

Scarlet Eligible: Yes

Bilirubin, Infant, Total/Direct, Serum

0046-3

Alternate Name: Total Bilirubin, Infant, Serum

Methodology: Colorimetry

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Turn Around Time: 1 days

CPT Code: 82247 (1), 82248 (1)

Scarlet Eligible: No

Immunoglobulin E (IgE), Serum/Plasma

0996-9

Alternate Name: Total IgE

Methodology: Immunospectrometric

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82785 (1)

Scarlet Eligible: Yes

Total IgG and subclasses IgG 1, 2, 3 & 4

H991-3

Alternate Name: Total IgG/Subclasses

Methodology: Immunospectrometric

Preferred Specimen: 1 mL SST Tube

Storage Instruction: Refrigerate

Turn Around Time: 3-4 days

CPT Code: 82784 (1), 82787 (4)

Scarlet Eligible: Yes

Protein, Total, Serum

0135-4

Alternate Name: Total Protein, Serum

Methodology: Colorimetry

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84155 (1)

Scarlet Eligible: Yes

Testosterone, Total/Free (LC/MS)

K630-4

Alternate Name: Total/Free Testosterone (LC/MS)

Collection Instruction: 1 SST and a Red must be submitted.

Turn Around Time: 3 days

CPT Code: 84403 (1), 84270 (1), 99999 (1)

Scarlet Eligible: Yes

Prealbumin

1006-6

Alternate Name: Transthyretin, Serum

Methodology: Immunoturbidimetric

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 84134 (1)

Scarlet Eligible: Yes

Tree Mix (tx1), IgE

2178-2

Alternate Name: Tree Mix

Methodology: Immunocap

Preferred Specimen: 5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86005 (1)

Scarlet Eligible: Yes

Trees Screen (6)

0742-7

Alternate Name: Trees Screen (6)

Preferred Specimen: 5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82785 (1), 86003 (6)

Scarlet Eligible: Yes

Serodia TPPA (Antibody to T. pallidum)

0654-4

Alternate Name: Treponema Pallidum Particle Agglutination, TPPA

Methodology: Agglutination

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Turn Around Time: 3 days

CPT Code: 86780 (1)

Scarlet Eligible: Yes

Liaison CIA (Antibody to T. pallidum)

0334-3

Alternate Name: Treponema pallidum Ab,Treponema Pallidum Ab Chemiluminescence Im

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86780 (1)

Scarlet Eligible: Yes

Lipase, Serum

0521-5

Alternate Name: Triacylglycerol Acylhydrolase, Serum

Methodology: Enzyme Immunoassay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 83690 (1)

Scarlet Eligible: Yes

Trichomonas By Multiplex PCR

6342-0

Alternate Name: Trichomonas By Mult PCR

Methodology: Polymerase Chain Reaction

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 87661 (1)

Scarlet Eligible: No

Trichomonas by MultiPlex PCR - swab

J836-8

Alternate Name: Trichomonas PCR Swab

Methodology: Multiplex PCR

Preferred Specimen: Swab Aptima Orange label

Collection Instruction: Collect specimen-break swab shaft into transport tube, then label with name and source

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87661 (1)

Scarlet Eligible: No

Allergen Trichophyton Rubrum (m205), IgE

A673-6

Alternate Name: Trichophyton Rubrum

Methodology: Immunocap

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Triglycerides

0155-2

Alternate Name: Triglycerides

Methodology: Enzyme Immunoassay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84478 (1)

Scarlet Eligible: Yes

Allergen Trout (f204), IgE

3643-4

Alternate Name: Trout**Methodology:** Immunocap**Preferred Specimen:** 2 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 3 days**CPT Code:** 86003 (1)**Scarlet Eligible:** Yes**Allergen Tuna (f40), IgE**

0698-1

Alternate Name: Tuna**Methodology:** Immunocap**Preferred Specimen:** 1 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 3 days**CPT Code:** 86003 (1)**Scarlet Eligible:** Yes**Allergen Turkey Feathers (e89), IgE**

1426-6

Alternate Name: Turkey Feathers**Methodology:** Immunocap**Preferred Specimen:** 3 mL SST Tube**Collection Instruction:** Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes**Storage Instruction:** Refrigerate**Turn Around Time:** 3 days**CPT Code:** 86003 (1)**Scarlet Eligible:** Yes

Allergen Turkey Meat (f284), IgE

1149-4

Alternate Name: Turkey Meat

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Acetaminophen, Serum

0622-1

Alternate Name: Tylenol

Methodology: Enzyme Immunoassay

Preferred Specimen: Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80143 (1)

Scarlet Eligible: Yes

Urinalysis w/Reflex to Culture

B518-2

Alternate Name: UA,URINE,REFLEX,

Methodology: Colorimetry

Preferred Specimen: 12 Urine Urinalysis Tube - Yellow,1 mL Urine Tube - Grey Top

Collection Instruction: In order to obtain reflex both specimens must be received Borincon or UGY and a UA tube, If only 1 specimen is received only the appropriate test will be performed.

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 81001 (1)

Scarlet Eligible: Yes

H.Pylori Urea Breath Test

6236-4

Alternate Name: UBIT-R

Methodology: Spectrophotometry

Preferred Specimen: Breath Bag, Blue/Pink

Collection Instruction: REFER TO KIT INSTRUCTIONS. HUMAN breath: Paired breath samples (pre and post) collected in BreathTek(tm) UBT Collection Kit bags and must be submitted together.

Storage Instruction: Room Temp

Turn Around Time: 7 days

CPT Code: 83013 (1)

Scarlet Eligible: No

Protein Electrophoresis, 24Hr Urine

Q069-7

Alternate Name: UPEP, 24 HR

Methodology: Electrophoresis

Preferred Specimen: Urine Cup 24 hour

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 81050 (1), 84156 (1), 84166 (1)

Scarlet Eligible: No

Urea Clearance, 24hr Urine

1698-0

Alternate Name: UREA CLEARANCE

Methodology: Urease

Preferred Specimen: Urine Container - 24hr,SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84545 (1)

Scarlet Eligible: No

Urea, Urine 24hr

0358-2

Alternate Name: UREA,URINE 24 HR

Preferred Specimen: 10 mL Urine Container - 24hr

Collection Instruction: Discard first morning sample, then collect all urines after for 24 hrs, on completion label with collection time and patient's name; record total volume on Requisition

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 81050 (1), 84540 (1)

Scarlet Eligible: No

Urea, Urine, Timed

2083-4

Alternate Name: UREA,URINE TIMED

Methodology: Urease

Preferred Specimen: 5 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84540 (1)

Scarlet Eligible: Yes

Ureaplasma By Multiplex PCR

6339-6

Alternate Name: UREAPLASMA BY MULTIPLEX PCR

Methodology: Polymerase Chain Reaction

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 87798 (1)

Scarlet Eligible: No

Uric Acid, Synovial

0826-8

Alternate Name: URIC ACID, SYNOVIAL

Methodology: Colorimetry

Preferred Specimen: 1 mL Fluid-source required

Collection Instruction: Submit fluid in sterile cup

Storage Instruction: Refrigerate

Turn Around Time: 9 days

CPT Code: 84560 (1)

Scarlet Eligible: No

Uric Acid, Urine, Timed Spec.

1663-4

Alternate Name: URIC ACID, URINE TIMED

Methodology: Colorimetry

Preferred Specimen: 3 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84560 (1)

Scarlet Eligible: Yes

Urinalysis, Routine

0159-4

Alternate Name: URINALYSIS (UA)

Methodology: Microscopic Examination

Preferred Specimen: 12 mL Urine Urinalysis Tube - Yellow

Collection Instruction: Collect urine in sterile cup, transfer 10 ml into tube, label with patient's name, keep refrigerated

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 81001 (1)

Scarlet Eligible: Yes

Cytopathology, Urine (Thinprep)

5254-8

Alternate Name: URINE CYTOPATHOLOGY

Preferred Specimen: 50 mL ThinPrep Cytolyt Solution

Collection Instruction: Add urine to equal part of 30ml thinprep cytolyte solution

Turn Around Time: 4 days

CPT Code: 88112 (1)

Scarlet Eligible: No

URINE FREE KAPPA+LAMBDA

L592-5

Alternate Name: URINE FREE KAPPA+LAMBDA

Methodology: Immunoturbidimetric

Preferred Specimen: 10.0 mL Urine Cup

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 83521 (2), 99999 (1)

Scarlet Eligible: No

Urogenital Mycoplasmas

P114-2

Alternate Name: UROGENITAL MYCOPLASMAS

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Turn Around Time: 5 days

CPT Code: 87798 (2)

Scarlet Eligible: No

U.R.R. PANEL

1829-1

Alternate Name: URR Panel

Preferred Specimen: SST Tube, Microtainer - Pediatric SST

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84520 (2)

Scarlet Eligible: No

Usher Syndrome Type 1F

A649-6

Alternate Name: USHER SYNDROME TYP 1F

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81400 (1)

Scarlet Eligible: Yes

Uric Acid, Serum

0157-8

Alternate Name: Urate

Methodology: Enzyme Immunoassay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84550 (1)

Scarlet Eligible: Yes

Uric Acid, Urine, 24 Hr.

0158-6

Alternate Name: Urate, Urine

Preferred Specimen: 5 mL Urine Container - 24hr

Collection Instruction: Discard first morning sample, then collect all urines after for 24 hrs, on completion label with collection time and patient's name; record total volume on Requisition

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 81050 (1), 84560 (1)

Scarlet Eligible: No

Urethral Smear

1500-8

Alternate Name: Urethral Smear

Methodology: Bacterial Culture

Preferred Specimen: Slides

Collection Instruction: Label slides with patient name, place in slide holder for transport

Turn Around Time: 2 days

CPT Code: 88160 (1)

Scarlet Eligible: No

Creatinine Urine, Timed or Random

1661-8

Alternate Name: Urine Creatinine

Methodology: Colorimetry

Preferred Specimen: 10 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82570 (1)

Scarlet Eligible: Yes

Urine Culture

0080-2

Alternate Name: Urine Culture

Methodology: Bacterial Culture

Preferred Specimen: 5-10 mL Urine Tube - Grey Top

Collection Instruction: Perform a clean catch process, collect urine in a sterile container, transfer 10-15ml into a Grey Vacutainer (UGY), label with patients name.

Storage Instruction: Refrigerate

Turn Around Time: 4 days

CPT Code: 87086 (1)

Scarlet Eligible: Yes

Urine Culture, Prenatal, W/GBS

TH76-3

Alternate Name: Urine Culture, Prenatal

Methodology: Bacterial Culture

Preferred Specimen: 5-10 mL Urine Tube - Grey Top

Storage Instruction: Refrigerate

Turn Around Time: 4 days

CPT Code: 87086 (1)

Scarlet Eligible: Yes

Urine Cytology Reflex Atypical or above to UroVysion

B415

Alternate Name: Urine Cytology Reflex Atypical or above to UroVysion

Scarlet Eligible: No

Urine Cytology, if Atypical, reflex to hTert

TG28-6

Alternate Name: Urine Cytology rflx hTert

Preferred Specimen: ThinPrep Cytolyt Solution

Storage Instruction: Room Temp

Turn Around Time: 4-6 days

CPT Code: 88112 (1)

Scarlet Eligible: No

Urine Cytopathology, TC (Alcohol)

5554

Alternate Name: Urine Cytopathology, TC (Alcohol)

Scarlet Eligible: No

Comprehensive Urine Pathology

5324-9

Alternate Name: UroVysion and Urine Cytology

Methodology: FISH

Preferred Specimen: 50 mL Voided Urine - Urocyte Collection Kit

Collection Instruction: Collect voided urine in Thinprep Cytolyte solution

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88112 (1), 88121 (4)

Scarlet Eligible: No

Comprehensive Urine Pathology (Technical Only)

5253-0

Alternate Name: UroVysion and Urine Cytology- TC Only

Methodology: FISH

Preferred Specimen: 50 mL Voided Urine - Urocyte Collection Kit

Collection Instruction: Collect voided urine in Thinprep Cytolyte solution

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88112 (1), 88121 (1)

Scarlet Eligible: No

UroVysion (Tech Only)

5249-8

Alternate Name: UroVysion for Bladder Cancer: Technical Component Only.

Methodology: FISH

Preferred Specimen: 50 mL Voided Urine - Urocyte Collection Kit

Collection Instruction: Follow kit instructions

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88121 (1)

Scarlet Eligible: No

Bladder Washing, FISH

J371-6

Alternate Name: Urovysion Bladder Wash

Methodology: FISH

Preferred Specimen: 50 mL BLADDER WASHINGS

Collection Instruction: Perform bladder-washing collection (50 mL) at the physicians office in a sterile 100 mL urine specimen container. Mix bladder washing with CytoLyt OR ethanol 50% in a 1:1 solution. Alternatively, mix bladder washing with preservative Carbowax (2% polyethylene glycol in 50% ethanol) 2:1(v:v). Stability is 72hrs, however the lab will not reject sample and will determine when received.

Storage Instruction: Room Temp

Turn Around Time: 3-4 days

CPT Code: 88120 (1)

Scarlet Eligible: No

Usher Syndrome Type 3

A647-0

Alternate Name: Usher Syndrome Type 3

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81400 (1)

Scarlet Eligible: Yes

Vaginitis/Vaginosis W/O PAP W/O HPV Expanded

F371-0

Alternate Name: VAG/VAG NO PAP/HPV EXPD

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Turn Around Time: 7 days

CPT Code: 87481 (2), 87512 (1), 87661 (1), 87798 (2), 87799 (2)

Scarlet Eligible: No

Vaginal Biopsy

5889-1

Alternate Name: VAGINAL BIOPSY

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Turn Around Time: 1 days

CPT Code: 88305 (1)

Scarlet Eligible: No

Vancomycin Peak, Serum

0575-1

Alternate Name: VANCOMYCIN, PEAK

Methodology: Enzyme Immunoassay

Preferred Specimen: Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80202 (1)

Scarlet Eligible: Yes

Vancomycin (Random)

A614-0

Alternate Name: VANCOMYCIN, RANDOM

Methodology: Enzyme Immunoassay

Preferred Specimen: Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80202 (1)

Scarlet Eligible: Yes

Vancomycin, Trough, Serum

1589-1

Alternate Name: VANCOMYCIN, TROUGH

Methodology: Fluor. Polarization Immunoassay

Preferred Specimen: Red Top

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80202 (1)

Scarlet Eligible: Yes

Venous Prothrom. Risk Factor

5875-0

Alternate Name: VENOUS PROTHROM.RISK FACTOR

Preferred Specimen: Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA and patient's name.

Turn Around Time: 4 days

CPT Code: 85240 (1), 85250 (1), 85270 (1), 85384 (1), 85610 (1), 85730 (1)

Scarlet Eligible: No

Venous Prothrombic Risk

5989-9

Alternate Name: VENOUS PROTHROMBOTIC RISK

Preferred Specimen: Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA and patient's name.

Turn Around Time: 4 days

CPT Code: 85240 (1), 85250 (1), 85270 (1), 85379 (1), 85384 (1), 85610 (1), 85730 (1)

Scarlet Eligible: No

Vitamin B12/Folate

0287-3

Alternate Name: VIT B12/FOLATE

Methodology: Chemiluminescence

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 82607 (1), 82746 (1)

Scarlet Eligible: Yes

Vitamin B12 Binding Capacity

0906-8

Alternate Name: VITAM.B12 BIND.CAP.

Methodology: Radio Immuno Assay

Preferred Specimen: 1 mL Aliquot Tube-Serum

Collection Instruction: Fill tube, invert gently 2-3 times, let stand for 20 minutes, spin for 10-15 minutes label with patient name.

Storage Instruction: Refrigerate

Turn Around Time: 14 days

CPT Code: 82608 (1)

Scarlet Eligible: Yes

Vitamin A + Beta Carotene

1332-6

Alternate Name: VITAMIN A + BETA CAROTENE

Preferred Specimen: 4 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82380 (1), 84590 (1)

Scarlet Eligible: Yes

Vitamin Panel

2861-3

Alternate Name: VITAMIN PANEL

Collection Instruction: All specimens must be received frozen and labeled accordingly. One of the serum from SST pour off's must be received in an Amber tube or vial wrapped in foil.

Turn Around Time: 3 days

CPT Code: 82380 (1), 84207 (1), 84252 (1), 84255 (1), 84425 (1), 84446 (1), 84590 (1)

Scarlet Eligible: Yes

Vitamin B1, Total Thiamine

0293-1

Alternate Name: VITB1,THIAMINE,TOTAL

Methodology: High Pressure Liquid Chromatography

Preferred Specimen: 2 mL Amber Serum

Collection Instruction: Draw SST or RED, let clot for 8-10 minutes. Centrifuge and pour into either an aliquot tube (must be covered in foil) or poured into an amber tube and label with patient name and specimen type (serum).

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 84425 (1)

Scarlet Eligible: Yes

Von Willebrand Screen- Basic

5825-5

Alternate Name: VON-WILLEBRAND SCREEN - BASIC

Preferred Specimen: Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 9 days

CPT Code: 85240 (1), 85245 (1), 85246 (1), 85610 (1), 85730 (1)

Scarlet Eligible: No

Vulvar Biopsy

5886-7

Alternate Name: VULVAR BIOPSY

Methodology: Histology

Preferred Specimen: Biopsy In Formalin

Collection Instruction: Submit in 10% Neutral Buffered Formalin, varied jar sizes.

Storage Instruction: Refrigerate

Turn Around Time: 2-3 days

CPT Code: 88305 (1)

Scarlet Eligible: No

Varicella-Zoster Virus

J839-2

Alternate Name: VZV

Methodology: Helicase Dependent Amplification

Preferred Specimen: 2 mL Swab-Viral Culturette

Collection Instruction: Collect sample, place in holder, label with patient name. NOT TO BE USED FOR BACTERIAL TRANSPORT**MUST SWAB COLLECTED FROM UNCAPPING THE LESION**. Frozen samples are stable for 7days prior to processing.

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 87798 (1)

Scarlet Eligible: No

Varicella Zoster Virus by IHC

TF85-8

Alternate Name: VZV by IHC

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Storage Instruction: Room Temp

Turn Around Time: 24-48 days

CPT Code: 88342 (1)

Scarlet Eligible: No

Varicella Zoster Virus by IHC Tech Only

TF88-2

Alternate Name: VZV by IHC Tech Only

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Storage Instruction: Room Temp

Turn Around Time: 24-48 days

CPT Code: 88342 (1)

Scarlet Eligible: No

Benzodiazepines, Urine, Screen (200)

0924-1

Alternate Name: Valium; Librium; Dalmane; Tranxene; Klonopin; Ativan ;Serax;Centrax; Restoril; Xanax; Halcion;Versed; Doral

Methodology: Enzyme Immunoassay

Preferred Specimen: 5 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 80307 (1)

Scarlet Eligible: Yes

Allergen Vanilla (rf234), IgE

2347-3

Alternate Name: Vanilla

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Varicella Zoster Ab IgG W/Reflex

T008-2

Alternate Name: Varicella Zos(IgG)Rfx

Methodology: Chemiluminescence

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86787 (1)

Scarlet Eligible: Yes

Allergen Profile, Insects (6)

2469-5

Alternate Name: Venom/Insect Profile (5)

Preferred Specimen: 5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (6), 82785 (1)

Scarlet Eligible: Yes

Vitamin C Dilution

6193-7

Alternate Name: Vitamin C Dilution

Methodology: High Pressure Liquid Chromatography

Preferred Specimen: 3 mL Amber Serum

Collection Instruction: Draw SST or RED, let clot for 8-10 minutes. Centrifuge and pour into either an aliquot tube (must be covered in foil) or poured into an amber tube and label with patient name and specimen type (serum). FREEZE.

Storage Instruction: STRICT Frozen

Turn Around Time: 10 days

CPT Code: 82180 (1)

Scarlet Eligible: No

Vit. D, 25-Hydroxy Tot.,D2+D3

8573-8

Alternate Name: Vitamin D3

Methodology: Liquid Chromatography Mass Spectrometry

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 82306 (1)

Scarlet Eligible: Yes

Ristocetin Cofactor ACT.

5805-7

Alternate Name: Von Willebrand Factor Activity

Methodology: Optical Platelet Aggregometry

Preferred Specimen: 1.5 mL Aliquot Plasma-Citrated-Blue Top

Collection Instruction: Collect full blue-top (citrated) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.

Storage Instruction: STRICT Frozen

Turn Around Time: 5 days

CPT Code: 85245 (1)

Scarlet Eligible: No

Vulvovaginal Candidiasis Panel by PCR

K853-2

Alternate Name: Vulvovaginal Candidiasis Pnl

Methodology: Polymerase Chain Reaction

Preferred Specimen: Swab-Aptima Genprobe

Collection Instruction: Collect specimen-break swab shaft into transport tube, then label with name and source

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 87481 (1)

Scarlet Eligible: No

Walker Warburg Syndrome

A648-8

Alternate Name: WALKER WARBURG

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81479 (1)

Scarlet Eligible: Yes

AAMI Water Analysis

1795-4

Alternate Name: WATER ANALYSIS

Preferred Specimen: 10 mL Cup Plain

Collection Instruction: Place sample in sterile cup label with patient's name and source

Storage Instruction: Refrigerate

Turn Around Time: 10 days

CPT Code: 99999 (1)

Scarlet Eligible: No

WBC

1497-7

Alternate Name: WBC

Methodology: Flow Cytometry

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 85048 (1)

Scarlet Eligible: Yes

Bioavailable Testosterone

2087-5

Alternate Name: WEAKLY BINDING TESTOSTERONE

Preferred Specimen: SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 82040 (1), 84270 (1), 84403 (1)

Scarlet Eligible: Yes

Wilson Disease Carrier Screen

B164-5

Alternate Name: WILSON DISEASE CARRIER

Methodology: Genotyping by Next Generation Sequencing

Preferred Specimen: 1 mL Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 10-14 days

CPT Code: 81479 (1)

Scarlet Eligible: Yes

Female High Risk Sexual Health Panel

F238-1

Alternate Name: WOMENS H/R SEX PNL

Collection Instruction: Submit 1 SST, 1 Viral Media, 1 Urine Aptima

Turn Around Time: 10 days

CPT Code: 86706 (1), 86803 (1), 87109 (2), 87389 (1), 87491 (1), 87591 (1), 87661 (1)

Scarlet Eligible: No

Wound Culture (Aerobic/Anaerobic) W/ Gram Stain (Panel)

0082-8

Alternate Name: WOUND CULTURE PANEL

Methodology: Bacterial Culture

Preferred Specimen: Swab-E

Collection Instruction: Collect specimen with swab then place swab into transport carrier; Label with name and source

Turn Around Time: 5 days

CPT Code: 87070 (1), 87075 (1), 87205 (1)

Scarlet Eligible: No

Allergen Walnut Tree (t10), IgE

0686-6

Alternate Name: Walnut Tree

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 2 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Walnuts (f256), IgE

1479-5

Alternate Name: Walnuts

Methodology: Immunocap

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Watermelon (f329), IgE

1941-4

Alternate Name: Watermelon

Methodology: Immunocap

Preferred Specimen: 3.0 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Weed Mix (wx1), IgE

2177-4

Alternate Name: Weed Mix

Methodology: Immunocap

Preferred Specimen: 5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86005 (1)

Scarlet Eligible: Yes

Weeds Screen (7)

0746-8

Alternate Name: Weeds Screen (7)

Preferred Specimen: 5 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 82785 (1), 86003 (7)

Scarlet Eligible: Yes

Allergen Western Ragweed (w2), IgE

3802-6

Alternate Name: West.Ragweed,Ige (W2)

Methodology: Immunocap

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Wheat (f4), IgE

0673-4

Alternate Name: Wheat

Methodology: Immunocap

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Whey (f236), IgE

1620-4

Alternate Name: Whey

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen White Ash (t15), IgE

1306-0

Alternate Name: White Ash

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen White Bean (f15), IgE

1241-9

Alternate Name: White Bean

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen White Pine (t16), IgE

3399-3

Alternate Name: White Pine (T16)

Methodology: Immunocap

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Whiteface Hornet (i2), IgE

0907-6

Alternate Name: Whiteface Hornet

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Willow (t12), IgE

1252-6

Alternate Name: Willow

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Wormwood (w5), IgE

1307-8

Alternate Name: Wormwood

Methodology: Immunocap

Preferred Specimen: 3.0 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Culture, Wound, Superficial

B656-0

Alternate Name: Wound Aerobic

Methodology: Bacterial Culture

Preferred Specimen: Swab-E

Collection Instruction: Collect specimen with swab then place swab into transport carrier; Label with name and source

Storage Instruction: Refrigerate

Turn Around Time: 5 days

CPT Code: 87070 (1)

Scarlet Eligible: No

XCE 12 by FISH

B513

Alternate Name: XCE 12 by FISH

Scarlet Eligible: No

XL RARA (17q21) BY FISH

J157-9

Alternate Name: XL RARA (17q21) BY FISH

Methodology: FISH

Preferred Specimen: Bone Marrow - Green Top

Collection Instruction: Place bone marrow into tube, label with patient's name

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88374 (1)

Scarlet Eligible: No

XL RARA (17q21) BY FISH - TC (NON-NY)

J160-3

Alternate Name: XL RARA (17q21) FISH-TC

Methodology: FISH

Preferred Specimen: Bone Marrow - Green Top

Collection Instruction: Place bone marrow into tube, label with patient's name

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88374 (1)

Scarlet Eligible: No

Xylose Absorption, Urine

0530-6

Alternate Name: XYLOSE ABSORPTION,UR

Methodology: Colorimetry

Preferred Specimen: 5 mL Urine Cup

Collection Instruction: Collect from a 5 hr collection transferred into a gray top tube. Record date and time of dose of xylose and the 5 hr total volume of urine collection.

Storage Instruction: Refrigerate

Turn Around Time: 8 days

CPT Code: 84620 (1)

Scarlet Eligible: No

Xylose, (Baseline) Whole Blood

1184-1

Alternate Name: Xylose, Whole Blood

Methodology: Colorimetry

Preferred Specimen: 4.0 mL Grey top- Sodium Fluoride

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 13 days

CPT Code: 84620 (1)

Scarlet Eligible: No

Allergen Yellow Hornet (i5), IgE

0885-4

Alternate Name: Yellow Hornet

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Yellow Jacket (i3), IgE

0884-7

Alternate Name: Yellow Jacket

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

ZAP-70 BY IHC

5410-6

Alternate Name: ZAP-70, IHC

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88342 (1)

Scarlet Eligible: No

Zika Molecular Serum+UrineK559-5

Alternate Name: Zika Molecular Serum+Urine**Turn Around Time:** 2 days**CPT Code:** 87798 (2)**Scarlet Eligible:** Yes

Zika Serum RNA PCR w/reflex to IgMJ405-2

Alternate Name: Zika RNA PCR w/RFX IgM**Methodology:** Transcription Mediated Amplification**Preferred Specimen:** 1 mL Aliquot Serum for Zika testing**Collection Instruction:** Fill Serum Separator (SST) tube completely, invert gently 2-3 times, let stand for 20 minutes, spin for 10-15 minutes. Transfer serum into the red plastic aliquot tube. Follow specimen labeling instructions provided.**Storage Instruction:** Refrigerate**Turn Around Time:** 2 days**CPT Code:** 87662 (1)**Scarlet Eligible:** Yes

Zika Serum RNA PCRJ406-0

Alternate Name: Zika Serum RNA PCR**Methodology:** Transcription Mediated Amplification**Preferred Specimen:** 1 mL Aliquot Serum for Zika testing**Collection Instruction:** Fill Serum Separator (SST) tube completely, invert gently 2-3 times, let stand for 20 minutes, spin for 10-15 minutes. Transfer serum into the red plastic aliquot tube. Follow specimen labeling instructions provided.**Storage Instruction:** Refrigerate**Turn Around Time:** 2 days**CPT Code:** 87662 (1)**Scarlet Eligible:** Yes

Anti-Zika Virus Elisa IgM Serum

J497-9

Alternate Name: Zika Serum Virus IgM

Methodology: Enzyme Immunoassay

Preferred Specimen: 1 mL Aliquot Serum for Zika testing

Collection Instruction: Fill Serum Separator (SST) tube completely, invert gently 2-3 times, let stand for 20 minutes, spin for 10-15 minutes. Transfer serum into the red plastic aliquot tube. Follow specimen labeling instructions provided.

Storage Instruction: Refrigerate

Turn Around Time: 7 days

CPT Code: 86790 (1)

Scarlet Eligible: Yes

Zika Urine RNA PCR

J407-8

Alternate Name: Zika Urine RNA PCR

Methodology: Transcription Mediated Amplification

Preferred Specimen: Urine Cup for Zika

Collection Instruction: Collect urine in sterile green container or the Aptima Urine Transport. Follow specimen labeling instructions.

Turn Around Time: 2 days

CPT Code: 87662 (1)

Scarlet Eligible: Yes

Zika Molecular Serum w/rfx and Urine

K583-5

Alternate Name: ZikaSerum w/rfx and Urine

Turn Around Time: 2 days

CPT Code: 87798 (2)

Scarlet Eligible: Yes

Zinc Protoporphyrin (ZPP),Whole Blood, Adult

6289-3

Alternate Name: Zinc Protoporphyrin (ZPP),Whole Blood,ADULT

Methodology: Fluorescence

Preferred Specimen: Tan Top

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 84202 (1)

Scarlet Eligible: Yes

Zinc Protoporphyrin (ZPP),Whole Blood, Child

0407-7

Alternate Name: Zinc Protoporphyrin (ZPP),Whole Blood,CHILD

Methodology: Fluorescence

Preferred Specimen: Lavender top- EDTA

Collection Instruction: Fill lavender-top (EDTA) tube completely, invert 8-10 times. DO NOT SHAKE TUBE!

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 84202 (1)

Scarlet Eligible: No

Zinc, RBC

1336-7

Alternate Name: Zn, RBC

Methodology: Inductively Coupled Plasma Mass Spectrometry

Preferred Specimen: Blue-Royal Blue Top K2 EDTA Plasma(*LAVENDER LINE)

Collection Instruction: Fill tube, invert gently 5-6 times, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 84630 (1)

Scarlet Eligible: No

Zinc, Serum

0285-7

Alternate Name: Zn, Serum

Methodology: Inductively Coupled Plasma Mass Spectrometry

Preferred Specimen: Blue-Royal Blue Trace Element Serum (**RED LINE**)

Collection Instruction: Fill tube, label with patient name

Storage Instruction: Refrigerate

Turn Around Time: 7 days

CPT Code: 84630 (1)

Scarlet Eligible: Yes

Zinc, Urine, Random

2440-6

Alternate Name: Zn, Urine, Random

Methodology: Inductively Coupled Plasma Mass Spectrometry

Preferred Specimen: 7 mL Urine Cup

Collection Instruction: Collect urine in cup, label with patients name

Storage Instruction: Refrigerate

Turn Around Time: 8 days

CPT Code: 82570 (1), 84630 (1)

Scarlet Eligible: Yes

Hepatitis Be Antibody

0540-5

Alternate Name: anti-HBe; HBeAb

Methodology: Enzyme Linked Immunoabsorbance

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 86707 (1)

Scarlet Eligible: Yes

Allergen Cashew (rf202), IgE

1412-6

Alternate Name: cashew nut

Methodology: Immunocap

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Cherries (f242), IgE

1583-4

Alternate Name: cherry

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Cotton-Wood (t14), IgE

0891-2

Alternate Name: cottonwood

Methodology: Immunocap

Preferred Specimen: 3 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

Allergen Milk (f2), IgE

0672-6

Alternate Name: cows milk; cow's milk

Methodology: Immunocap

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86003 (1)

Scarlet Eligible: Yes

dsDNA Antibody IgG, Crithidia IFA

3825-7

Alternate Name: dsDNA(crithidia)Ab IgG

Methodology: Indirect Fluorescence Assay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 86255 (1)

Scarlet Eligible: Yes

HCG, Quantitative, Serum

0327-7

Alternate Name: hCG, Quantitative, Serum; BHCG., Quant.

Methodology: Electrochemiluminescence Immunoassay

Preferred Specimen: 2 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84702 (1)

Scarlet Eligible: Yes

Hcg, Quantitative, Tumor Marker

1201-3

Alternate Name: hCG, Quantitative, Tumor Marker

Methodology: Electrochemiluminescence Immunoassay

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 1 days

CPT Code: 84702 (1)

Scarlet Eligible: Yes

Human Telomerase Reverse Transcriptase - Tech Only

TF87-4

Alternate Name: hTert - Tech Only

Methodology: Immunohistochemistry

Preferred Specimen: ThinPrep Cytolyt Solution

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 88342 (1)

Scarlet Eligible: No

p16 by IHC

J356-7

Alternate Name: p16 by IHC

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Storage Instruction: Room Temp

Turn Around Time: 1 days

CPT Code: 88342 (1)

Scarlet Eligible: No

17P13 BY FISH-TC ONLY

A874-0

Alternate Name: p53 by FISH, TC Only

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube. Also acceptable: Whole blood in green (sodium heparin) tube.

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88367 (1)

Scarlet Eligible: No

p53 by IHC w/digital imaging

5153-2

Alternate Name: p53-TUMOR SUPP.GENE

Methodology: Immunohistochemistry

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in a block or slide form from client with surgical number imprint.

Turn Around Time: 1 days

CPT Code: 88361 (1)

Scarlet Eligible: No

p57 by IHC

J546-3

Alternate Name: p57 by IHC

Methodology: Microscopic Examination

Preferred Specimen: Formalin-fixed, Paraffin-embedded Tissue

Collection Instruction: This comes in block form from client with surgical number imprint.

Storage Instruction: Room Temp

Turn Around Time: 1 days

CPT Code: 88342 (1)

Scarlet Eligible: No

Reverse Triiodothyronine Serum/Plasma (rT3)

T397-9

Alternate Name: rT3

Methodology: Liquid Chromatography Tandem Mass Spectrometry

Preferred Specimen: 1 mL SST Tube

Collection Instruction: Fill tube, invert gently 5 times, label with patient name, let stand for minimum of 30 minutes, maximum of 1 hr, spin for 10-15 minutes

Storage Instruction: Refrigerate

Turn Around Time: 3 days

CPT Code: 84482 (1)

Scarlet Eligible: Yes

Pap (Reflex to HPV DNA Genotyping 16, 18 when ASC-US or >)

6196-0

Alternate Name: rfx HPV DNA when ASCUS or >

Preferred Specimen: ThinPrep Vial

Collection Instruction: Collect specimen, swish in Thin prep vial, label with pt name.

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 88175 (1)

Scarlet Eligible: No

Lymphoma - Mantle Cell (MCL), BCL1/IGH by FISH

5026-0

Alternate Name: t(11;14), CCND1, CYCLIN D1

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube. Also acceptable: Whole blood in green (sodium heparin) tube

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88374 (1)

Scarlet Eligible: No

Lymphoma - Follicular, BCL2/IGH by FISH

5270-4

Alternate Name: t(14;18)

Methodology: FISH

Preferred Specimen: 2 mL Bone Marrow - Green Top

Collection Instruction: Non-diluted bone marrow aspirate collected and transferred into sodium heparin (green top) tube. Also acceptable: Whole blood in green (sodium heparin) tube.

Storage Instruction: Refrigerate

Turn Around Time: 3-5 days

CPT Code: 88377 (1)

Scarlet Eligible: No

Lymphoma - Follicular, BCL2/IGH by PCR

5271-2

Alternate Name: t(14;18), Follicular Lymphoma

Methodology: Polymerase Chain Reaction

Preferred Specimen: 5 Bone Marrow - Lavender Top

Collection Instruction: Place bone marrow into tube, label with patient's name. Send in bone marrow kit.

Storage Instruction: Refrigerate

Turn Around Time: 7 days

CPT Code: 81402 (1)

Scarlet Eligible: No

t(1;19) (q23.3: p13.3) FISH

B611-5

Alternate Name: t(1;19) (q23.3: p13.3)

Methodology: FISH

Preferred Specimen: 1 mL Bone Marrow - Green Top

Collection Instruction: Place bone marrow into tube, label with patient's name

Storage Instruction: Room Temp

Turn Around Time: 5 days

CPT Code: 88271 (2), 88275 (1)

Scarlet Eligible: No

Von Willebrand Factor, Antigen

5815-6

Alternate Name: vWF Antigen**Methodology:** Immunospectrometric**Preferred Specimen:** 1 mL Aliquot Plasma-Citrate-Blue Top**Collection Instruction:** Collect full blue-top (citrate) tube. Centrifuge Blue-top at high speed for 15 minutes. Transfer plasma into plastic transfer tube. Label with date/time collected, PLASMA, patient's name and freeze.**Storage Instruction:** STRICT Frozen**Turn Around Time:** 9 days**CPT Code:** 85246 (1)**Scarlet Eligible:** No

BioReference®

Index by Test Code

Index by Test Name

XCE 8 BY FISH - B514-1	B1
-5/5Q- BY FISH-TC ONLY - A869-0	B1
D5S2064 FOR -5/5q- BY FISH - B507-5	B1
-7/7Q- BY FISH-TC ONLY - A868-2	B2
D7S486/XCE 7 FOR -7/7q- BY FISH - B506-7	B2
GLYCOMARK - J114-0	B2
10-Color Flow Cytometry Test Codes -	B1
10-Color Screening (Short) Panel - TK69	B3
10-Color Screening (Short) Panel TC Only - TK81	B3
ATM FOR 11q22.3 BY FISH - B512-5	B3
13Q14 by FISH TC - A872	B3
D13S319 (13q14.3) BY FISH - B510-9	B4
Estradiol, Serum - 0516-5	B4
Hydroxyprogesterone(17), Alpha Serum/Plasma - 0997-7	B4
P53 by FISH - A314-7	B5
1P/19Q,D(1;19)DELETION BY FISH - A311-3	B5
1q21/CKS1B by FISH - B509	B5
20Q12 BY FISH-TC ONLY - A871-6	B6
D20C108 FOR DEL (20q) BY FISH - B508-3	B6
Osmolality, Urine 24hr. - 0988-6	B6
4K PSA rflx >3.0 - J264-3	B7
4K PSA Total + Free Reflex >3.0 to 4Kscore - K135-4	B7
4K Prostate Score (Serum) - J148-8	B7
5-Flucytosine, S - 6180-4	B8
A. Fumigatus Mold (m3), IgE - 0694-0	B8
A. Vaginalis by RT-PCR - J212-2	B8
Abetalipoproteinemia - B148-8	B9
Onkosight NGS ABL1 Sequencing - J132-2	B9
ABO and Rh Blood Typing - 0156-0	B9
Absolute Neutrophil Count - 3617-8	B10
Cardiolipin Antibodies, IgG, IgM, IgA - 0431-7	B10
Adrenocorticotrophic Hormone (ACTH) - 0507-4	B10
Acute Hepatitis Panel - 3283-9	B11
Acute Phase Reactant Profile - 0584-3	B11
Adenocarcinoma vs. Mesothelioma IHC Panel - 5181-3	B11
GP/AE1/AE3 (Stain Only - 5108-6	B12
Maternal Serum Alpha-Fetoprotein Screen (MSAFP) - ONTD Only - 0274-1	B12

Amniotic fluid AFP - 2122-0.....	B12
Alpha Fetoprotein, Tumor Marker (AFP) - 0025-7	B13
AFP/HCG Tumor Marker - 4514-6	B13
Ashkenazi Jewish Expanded Profile - J293-2.....	B13
Ashkenazi Jewish DNA Panel (NON NY) - A282-6.....	B14
Post Dialysis Albumin - 2256-6.....	B14
Alcohol Biomarkers (ETG), Qual. (500) - 6284-4.....	B14
Aldosterone, Serum - 0302-0.....	B15
ALK (EML4/ALK) BY IHC - B565-3.....	B15
ALK (EML4/ALK) by IHC TC - B603	B15
ALK BY FISH (TC) - B236-1.....	B16
ALK BY FISH, IF NEG REFLEX TO ROS1 BY FISH - J015-9	B16
EML4-ALK by FISH (FDA): Lung - A241-2.....	B16
ALK-2P23: Anaplastic Lymphoma - B684-2.....	B17
ALL Prognosis Panel by FISH - P260-3.....	B17
Allergen Cedar Red (t212), IgE - 6305-7.....	B17
Allergen Chestnut (t206), IgE - 2420-8.....	B18
Allergen Herring (f205), IgE - 6267-9.....	B18
Allergen White Mulberry (T70), IgE - 3835-6.....	B18
Allergy Outdoor Panel - P267-8.....	B19
Alkaline Phosphatase, Serum - 0185-9	B19
Alpha-2-Macroglobulin - 3882-8.....	B19
Alpha-Globin mut. Fetus - B669-3.....	B20
Aluminum, Water, Dialysis - 1817-6	B20
Mitochondrial Antibody (AMA) with titer - 0939-9.....	B20
AML Diagnostic Panel by FISH - P261-1	B21
AML-M4Eo by FISH TC - J380.....	B21
AML-M2 (AML/ETO) by Fish- TC ONLY - B891-3	B21
Amylase,Pleural Fluid - 6218-2.....	B22
Amylase, Urine, Random - 0037-2	B22
Antinuclear Antibody Evaluation (ANA) - 0594-2	B22
ANA with Reflex to IFA - A233-9	B23
Antinuclear Ab. (ANA) W/Reflex - 0038-0	B23
Anal GC/CT ThinPrep (Non-NY) - K987-8.....	B23
Anemia I Profile - 0998-5.....	B24
Anemia II Profile - 0999-3.....	B24
Anemia III Profile - 2073-5.....	B24
Anemia Profile, Expanded - 9893-9	B24

Anthrax Culture Screen - 3573-3.....	B25
Anti-Annexin V. Antibodies, IgG, IgM - 5706-7	B25
Beta-2 Glycoprotein 1 Antibody, IgA - 3905-7	B25
Beta-2 Glycoprotein 1 Antibody, IgG - 3904-0.....	B26
Beta-2 Glycoprotein 1 Antibody, IgM - 3906-5.....	B26
Anti CMV Avidity, IgG - J057-1	B26
Anti-Epidermal Antibody, S - 0184-2	B27
Anti-Mullerian Hormone - 3707-7	B27
Parietal Cell Antibody - 0263-4	B27
Antiphosphatidylserine Antibodies, IgG, IgM, Special Coagulation - 5790-1	B28
Smith Antibody - 0851-6.....	B28
Antithrombin III Activity (AT III) - 5714-1	B28
Antithrombin III Antigen (ATIII) - 5715-8	B29
Arsenic, Blood - J019-1.....	B29
Arthritis Profile - 0325-1	B29
Arthritis Profile - Expanded - 9892-1	B30
Ashkenazi Jewish 11 - P658-8	B30
Ashkenazi Jewish 18 - P659-6	B30
Ashkenazi Jewish 25 - Q161-2.....	B31
Ashkenazi Jewish DNA Carrier - 3267-2	B31
Ashkenazi Jewish DNA (NY patients only) - P223-1	B31
Antistreptolysin O (ASO) Titer - 0042-2	B32
Aspergillus Antibodies, Immunodiffusion, Serum - 1741-8	B32
Salicylates, Serum - 0526-4.....	B32
Atenolol Level - 2392-9	B33
Atherosclerotic Risk Expanded - 6802-3	B33
Antithrombin III Anti-Xa Assay - 5723-2	B33
11Q22.3 BY FISH-TC ONLY - A879-9.....	B34
Automated UroVysion By FISH - 5030-2	B34
AV Panel w/o Lacto Reuteri and Rhamnosis (Non NY) - L767-3	B34
Abacavir Genotype HLA B57 - J096-9	B35
Abnormal APTT Assessment - 5851-1	B35
Abnormal PT and APTT Assessment - 5853-7	B35
Eosinophil Count, Absolute, Blood - 0381-4.....	B36
Activation Markers - 5811-5.....	B36
Acute 10-Color Leukemia/Lymphoma Panel - TH61	B36
Acute 10-Color Leukemia/Lymphoma Panel TC Only - TH64	B36
Acute Leukemia + Myeloid Disorders Analysis (Technical Only) and StormPath BMM Report	

Template - R735-2	B37
Acute Leukemia/MDS IHC Panel - 5116-9	B37
AML - M2 (AML/ETO) (RUNX1/RUNX1T1) by FISH - 5024-5	B37
AML - M2 (AML/ETO) by RT-PCR - 5034-4	B38
AML - M4Eo by FISH - 5025-2	B38
AML - M3 (PML/RARA) by RT-PCR - 5261-3	B38
AML - M3 (PML/RARA) by FISH - 5260-5	B39
Acylglycines, Quant. Urine - 6310-7	B39
OnkoSight Advanced CNS NGS Panel - TH57-3	B39
OnkoSight Advanced Gastrointestinal Stromal Tumor NGS Panel - TH58-1	B40
OnkoSight Advanced GYN Tumor NGS - TH53-2	B40
OnkoSight Advanced Lung Cancer NGS Panel - TH49-0	B40
OnkoSight Advanced Melanoma NGS Panel - TH59-9	B41
OnkoSight Advanced Prostate NGS Panel - TH48-2	B41
OnkoSight Advanced Solid Tumor NGS Panel - TJ16-5	B41
Aerobic Vaginitis Profile (Non NY) - R732-9	B41
Aluminum, Serum - 0378-0	B42
Albumin, Serum - 0033-1	B42
Alcohol, Methyl - A473-1	B42
Allergen Lime (RF306), IgE - 3522-0	B43
Allergen Yeast Bakers, IgE - 0813-6	B43
Digoxin, Serum - 0083-6	B43
Allergen Almond Nut (f20), IgE - 0815-1	B44
Allergen Alpha Lactalbumin (f76), IgE - 1406-8	B44
Vitamin E, Serum - 0447-3	B44
Triple Screen AFP3 (AFP, UE3, HCG) - 0825-0	B45
Alport Syndrome, Autosomal Recessive - J316-1	B45
Allergen Alternaria Tenuis Mold (m6), IgE - 0660-1	B45
Theophylline, Serum - 0380-6	B46
Barbiturates, Urine, Screen (200) - 0921-7	B46
Amylase, Serum - 0036-4	B46
Anal CT on ThinPrep (Non NY) - J792-3	B47
Anal GC on ThinPrep (Non NY) - J791-5	B47
Allergen Anchovy (rf313), IgE - 1463-9	B47
PCP (6-Phencyclidine), Urine, Screen (25) - 0926-6	B48
Animal Mix (ex1), IgE - 2180-8	B48
Animals Screen(6) - 0902-7	B48
Endomysial Antibody (IgA) Screen with Titer - 1206-2	B49

Anti-Jo-1, IgG - T974-5	B49
Double-Stranded DNA (dsDNA) Antibody - 0364-0	B49
Antibody Screen (Ind.Coombs) - 0039-8	B50
Apolipoprotein A-1 (Apo A-1) - 0451-5	B50
Apolipoprotein A-1 and B (Apo A1/B) - 0895-3.....	B50
Apolipoprotein B (Apo B) - 0457-2	B51
Allergen Apple (f49), IgE - 1134-6	B51
Allergen Apricot (f237), IgE - 2344-0	B51
Vaginitis/Vaginosis by RT-PCR, Aptima Tube - L050-4.....	B52
Arthrogryposis, Mental Retardation, and Seizures - J297-3.....	B52
Arsenic, Urine, 24 Hours - 0427-5	B52
Vitamin C - 0291-5.....	B53
Lorazepam (Ativan)), Urine - 5924-6.....	B53
Allergen Australian Pine, IgE (t73) - 3837-2	B53
ISH Single Antibody Stain Only - 5289-4	B54
Allergen Avocado (f96), IgE - 1945-5	B54
B. fragilis by RT PCR - J216-3	B54
Beta-2-Microglobulin, Serum - 0262-6	B55
Beta-2-Microglobulin, Urine, Random - 1754-1	B55
Vaginitis/Vaginosis, DNA - BD Affirm - 5093-0	B55
Basic Antiphospholipid/Lupus - 5994-9	B56
Basic Food Allergy Panel - 4633-4.....	B56
BCLC (3Q27) by FISH TC - B348	B56
BCR/ABL by FISH TC t(9:22) - A880.....	B56
BCR/ABL Quantification by RT-PCR - 5858-6	B57
Vaginitis/Vaginosis by DNA probe, Thin Prep - P576-2.....	B57
Bean Profile - Q032-5	B57
Benzodiazepines, Urine, Screen (100) - 3233-4	B58
Beta Carotene, Total - 0512-4.....	B58
Beta-2-Glycoprotein I (G,A,M) - 2675-7.....	B58
Beta-Galactosidase, Leukocytes - 3514-7	B59
Biopsy Level II, PC ONLY - 5233-2	B59
Biopsy Level III, PC ONLY - 5234-0	B59
Biopsy Level IV, PC ONLY - 5235-7	B60
Biopsy Level V, PC ONLY - 5236-5	B60
Biopsy Pod 2-@ _ - 5002-1.....	B60
Biopsy Pod 3-@ _ - 5003-9.....	B61
Biopsy Pod -@ _ - 5000-5.....	B61

Biopsy, general-@ (AREA) - 4000-6.....	B61
Biopsy Gen 2-@ (AREA) - 4007-1	B62
Biopsy Gen 3@AREA - 4003-0.....	B62
Biopsy Gen 4-@ (AREA) - 4004-8	B62
Biopsy Gen 5-@ (AREA) - 4005-5	B63
Biopsy Gen 6-@ (AREA) - 4006-3	B63
Biopsy, GI (1 Jar) - 5631-7	B63
Biopsy, GI (10 JARS) - A296-6.....	B64
Biopsy, GI (2 Jars) - 5632-5	B64
Biopsy, GI (3 Jars) - 5633-3	B64
Biopsy, GI (4 Jars) - 5634-1	B65
Biopsy, GI (5 Jars) - 5635-8	B65
Biopsy, GI (6 Jars) - 5636-6	B65
Biopsy, GI (7 JARS) - A293-3.....	B66
Biopsy, GI (8 JARS) - A294-1	B66
Biopsy, GI (9 JARS) - A295-8.....	B66
Biopsy, Level VI Global - 5237-3	B67
Liver Biopsy - 5009-6	B67
Bleeding Diatheses Evaluation Profile I - 5850-3.....	B67
Bleeding Diathesis Evaluation Profile II - 5997-2	B68
Bloom Syndrome Carrier Test - A658-7	B68
Bone Marrow Morph.- Tech Only - B846-7.....	B68
Bone Marrow Morph. Clot Only- Tech Only - B843-4.....	B69
Bone Marrow Morph. Core Only- Tech Only - B845-9	B69
Bone Marrow Morph. Smear Only- Tech Only - B844-2	B69
Basic Metabolic Profile (BMP) - 2555-1	B70
B-Type Natriuretic Peptide (BNP) - 3553-5	B70
Bone Marrow Morphology - 5199-5	B70
Bone Marrow Analysis - Clot Only - 5202-7	B71
Bone Marrow Morphology - Core Only - 5200-1	B71
Bone Marrow Slide Consult - 5207-6.....	B71
X/Y Bone Marrow Transplant Monitoring - 5066-6	B72
Bone Marrow Analysis - Smears - 5211-8.....	B72
Bone/Tissue Culture - P210-8.....	B72
BRAF, if neg reflex to MLH1 Hypermethylation- Non NY - J107-4	B73
BRAF, if neg reflex to MLH1 Hypermethylation NY Only - J463-1	B73
BRAF Sequencing- Non NY - J021-7.....	B73
BRAF Sequencing NY Only - J456-5.....	B74

OnkoSight BRAF Sequencing - BM/PB - T979-4.....	B74
BRAF/EGFR NGS if EGFR neg rfx to ALK FISH if neg rfx ROS1 - TB34-4.....	B74
BRCA1/2 Sequencing and Del/Dup Analysis - B362-5.....	B75
Breast Carcinoma/Lymph Nodes/MM Metastasis - 5170-6.....	B75
Breast Carcinoma, R/O Microinvasion BY IHC-TC - A952-4.....	B75
BREAST SMEAR (slide) - 1301-1.....	B76
Blood Urea Nitrogen (BUN), Serum - 0049-7.....	B76
BUN, Post Dialysis - 1679-0.....	B76
BUN + Creatinine with Ratio - 4030-3.....	B77
Buprenorphine (Suboxone), Urine, Screen (5) - 3956-0.....	B77
Bacterial Vaginosis and Aerobic Vaginitis (Non NY) - L571-9.....	B77
BV PANEL BY RT-PCR - H870-9.....	B77
BV Panel by RT-PCR - K835-9.....	B78
Bacterial Vag. W/RFX TP-NY. - H575-4.....	B78
BV If Negative Reflex to Aerobic Vaginitis (Non NY) - L570-1.....	B78
BVAB2 by RT-PCR - J215-5.....	B78
Allergen Bahia Grass (g17), IgE - 1825-9.....	B79
Allergen Bamboo Shoot (f51), IgE - 3642-6.....	B79
Allergen Banana (f92), IgE - 0863-1.....	B79
Phenobarbital, Serum - 0289-9.....	B80
Bardet-Biedl Syndrome 2 - J317-9.....	B80
Allergen Barley (f6), IgE - 0816-9.....	B80
Lipid Screen (Basic Lipid Profile) - 0009-1.....	B81
Allergen Basil (rf269), IgE - 1461-3.....	B81
Allergen Beech Tree, IgE (t5) - 0948-0.....	B81
Allergen Beef (f27), IgE - 0861-5.....	B82
Allergen Beetroot (f319), IgE - 3803-4.....	B82
Immunotyping, Urine - 1644-4.....	B82
Amphetamines, Urine, Screen (1000) - 0920-9.....	B83
Allergen Bermuda Grass (g2), IgE - 1198-1.....	B83
Allergen Betalactoglobulin (f77), IgE - 1268-2.....	B83
Pregnancy, Urine - 0133-9.....	B84
Procainamide, Ser.(+napa #0525) - 1846-5.....	B84
Special Stain, Halls Bile - TC - T409-2.....	B84
Bilirubin, Direct, Serum - 0044-8.....	B85
Bilirubin, Total + Direct, Serum - 0045-5.....	B85
Bilirubin, Total, Serum - 0043-0.....	B85
Quinidine, Serum - 0140-4.....	B86

Allergen Black Pepper (f250), IgE - 1422-5.....	B86
Allergen Blackberry (rf211), IgE - 1486-0.....	B86
GP/Non-Automated Urovysion BY FISH - 6274-5	B87
Allergen Blomia Tropicalis (d201), IgE - B235-3	B87
Culture, Blood (spec.#2) - 1350-8	B87
Culture, Blood (spec.#3) - 1351-6	B88
Culture, Blood - 0341-8	B88
Malaria/Babesia/Other Blood Parasites - J897-0	B88
Occult Blood,Stool (1st Specimen) - 0122-2	B89
Occult Blood,Stool (3rd Spec.) - 0323-6.....	B89
Occult Blood,Stool (2nd Spec.) - 0319-4.....	B89
Allergen Blue Mussel (f37), IgE - 0697-3.....	B90
Allergen Blueberry (rf228), IgE - 1487-8	B90
Culture, Body Fluid W/Gram Stain - 0556-1	B90
Ostase - Alkaline Phosphatase Bone Specific - T447-2	B91
Allergen Brazil Nut (f18), IgE - 1405-0.....	B91
Breast Basal Layer Triple Stain by IHC Global - TA18-9	B91
ER / PR / HER2 (Tech Only)w/digital imaging - 5127-6.....	B92
Breast Cancer Management Panel - J055-5.....	B92
Allergen Broccoli (f260), IgE - 1229-4.....	B92
Allergen Brome Grass (G11), IgE - 2394-5.....	B93
Breast Basal Layer Triple Stain by IHC Tech Only - TA17-1	B93
Allergen Brussel Sprouts (rf217), IgE - 1814-3.....	B93
Allergen Buckwheat (f11), IgE - 1629-5.....	B94
Clostridium difficile Toxin/GDH antigen with Reflex to PCR - 0259-2.....	B94
C-Peptide, Serum - 0578-5.....	B94
C-Reactive Protein (hs-CRP), High Sensitivity - 3320-9	B95
C-Reactive Protein (CRP), Serum - 0069-5.....	B95
Candida glabrata by PCR - T783-0	B95
Candida dubliniensis by PCR - T786-3.....	B96
Candida krusei - T449-8	B96
Candida parapsilosis - T536-2	B96
Candida tropicalis by PCR - T745-9	B97
Chlamydia Trachomatis, Urine, rRNA - 6368-5	B97
Complement C3, Serum - 0532-2.....	B97
Complement C4, Serum - 0533-0.....	B98
CA 15-3, Serum - 2130-3.....	B98
CA 19-9, Serum - 0535-5.....	B98

CA 27.29, Serum - 0823-5.....	B99
CA 125, Serum - 0536-3.....	B99
Cadmium, Blood - J018-3	B99
Vitamin D, 25-Hydroxy, Serum - 0286-5.....	B100
CALL RESULT TO DOCTOR - 9990-3.....	B100
Canavan Disease Test - A657-9	B100
Vulvovaginal Candidiasis Panel by PCR w/ reflex to species - B976-2.....	B101
Breast Carcinoma, Rule Out Microinvasion, with Interp. - 5176-3.....	B101
Cathepsin D by IHC - 5151-6.....	B101
CCP (CYCLIC CITRULLINATED PEPTIDE) IGG, SERUM - 3747-3.....	B102
T-NK CD57 Flow Cytometry Analysis - 5221-7.....	B102
Clostridium Difficile Toxins A+B Genes By DNA Amplification - B137-1.....	B102
Carcinoembryonic Antigen (CEA) - 0055-4.....	B103
CEA, Peritoneal Fluid - 6169-7	B103
AML - CEBPA Mutational Analysis - 8860-9.....	B103
Celiac Disease Antibody Evaluation - 1186-6	B104
Celiac Disease Eval + HLA Genetics - K002-6.....	B104
Cell Count, Fluid - 0056-2	B104
Centromere Autoantibodies - 0587-6	B105
Cervicitis W/O PAP W/O HPV PLUS - Q265-1	B105
CF/Fragile X/SMA - 2115-4.....	B105
Cystic Fibrosis (CFTR) Intron 9 Poly-T Analysis (5T/7T/9T) - J520-8.....	B106
Chem 12 Profile - 0004-2.....	B106
Chem 7 Panel - 0768-2	B106
Chlamydia Trachomatis by Multiplex PCR - 6337-0	B107
Chloramphenicol - 3792-9	B107
Cholinesterase, Serum - 0545-4.....	B107
Chromosome Analysis- Breakage, Fanconi Anemia - B534-9.....	B108
Chromosome Analysis w/Reflex Microarray - 5394-2.....	B108
CPK Isoenzymes - 0558-7	B108
CK-903 for Prostate Cancer - 5434-6.....	B109
CPK (CK), Total - 0068-7	B109
Creatine Kinase-MB Fraction(CKMB) - 2922-3.....	B109
1Q21 BY FISH-TC ONLY - A876-5.....	B110
OnkoSight Advanced Chronic Lymphoid Neoplasm NGS Panel - TH55-7.....	B110
CLL FISH Panel - TC Only - P962-4.....	B110
Minimal Residual Disease for CLL/SLL - 5155-7.....	B111
CMV Ab.,IGG *POS. Reflex To IGM - 2377-0	B111

Cytomegalovirus (CMV) by Multiplex PCR (NON-NY) - 6353-7	B111
Cytomegalovirus (CMV), DNA, Ultraquant - 1161-9.....	B112
Carbon Dioxide, Serum (CO2) - 0052-1	B112
COA/BMP - K935-7	B112
AtheroSclerotic Risk (Coag) - 5874-3	B113
Coagulation Profile - 0994-4	B113
Cocaine, Urine, Screen (100) - 1478-7	B113
Cocaine, Urine, Screen (300) - 0923-3	B114
Collagen Binding Assay - 5744-8.....	B114
Colon Biopsy (2 Specimens) - 5652-3	B114
Comprehensive Heart Health Baseline - F306-6.....	B115
Comprehensive Heart Health Monitoring - F307-4.....	B115
Comprehensive Metabolic Panel - 3427-2.....	B115
Comp. Metabolic/GGTP - 9150-4.....	B116
Comprehensive Slide Consult - 5256-3	B116
Cone Biopsy/LEEP - 5672-1	B116
Congo Red Amyloid Stain - 5011-2	B117
Coronary Risk Profile II - 0993-6.....	B117
Cortisol, AM, baseline - 0066-1.....	B117
Novel Coronavirus COVID-19 Nasal/Nasopharynx - TH68-0.....	B118
Novel Coronavirus COVID-19 Oropharynx - TH69-8	B118
Novel Coronavirus COVID-19 Pooled N/NP/OP - TH71-4	B118
Carnitine Palmitoyltransferase Def. Type 2 Carrier Screening - J315-3.....	B119
OnkoSight Advanced Colorectal Cancer NGS Panel - TH50-8.....	B119
Post Dialysis Creatinine - 1687-3	B119
PDF Creatinine (Timed Test code) - 2063-6	B120
Creatinine, Arterial - 4534-4	B120
Creatinine Clearance (24Hr. Urine) - 0072-9	B120
Creatinine, Systemic - 4535-1	B121
Creatinine, Urine, 24 HR. - 0071-1	B121
Cerebrospinal Fluid (CSF) Cytology - 5112-8	B121
Chlamydia (CT) PCR Urine - T407-6.....	B122
Chlamydia, Anal Aptima - A114-1	B122
Chlamydia, Oral Aptima - A113-3.....	B123
Chlamydia, swab, TMA, Aptima - 1004-1	B123
High Risk STI CT/GC - 3851-3	B123
Chlamydia/N. Gonorrhoea CT/GC PCR Urine - L344-1	B124
Chlamydia + Gonorrhoea, Urine, Aptima - 0257-6	B124

Chlamydia + Gonorrhoeae, Anal, Swab, Aptima - P011-0	B124
Chlamydia, Gonorrhea, and Trichomonas rRNA,TMA Aptima Urine - H001-1	B125
High Risk STI CT/GC/Trich/M.Gen. - P006-0	B125
Culture, Catheter Tip - J010-0	B125
Culture, Mouth - 0960-5	B125
Culture, Stool (2nd Specimen) - 0226-1	B126
Culture Stool (3rd Specimen) - 0225-3	B126
Culture,ID Only (Any Source) - 1043-9	B126
Cyclospora ID, Modified Acid Fast - 1468-8	B127
PAP Smear, 2 Slides, Autopap - 1102-3	B127
Cytology, Respiratory System - 1402-7	B127
Cytology, Thyroid By FNA - 4300-0	B127
Cytology, Urine With Interpretation - 2100-6	B128
Cytopatho. Urine-PC (Alcohol) - 5314-0	B128
Cytopathology Urine-TC (Thinprep) - 5218-3	B128
Anal Cytopathology, ThinPrep - 6626-6	B129
Calcium, Ionized, Serum - 0051-3	B129
Calcium, Serum - 0050-5	B129
Calcium, Urine, 24 Hr. - 0359-0	B130
Calcium, Urine, Timed Specimen - 1660-0	B130
Vitamin D, 1, 25-Dihydroxy, Serum - 0288-1	B130
Chromosome Analysis (Oncology)- Tissue - J913-5	B131
Allergen C.Albicans (m5), IgE - 0892-0	B131
Candida Species by M. PCR - Swab - L416-7	B131
Allergen Cardamon (rf267), IgE - 2308-5	B132
Homocysteine Level - 1311-0	B132
Allergen Carrots (f31), IgE - 0934-0	B132
Allergen Cat Dander (e1), IgE - 0665-0	B133
Urine Culture, Catheterized - 2603-9	B133
Allergen Cauliflower (rf291), IgE - 1815-0	B133
Allergen Celery (f85), IgE - 1380-5	B134
Cereal Mix (fx3), IgE - 2185-7	B134
Culture, Cerebrospinal Fluid - 0966-2	B134
Protein Electrophoresis, CSF - 3946-1	B135
Cervical Biopsy - 5660-6	B135
Allergen Cheddar Cheese (f81), IgE - 0850-8	B135
Allergen Cheesemolds (f82), IgE - 1249-2	B136
Chem 29 Profile - 0165-1	B136

CH24/HDL,CBC/D/PLT - 0855-7	B136
Allergen Chick Pea (f309), IgE - 1163-5.....	B137
Allergen Chicken Feathers (e85), IgE - 2570-0.....	B137
Allergen Chicken Meat (f83), IgE - 0817-7	B137
Childhood Allergy Profile - K726-0.....	B138
Allergen Chili Pepper (rf279), IgE - 1423-3.....	B138
Chlamydia NAAT - 3852-1	B138
Chlamydia Trachomatis MultiPlex PCR - swab - J835-0	B139
High Risk-STI CT/GC/Trich-TP/SP - 6116-8	B139
Chlamydia + Gonorrhea ,Oral, Swab, Aptima - P010-2.....	B139
Chlamydia and Gonorrhea, rRNA,TMA Aptima Swab - 1005-8.....	B140
Allergen Chocolate Cocoa (f93), IgE - 1007-4.....	B140
Cholesterol - 0058-8.....	B140
Chromium, Plasma/Serum - 0277-4.....	B141
Chromogranin A - TA34-6	B141
Karyotype CVS - A587-8.....	B141
Chromosome Analysis, POC, Tissue - 1053-8.....	B142
Chromosome Analysis,Amnio.FLD(IG) - 2136-0	B142
Chromosome Analysis/Karyotype, Blood - 0559-5.....	B142
Chronic 10-Color Lymphoproliferative Panel - TH66	B143
Chronic 10-Color Lymphoproliferative Panel TC Only - TH65	B143
Chronic Fatigue Syndrome - F100-3	B143
Lymphoproliferative Disorders Analysis (TC Only) and StormPath BMM ReportTemplate - R734-5.....	B143
Hepatitis Comprehensive II - 2302-8.....	B143
Lithium, Serum - 0119-8.....	B144
Cyclosporine, Whole Blood - 0022-4	B144
Allergen Coriander (rf317), IgE - 2361-4.....	B144
Allergen Cinnamon (rf220), IgE - 2345-7	B145
Chloride, Serum - 0057-0	B145
Chloride, Urine, Random - 1794-7	B145
Chloride, Urine 24HR. - 0439-0.....	B146
Allergen C. Herbarum (m2), IgE - 0661-9	B146
Allergen Clam (f207), IgE - 1169-2.....	B146
Drug Confirmation, Clonazepam, Urine (LC/MS) - L204-7.....	B147
Allergen Clove (rf268), IgE - 2350-7.....	B147
CoEnzyme Q10 Deficiency, Primary, 7 - J441-7	B147
Allergen Cocklebur (w13), IgE - 1939-8.....	B148

Allergen Cockroach (i6), IgE - 0716-1	B148
Allergen Coconut (f36), IgE - 1318-5	B148
Allergen Codfish (f3), IgE - 0840-9	B149
Coenzyme Q10 - J322-9	B149
Allergen Coffee (rf221), IgE - 1133-8	B149
Cold Agglutinins Titer - 0063-8	B150
B-CrossLaps/Serum - T421-7	B150
Colorectal Cancer Panel - B274	B150
Allergen Common Birch (t3), IgE - 0685-8	B151
Allergen Common Ragweed (w1) ,IgE - 0717-9	B151
Allergen Common Reed (G7) IgE+ - 2395-2	B151
Comprehensive Thyroid-Expanded - 9891-3	B152
CBC w/Diff, Platelet Ct. - 0053-9	B152
CBC w/o DIFF (Hemogram)/PLAT. CT - 0034-9	B152
Cortisol, Serum AM and PM - 0374-9	B153
Cortisol, Random - 0900-1	B153
Cortisol, PM - 0067-9	B153
Comprehensive Common Cancer Panel - B275-9	B154
Congenital Amegakaryocytic Thrombocytopenia (CAMT) - J306-2	B154
Conventional Pap Slide(s) - 1100-7	B155
Allergen Corn (maize)(f8), IgE - 0821-9	B155
Allergen Cotton-Fiber (o1), IgE - 1436-5	B155
Warfarin, Pharmacogenomic Assay - 6261-2	B156
Allergen Cow Dander (e4), IgE - 0692-4	B156
Allergen Crab (f23), IgE - 1168-4	B156
Allergen Cranberry (rf341), IgE - 3796-0	B157
Creatinine, Serum - 0070-3	B157
Cryptosporidium EIA - 0348-3	B157
Cryptosporidium ID, Modified Acid Fast - J745-1	B158
Crystal Identification, Body Fluid - 0579-3	B158
Copper, Plasma/Serum - 0479-6	B158
Allergen Cucumber (f244), IgE - 1664-2	B159
Allergen Cultivated Rye (G12), IgE - 2406-7	B159
Allergen Cultivated Wheat (g15), IgE - 3221-9	B159
Culture, Actinomyces - J011-8	B160
Medical Implants/Devices Culture - TG39-3	B160
Culture, Anal (Beta strep. Group B) - 3773-9	B160
Culture, Environmental - 0819-3	B161

Culture, Tracheal Airway - K909-2.....	B161
Allergen Curry (rf281), IgE - 2356-4.....	B161
Allergen Curvularia Lunata (m16), IgE - 2476-0	B162
Cystatin C - J013-4	B162
Cystic Fibrosis 40- NGS - J435-9.....	B162
Cytology, Bladder Washing (50% Alcohol) - J529-9	B163
Cytology, Bladder Washing (ThinPrep) - J534-9	B163
Oncology Chromosome Analysis- BM/PB - 5250-6	B163
CMV IGM - 0461-4	B164
CMV Ab. (IGG) - 0400-2	B164
Cytopath Conc Tech - TC only - 5210.....	B164
Cytopathology Non-Gyn TC only - B003.....	B164
D-Dimer - 5743-0.....	B165
Allergen D. Pteronyssinus (d1), IgE - 0715-3.....	B165
13Q14.3 BY FISH-TC ONLY - A878-1.....	B165
Decalcification (Nails/Bone) - 3188-0	B166
Dehydroepiandrosterone (DHEA) - 0178-4	B166
DHEA Sulfate - 0406-9	B166
Dihydrotestosterone - J892-1	B167
Diabetes Profile - 0391-3.....	B167
Diabetic Screen - 0222-0	B167
DIC Level 2 - 5872-7	B168
DIC Assessment Panel Level-1 - 5839-6	B168
DIC Assessment - 5992-3	B168
Dihydropyrimidine Dehydrogenase - A644-7	B169
DLBCL FISH PANEL: MYC BA, BCL2, BCL6 - L218-7	B169
Lymphoma - Diffuse large B-Cell, BCL6 by FISH - 5028-6.....	B169
DNA Ploidy for Molar Pregnancy - 5555	B170
DNA Ploidy for Molar Pregnancy + P57 - K735	B170
Dihydropyrimidine Dehydrogenase (DPD) Gene Mutation Analysis - 6285-1	B170
Drug Abuse Screen, Urine. (8 drug) - 0628-8.....	B170
Allergen Casein (f78), IgE - 1267-4	B171
Allergen Dandelion (w8), IgE - 1363-1	B171
Deafness, Autosomal Recessive 1A (GJB2) - J308-8	B171
Gliadin Antibodies, IgG, IgA - 0608-0	B172
Meperidine, Urine, Screen, Qual. (200) - 8545-6.....	B172
Valproic Acid - 0542-1	B172
Allergen D. Farinae (Dust Mite) (d2), IgE - 1372-2	B173

Culture, Dialysate - 0857-3.....	B173
Culture, Dialysis Water - 0856-5.....	B173
Allergen Dill (f277), IgE - 3605-3	B174
DRVVT Complete - 5730-7	B174
Coombs Test, Direct - 0064-6	B174
Sickle Cell Screen - 0366-5	B175
Allergen Dog Dander (e5), IgE - 0664-3	B175
Allergen Duck Feathers (e86), IgE - 0846-6	B175
Dust Mix (hx2), IgE - 2182-4	B176
Dyskeratosis Congenita, Autosomal Recessive 5 - J299-9	B176
Escherichia coli (Non NY) - J968-9	B176
Enterococcus faecalis (Non NY) - J967-1	B177
EDDP (Methadone Metabolite), Urine, Screen, Qual. (300) - 6286-9	B177
EGFR, if neg reflex to ALK- Non Ny - J105-8	B177
EGFR, if neg reflex to ALK NY Only - J461-5	B178
EGFR, if neg reflex to ALK - TC only Non NY - J202-3	B178
EGFR, if neg reflex to ALK-TC only- NY Only - J466-4	B178
EGFR, if neg reflex to ALK, if neg reflex to ROS1- Non NY - J106-6.....	B179
EGFR, if neg reflex to ALK, if neg reflex to ROS1 NY Only - J460-7	B179
EGFR by FISH - A240-4	B179
EGFR Sequencing NY Only - J454-0	B180
EGFR and KRAS, if both neg, reflex ALK, if neg reflex ROS1 Non NY - H579-6	B180
EGFR by NGS, ALK, ROS1, and MET by FISH- Non NY - H581-2.....	B180
EGFR by NGS, ALK, ROS1, and MET by FISH (Ny Only) - K335-0.....	B180
EGFR and KRAS, if both neg, reflex ALK, if neg reflex ROS 1 NY Only - J465-6.....	B181
Electrolytes (24 hr urine) - 8017-6.....	B181
Electrolytes, Urine, Random - 0989-4.....	B181
Encainide + Metabolite Level - 0410-1	B182
Endo-Cervical Biopsy - 5662-2.....	B182
Endometrial Biopsy - 5664-8	B182
Enterococcus/VRE Screen-@ - 6108-5	B183
Enterovirus RNA, Qualitative - 3788-7.....	B183
Erythropoietin, Serum - 0183-4.....	B183
Epstein-Barr Virus By PCR - A315-4	B184
ER,PR,KI-67,P53,HER2 FISH, TC - A278-4	B184
ER / PR / Ki-67/ HER2 (Tech Only) - 5433-8	B184
ER/PR (Tech Only) - 5425-4	B185
ER/PR - w/Interpretation w/digital imaging - 5198-7.....	B185

ER/ PR/ DNA/ KI-67/ HER2 - Reflex to HER2 by FISH - 5408-0	B185
ER/PR /DNA/ Ki-67 - w/Interp w/digital imaging - 5160-7	B186
ER/PR/HER2 - w/Interpretation w/digital imaging - 5161-5	B186
ER/PR/HER2w/digital imaging- Rflex to HER2 by FISH - 5406-4	B186
ER/ PR/ KI-67/ HER2 - Reflex to HER2 by FISH - 5407-2	B187
ER/PR/ Ki-67/ HER2 - w/Interp w/digital imaging - 5163-1	B187
ER/PRKI-67/HER2/P53 BY IHC with Interp. w/digital imaging - A954-0	B187
ER/PRKI-67/HER2/P53 BY IHC,TC - A951-6	B188
ERG by IHC (EP111) - T411-8	B188
ERG by IHC (EP111) - Tech Only - T412-6	B188
Esophageal Biospy X1 - 5630-9	B189
Estradiol (E2), Ultrasensitive, LC/MS, Serum - J258-5	B189
Estrogen Panel - 2317-6	B189
Estrone (E1), Serum - J259-3	B190
Drug Confirmation, Ethanol, Urine (GC/HS) - 6239-8	B190
Expanded CF, Fragile X, SMA - P150-6	B190
Expanded Anti-Phos/Lupus - 5990-7	B191
ExtraTube(s) Drawn - B978-8	B191
Culture, Ear - 0087-7	B191
Culture, Ear (Left) - 0462-2	B192
Culture, Ear (Right) - 0634-6	B192
Egg Components IgE - H992-1	B192
Allergen Egg White (f1), IgE - 0671-8	B193
Egg White w/Reflex to Component Panel - M066-9	B193
Allergen Egg Yolk (f75), IgE - 0849-0	B193
Allergen Eggplant (rf262), IgE - 2265-7	B194
Ehlers-Danlos Syndrome, Type VIIC - J294-0	B194
Electrolytes - 0002-6	B194
Protein Electrophoresis - 0085-1	B195
Allergen Elm Tree (t8), IgE - 0659-3	B195
Allergen Eng. Plantain (w9), IgE - 0650-2	B195
Nasal Smear (for Eosinophils) - 0248-5	B196
Eosinophil Ct, Urine - 1183-3	B196
EBV Panel (Capsid/EBNA G+M) - 0572-8	B196
EBV, Early Antigen Ab - 0582-7	B197
EBV, Nuclear Antigen Ab, IgG - 0583-5	B197
EBV Capsid Ab, IgG - 0234-5	B197
EBV Capsid Ab, IgM - 0580-1	B198

ESR (Sed-Rate) - 0086-9	B198
Estriol, Unconjugated - 3116-1	B198
Estrogens, Fractionated - H906-1	B199
Drug Confirmation, Ethanol, Blood (GC/HS) - B784-0	B199
Alcohol (Ethanol), Blood - 0985-2.....	B199
Alcohol (Ethanol), Urine, Screen (20) - 0995-1	B200
Allergen Eucalyptus (t18), IgE - 1480-3	B200
EWSR1(22Q12) Gene Rearrang. FISH - B556-2	B200
Expanded Cystic Fibrosis Panel - J434-2	B201
Factor II Activity - 5724-0	B201
Factor IX Activity - 5769-5.....	B201
Factor V Activity - 5725-7	B202
Factor VII Activity - 5728-1	B202
Factor VIII Activity - 5729-9.....	B203
Factor X Activity - 5735-6.....	B203
Factor XI Activity - 5737-2.....	B204
Factor XII Activity - 5738-0.....	B204
Factor V Mutation (Leiden) - 5726-5	B204
Factor X Antigen - 5736-4.....	B205
Factor X Chromogenic - 5707-5.....	B205
Factor XIII Activity - 5739-8.....	B205
Familial Hyperinsulinism - A643-9	B206
Familial Dysautonomia - A697-5	B206
Fanconi Anemia, Type C - A656-1	B206
Fax Result To Doctor - 9993-7	B207
Glucose, Fasting - 0095-0	B207
Fecal Calprotectin - J920-0	B207
Fecal Lactoferrin - 3374-6	B208
Fentanyl, Urine, Screen, Qual. (2) - 6194-5	B208
Ferritin - 0088-5	B208
Fetal Fibronectin - 3576-6	B209
FGFR1 by FISH (Tech Only) - TG46-8	B209
FGFR2/IGH t(4;14) by FISH - A291	B209
Fibrinogen, Activity - 0089-3	B210
Fibrinolytic Evaluation - Standard - 5860-2	B210
Chronic Eosinophilic Leukemia: FIP1L1-PDGFR4 by FISH - 5182-1	B210
Lymphoma - Non-Hodgkin, BCL1, BCL2, BCL6 (FISH Panel) - 5273-8	B211
Multiple Myeloma FISH Panel - 5282-9	B211

FISH, FGFR1, 8p11-12 - B557-0	B211
OnkoSight FLT3 - A044-0	B212
Onkosight NGS FLT3 IF NEG. RFLX TO CEBPA+NPM1 NGS - J154-6	B212
Onkosight NGS FLT3 AND NPM1 Sequencing - J155-3.....	B212
Influenza A/B, RSV BY PCR - 1410-0	B213
Cytology, Alcohol-Based FNA - 5243-1	B213
Fine Needle Aspirate - 5242-3	B213
Folate, RBC - 0873-0.....	B214
Food Allerpanel (4 mix) - 2189-9	B214
Food Panel II - 2797-9.....	B214
FOXO1 (FKHR) (13q14) GENE BY FISH - B665-1	B215
Fractional Excretion Of Sodium - P791-7	B215
Fragile X - 2383-8.....	B215
FSH (Follicle Stimulating Hormone) - 0092-7.....	B216
Factor XI carrier screening (Hemophilia C) - J310-4.....	B216
Allergen False Ragweed (w4), IgE - 2407-5.....	B216
Familial Hypercholesterolemia LDLR - B151-2.....	B217
Familial Mediterranean Fever - J307-0	B217
Feather Mix (ex71), IgE - 1300-3.....	B217
Occult Blood, Stool, Immunochem (FIT) - 6182-0	B218
Female STI Urine Profile - P863-4	B218
Feulgen Stain - Tech Only - T983-6	B218
Fibrotest/Actitest - H510-1	B219
Allergen Fig (f328), IgE - 1937-2	B219
Allergen Fire Ant (i70), IgE - 0904-3.....	B219
Combined First (NT, Total BHCG, PAPP-A) - 6246-3.....	B220
Combined First Plus (NT, Total BHCG, Inhibin-A, PAPP-A) - 8622-3	B220
Folate (Folic Acid) - 0090-1	B220
Food Allergy Profile w/ rfx to Milk and Peanut - K361-6.....	B221
Food Allergy Profile - K426-7	B221
Foods I Screen (fish/6) - 0743-5.....	B221
Foods II Screen (10) - 0745-0	B222
Allergen Formaldehyde (k80), IgE - 2106-3.....	B222
Catecholamines, Fractionation, Plasma - 0214-7	B223
Thyroxine, Free (FT4) - 0091-9	B223
GGTP (Gamma-GT) - 0093-5	B224
G. vaginalis by RT-PCR - J211-4	B224
G.T.T. 2 hr (includ. 30 min) - Q605-8.....	B224

Galactosemia Carrier Screen - B155-3	B225
Gastric Biopsy (1 Specimen) - 5640-8.....	B225
Gastric Washings/Brushings - 5643-2	B225
GC MultiPlex PCR - swab - J834-3	B226
GC, Anal Aptima - A112-5.....	B226
GC, Oral Aptima - A111-7	B226
GC,rRNA,TMA Aptima Swab - 1003-3	B227
GC/CT SWAB BY PCR - L583-4.....	B227
General Inhibitor Screen - 5800-8.....	B227
General Thrombophilia Assessment - 6800-7.....	B228
Culture, Genital - 0299-8	B228
Gentamicin, Peak - 0730-2	B228
Gentamicin, Peak + Trough - 9081-1	B229
Gentamicin (Random) - A616-5	B229
GENTAMICIN, TROUGH - 1930-7	B229
Germ Cell Tumor IHC Panel - 5194-6	B230
Giardia Antigen EIA (2ND SPECIMEN) - 3898-4	B230
Giardia Antigen EIA (3RD SPECIMEN) - 3899-2.....	B230
GIST Profile (IHC) - 5175-5.....	B231
GIST Profile by IHC TC - J687	B231
Gleason Score 3+4, Reflex To PTEN FISH - J143-9.....	B231
Gliadin IGA, ABS. - 5868-5	B231
Gliadin IgG, ABS - 3822-4.....	B232
GLU.TOL.3HR.(Pregnancy)100GM - 1015-7	B232
Glucose, Post Dialysis - 3110-4.....	B232
Glucose Tolerance 2 Hour - 8979-7	B233
Glucose Tolerance 3 HRS. - 0097-6.....	B233
Glucose Tolerance 4 HRS. - 0098-4.....	B233
Glucose Tolerance 5 HRS. - 0099-2.....	B233
Glucose, Urine, Random - 5042-7	B234
Glucose, 2HR. PP W/Glucola - 3157-5.....	B234
Glucose, Urine 24HR Quantitative - 0096-8.....	B234
Glycogen Storage, Type IA (GSDIA) - A651-2	B235
Gonadotropin I - 2318-4.....	B235
GP/Bile Duct Brushing - 5119-3.....	B235
Bone Marrow Smear Interpretation - 5285-2	B236
ISH Single Antibody, W/Interp. - 5286-0.....	B236
GP/Level III Histo (2 Slides) - 5159-9	B236

GP/Level III Histo. TC (1 Slide) - 5158-1	B237
Gram Stain - 0103-2	B237
Group A Strep (Non NY) - J965-5	B237
Group B Strep (Non NY) - J966-3	B238
Growth Hormone - 0510-8.....	B238
GTT,GEST.,2HR.(ADA) PANEL - P782-6.....	B238
GTT,GEST,3HR.(ACOG,CARP/COUS) - 6794-2.....	B239
Allergen Garlic (f47), IgE - 0965-4	B239
Gastric Distress - F101-1	B239
Gastrin, Serum - 0508-2	B240
Gaucher Disease Carrier Screening - A655-3	B240
Solid Tumor Gene Fusion Assay - J355-9	B240
General Health Panel - 9149-6	B241
Chlamydia Swab by PCR - T832-5	B241
Genital Cult. R/O Yeast/Gard. - 1295-5.....	B242
Gonorrhea Swab by PCR - T833-3.....	B242
Culture, Genital (Beta strep. Group B) - 3400-9.....	B243
Allergen Giant Ragweed (w3), IgE - 1197-3	B243
Giardia Antigen (EIA) - 1396-1.....	B243
Allergen Ginger (rf270), IgE - 2346-5.....	B244
Glioma: IDH1/IDH2 by NGS - J635-4	B244
Glucose, 1hr, PP (Preg) w/Glucola - 0368-1	B244
Glucose - 1976-0.....	B245
Glucose 2hr PP - 0369-9	B245
Glucose, 1hr PP - 0171-9	B245
Allergen Gluten (f79), IgE - 0935-7	B246
Fructosamine, Serum - 0485-3.....	B246
Glycemic Control Profile - H959-0.....	B246
Hemoglobin A1c (GlycoHgb) - 0102-4.....	B247
Allergen Glycyphagus Domesticus (d73), IgE - 3813-3.....	B247
Allergen Goat Milk (rf300), IgE - 2263-2.....	B247
Allergen Golden Rod (w12), IgE - 1364-9.....	B248
Culture, Neisseria gonorrhoeae (GC), Anal - 1353-2	B248
Culture, Neisseria gonorrhoeae (GC), Throat - 1352-4	B248
Culture, Neisseria gonorrhoeae (GC) - 0075-2.....	B249
Gonorrhea NAAT - 3853-9.....	B249
GC Smear - 0094-3	B249
Allergen Goose Feathers (e70), IgE - 0847-4.....	B250

Allergen Grape (f259), IgE - 1225-2	B250
Allergen Grapefruit (f209), IgE - 1349-0.....	B250
Grass Mix (gx1), IgE - 1841-6.....	B251
Grasses Screen (6) - 0741-9	B251
Allergen Green (Bell) Pepper (rf263), IgE - 1243-5.....	B251
Allergen String Bean (rf315), IgE - 1230-2.....	B252
Allergen Green Pea (f12), IgE - 1228-6.....	B252
Allergen Grey Alder (t2), IgE - 0561-1	B252
Culture, Genital - Beta strep. Group B (Penicillin Allergic) - A265-1	B253
Allergen Guinea Pig (e6), IgE - 3566-7	B253
H. Pylori BY IHC Tech Only - J159-5.....	B253
H. Pylori Antibody Panel (G,A,M) - 0953-0	B254
H. Pylori Antigen, Stool - 2441-4	B254
H. PYLORI IgM - 7736-2.....	B254
H. Pylori,IGA (For Panel #0953) - 1766-5.....	B255
H. Pylori,IGG (For Panel #0953) - 1765-7	B255
H.Pylori (Stain Ident.) - 5645-7	B255
H. Pyloril BY IHC W/ Interp - J158-7.....	B256
H. Pylori by IHC. If+ rflx to H. pylori PCR with rflx to H. pylori Clarithromycin resistance detecti - L282-3	B256
Haptoglobin, Serum - 0514-0	B256
Hepatitis B Core Antibody, Total - 0108-1	B257
Hepatitis B Core Antibody, IgM - 0206-3.....	B257
Hepatitis B Surface Antibody (Qual/Immunity) - 0107-3	B257
Hepatitis B Surface Antigen - 0106-5.....	B258
HCG Quant.,Progesterone - 2379-6.....	B258
Hematocrit - 0019-0.....	B258
HCV Viral Load reflex to NS3/NS5 (Non-NY) - K943-1	B259
Hepatitis C RNA, Quant., PCR - 3376-1	B259
HDL Cholesterol - 0059-6	B259
Human Epididymal Protein 4 - J032-4	B260
Heart Health Monitoring - Q472-3.....	B260
Heart Health Baseline - Q471-5	B260
Hemoglobin + Hematocrit - 0104-0.....	B260
Hepatitis A Antibody, IgM - 0538-9	B261
Hepatitis A Ab.(non-reflex) - 2126-1	B261
Hepatitis A Ab.,Total W/Reflex - 0105-7	B261
Hepatitis B Core W/Reflex To IgM - 3138-5	B262

Hepatitis B Surface Antigen, Quant. - 1608-9	B262
Hepatitis B Surface Antibody, Quant. - 1606-3	B262
Hepatitis Be Antigen - 0539-7	B263
Hepatitis C Genotype - 2161-8	B263
Heparin Level (ANTI-XA) - LMWH - 5761-2	B263
Hepatic Function Panel - 3422-3	B264
Hepatitis B Recovery - 2075-0	B264
Hepatitis/ Comprehensive - 0992-8	B264
Hepatocellular vs. Cholangiocarcinoma vs. Metastasis IHC Pnl - 5193-8	B265
HER2 by IHC (Tech Only)w/digital imaging - 5405-6	B265
HER2 by IHC - reflex +2 to HER2 by FISH - 5428-8	B265
HER-2, IHC - A301-4	B266
HER2 by IHC w/digital imaging - 5171-4	B266
HER-2/neu,FISH (Manual) - 5247-2	B266
HER2 by FISH with Interpretation - 5262-1	B267
HER2 FISH Global Gastric+GEJ - A427-7	B267
HER2 FISH-Reflex Equivocal to HER2 IHC - B812-9	B267
Automated HER2 FISH TC Only - 5259-7	B268
HER2 BY IHC GI/ESO Carcinoma - A972-2	B268
HER2 IHC Gastric/GEJ, TC Only w/digital imaging - B269-2	B268
Hereditary Hemochromatosis DNA Mutation Analysis (NY only) - 2696-3	B269
Hereditary Hemochromatosis DNA Mutation Analysis (Non-NY) - 3420-7	B269
Histo Chemical Micro Organism TCx2 - B671-9	B269
Histo Chemical Micro Organism x3 W/IN - B673-5	B270
Histo Chemical Micro Organism TC - B672-7	B270
Histo Chemical Micro Organism x2 W/IN - B674-3	B270
Histo Chemical Micro Organism W/INTERP - B675-0	B271
Histo Chemical Micro Organism TCx3 - B670-1	B271
Histology Processing - 5416-3	B271
HIT Profile-Level 1 - 5871-9	B271
HIV AG/AB 4th Generation - B688-3	B272
HIV-1,RNA Viral Reflex Genotype - 6364-4	B272
HIV-1, RNA, Ultra/PCR, Viral Load - 1010-8	B272
HLA Typing for Celiac Disease - J097-7	B273
HOLD BLOCK FOR ADDT'L TESTING - B700-6	B273
Hold Cells For Further Testing - B232-0	B273
HPV DNA 16, 18, NON 16/18 - A531-6	B273
HPV High Risk by ISH - 8714-8	B274

HPV High by ISH Stain ONLY - 8689-2	B274
HPV High/Low Risk by ISH - 3412-4	B274
HPV High/Low Risk by ISH (Stain Only) - 8696-7	B275
HPV PRIMARY SCREENING >=25 YRS, WHEN + RFX TO PAP - H852-7	B275
HPV mRNA screen, if pos, Reflex to 16, 18/45 - B342-7	B275
HPV DNA, High Risk Anal - TP - J025-8	B276
HPV, High Risk - 3910	B276
HSV 1+2 Reflex to VZV - L009-0	B276
Herpes Virus I/II DNA - H760-2	B277
HSV 1+ 2/ VZV Assay - L008-2	B277
HERPES SIMPLEX VIRUS I DNA - J164-5	B277
HERPES SIMPLEX VIRUS II DNA - J165-2	B278
Herpes I/II Ab. Screen, IgG - 0109-9	B278
HSV/LGV Lesion Profile (NON-NY) - P764-4	B278
HSV 1+2 (NON NY) - L010-8	B279
Hyperactive Coag. Panel - 5993-1	B279
Hypersensitivity Pneumonitis - 1366-4	B279
Hyperthyroid Panel - 0403-6	B280
Allergen Halibut (rf303), IgE - 1485-2	B280
Allergen Hazelnut (f17), IgE - 1404-3	B280
Hemoglobin - 1499-3	B281
Hemoglobin A1c with reflex to Glycomark - J218-9	B281
Allergen Helminthosporium (m8), IgE - 1894-5	B281
Complement, Total (CH50), Serum - 0534-8	B282
Hepatitis B Panel - 0197-4	B282
Hepatitis B Virus, DNA, Quantitative, RT-PCR - 3389-4	B282
Hepatitis C Antibody W/ Reflex RT PCR - B125-6	B283
Hepatitis Screen - 0330-1	B283
Hepatitis Management Profile - 0984-5	B283
HER2 FISH TC-Gastric+GEJ - A428-5	B284
Hereditary MDS/Leukemia Panel - T830-9	B284
Pancreatic Cancer Panel - B343-5	B284
Hermansky-Pudlak Syndrome 3 - B150-4	B285
Opiates, Urine, Screen (300) - 0433-3	B285
Herpes I AB., IgG - 0110-7	B285
Herpes II Ab., IgG - 0570-2	B286
Herpes Simplex II by Multiplex (Non NY) - 6352-9	B286
Herpes Simplex I by Multiplex (Non NY) - 6351-1	B286

Varicella Zoster Ab., IgG - 0597-5	B287
Histology Processing 8 H+E - J858-2	B287
Histology Processing (GI) - J862-4	B287
Histology Processing 10 H+E - J861-6	B288
Histology Processing 6 H+E - J857-4	B288
Histology Processing 7 H+E - J859-0	B288
Histology Processing 9 H+E - J860-8	B289
Histone Ab, IgG - T841-6.....	B289
Allergen Honey Bee (i1), IgE - 0883-9	B289
Allergen Hops (RF324), IgE - A333-7	B290
Allergen Horse Dander (e3), IgE - 0691-6.....	B290
Allergen Horsefly (i204), IgE - 3556-8.....	B290
Allergen House Dust (Greer) (h1), IgE - 1631-1	B291
Allergen House Dust H/S (h2), IgE - 0662-7	B291
House Dust Screen (5) - 0910-0.....	B291
HLA-B27 - 0375-6.....	B292
Prothrombin Fragment 1+2 - 5797-6	B292
IBD Differentiation Panel - H583-8	B292
IGF Binding Protein-3 - 1316-9	B293
Insulin-Like Growth Factor I - 1072-8	B293
Lymphoma - B-Cell, IGH by FISH - 5032-8.....	B293
IGH/BCL2 t(14;18) by FISH TC - b350	B294
IGH/MAF t(14;16) by FISH - B564.....	B294
FISH, IGH/MAF, t(14;16)- TC - B760-0.....	B294
FISH, IGH/MAFB, t(14;20) - B545-5	B294
IGH/MAFB T(14;20) - TC ONLY - TB61-7.....	B295
QFT - TB GOLD PLUS (QUANTIFERON) - T814-3.....	B295
IHC 1 Antibody Interp. - B024-1	B295
IHC 1 Antibody/Stain Only - 5144-1	B296
IHC 2 Antibody/Stain Only - 5145-8	B296
IHC 3 Antibody/Stain Only - 5146-6	B296
IHC 4 Antibody/Stain Only - 5147-4	B297
IHC 5 Antibody/Stain Only - 5148-2	B297
IHC Antibody (2) With Interpretation - 5137-5	B297
IHC Antibody (3) With Interpretation - 5138-3	B298
IHC Antibody (4) With Interpretation - 5139-1	B298
IHC Antibody (5) With Interpretation - 5141-7	B298
IHC Single Antibody With Interpretation - 5136-7	B299

Interleukin 6 - J785-7	B299
Imaged IHC (1) W/Interp. - 5308-2	B299
Imaged IHC, 1 Slide TC Only - 5842-0	B300
Imaged IHC, 2 Slides TC Only - 5838-8	B300
Imaged IHC, 2 Slides W/Interp. - 5845-3	B300
Imaged IHC, 3 Slides TC Only - 5837-0	B301
Imaged IHC, 3 Slides W/Interp. - 5844-6	B301
Imaged IHC, 4 Slides TC Only - 5834-7	B301
Imaged IHC, 4 Slides W/Interp. - 5843-8	B302
Immunoglobulins, Serum - 0520-7	B302
Indoor Allergy Panel - P268-6	B302
Inherigen - A662-9.....	B303
Inherigen Plus - A659-5	B303
Inhibitor Titer - 5757-0	B303
Prenatal Rapid Aneuploidy FISH - 3582-4.....	B304
Insulin Tolerance (3 Hour) - 0024-0	B304
Insulin, 1st HR. Specimen - 1649-3	B304
Insulin, Fasting - 1648-5	B305
Interferon-ALPHA, EIA - 3997-4.....	B305
Interleukin 28B POLYMORPHISM (NON-NY) - A304-8.....	B305
Interleukin 5 - J223-9.....	B306
Iodine, Serum/Plasma - J076-1	B306
Iodine, Urine Panel - H474-0	B306
Iron + TIBC - 0250-1	B307
Iron, Liver Studies, Quant. - 5099-7	B307
ISH Single Antibody -PC Only - B960-6.....	B307
Isospora Belli ID, Modified Acid Fast - 2670-8.....	B308
Accutype Ribavirin - B832-7	B308
Immunoglobulin A, (IgA) Serum - 0211-3.....	B308
Immunoglobulin G (IgG), Serum - 0541-3	B309
Lymphoma - B-Cell, IGH by PCR - 5278-7.....	B309
Immunoglobulin M, (IgM) Serum - 0549-6.....	B309
CLL Prognosis: IgVH Mutation Analysis by PCR - 5223-3	B310
Immunotyping, Serum - 0413-5	B310
Insulin, Total/Random - 0113-1	B310
InheriGen Plus w/o Fragile X - F021-1	B311
InheriGenTx - B735-2.....	B311
Inhibin-A - TF79-1	B311

Phosphorus, Inorganic - 0127-1	B312
Integrated Part 1 (NT, PAPP-A) - 6375-0.....	B312
Integrated Part 2 (AFP, UE3, HCG, INHIBIN-A) - 6376-8	B312
Intrinsic Factor AB - T835-8.....	B313
Iron, Serum - 0114-9	B313
Iron Reflex to TIBC - M064-4.....	B313
Transferrin - 0192-5	B314
Iron/Ferritin Panel - 1021-5	B314
Allergen Isocyanate MD1 (k76), IgE - 2104-8	B314
Allergen Isocyanate TD1 (k75), IgE - 2103-0	B315
Allergen Ispaghula (Psyllium) (k72), IgE - 3968-5	B315
JAK2 (including V617F and exon 12), MPL and CALR if negative reflex to MPN Panel - J632.....	B315
JAK2 (including V617F and exon 12), MPL, and CALR if negative, reflex to MPN Panel - J632-1	B316
JCV Ab w/ Index w/ RFX - B542-2	B316
Joubert Syndrome - A645-4	B316
Allergen Johnson Grass (g10), IgE - 1824-2.....	B317
Allergen June/Ken.Blue Grass (g8), IgE - 0648-6.....	B317
Potassium, Serum - 0129-7.....	B317
Potassium, Urine 24 Hr. - 0301-2.....	B318
Cannabinoid, synthetic (K2,Spice), Urine, Screen, group 1, Qual.(20) - A292-5.....	B318
Cannabinoid, synthetic (K2,Spice), Urine, Screen, group 2, Qual.(10) - J090-2.....	B318
Kappa Free Light Chain Serum - 6241-4	B319
Free Kappa + Lambda Light Chain (serum) - 3893-5.....	B319
Keppra Level (Levetiracetam) - 3465-2.....	B319
KI-67 Imaged By IHC-Tech Only - A346-9	B320
KIT Sequencing (Solid Tumors) Non NY - J020-9	B320
Onkosight NGS KIT Sequencing (Non-NY) - J138-9	B320
KIT and PDGFRA Sequencing NY Only - J464-9	B321
KIT and PDGFRA Sequencing- Non Ny - H578-8	B321
KIT Sequencing (Solid Tumors) NY Only - J453-2	B321
KRAS Sequencing- Non NY - J024-1.....	B322
KRAS, if neg reflex to BRAF NY Only - J462-3	B322
KRAS, if neg reflex to BRAF- Non NY - J104-1	B322
KRAS Sequencing NY Only - J458-1	B323
KRAS by NGS, if negative reflex to BRAF; + MSI by PCR, MSI-High, reflex to MLH1 - H580-4.....	B323
KRAS/NGS if neg rflx to BRAF/ MSI-PCR, if MSI-High, rflx MLH(NY) - J452-4.....	B323
Peritoneal KT/V, Measured - 6238-0.....	B324

Cabbage (f216), IgE - 1244-3	B324
Ketones (Acetoacetic), Serum - 0028-1	B324
Ki-67 by IHC - 5152-4	B325
Allergen Kidney Bean (rf287), IgE - 3122-9	B325
Kidney Profile - L325-0	B325
Allergen Kiwi (f84), IgE - 0842-5	B326
Lupus-Hexagonal Phase - 5749-7	B326
Lactobacillus (Species) by PCR - J213-0	B326
Lactose Tolerance Test - 1178-3	B327
LAL (Limulus Amebocyte Lysate) - 2475-2	B327
Lambda Free Light Chain Serum - 6240-6	B327
Liver Cytosol (LC-1) Autoantibodies - B608-1	B328
LDL Direct - 2194-9	B328
Lead, Blood + E.P. - 0588-4	B328
Legionella Culture - 3874-5	B329
Leptin (ELISA) - J220-5	B329
Level 3 Surgical Pathology - 5695-2	B329
Level 4 Surgical Pathology - 5696-0	B330
Level 5 Surgical Pathology - 5697-8	B330
Level II Histo,TC - A229-7	B330
LEVEL II-PC SURGICAL PATHOLOGY - B012-6	B331
Level III Histo. (3 Slides) - 5164-9	B331
Level III Histo. (4 Slides) - 5172-2	B331
Level III Histo.(5 Slides) - 5173-0	B332
LEVEL III-PC SURGICAL PATH. - B013-4	B332
Level IV Histo, TC - A231-3	B332
LEVEL IV-PC SURGICAL PATHOLOGY - B014-2	B333
Level V Histo,TC - A228-9	B333
LEVEL V-PC SURGICAL PATHOLOGY - B015-9	B333
Level VI Histo,TC - A230-5	B334
LEVEL VI-PC SURGICAL PATHOLOGY - B016-7	B334
Lymphogranuloma Venereum (LGV) (Non-NY) - 6356-0	B334
LH, FSH, Estradiol - 0982-9	B335
Luteinizing Hormone (LH) - 0342-6	B335
Lipase, Peritoneal Fluid - 6167-1	B335
Lipoprotein Electrophoresis w/Lipids - 0551-2	B336
Liver Panel - 1074-4	B336
Lupus Antibody Panel (SLE) - 4291-1	B336

Lyme AB West-Blot IgG/IgM - 4033-7	B337
Lyme Disease AB RFLX IGM TO WB G/M - 2297-0	B337
Lyme Ab.,Total (No Reflex) - 2382-0	B337
Lyme Disease Antibody, Serum - 0568-6	B338
Lyme Blot, IgG W/Bands - 6272-9	B338
Lyme Blot, IgM W/Bands - 6273-7	B338
Sjogren's Abs (SS-A, SS-B) - 0668-4	B339
LabFinder Report - T420-9	B339
Lactate Dehydrogenase (LDH) - 0117-2	B339
Lactobacillus reuteri/rhamnosus (Non NY) - J964-8.....	B340
Lamb/Mutton (f88) , IgE - 0862-3	B340
Allergen Latex (k82), IgE - 1371-4.....	B340
Allergen Lemon (f208), IgE - 1224-5.....	B341
Allergen Lentil (f235), IgE - 3521-2	B341
Allergen Lettuce (f215), IgE - 1409-2.....	B341
Syphilis Reverse Algorithm - J275-9	B342
Lipid Screen W/Reflex to Direct LDL - M063-6	B342
Liver Profile - 0012-5.....	B342
Allergen Lobster (f80), IgE - 0841-7.....	B343
Lipoprotein-a (Lp-a) - 1639-4	B343
Lp-Pla2 (PLAC) - 6155-6	B343
Lupus Anticoagulant Screen - 5817-2	B344
Lyme Ab., IgM - 1615-4.....	B344
Lynch/Colorectal High Risk Panel - B522	B344
Mycoplasma Genitalium by MultiPlex PCR- swab - J837-6	B345
Multiple Myeloma FISH Panel - TC Only - P961-6	B345
Magnesium, Urine, Random - 3135-1	B345
Male STI Urine Profile (GP) - P310-6	B346
Maple Syrup Urine Disease(Type 1A & 1B) - A650-4	B346
MDM2 Gene Amplification by FISH - B544-8.....	B346
Ecstasy (MDMA), Urine, Screen, Qual. (300) - 3174-0.....	B347
Megasphaera Type 1 by RT PCR - J214-8	B347
Metastatic Melanoma - Lymph Node - 5126-8.....	B347
Menopause/Perimenopause Risk Management Profile - P244-7	B348
Male High Risk Sexual Health Panel - F239-9.....	B348
MET BY FISH-Global - B335-1	B348
MET BY FISH-Tech Only - B336-9.....	B349
MET by IHC - A612-4	B349

Methadone, Urine, Screen (300) - 0925-8.....	B349
MGMT Promoter Region Methylation - J876-4.....	B350
Microalbumin, Urine, 24Hr. - 3171-6.....	B350
Microsporidium ID, Stool - 2673-2.....	B350
Milk Component Panel IgE - H191-0.....	B351
Mixing Studies-APTT - 5830-5.....	B351
Mixing Studies-PT - 5840-4.....	B351
MLH-1 Promo. Methylation - A313-9.....	B352
AML-M4 and M5 by Fish- TC ONLY - B889-7.....	B352
AML - M4 + M5 by FISH - 5427-0.....	B352
Measles,Mumps,Rubella - 2565-0.....	B353
Mobiluncus (Species) by RT PCR - J217-1.....	B353
6-Monoacetyl Morphine (6-MAM), Urine, Screen (10) - 3911-5.....	B353
Rule Out Mosaicism Chromosome Analysis, Blood - T982-8.....	B354
MPN: JAK2/MPL/CALR by NGS - J633-9.....	B354
Myeloperoxidase Ab. - 2192-3.....	B354
Minimal Residual Disease for CLL/SLL- Tech Only - J496-1.....	B355
MRSA Screen - 6119-2.....	B355
Mismatch Repair(MMR) by IHC w/Interp - A943-3.....	B355
Mismatch Repair(MMR) by IHC, Tech Only - A944-1.....	B356
MSI: if MSI-H reflex to BRAF, if BRAF neg. reflex to MLH1 - J761-8.....	B356
Microsatellite Instability (MSI) by PCR - 3371-2.....	B356
MUC4 by IHC - TF93-2.....	B357
MUC4 by IHC Tech Only - TF92-4.....	B357
Mucopolidosis Type IV - A654-6.....	B357
Allergen Mucor Racemosus IgG - B770-9.....	B358
Mumps Virus Antibody, IgG - 0316-0.....	B358
MYC-BA, 8q24 Rearrangement FISH- TECH ONLY - J144-7.....	B358
MYC-BA, 8q24 Rearrangement FISH - J141-3.....	B359
MYC/IGH t(8;14) by FISH TC - B439.....	B359
Lymphoma - Burkitt, MYC/IGH by FISH - 5027-8.....	B359
Mycoplasma Culture, Blood - 3247-4.....	B360
Mycoplasma Genitalium By Multiplex PCR - 6340-4.....	B360
Mycoplasma Hominis By Multiplex PCR (NON-NY) - 6341-2.....	B360
Mycoplasma pneumoniae IgG/IgM Antibody, IFA (CSF) (NON-NY) - 3442-1.....	B361
Mycoplasma/Ureaplasma Culture - 2523-9.....	B361
Myoglobin, Serum - 0789-8.....	B361
Allergen Macadamia Nuts (RF345), IgE - 3644-2.....	B362

Allergen Malt (f90), IgE - 1630-3	B362
Allergen Mandarin (rf302), IgE - 3520-4	B362
Allergen Mango (f91), IgE - 3456-1	B363
Manual Diff/Platelet only - 1340-9	B363
Citrated Platelet Count - J661-0.....	B363
Allergen Maple (Box Elder) (t1), IgE - 0684-1	B364
Allergen Marguerite/Ox-eye Daisy (w7), IgE - 3879-4.....	B364
Cannabinoids (Marijuana), Urine, Screen (100) - 0922-5.....	B364
Allergen Marjoram (f274), IgE - 6206-7	B365
Allergen Marsh Elder (w16), IgE - 1481-1	B365
Allergen Meadow Fescue (g4), IgE - 0681-7	B365
Rubeola/Measles IgG W/Reflex - T007-4.....	B366
Rubeola/Measles IgG Antibody - 0567-8.....	B366
Allergen Melaleuca Tree (t212), IgE - 2642-7	B366
Allergen Melon (Cantaloupe) (f87), IgE - 1227-8	B367
Meprobamate - 1936-4	B367
Mesoridazine - 3709-3.....	B367
Vitamin B12 - 0160-2	B368
Methylmalonic Acid, Serum - 0349-1.....	B368
Magnesium, Serum - 0120-6	B368
Magnesium, Urine, 24Hr. - 1171-8.....	B369
Magnesium, RBC - 1383-9.....	B369
Microalbumin, Urine, Random - 0228-7.....	B369
MSI by PCR if MSI-H Reflex to MLH1 - P264-5	B370
Mitochondrial Complex I Deficiency - J305-4	B370
Mold Mix (mx1), IgE - 2181-6	B370
Molds Screen (6) - 0901-9	B371
Skin Biopsy, (1 Jar) - 5602-8.....	B371
Mono Screen - 0121-4	B371
Allergen Mosquito (i71), IgE - 2179-0.....	B372
Allergen Mountain Juniper (T6), IgE - 6334-7.....	B372
Allergen Mouse Epithelium (e71), IgE - 1359-9	B372
Allergen Mouse U. Protein (e72), IgE - 1360-7.....	B373
Allergen Mucor Racemosus (m4), IgE - 1315-1	B373
Allergen Mugwort (w6), IgE - 0718-7.....	B373
Multiple Double Stain by IHC Global - B962.....	B374
Multiple Double Stain by IHC TC Only - B963.....	B374
Multiple Sulfatase Deficiency - J296-5	B374

Mumps Virus Ab IgG W/Reflex - T009-0	B374
Allergen Mushroom (rf212), IgE - 1462-1	B375
Allergen Mustard (f89), IgE - 1823-4	B375
MDS FISH Panel - 5281-1	B375
MDS FISH Panel- TC Only - P960-8.....	B376
Primidone (incl. Phenobarb) - 0464-8	B376
N. Gonorrhea, Urine, rRNA - 6369-3.....	B376
Infant Tot.Bili (tot/dir=#0046) - 1900-0	B377
Chlamydia Direct Antigen - 0782-3.....	B377
Infant Dir.Bili (tot/dir=#0046) - 1901-8.....	B377
Nefazodone Level - 3147-6	B378
Allergen Nettle (w20), IgE - 6308-1.....	B378
Neuroendocrine Neoplasm IHC Panel - 5195-3.....	B378
Neutrophil Function Panel - B124-9	B379
N. Gonorrhea/NG (GC) PCR Urine - T408-4	B379
Ammonia, Plasma (As Nitrogen) - 0035-6.....	B380
Niemann-Pick Type A+B Carrier Test - A653-8.....	B380
Nipple Discharge Smear - 1303-7	B380
ClariTest Core Non-Invasive Prenatal Screen (NON-NY) - TH18-5	B381
ClariTest Core Non-Invasive Prenatal Screen with 22q11.2 Microdeletion (DiGeorge Syndrome)(NON-NY) - TH19-3.....	B381
Hodgkin Lymphoma IHC Panel - 5101-1	B381
Northern Allerpanel (3 mix) - 2187-3.....	B382
Onkosight NGS NPM1 Sequencing - J137-1	B382
NRAS Sequencing- Non NY - J023-3	B382
NRAS Sequencing NY Only - J455-7.....	B383
Sodium, Serum - 0148-7.....	B383
Sodium, Urine, 24 Hr. - 0308-7	B383
Sodium, Urine, Random - 1793-9.....	B384
Naproxen (Aleve,Anaprox,Naprosyn) - 6223-2	B384
Culture, Nasopharyngeal/Nasal - 0074-5.....	B384
Nasal Smear: Complete WBC w/ Diff - 4009-7	B385
GC By Multiplex PCR - 6338-8.....	B385
Nemaline Myopathy 2 - A646-2.....	B385
New PD-L1 -	B1
Cytology, Thinprep Fluids - 5244-9	B386
Cytology, Alcohol-Based Fluids - 5245-6	B386
Non-Gyn Other (Slides) - Interp Only - J949-9	B386

Allergen Profile, Northeast - 0740-1	B387
Allergen Nut Mix (fx1), IgE - 2183-2	B387
Allergen Nutmeg (rf282), IgE - 2352-3	B387
Ova and Parasite (O&P) Panel, Stool W/ Trichrome(3 specimens) - 2326-7	B388
Ova and Parasite (O&P), Sputum - J002-7.....	B388
Ova and Parasite (O&P), Urine - 2534-6.....	B388
Obstetric/Prenatal I - 0008-3	B389
Obstetric/Prenatal II - 0981-1	B389
OnkoSight Lung Cancer Panel, ALK and ROS1 by FISH - H125-8.....	B389
Drug Confirmation, Opiates, Urine (GC/MS) - 2596-5	B390
Opiates, Urine, Screen (200) - 6258-8.....	B390
Allergen Oak Tree (t7), IgE - 0652-8	B390
Oat (f7), IgE - 0814-4	B391
Culture, Eye - 0076-0.....	B391
Culture, Eye (Left) - 0565-2.....	B391
Culture, Eye (Right) - 0633-8	B392
Olive Tree (t9), IgE - 1251-8.....	B392
OneCheck Hematopathology - 5500-4.....	B392
Allergen Onion (f48), IgE - 0864-9	B393
OnkoSight AML Panel Sequencing - B894-7	B393
OnkoSight Colorectal Cancer Panel - B822-8	B393
OnkoSight Glioma Panel - J634-7	B394
OnkoSight Lung Cancer Panel - B819-4	B394
OnkoSight Lung Comprehensive - K924-1	B394
OnkoSight Lung Panel, if EGFR- neg. rflx to (J015) ALK by FISH, if neg. rflx to ROS1 by FISH - L113-0	B395
OnkoSight MDS Panel Sequencing - B893-9	B395
OnkoSight MPN Panel Sequencing - B818-6	B395
OnkoSight MYD88 Sequencing (Non-NY) - J130-6	B396
OnkoSight Melanoma Panel - B820-2	B396
OnkoSight Myeloid Disorder Panel - B823-6.....	B396
OnkoSight Solid Tumor Panel - B821-0.....	B397
Opiates, Serum Screen - J865-7.....	B397
Allergen Orange (f33), IgE - 0835-9	B397
Allergen Orchard Grass (g3), IgE - 0649-4.....	B398
Allergen Oregano (rf283), IgE - 2351-5.....	B398
Ova and Parasite (O&P), Stool W/ Trichrome (2nd specimen) - 0355-8	B398
Ova and Parasite (O&P), Stool W/ Trichrome (3rd specimen) - 0227-9	B399

Allergen, Ovalbumin (F232), IgE - J281-7	B399
Allergen Ovomuroid (f233), IgE - J280-9.....	B399
Oxycodone, Urine, Screen (300) - A595-1	B400
Allergen Oyster (f290), IgE - 3121-1	B400
Vitamin B6 - 0296-4.....	B400
Neutrophil Cytoplasmic AB/ANCA - 0636-1	B401
Prostate Specific Acid Phos (P.A.P.) - 0030-7	B401
PTT (Partial Thromboplastin Time) - 0139-6.....	B401
PAI 4G/5G Polymorphism - A269-3.....	B402
Chronic Lymphocytic Leuk.(CLL): FISH Pnl(+12, 13q, 17p, 11q) - 5280-3.....	B402
PANEL L202 - L202-1	B402
PAP (RFLX HPV HIGH ASC-US or >) - 6630-8	B403
PAP (RFLX HPV HIGH ASC-US) - 6629-0	B403
PAP + CT/GC - 9241-1.....	B403
PAP + HPV HI - 9102-5	B404
PAP + HPV HI + CT/GC - 9235-3	B404
PAP RFLX HPV mRNA ASCUS,RFX 16/18/45 w/CT/GC - A253-7	B404
PAP Reflex to HPV mRNA when ASCUS or greater,RFX 16,18/45 - A254-5	B405
PAP W/O Image Reflex HR HPV - 5240-7.....	B405
PAP(RFX HPV Plus) - A489-7	B405
PAP+CT+GC(RF HPV+ LO ASC OR >) - P383-3	B405
PAP+CT+GC(RFX HPV PLUS) - P374-2	B406
PAP + HPV Plus + CT + GC - P372-6.....	B406
PAP, Liquid-Based - 1962-0.....	B406
StormPath PAP, Liquid based - B352-6	B406
PAP/CT/GC (Reflex HPV PLUS when ASC-US) - 6431-1	B407
PAP/CT/GC (Reflex HPV PLUS when ASC-US or UNSAT) - L327-6	B407
PAP RFX HPV Plus mRNA ASCUS>,RFX 16/18/45 w/CT/GC - A255-2.....	B407
PAP/CT/GC/ HPV mRNA,RFX 16/18/45 GENOTYPE - A252-9.....	B407
PAP/HPVHI/GC RFX>GENCERV-P+HPV - P699-2.....	B408
PAP/PSA - 0956-3.....	B408
PAP/HPV+/GC RFX>GENCERV-PAP+HP (NON-NY) - P701-6	B408
Parainfluenza Antibody - 0946-4.....	B409
Path Cons On Ref.BLK/SLD GLOB. - 5283-7	B409
Pathology Slide Consultation - 5111-0	B409
PBG Deaminase, WB - 3245-8	B410
Procalcitonin - T969-5.....	B410
PD-1 BY IHC - B994-5	B410

PD-L1 (SP142) , TECENTRIQ. - T413-4	B411
PD-L1 22C3 for Keytruda Other Tumor Types - TK85	B411
PD-L1 IHC, Keytruda for Cervical Cancer W/INTERP. - T988-5	B411
PD-L1 (Clone E1L3N) IHC w/Interp. - J257-7	B412
Keytruda PD-L1 Head & Neck Squamous Carcinoma (HNSCC) - TH63-1	B412
PD-L1 IHC 22C3 pharmDx, Keytruda W/INTERP. - B993-7	B412
PD-L1 IHC 28-8 pharmDx (OPDIVO) W/INTERP. - J221-3.....	B413
PD-L1 (Clone E1L3N) IHC Tech Only - TB35-1	B413
PD-L1 IHC , Keytruda for Urothelial Carcinoma W/INTERP. - TA31-2.....	B413
PDGFRA Sequencing- Non NY - J045-6	B414
PDGFRA Sequencing - J459-9	B414
PDGFRB by FISH (Global)- Non NY - T781-4.....	B414
PDGFRB by FISH (Tech Only) Non NY - T782-2	B415
Keytruda PD-L1 Esophageal Squamous Cell Carcinoma by IHC - TH62-3.....	B415
PDL1 (22C3 Keytruda) IHC + EGFR NGS + ALK FISH + ROS1 FISH - L393-8.....	B415
Pediatric Arthritis Profile - 0596-7	B416
Pediatric Blood Culture #1 - B331-0	B416
Pediatric Blood Culture #2 - B332-8	B416
Pediatric Blood Culture #3 - B333-6	B417
Pediatric Profile I - 0980-3.....	B417
Pediatric Profile II - 0979-5.....	B417
Perennial Allergypanel (3 mix) - 2188-1	B418
PET II (Dialysis Fluid) - 4447-9	B418
Phenosense Entry (incl.Fuzeon) - 6235-6.....	B418
Phlebotomy Charge - 0305-3.....	B419
Phosphatidyl Ethanolamine - 5722-4	B419
Phosphorus, Urine, Random - 3134-4.....	B419
Phosphates, Urine 24 Hr. - 0411-9.....	B419
PID/Infertility W/O PAP Expanded - F368-6	B420
Pinworm Tape Test - 0204-8.....	B420
Pituitary I - 2319-2	B420
Plasma Cell Dyscrasia/Multiple Myeloma IHC Panel - 5579-8	B421
Plasma Cell Analysis by Flow Cytometry (Tech Only) and StormPath BMM Report Template - R736-0.....	B421
Plasminogen Quant. - 5775-2	B421
Platelet Antibody Identification - 5776-0	B422
PNH by Flow Cytometry - 5564-0	B422
Platelet Neutralizing Procedure - 5784-4.....	B422

Polycyst.Ovary Synd. Treat. Panel - P259-5	B423
Polycyst.Ovary Synd. Conf.Panel - P258-7	B423
Post Dialysis Potassium - 1688-1	B423
Potassium, Urine, Timed Specimen - 1659-2	B424
Proteinase-3 AutoAbs - 2663-3	B424
15q Del Characterization, FISH - B796-4	B424
Preg.+Throm.Assessment - 6801-5.....	B425
Pregnenolone, Serum - J150-4	B425
Prenatal + HIV - 2765-6.....	B425
Procainamide/NAPA - 0525-6.....	B426
PROCESS AND HOLD FOR FISH - B664-4	B426
Progesterone, Serum - 0335-0	B426
Prolactin, Serum - 0134-7	B427
Propoxyphene, Urine, Screen (300) - 0963-9	B427
Prostate Cancer Triple Stain By IHC (PC) - B964-8.....	B427
Prostate Cancer/ Lymph Nodes - 5122-7.....	B428
Prostate Core Biopsy (1 Spec.) - 5611-9.....	B428
Prostate Core Biopsy (8 Spec.) - 5618-4.....	B428
Prostate Core Biopsy (4 Spec.) - 5614-3.....	B429
Prostate Core Biopsy (3 Spec.) - 5613-5.....	B429
Prostate Core Biopsy (20 Spec.) - 5174-8.....	B429
Prostate Core Biopsy (2 Spec.) - 5612-7	B430
Prostate Core Biopsy (7 Spec.) - 5617-6.....	B430
Prostate Core Biopsy (5 Spec.) - 5615-0.....	B430
Prostate Core Biopsy (6 Spec.) - 5616-8.....	B431
Prostate Core Biopsy (12 Spec.) - 5622-6.....	B431
Activated Protein C Resistance Modified - 5705-9.....	B431
Protein C Antigen - 5786-9.....	B432
Activated Protein C Resistance - 5704-2	B432
Protein S Activity - 5787-7	B433
Protein S Antigen Total - 5789-3	B433
Protein S Antigen, Free - 5788-5	B434
Protein, Urine, 24 Hr. - 0136-2	B434
Protein/Creatinine w/Ratio - 3329-0.....	B434
Prothrombin ABS IGG - 5820-6	B435
COAG/Prothrombin Time - 5798-4	B435
PSA Free + Total - 2088-3	B435
PSA Reflex to Free PSA - M067-7.....	B436

PT + PTT - 2237-6.....	B436
Prothrombin G20210A Mutation - 5795-0.....	B436
Prothrombin Time/INR (PT) - 0137-0.....	B437
PTEN BY FISH - B854-1.....	B437
PTEN BY IHC With Interpretation - A985-4.....	B437
PTH, Intact w/ Calcium - 1689-9.....	B438
PTH, Intact (use 1689 if with Calcium) - 0598-3.....	B438
PTH, Intact, FNA - J047-2.....	B438
Part.Thromb.Time(PTT)-LA - 5771-1.....	B439
Panel A (Inhalants) - 1632-9.....	B439
Pap + Pap Dependent HPV - B975-4.....	B439
Pap + Pap Dependent HPV + CT/GC - P734-7.....	B440
Pap (Reflex to HPV DNA Genotyping 16,18 when ASC-US or UNSAT) - L326-8.....	B440
Pap (Reflex to HPV DNA Genotyping 16,18 when ASC-US) - 6195-2.....	B440
Pap (Reflex to HPV mRNA Genotyping 16, 18/45 when ASC-US) - A251-1.....	B441
Pap and HPV DNA Genotyping 16,18 - P079-7.....	B441
Pap and HPV mRNA Genotyping 16, 18/45 - A250-3.....	B441
PAP/CT/GC(RFLX HPV PLUS ASC-US or >) - 9985-3.....	B442
Allergen Papaya (f293), IgE - 1942-2.....	B442
Allergen Paper Wasp (i4), IgE - 0886-2.....	B442
Allergen Paprika (f218),IgE - 1424-1.....	B443
Ova and Parasite (O&P), Stool W/ Trichrome (1st specimen) - 0377-2.....	B443
Allergen Parsley (f86), IgE - 1355-7.....	B443
Parvovirus B19 Ab, IgG - T790-5.....	B444
Parvovirus B19 Ab, IgM - T791-3.....	B444
Parvovirus B19 Ab (IgG/IgM) - L399-5.....	B444
Allergen Passion Fruit (rf294), IgE - 3480-1.....	B445
Lead, Blood (Adult) - 1438-1.....	B445
Lead, Blood (Child) - 0398-8.....	B445
Allergen Peach (f95), IgE - 1226-0.....	B446
Allergen Peanut (f13),IgE - 0674-2.....	B446
Allergen Peanut Component Ara (h1), IgE - J284-1.....	B446
Peanut Component Ara (h2), IgE - J283-3.....	B447
Peanut Component Ara (h3), IgE - J285-8.....	B447
Peanut Component Ara (h6), IgE - TG35-1.....	B447
Peanut Component Ara (h8), IgE - J287-4.....	B448
Peanut Component Ara (h9), IgE - J286-6.....	B448
Peanut Components IgE - H993-9.....	B448

Peanut w/Reflex to Component Panel - M065-1	B449
Allergen Pear (f94), IgE - 0930-8	B449
Allergen Pecan Nuts (f201), IgE - 1592-5	B449
Allergen Pecan/Hickory Tree (t22), IgE - 1938-0	B450
Pediatric Mix (fx5), IgE - 2186-5	B450
Allergen Penicillium Notatum Mold (m1), IgE - 2167-5	B450
Allergen Perennial Rye Grass (g5), IgE - 0682-5	B451
Perimenopausal Hormonal Profile - 8431-9	B451
Peroxisome Biogenesis Disorder 5A (Zellweger) - J303-9	B451
Culture, Throat, Routine - 0078-6	B452
Phenylketonuria (PKU) Carrier Screening - J304-7	B452
Dilantin - 0084-4	B452
BCR/ABL: t(9;22) by FISH - 5265-4	B453
Allergen Phoma Betae (m13), IgE - 2477-8	B453
Phosphoglycerate dehydrogenase deficiency - J302-1	B453
Allergen Pignolia (Pine) Nuts (t213), IgE - 3638-4	B454
Allergens Pineapple (f210), IgE - 1389-6	B454
Allergen Pistachio (rf203), IgE - 1195-7	B454
Plasma 10-Color (Short) Panel - TK94	B455
Plasma 10-Color (Short) Panel TC Only - TK95	B455
Plasma 10-Color Myeloma/Monoclonal Gammopathy - TK97	B455
Plasma 10-Color Myeloma/Monoclonal Gammopathy TC Only - TK96	B455
Glucose, Non-Fasting - 0223-8	B455
Heparin Platelet Antibody - 5746-3	B456
Allergen Plum (f255), IgE - 1584-2	B456
Pneumococcal Antibody Panel (14 Serotypes) - 1250-0	B456
Polycystic Kidney Disease, Autosomal Recessive (ARPKD) - J301-3	B457
Allergen Poplar White IgE - 3704-4	B457
Allergen Poppy Seed (rf224), IgE - 1031-4	B457
Allergen Pork (f26), IgE - 0836-7	B458
Potassium, RBC - J571-1	B458
Allergen Potato (f35), IgE - 0844-1	B458
Obstetric (AMA) Panel - 0010-9	B459
Ciliary Dyskinesia, Primary 1 - J311-2	B459
Allergen Privet Tree (T210), IgE - A321-2	B459
Process and Hold for Cytogenetics - J655-2	B460
Tacrolimus (FK506) - 1871-3	B460
Prostate Cancer (Stain Only) - 5104-5	B460

Hereditary Prostate Cancer Panel - J665-1	B461
Prostate Carcinoma IHC Panel (CK903, p63, P504S) - 5097-1	B461
PSA Total - 0190-9	B462
Protein C Activity - 5785-1	B462
Protein Electrophoresis, Urine, Random - 0404-4.....	B462
Protein,Urine,Timed or Random - 1662-6	B463
Allergen Pumpkin (f225), IgE - 1488-6	B463
Allergen Pumpkin Seed (rf226), IgE - 1484-5	B463
Pyruvic Acid (Pyruvate), Blood - 3767-1	B464
Methaqualone, Urine, Screen (300) - 3239-1	B464
Alpha-Fetoprotein (AFP 4) includes HCG, UE3, Inhibin-A - 3158-3	B464
Pregnancy, Serum (Qualitative) - T824-2	B465
Allergen R. Nigricans (m11), IgE - 1373-0	B465
Rf (Rheumatoid Arthritis) Titer - 0796-3.....	B465
Allergen Red Raspberry (rf343), IgE - 3508-9.....	B466
RB1 for 13q14 by FISH - B511	B466
RBC - 1498-5.....	B466
Reactive Hyperplasia vs. Lymphoma (NHL) IHC Panel - 5103-7	B467
Renal Function (Medicare) - 3284-7.....	B467
Reptilase Time - 5799-2.....	B467
Rheumatoid Factor IgA - TA08-0	B468
Rheumatoid Factor IgM - TA09-8	B468
ROMA - L012-4.....	B468
ROS-1 BY FISH-TECH Only - B337-7	B469
ROS1 BY FISH Global - B334-4	B469
POST TREATMENT SYPHILIS (RPR) - R733-7	B469
RPR Serology - 0142-0	B470
Resp. Syncytial Virus - 1397-9.....	B470
Rubella Ab, IgG, Qual, Immune Status - 3151-8.....	B470
Sirolimus, Blood - 3481-9	B471
Allergen Redtop (Bent) Grass (g9), IgE - 1199-9	B471
Rest of common cancer panel after Lynch Syndrome Panel B985 - B984-6	B471
Rest of Comprehensive Cancer Panel - B363-3	B472
Rest of Comprehensive Cancer Panel (after Prostate Panel) - J664-4	B472
Reflex to Remainder of Common Cancer Panel (after J055) - J054-8.....	B473
Region 1 Allergy Profile (North East) - 6823-9	B473
Region 10 Allergy Profile (TX,OK) - 6827-0	B473
Region 11 Allergy Profile (AZ, ID, NM, WY, CO, UT) - K908-4	B473

Region 13 Allergy Profile (S. CA Coast) - 6820-5.....	B474
Region 14 Allergy Profile (Central CA) - 6821-3.....	B474
Region 17 Allergy Profile (North West) - 6822-1.....	B474
Region 2 Allergy Profile (Mid Atlantic) - 6824-7.....	B474
Region 3 Allergy Profile (South East) - 6825-4.....	B475
Region 4 Allergy Profile (S.Florida) - 6826-2.....	B475
Region 5 Allergy Profile (IN, OH, TN, WV, KY) - 6829-6.....	B475
Renal Cancer Panel - B394.....	B475
Renal Cancer Panel - B394.....	B475
Renal Cancer Panel - B394.....	B475
Rest of Demyelinating CMT Panel (after PMP22 Del/Dup) - J795-6.....	B476
Reticulocyte Count - 0141-2.....	B476
RET-HE - 5848-7.....	B476
Retinosis Pigmentosa 59 - J313-8.....	B477
Vitamin A (Retinol), Serum - 0337-6.....	B477
Reflex to GenCerv - A780-9.....	B477
Rheumatoid Factor IgG, IgA, IgM - L626-1.....	B478
Vitamin B2 - 3325-8.....	B478
Allergen Rice (f9), IgE - 0843-3.....	B478
Allergen Rough Careless Pigweed (w14), IgE - 1627-9.....	B479
Rubella Ab IgG W/Reflex - T385-4.....	B479
RUBELLA IGM - J900-2.....	B479
Rubella Ab, IgG, Quant - 0973-8.....	B480
Allergen Russian Thistle (w11), IgE - 0732-8.....	B480
Allergen Rye (food) (f5), IgE - 0822-7.....	B480
Hemoglobin Fractionation, Capillary Electrophoresis - 0216-2.....	B481
Staphylococcus aureus (Non NY) - J970-5.....	B481
S. Cerevisiae Antibodies (IGG/IGA) - H940-0.....	B481
Allergen Sardine, IgG (NON-NY) - B424-3.....	B482
Scleroderma Profile - 0593-4.....	B482
Selenium Serum/Plasma - 0232-9.....	B482
Send Copy Of Results To: - 9991-1.....	B483
Sequential Screen Part 1 (NT, Total BHCG, PAPP-A) - 8602-5.....	B483
Sequential Screen Part 2 (AFP, UE3, HCG, Inhibin-A) - 8601-7.....	B483
Serum Integrated Part 1 (PAPP-A) - A335-2.....	B484
Serum Integrated Part 2 (AFP, UE3, HCG, Inhibin-A) - A334-5.....	B484
Aspartate Aminotransferase (AST/SGOT) - 0146-1.....	B484
SGPT (ALT) - 0147-9.....	B485

Sex-Hormone Binding Globulin - 0658-5	B485
Silicate, Water - 6219-0.....	B485
Sinus Culture W/ Gram Stain - F245-6.....	B486
Skin Biopsy, (2 Jars) - 5603-6.....	B486
Skin Biopsy, (3 Jars) - 5604-4.....	B486
Skin Biopsy, (4 Jars) - 5605-1	B487
Skin Biopsy, (5 Jars) - 5606-9.....	B487
Skin Biopsy, (6 Jars) - 5607-7.....	B487
Spinal Muscular Atrophy (SMA) Carrier (SMN1 Dosage) - 6380-0.....	B488
Chem 24 (SMAC) - 0005-9	B488
Small Round Cell Tumor IHC Panel - 5177-1	B488
Smith-Lemli-Opitz Synd.Prenat (Non-NY) - B139-7	B489
SOX-10 BY IHC W/Interp - J147-0.....	B489
SOX-10 BY IHC TECH Only - J146-2.....	B489
Special Stain Enzyme Constituents - PC - B756-8	B490
Special Stain Enzyme Constituents -TC - B758-4	B490
Special Stain Enzyme Constituents - Global - B757-6	B490
Specific Factor Inhibitor - 5818-0	B491
Sperm Count (Post Vasectomy) - 0145-3	B491
Spindle Cell Neoplasm IHC Panel - 5191-2.....	B491
Spore Strips - 0402-8	B492
Lung - Adenocarcinoma vs. Squamous cell carcinoma - A130-7.....	B492
Synovial Sarcoma 18q11.2 FISH - J331-0.....	B492
Allergen StemphyliumBotryosum(m10),IgE - 3828-1	B493
STI Lesion (NON-NY) - Q301-4	B493
STI Serology - P216-5	B493
Stormpath PAP + HPV HR W/Rfx to HPV 16 and 18/45 - Q723-9.....	B494
Stormpath PAP+CT/GC+HPV HR w/RFX to HPV 16 and 18/45. - Q724-7	B494
Strep Group B (Amplified Probe) - 0202-2.....	B494
Pneumococcal Antibody Panel (23 Serotype) - J048-0	B495
Surgical Pathology,Level IV - 5619-2.....	B495
SURGICAL PATHOLOGY LEVEL 2 - 5224-1	B495
Allergen Salmon (f41), IgE - 0714-6.....	B496
Beta-Hexosaminidase B Enzyme Analysis (Sandhoff Disease) - B761-8	B496
Allergen Sardine (Pilchard) (f308), IgE - 6266-1	B496
Allergen Scallop (f338), IgE - 2543-7	B497
Hypercoagulation - 2321-8	B497
Rast/Seafood Mix(fx2) IgE + - 1843-2.....	B497

Culture, Seminal Fluid - 1827-5	B498
Serum Drug Screen - L063-7.....	B498
Allergen Sesame Seed (f10), IgE - 0936-5.....	B498
Allergen Sheep Epithelium (e81), IgE - 1948-9	B499
Allergen Sheep Sorrel (w18), IgE - 0894-6	B499
Allergen Shrimp (f24),IgE - 0696-5	B499
Allergen Silk (k74), IgE - 1253-4	B500
Extractable Nuclear Antigen (ENA) Antibodies (Sm/RnP) - 0318-6.....	B500
Smith-Lemli-Opitz Syndrome - J314-6	B500
Smooth Muscle Antibody (ASM) - 0419-2	B501
Allergen Soybean (f14), IgE - 0680-9	B501
Spastic Tetraplegia, Thin Corpus Callosum and Progressive Microcephaly - J298-1	B501
Allergen Trichoderma Viride (m15), IgE - 3829-9	B502
Culture, Sputum W/ Gram Stain - 0382-2	B502
Culture, Stool (1st Specimen) - 0077-8.....	B502
Allergen Storage Mite (Acarus) (d70), IgE - 1605-5	B503
StormPath PAP(RFLX to HPV mRNA GENOTYPING 16,18/45 w/ASC-US) - TB33-6	B503
StormPath PAP, liquid based + HPV DNA Genotyping 16/18 - L079-3	B503
StormPath PAP, liquid based + HPV DNA + CT/GC - L080-1	B504
Stormpath Pap/GC/CT reflex HPV when ASC-US - Q725	B504
Allergen Strawberry (f44), IgE - 0931-6.....	B504
Strep. Group A Screen DNA - 0079-4.....	B504
Allergen Sunflower Seed (k84), IgE - 1414-2.....	B505
High Risk STI- CT/GC- Aptima Swab - L062-9.....	B505
High Risk STI-CT/GC/Trich-Aptima Swab - L047-0.....	B505
High Risk STI- CT/GC/Trich/M.Gen-Aptima Swab - L048-8.....	B506
Aptima Swab Cervicitis (NON-NY) - L049-6.....	B506
Allergen Sweet Potato (f54), IgE - 1231-0.....	B506
Allergen Sweet Vernal (g1), IgE - 1196-5	B507
Allergen Sycamore (t11), IgE - 1305-2.....	B507
C-KIT (D816V) - 5179-7	B507
Immune Deficiency Panel CD4/CD8 w/CBC - 1763-2	B508
T(11;14)+11 BY FISH-TC Only - A873-2.....	B508
T(4;14) BY FISH-TC Only - A875-7	B508
T.Pallidum By Multiplex PCR (NON-NY) - 6354-5	B509
T3, Total - 0150-3.....	B509
T3 Uptake - 0152-9.....	B509
T3 Free, Serum - 0271-7	B510

T4 (Thyroxine) Total - 0151-1	B510
EGFR Sequencing- Non NY - J022-5.....	B510
TAO Brush Endometrial Cytology - B685-9.....	B511
Tay-Sachs Disease Carrier Test - A652-0.....	B511
Osteopetrosis, Autosomal Recessive 1 - J295-7.....	B511
Lymphoma - T-Cell, TCR Gamma by PCR - 5031-0.....	B512
TEL/AML1: t(12;21) by FISH - 5038-5.....	B512
Testosterone Free + Total - 7301-5.....	B512
Testosterone, Serum - 0379-8.....	B513
Cannabinoids (Marijuana), Urine, Screen (50) - 3971-9.....	B513
Cannabinoids (Marijuana), Urine, Screen (20) - 1576-8.....	B513
Thrombin Time - 5807-3.....	B514
Thrombin/Antithrombin - 5808-1	B514
Thyroglobulin, Serum - 0577-7.....	B514
Thyroid Comprehensive Panel - 0867-2.....	B515
Thyroid PNL-T3U,T4,TSH - 0007-5.....	B515
Thyroid Screen - 0006-7	B515
Thyroiditis - 2322-6	B516
Tissue Path.(POC) - 3001-5.....	B516
Tissue Pathology, Gross Only - 4002-2.....	B516
Tobramycin, Peak - 1235-1	B517
Tobramycin (Random) - A615-7	B517
Tobramycin, Trough - 2132-9	B517
Toxoplasma Antibodies, (IgG, IgM) - 2721-9.....	B518
Toxoplasma IgG W/Reflex - 2376-2.....	B518
Toxoplasma Ab. (IgG), EIA - 0552-0.....	B518
Toxoplasma Ab. (IgM) - 0877-1	B519
TP53 Sequencing - J044-9.....	B519
OnkoSight NGS TP53 Sequencing for Solid Tumors - T980-2	B519
Tramadol, Urine, Screen, Qual. (200) - 5917-0.....	B520
Transglutaminase, IgA Abs. - 3472-8	B520
Transglutaminase, IGG ABS. - 5869-3.....	B520
Trichomonas Vaginalis, Probe rRNA - A861-7	B521
Trimipramine, Serum - 0613-0	B521
Trisomy 12 by FISH-TC Only - A877-3	B521
Trisomy 8 By FISH -TC Only - A870-8	B522
Trofile DNA - A205-7	B522
Troponin I - 2163-4.....	B522

TSH W/ Reflex To Free T4 - A518-3	B523
Thyroid Stimulating Hormone (TSH) - 0153-7	B523
TSH (Pregnancy) - J265-0	B523
Thyroid Stimulating Immunoglobulin - J510-9	B524
Transglutaminase IGG/A,AB, Evaluation - 3473-6	B524
Tumor of Unknown Primary IHC Panel - 5190-4	B524
Tyrosinemia, Type 1 - B149-6	B525
Tysabri Antibodies - 6293-5	B525
Hexosaminidase Leukocytes (Tay-Sachs Enzyme) - 0865-6	B525
Allergen Tea (rf222), IgE - 1390-4	B526
Carbamazepine - 0154-5	B526
Total Testosterone by LC/MS/MS - J451-6	B526
Throat Culture, R/O Group A Strep. ONLY - TJ06-6	B527
Culture, Fungal (Yeast) - A502-7	B527
Calcitonin, Serum - 0196-6.....	B527
Anti-Thyroglobulin Ab. (ATG) - 0041-4	B528
Anti-Microsomal Ab(Thyroid Peroxidase)TPO - 0412-7	B528
Tick Identification - J291-6.....	B528
Allergen Tilapia (RF414), IgE - 8963-1	B529
Allergen Timothy Grass (g6), IgE - 0683-3	B529
Tissue Pathology One Check - T793	B529
Allergen Tomato (f25), IgE - 0837-5	B530
SCL70 (Scleroderma Ab.), Serum - 0315-2.....	B530
TORCH Panel (IgG/IgM) - 1297-1	B530
Bilirubin, Infant, Total/Direct, Serum - 0046-3	B531
Immunoglobulin E (IgE), Serum/Plasma - 0996-9	B531
Total IgG and subclasses IgG 1, 2, 3 & 4 - H991-3	B531
Protein, Total, Serum - 0135-4	B532
Testosterone, Total/Free (LC/MS) - K630-4.....	B532
Prealbumin - 1006-6	B532
Tree Mix (tx1), IgE - 2178-2	B533
Trees Screen (6) - 0742-7.....	B533
Serodia TPPA (Antibody to T. pallidum) - 0654-4	B533
Liaison CIA (Antibody to T. pallidum) - 0334-3	B534
Lipase, Serum - 0521-5.....	B534
Trichomonas By Multiplex PCR - 6342-0.....	B534
Trichomonas by MultiPlex PCR - swab - J836-8.....	B535
Allergen Trichophyton Rubrum (m205), IgE - A673-6	B535

Triglycerides - 0155-2	B535
Allergen Trout (f204), IgE - 3643-4	B536
Allergen Tuna (f40), IgE - 0698-1	B536
Allergen Turkey Feathers (e89), IgE - 1426-6	B536
Allergen Turkey Meat (f284), IgE - 1149-4	B537
Acetaminophen, Serum - 0622-1	B537
Urinalysis w/Reflex to Culture - B518-2	B537
H.Pylori Urea Breath Test - 6236-4	B538
Protein Electrophoresis, 24Hr Urine - Q069-7	B538
Urea Clearance, 24hr Urine - 1698-0	B538
Urea, Urine 24hr - 0358-2	B539
Urea, Urine, Timed - 2083-4	B539
Ureaplasma By Multiplex PCR - 6339-6	B539
Uric Acid, Synovial - 0826-8	B540
Uric Acid, Urine, Timed Spec. - 1663-4	B540
Urinalysis, Routine - 0159-4	B540
Cytopathology, Urine (Thinprep) - 5254-8	B541
URINE FREE KAPPA+LAMBDA - L592-5	B541
Urogenital Mycoplasmas - P114-2	B541
U.R.R. PANEL - 1829-1	B542
Usher Syndrome Type 1F - A649-6	B542
Uric Acid, Serum - 0157-8	B542
Uric Acid, Urine, 24 Hr. - 0158-6	B543
Urethral Smear - 1500-8	B543
Creatinine Urine, Timed or Random - 1661-8	B543
Urine Culture - 0080-2	B544
Urine Culture, Prenatal, W/GBS - TH76-3	B544
Urine Cytology Reflex Atypical or above to UroVysion - B415	B544
Urine Cytology, if Atypical, reflex to hTert - TG28-6	B544
Urine Cytopathology, TC (Alcohol) - 5554	B545
Comprehensive Urine Pathology - 5324-9	B545
Comprehensive Urine Pathology (Technical Only) - 5253-0	B545
UroVysion (Tech Only) - 5249-8	B545
Bladder Washing, FISH - J371-6	B546
Usher Syndrome Type 3 - A647-0	B546
Vaginitis/Vaginosis W/O PAP W/O HPV Expanded - F371-0	B546
Vaginal Biopsy - 5889-1	B547
Vancomycin Peak, Serum - 0575-1	B547

Vancomycin (Random) - A614-0	B547
Vancomycin, Trough, Serum - 1589-1	B548
Venous Prothrom. Risk Factor - 5875-0	B548
Venous Prothrombic Risk - 5989-9	B548
Vitamin B12/Folate - 0287-3	B549
Vitamin B12 Binding Capacity - 0906-8	B549
Vitamin A + Beta Carotene - 1332-6	B549
Vitamin Panel - 2861-3	B550
Vitamin B1, Total Thiamine - 0293-1	B550
Von Willebrand Screen- Basic - 5825-5	B550
Vulvar Biopsy - 5886-7	B551
Varicella-Zoster Virus - J839-2.....	B551
Varicella Zoster Virus by IHC - TF85-8	B551
Varicella Zoster Virus by IHC Tech Only - TF88-2	B552
Benzodiazepines, Urine, Screen (200) - 0924-1	B552
Allergen Vanilla (rf234), IgE - 2347-3.....	B552
Varicella Zoster Ab IgG W/Reflex - T008-2	B553
Allergen Profile, Insects (6) - 2469-5	B553
Vitamin C Dilution - 6193-7	B553
Vit. D, 25-Hydroxy Tot.,D2+D3 - 8573-8	B554
Ristocetin Cofactor ACT. - 5805-7	B554
Vulvovaginal Candidiasis Panel by PCR - K853-2	B554
Walker Warburg Syndrome - A648-8.....	B555
AAMI Water Analysis - 1795-4	B555
WBC - 1497-7	B555
Bioavailable Testosterone - 2087-5	B556
Wilson Disease Carrier Screen - B164-5.....	B556
Female High Risk Sexual Health Panel - F238-1	B556
Wound Culture (Aerobic/Anaerobic) W/ Gram Stain (Panel) - 0082-8	B557
Allergen Walnut Tree (t10), IgE - 0686-6.....	B557
Allergen Walnuts (f256), IgE - 1479-5.....	B557
Allergen Watermelon (f329), IgE - 1941-4	B558
Weed Mix (wx1), IgE - 2177-4.....	B558
Weeds Screen (7) - 0746-8	B558
Allergen Western Ragweed (w2), IgE - 3802-6.....	B559
Allergen Wheat (f4), IgE - 0673-4.....	B559
Allergen Whey (f236), IgE - 1620-4	B559
Allergen White Ash (t15), IgE - 1306-0	B560

Allergen White Bean (f15), IgE - 1241-9	B560
Allergen White Pine (t16), IgE - 3399-3	B560
Allergen Whiteface Hornet (i2), IgE - 0907-6	B561
Allergen Willow (t12), IgE - 1252-6	B561
Allergen Wormwood (w5), IgE - 1307-8	B561
Culture, Wound, Superficial - B656-0	B562
XCE 12 by FISH - B513	B562
XL RARA (17q21) BY FISH - J157-9	B562
XL RARA (17q21) BY FISH - TC (NON-NY) - J160-3	B563
Xylose Absorption, Urine - 0530-6	B563
Xylose, (Baseline) Whole Blood - 1184-1	B563
Allergen Yellow Hornet (i5), IgE - 0885-4	B564
Allergen Yellow Jacket (i3), IgE - 0884-7	B564
ZAP-70 BY IHC - 5410-6	B564
Zika Molecular Serum+Urine - K559-5	B565
Zika Serum RNA PCR w/reflex to IgM - J405-2	B565
Zika Serum RNA PCR - J406-0	B565
Anti-Zika Virus Elisa IgM Serum - J497-9	B566
Zika Urine RNA PCR - J407-8	B566
Zika Molecular Serum w/rfx and Urine - K583-5	B566
Zinc Protoporphyrin (ZPP), Whole Blood, Adult - 6289-3	B567
Zinc Protoporphyrin (ZPP), Whole Blood, Child - 0407-7	B567
Zinc, RBC - 1336-7	B567
Zinc, Serum - 0285-7	B568
Zinc, Urine, Random - 2440-6	B568
Hepatitis Be Antibody - 0540-5	B568
Allergen Cashew (rf202), IgE - 1412-6	B569
Allergen Cherries (f242), IgE - 1583-4	B569
Allergen Cotton-Wood (t14), IgE - 0891-2	B569
Allergen Milk (f2), IgE - 0672-6	B570
dsDNA Antibody IgG, Crithidia IFA - 3825-7	B570
HCG, Quantitative, Serum - 0327-7	B570
Hcg, Quantitative, Tumor Marker - 1201-3	B571
Human Telomerase Reverse Transcriptase - Tech Only - TF87-4	B571
p16 by IHC - J356-7	B571
17P13 BY FISH-TC ONLY - A874-0	B572
p53 by IHC w/digital imaging - 5153-2	B572
p57 by IHC - J546-3	B572

Reverse Triiodothyronine Serum/Plasma (rT3) - T397-9	B573
Pap (Reflex to HPV DNA Genotyping 16, 18 when ASC-US or >) - 6196-0	B573
Lymphoma - Mantle Cell (MCL), BCL1/IGH by FISH - 5026-0.....	B573
Lymphoma - Follicular, BCL2/IGH by FISH - 5270-4	B574
Lymphoma - Follicular, BCL2/IGH by PCR - 5271-2	B574
t(1;19) (q23.3: p13.3) FISH - B611-5.....	B574
Von Willebrand Factor, Antigen - 5815-6	B575

Index by Test Code

0002-6 - Electrolytes B194

0004-2 - Chem 12 Profile..... B106

0005-9 - Chem 24 (SMAC) B488

0006-7 - Thyroid Screen B515

0007-5 - Thyroid PNL-T3U,T4,TSH..... B515

0008-3 - Obstetric/Prenatal I B389

0009-1 - Lipid Screen (Basic Lipid Profile)..... B81

0010-9 - Obstetric (AMA) Panel B459

0012-5 - Liver Profile..... B342

0019-0 - Hematocrit..... B258

0022-4 - Cyclosporine, Whole Blood B144

0024-0 - Insulin Tolerance (3 Hour) B304

0025-7 - Alpha Fetoprotein, Tumor Marker (AFP) B13

0028-1 - Ketones (Acetoacetic), Serum B324

0030-7 - Prostate Specific Acid Phos (P.A.P.)..... B401

0031-5 - Acid Phosphatase, Total, Serum..... B0

0033-1 - Albumin, Serum B42

0034-9 - CBC w/o DIFF (Hemogram)/PLAT. CT B152

0035-6 - Ammonia, Plasma (As Nitrogen)..... B380

0036-4 - Amylase, Serum B46

0037-2 - Amylase, Urine, Random B22

0038-0 - Antinuclear Ab. (ANA) W/Reflex B23

0039-8 - Antibody Screen (Ind.Coombs) B50

0041-4 - Anti-Thyroglobulin Ab. (ATG)..... B528

0042-2 - Antistreptolysin O (ASO) Titer B32

0043-0 - Bilirubin, Total, Serum..... B85

0044-8 - Bilirubin, Direct, Serum..... B85

0045-5 - Bilirubin, Total + Direct, Serum..... B85

0046-3 - Bilirubin, Infant, Total/Direct, Serum B531

0047-1 - Coxsackie Virus Antibody, Group A Only..... B0

0049-7 - Blood Urea Nitrogen (BUN), Serum B76

0050-5 - Calcium, Serum..... B129

0051-3 - Calcium, Ionized, Serum B129

0052-1 - Carbon Dioxide, Serum (CO2)..... B112

0053-9 - CBC w/Diff, Platelet Ct..... B152

0055-4 - Carcinoembryonic Antigen (CEA)..... B103

0056-2 - Cell Count, Fluid B104

0057-0 - Chloride, Serum	B145
0058-8 - Cholesterol.....	B140
0059-6 - HDL Cholesterol.....	B259
0063-8 - Cold Agglutinins Titer	B150
0064-6 - Coombs Test, Direct	B174
0066-1 - Cortisol, AM, baseline.....	B117
0067-9 - Cortisol, PM	B153
0068-7 - CPK (CK), Total	B109
0069-5 - C-Reactive Protein (CRP), Serum.....	B95
0070-3 - Creatinine, Serum.....	B157
0071-1 - Creatinine, Urine, 24 HR.	B121
0072-9 - Creatinine Clearance (24Hr. Urine)	B120
0074-5 - Culture, Nasopharyngeal/Nasal.....	B384
0075-2 - Culture, Neisseria gonorrhoeae (GC).....	B249
0076-0 - Culture, Eye.....	B391
0077-8 - Culture, Stool (1st Specimen).....	B502
0078-6 - Culture, Throat, Routine	B452
0079-4 - Strep. Group A Screen DNA.....	B504
0080-2 - Urine Culture	B544
0082-8 - Wound Culture (Aerobic/Anaerobic) W/ Gram Stain (Panel)	B557
0083-6 - Digoxin, Serum.....	B43
0084-4 - Dilantin	B452
0085-1 - Protein Electrophoresis	B195
0086-9 - ESR (Sed-Rate)	B198
0087-7 - Culture, Ear	B191
0088-5 - Ferritin	B208
0089-3 - Fibrinogen, Activity	B210
0090-1 - Folate (Folic Acid)	B220
0091-9 - Thyroxine, Free (FT4).....	B223
0092-7 - FSH (Follicle Stimulating Hormone).....	B216
0093-5 - GGTP (Gamma-GT)	B224
0094-3 - GC Smear	B249
0095-0 - Glucose, Fasting.....	B207
0096-8 - Glucose, Urine 24HR Quantitative.....	B234
0097-6 - Glucose Tolerance 3 HRS.....	B233
0098-4 - Glucose Tolerance 4 HRS.....	B233
0099-2 - Glucose Tolerance 5 HRS.....	B233
0102-4 - Hemoglobin A1c (Glycohb).....	B247

0103-2 - Gram Stain	B237
0104-0 - Hemoglobin + Hematocrit	B260
0105-7 - Hepatitis A Ab.,Total W/Reflex	B261
0106-5 - Hepatitis B Surface Antigen.....	B258
0107-3 - Hepatitis B Surface Antibody (Qual/Immunity)	B257
0108-1 - Hepatitis B Core Antibody, Total	B257
0109-9 - Herpes I/II Ab. Screen, IgG.....	B278
0110-7 - Herpes I AB., IgG	B285
0113-1 - Insulin, Total/Random	B310
0114-9 - Iron, Serum.....	B313
0117-2 - Lactate Dehydrogenase (LDH)	B339
0119-8 - Lithium, Serum.....	B144
0120-6 - Magnesium, Serum	B368
0121-4 - Mono Screen	B371
0122-2 - Occult Blood,Stool (1st Specimen)	B89
0127-1 - Phosphorus, Inorganic.....	B312
0129-7 - Potassium, Serum.....	B317
0133-9 - Pregnancy, Urine	B84
0134-7 - Prolactin, Serum	B427
0135-4 - Protein, Total, Serum	B532
0136-2 - Protein, Urine, 24 Hr.	B434
0137-0 - Prothrombin Time/INR (PT).....	B437
0139-6 - PTT (Partial Thromboplastin Time).....	B401
0140-4 - Quinidine, Serum.....	B86
0141-2 - Reticulocyte Count.....	B476
0142-0 - RPR Serology	B470
0145-3 - Sperm Count (Post Vasectomy)	B491
0146-1 - Aspartate Aminotransferase (AST/SGOT).....	B484
0147-9 - SGPT (ALT).....	B485
0148-7 - Sodium, Serum.....	B383
0150-3 - T3, Total.....	B509
0151-1 - T4 (Thyroxine) Total.....	B510
0152-9 - T3 Uptake	B509
0153-7 - Thyroid Stimulating Hormone (TSH)	B523
0154-5 - Carbamazepine	B526
0155-2 - Triglycerides	B535
0156-0 - ABO and Rh Blood Typing.....	B9
0157-8 - Uric Acid, Serum.....	B542

0158-6 - Uric Acid, Urine, 24 Hr.....	B543
0159-4 - Urinalysis, Routine.....	B540
0160-2 - Vitamin B12.....	B368
0165-1 - Chem 29 Profile.....	B136
0168-5 - Androstenedione, LC/MS/MS.....	B0
0171-9 - Glucose, 1hr PP.....	B245
0178-4 - Dehydroepiandrosterone (DHEA).....	B166
0183-4 - Erythropoietin, Serum.....	B183
0184-2 - Anti-Epidermal Antibody, S.....	B27
0185-9 - Alkaline Phosphatase, Serum.....	B19
0190-9 - PSA Total.....	B462
0192-5 - Transferrin.....	B314
0196-6 - Calcitonin, Serum.....	B527
0197-4 - Hepatitis B Panel.....	B282
0202-2 - Strep Group B (Amplified Probe).....	B494
0204-8 - Pinworm Tape Test.....	B420
0206-3 - Hepatitis B Core Antibody, IgM.....	B257
0207-1 - Fluphenazine, Serum/Plasma.....	B0
0211-3 - Immunoglobulin A, (IgA) Serum.....	B308
0214-7 - Catecholamines, Fractionation, Plasma.....	B223
0216-2 - Hemoglobin Fractionation, Capillary Electrophoresis.....	B481
0222-0 - Diabetic Screen.....	B167
0223-8 - Glucose, Non-Fasting.....	B455
0225-3 - Culture Stool (3rd Specimen).....	B126
0226-1 - Culture, Stool (2nd Specimen).....	B126
0227-9 - Ova and Parasite (O&P), Stool W/ Trichrome (3rd specimen).....	B399
0228-7 - Microalbumin, Urine, Random.....	B369
0232-9 - Selenium Serum/Plasma.....	B482
0234-5 - EBV Capsid Ab, IgG.....	B197
0248-5 - Nasal Smear (for Eosinophils).....	B196
0250-1 - Iron + TIBC.....	B307
0257-6 - Chlamydia + Gonorrhea, Urine, Aptima.....	B124
0259-2 - Clostridium difficile Toxin/GDH antigen with Reflex to PCR.....	B94
0262-6 - Beta-2-Microglobulin, Serum.....	B55
0263-4 - Parietal Cell Antibody.....	B27
0271-7 - T3 Free, Serum.....	B510
0274-1 - Maternal Serum Alpha-Fetoprotein Screen (MSAFP) - ONTD Only.....	B12
0277-4 - Chromium, Plasma/Serum.....	B141

0285-7 - Zinc, Serum.....	B568
0286-5 - Vitamin D, 25-Hydroxy, Serum.....	B100
0287-3 - Vitamin B12/Folate	B549
0288-1 - Vitamin D, 1, 25-Dihydroxy, Serum	B130
0289-9 - Phenobarbital, Serum	B80
0291-5 - Vitamin C.....	B53
0293-1 - Vitamin B1, Total Thiamine	B550
0296-4 - Vitamin B6.....	B400
0299-8 - Culture, Genital	B228
0301-2 - Potassium, Urine 24 Hr.....	B318
0302-0 - Aldosterone, Serum.....	B15
0305-3 - Phlebotomy Charge.....	B419
0308-7 - Sodium, Urine, 24 Hr.....	B383
0309-5 - Aldolase,serum	B0
0315-2 - SCL70 (Scleroderma Ab.), Serum.....	B530
0316-0 - Mumps Virus Antibody, IgG.....	B358
0318-6 - Extractable Nuclear Antigen (ENA) Antibodies (Sm/RnP).....	B500
0319-4 - Occult Blood,Stool (2nd Spec.).....	B89
0323-6 - Occult Blood,Stool (3rd Spec.).....	B89
0325-1 - Arthritis Profile	B29
0327-7 - HCG, Quantitative, Serum	B570
0330-1 - Hepatitis Screen.....	B283
0334-3 - Liaison CIA (Antibody to T. pallidum)	B534
0335-0 - Progesterone, Serum	B426
0337-6 - Vitamin A (Retinol), Serum.....	B477
0341-8 - Culture, Blood	B88
0342-6 - Luteinizing Hormone (LH).....	B335
0348-3 - Cryptosporidium EIA	B157
0349-1 - Methylmalonic Acid, Serum.....	B368
0355-8 - Ova and Parasite (O&P), Stool W/ Trichrome (2nd specimen)	B398
0358-2 - Urea, Urine 24hr	B539
0359-0 - Calcium, Urine, 24 Hr.	B130
0364-0 - Double-Stranded DNA (dsDNA) Antibody	B49
0366-5 - Sickle Cell Screen	B175
0368-1 - Glucose, 1hr, PP (Preg) w/Glucola	B244
0369-9 - Glucose 2hr PP	B245
0374-9 - Cortisol, Serum AM and PM	B153
0375-6 - HLA-B27.....	B292

0377-2 - Ova and Parasite (O&P), Stool W/ Trichrome (1st specimen).....	B443
0378-0 - Aluminum, Serum.....	B42
0379-8 - Testosterone, Serum.....	B513
0380-6 - Theophylline, Serum	B46
0381-4 - Eosinophil Count, Absolute, Blood.....	B36
0382-2 - Culture, Sputum W/ Gram Stain	B502
0391-3 - Diabetes Profile.....	B167
0398-8 - Lead, Blood (Child).....	B445
0400-2 - CMV Ab. (IGG)	B164
0402-8 - Spore Strips	B492
0403-6 - Hyperthyroid Panel	B280
0404-4 - Protein Electrophoresis, Urine, Random.....	B462
0406-9 - DHEA Sulfate	B166
0407-7 - Zinc Protoporphyrin (ZPP),Whole Blood, Child.....	B567
0410-1 - Encainide + Metabolite Level	B182
0411-9 - Phosphates, Urine 24 Hr.....	B419
0412-7 - Anti-Microsomal Ab(Thyroid Peroxidase)TPO	B528
0413-5 - Immunotyping,Serum	B310
0419-2 - Smooth Muscle Antibody (ASM)	B501
0421-8 - Mercury, Urine 24HR.....	B0
0427-5 - Arsenic, Urine, 24 Hours.....	B52
0431-7 - Cardiolipin Antibodies, IgG, IgM, IgA.....	B10
0433-3 - Opiates, Urine, Screen (300).....	B285
0439-0 - Chloride, Urine 24HR.....	B146
0447-3 - Vitamin E, Serum.....	B44
0451-5 - Apolipoprotein A-1 (Apo A-1)	B50
0455-6 - Mycoplasma Pneumoniae IgG/IgM Abs.....	B0
0457-2 - Apolipoprotein B (Apo B)	B51
0461-4 - CMV IGM	B164
0462-2 - Culture, Ear (Left).....	B192
0464-8 - Primidone (incl. Phenobarb)	B376
0479-6 - Copper, Plasma/Serum.....	B158
0485-3 - Fructosamine, Serum	B246
0499-4 - Viscosity, Serum	B0
0507-4 - Adrenocorticotropic Hormone (ACTH)	B10
0508-2 - Gastrin, Serum	B240
0509-0 - Amitriptyline and Metabolite, Serum/Plasma	B0
0510-8 - Growth Hormone.....	B238

0512-4 - Beta Carotene, Total.....	B58
0513-2 - Catecholamines, Fractionation, Urine	B0
0514-0 - Haptoglobin, Serum	B256
0516-5 - Estradiol, Serum.....	B4
0520-7 - Immunoglobulins, Serum	B302
0521-5 - Lipase, Serum.....	B534
0525-6 - Procainamide/NAPA.....	B426
0526-4 - Salicylates, Serum	B32
0530-6 - Xylose Absorption, Urine	B563
0532-2 - Complement C3, Serum.....	B97
0533-0 - Complement C4, Serum.....	B98
0534-8 - Complement, Total (CH50), Serum.....	B282
0535-5 - CA 19-9, Serum.....	B98
0536-3 - CA 125, Serum.....	B99
0538-9 - Hepatitis A Antibody, IgM	B261
0539-7 - Hepatitis Be Antigen	B263
0540-5 - Hepatitis Be Antibody	B568
0541-3 - Immunoglobulin G (IgG), Serum	B309
0542-1 - Valproic Acid	B172
0545-4 - Cholinesterase, Serum.....	B107
0549-6 - Immunoglobulin M, (IgM) Serum.....	B309
0551-2 - Lipoprotein Electrophoresis w/Lipids	B336
0552-0 - Toxoplasma Ab. (IgG), EIA.....	B518
0556-1 - Culture, Body Fluid W/Gram Stain	B90
0558-7 - CPK Isoenzymes	B108
0559-5 - Chromosome Analysis/Karyotype, Blood.....	B142
0561-1 - Allergen Grey Alder (t2), IgE	B252
0565-2 - Culture, Eye (Left).....	B391
0567-8 - Rubeola/Measles IgG Antibody.....	B366
0568-6 - Lyme Disease Antibody, Serum	B338
0570-2 - Herpes II Ab., IgG	B286
0572-8 - EBV Panel (Capsid/EBNA G+M)	B196
0575-1 - Vancomycin Peak, Serum	B547
0577-7 - Thyroglobulin, Serum.....	B514
0578-5 - C-Peptide, Serum	B94
0579-3 - Crystal Identification, Body Fluid.....	B158
0580-1 - EBV Capsid Ab, IgM.....	B198
0582-7 - EBV, Early Antigen Ab	B197

0583-5 - EBV, Nuclear Antigen Ab, IgG	B197
0584-3 - Acute Phase Reactant Profile	B11
0587-6 - Centromere Autoantibodies	B105
0588-4 - Lead, Blood + E.P.	B328
0593-4 - Scleroderma Profile	B482
0594-2 - Antinuclear Antibody Evaluation (ANA)	B22
0596-7 - Pediatric Arthritis Profile	B416
0597-5 - Varicella Zoster Ab., IgG	B287
0598-3 - PTH, Intact (use 1689 if with Calcium)	B438
0608-0 - Gliadin Antibodies, IgG, IgA	B172
0613-0 - Trimipramine, Serum	B521
0622-1 - Acetaminophen, Serum	B537
0628-8 - Drug Abuse Screen, Urine. (8 drug).....	B170
0633-8 - Culture, Eye (Right)	B392
0634-6 - Culture, Ear (Right).....	B192
0636-1 - Neutrophil Cytoplasmic AB/ANCA	B401
0648-6 - Allergen June/Ken.Blue Grass (g8), IgE.....	B317
0649-4 - Allergen Orchard Grass (g3), IgE.....	B398
0650-2 - Allergen Eng. Plantain (w9), IgE.....	B195
0652-8 - Allergen Oak Tree (t7), IgE	B390
0654-4 - Serodia TPPA (Antibody to T. pallidum).....	B533
0658-5 - Sex-Hormone Binding Globulin	B485
0659-3 - Allergen Elm Tree (t8), IgE.....	B195
0660-1 - Allergen Alternaria Tenuis Mold (m6), IgE	B45
0661-9 - Allergen C. Herbarum (m2), IgE	B146
0662-7 - Allergen House Dust H/S (h2), IgE	B291
0664-3 - Allergen Dog Dander (e5), IgE	B175
0665-0 - Allergen Cat Dander (e1), IgE	B133
0668-4 - Sjogren's Abs (SS-A, SS-B)	B339
0671-8 - Allergen Egg White (f1), IgE.....	B193
0672-6 - Allergen Milk (f2), IgE	B570
0673-4 - Allergen Wheat (f4), IgE.....	B559
0674-2 - Allergen Peanut (f13),IgE.....	B446
0680-9 - Allergen Soybean (f14), IgE	B501
0681-7 - Allergen Meadow Fescue (g4), IgE	B365
0682-5 - Allergen Perennial Rye Grass (g5),IgE	B451
0683-3 - Allergen Timothy Grass (g6), IgE	B529
0684-1 - Allergen Maple (Box Elder) (t1), IgE	B364

0685-8 - Allergen Common Birch (t3), IgE.....	B151
0686-6 - Allergen Walnut Tree (t10), IgE.....	B557
0691-6 - Allergen Horse Dander (e3), IgE.....	B290
0692-4 - Allergen Cow Dander (e4), IgE.....	B156
0694-0 - A. Fumigatus Mold (m3), IgE.....	B8
0696-5 - Allergen Shrimp (f24),IgE	B499
0697-3 - Allergen Blue Mussel (f37), IgE.....	B90
0698-1 - Allergen Tuna (f40), IgE	B536
0714-6 - Allergen Salmon (f41), IgE.....	B496
0715-3 - Allergen D. Pteronyssinus (d1), IgE.....	B165
0716-1 - Allergen Cockroach (i6), IgE.....	B148
0717-9 - Allergen Common Ragweed (w1) ,IgE	B151
0718-7 - Allergen Mugwort (w6),IgE.....	B373
0730-2 - Gentamicin, Peak.....	B228
0732-8 - Allergen Russian Thistle (w11), IgE	B480
0740-1 - Allergen Profile, Northeast.....	B387
0741-9 - Grasses Screen (6)	B251
0742-7 - Trees Screen (6).....	B533
0743-5 - Foods I Screen (fish/6).....	B221
0745-0 - Foods II Screen (10)	B222
0746-8 - Weeds Screen (7)	B558
0768-2 - Chem 7 Panel	B106
0769-0 - Chlamydia and Chlamydophila AB Panel 3 (IgG, IgA, IgM).....	B0
0782-3 - Chlamydia Direct Antigen	B377
0789-8 - Myoglobin, Serum.....	B361
0796-3 - Rf (Rheumatoid Arthritis) Titer	B465
0808-6 - Lead, Urine, 24 Hr.....	B0
0813-6 - Allergen Yeast Bakers, IgE	B43
0814-4 - Oat (f7), IgE	B391
0815-1 - Allergen Almond Nut (f20), IgE.....	B44
0816-9 - Allergen Barley (f6), IgE	B80
0817-7 - Allergen Chicken Meat (f83), IgE	B137
0819-3 - Culture, Environmental.....	B161
0821-9 - Allergen Corn (maize)(f8), IgE.....	B155
0822-7 - Allergen Rye (food) (f5), IgE	B480
0823-5 - CA 27.29, Serum.....	B99
0825-0 - Triple Screen AFP3 (AFP, UE3, HCG).....	B45
0826-8 - Uric Acid, Synovial.....	B540

0830-0 - ALDOSTERONE, URINE, 24HR.	B0
0835-9 - Allergen Orange (f33), IgE	B397
0836-7 - Allergen Pork (f26), IgE	B458
0837-5 - Allergen Tomato (f25), IgE	B530
0840-9 - Allergen Codfish (f3), IgE	B149
0841-7 - Allergen Lobster (f80), IgE.....	B343
0842-5 - Allergen Kiwi (f84), IgE	B326
0843-3 - Allergen Rice (f9), IgE	B478
0844-1 - Allergen Potato (f35), IgE	B458
0846-6 - Allergen Duck Feathers (e86), IgE	B175
0847-4 - Allergen Goose Feathers (e70), IgE.....	B250
0849-0 - Allergen Egg Yolk (f75), IgE.....	B193
0850-8 - Allergen Cheddar Cheese (f81), IgE	B135
0851-6 - Smith Antibody.....	B28
0855-7 - CH24/HDL,CBC/D/PLT	B136
0856-5 - Culture, Dialysis Water	B173
0857-3 - Culture, Dialysate.....	B173
0861-5 - Allergen Beef (f27), IgE	B82
0862-3 - Lamb/Mutton (f88) , IgE	B340
0863-1 - Allergen Banana (f92), IgE.....	B79
0864-9 - Allergen Onion (f48), IgE	B393
0865-6 - Hexosaminidase Leukocytes (Tay-Sachs Enzyme)	B525
0867-2 - Thyroid Comprehensive Panel.....	B515
0873-0 - Folate, RBC.....	B214
0877-1 - Toxoplasma Ab. (IgM).....	B519
0883-9 - Allergen Honey Bee (i1), IgE	B289
0884-7 - Allergen Yellow Jacket (i3),IgE.....	B564
0885-4 - Allergen Yellow Hornet (i5), IgE.....	B564
0886-2 - Allergen Paper Wasp (i4), IgE	B442
0891-2 - Allergen Cotton-Wood (t14), IgE	B569
0892-0 - Allergen C.Albicans (m5), IgE.....	B131
0894-6 - Allergen Sheep Sorrel (w18), IgE.....	B499
0895-3 - Apolipoprotein A-1 and B (Apo A1/B).....	B50
0900-1 - Cortisol, Random	B153
0901-9 - Molds Screen (6)	B371
0902-7 - Animals Screen(6)	B48
0904-3 - Allergen Fire Ant (i70), IgE.....	B219
0906-8 - Vitamin B12 Binding Capacity	B549

0907-6 - Allergen Whiteface Hornet (i2), IgE	B561
0910-0 - House Dust Screen (5).....	B291
0920-9 - Amphetamines, Urine, Screen (1000)	B83
0921-7 - Barbiturates, Urine, Screen (200)	B46
0922-5 - Cannabinoids (Marijuana), Urine, Screen (100).....	B364
0923-3 - Cocaine, Urine, Screen (300)	B114
0924-1 - Benzodiazepines, Urine, Screen (200)	B552
0925-8 - Methadone, Urine, Screen (300).....	B349
0926-6 - PCP (6-Phencyclidine), Urine, Screen (25)	B48
0930-8 - Allergen Pear (f94), IgE	B449
0931-6 - Allergen Strawberry (f44), IgE.....	B504
0934-0 - Allergen Carrots (f31), IgE	B132
0935-7 - Allergen Gluten (f79), IgE	B246
0936-5 - Allergen Sesame Seed (f10), IgE.....	B498
0939-9 - Mitochondrial Antibody (AMA) with titer.....	B20
0946-4 - Parainfluenza Antibody.....	B409
0948-0 - Allergen Beech Tree, IgE (t5)	B81
0953-0 - H. Pylori Antibody Panel (G,A,M)	B254
0956-3 - PAP/PSA.....	B408
0960-5 - Culture, Mouth	B125
0963-9 - Propoxyphene, Urine, Screen (300)	B427
0965-4 - Allergen Garlic (f47), IgE.....	B239
0966-2 - Culture, Cerebrospinal Fluid	B134
0973-8 - Rubella Ab, IgG, Quant.....	B480
0978-7 - Plasma Renin Activity, LC/MS/MS	B0
0979-5 - Pediatric Profile II.....	B417
0980-3 - Pediatric Profile I.....	B417
0981-1 - Obstetric/Prenatal II	B389
0982-9 - LH, FSH, Estradiol.....	B335
0984-5 - Hepatitis Management Profile	B283
0985-2 - Alcohol (Ethanol), Blood.....	B199
0988-6 - Osmolality, Urine 24hr.....	B6
0989-4 - Electrolytes, Urine, Random.....	B181
0992-8 - Hepatitis/ Comprehensive.....	B264
0993-6 - Coronary Risk Profile II.....	B117
0994-4 - Coagulation Profile	B113
0995-1 - Alcohol (Ethanol), Urine, Screen (20).....	B200
0996-9 - Immunoglobulin E (IgE), Serum/Plasma	B531

0997-7 - Hydroxyprogesterone(17), Alpha Serum/Plasma	B4
0998-5 - Anemia I Profile.....	B24
0999-3 - Anemia II Profile.....	B24
1003-3 - GC,rRNA,TMA Aptima Swab	B227
1004-1 - Chlamydia, swab, TMA, Aptima	B123
1005-8 - Chlamydia and Gonorrhea, rRNA,TMA Aptima Swab.....	B140
1006-6 - Prealbumin	B532
1007-4 - Allergen Chocolate Cocoa (f93), IgE.....	B140
1010-8 - HIV-1, RNA, Ultra/PCR, Viral Load	B272
1015-7 - GLU.TOL.3HR.(Pregnancy)100GM	B232
1021-5 - Iron/Ferritin Panel	B314
1031-4 - Allergen Poppy Seed (rf224), IgE.....	B457
1043-9 - Culture,ID Only (Any Source).....	B126
1053-8 - Chromosome Analysis, POC, Tissue.....	B142
1072-8 - Insulin-Like Growth Factor I	B293
1074-4 - Liver Panel	B336
1100-7 - Conventional Pap Slide(s)	B155
1102-3 - PAP Smear, 2 Slides, Autopap.....	B127
1133-8 - Allergen Coffee (rf221), IgE	B149
1134-6 - Allergen Apple (f49), IgE	B51
1139-5 - CLOMIPRAMINE	B0
1149-4 - Allergen Turkey Meat (f284), IgE.....	B537
1161-9 - Cytomegalovirus (CMV), DNA, Ultraquant.....	B112
1163-5 - Allergen Chick Pea (f309), IgE	B137
1168-4 - Allergen Crab (f23), IgE.....	B156
1169-2 - Allergen Clam (f207), IgE	B146
1171-8 - Magnesium, Urine, 24Hr.....	B369
1178-3 - Lactose Tolerance Test	B327
1183-3 - Eosinophil Ct, Urine	B196
1184-1 - Xylose, (Baseline) Whole Blood.....	B563
1186-6 - Celiac Disease Antibody Evaluation	B104
1194-0 - HSV I/II DNA PCR, Quant.....	B0
1195-7 - Allergen Pistachio (rf203), IgE	B454
1196-5 - Allergen Sweet Vernal (g1), IgE	B507
1197-3 - Allergen Giant Ragweed (w3), IgE	B243
1198-1 - Allergen Bermuda Grass (g2), IgE	B83
1199-9 - Allergen Redtop (Bent) Grass (g9), IgE	B471
1201-3 - Hcg, Quantitative, Tumor Marker.....	B571

1206-2 - Endomysial Antibody (IgA) Screen with Titer	B49
1219-5 - Allergen Penicillium Chrysogenum/notatum IgG.....	B0
1224-5 - Allergen Lemon (f208), IgE.....	B341
1225-2 - Allergen Grape (f259), IgE	B250
1226-0 - Allergen Peach (f95), IgE	B446
1227-8 - Allergen Melon (Cantaloupe) (f87), IgE	B367
1228-6 - Allergen Green Pea (f12), IgE.....	B252
1229-4 - Allergen Broccoli (f260), IgE.....	B92
1230-2 - Allergen String Bean (rf315), IgE	B252
1231-0 - Allergen Sweet Potato (f54), IgE.....	B506
1235-1 - Tobramycin, Peak	B517
1241-9 - Allergen White Bean (f15), IgE.....	B560
1243-5 - Allergen Green (Bell) Pepper (rf263), IgE.....	B251
1244-3 - Cabbage (f216), IgE	B324
1249-2 - Allergen Cheesemolds (f82), IgE	B136
1250-0 - Pneumococcal Antibody Panel (14 Serotypes)	B456
1251-8 - Olive Tree (t9), IgE	B392
1252-6 - Allergen Willow (t12), IgE	B561
1253-4 - Allergen Silk (k74), IgE	B500
1267-4 - Allergen Casein (f78), IgE	B171
1268-2 - Allergen Betalactoglobulin (f77), IgE	B83
1295-5 - Genital Cult. R/O Yeast/Gard.....	B242
1297-1 - TORCH Panel (IgG/IgM)	B530
1300-3 - Feather Mix (ex71), IgE.....	B217
1301-1 - BREAST SMEAR (slide)	B76
1303-7 - Nipple Discharge Smear	B380
1305-2 - Allergen Sycamore (t11), IgE.....	B507
1306-0 - Allergen White Ash (t15), IgE	B560
1307-8 - Allergen Wormwood (w5), IgE.....	B561
1311-0 - Homocysteine Level	B132
1315-1 - Allergen Mucor Racemosus (m4), IgE	B373
1316-9 - IGF Binding Protein-3	B293
1318-5 - Allergen Coconut (f36), IgE	B148
1325-0 - Aluminum, 24-Hour Urine without creatinine.....	B0
1328-4 - COXSACKIE VIRUS Ab.(GRP. B)	B0
1332-6 - Vitamin A + Beta Carotene	B549
1336-7 - Zinc, RBC	B567
1340-9 - Manual Diff/Platelet only	B363

1349-0 - Allergen Grapefruit (f209), IgE.....	B250
1350-8 - Culture, Blood (spec.#2)	B87
1351-6 - Culture, Blood (spec.#3)	B88
1352-4 - Culture, Neisseria gonorrhoeae (GC), Throat	B248
1353-2 - Culture, Neisseria gonorrhoeae (GC), Anal	B248
1355-7 - Allergen Parsley (f86), IgE	B443
1359-9 - Allergen Mouse Epithelium (e71), IgE	B372
1360-7 - Allergen Mouse U. Protein (e72), IgE.....	B373
1363-1 - Allergen Dandelion (w8), IgE	B171
1364-9 - Allergen Golden Rod (w12), IgE.....	B248
1366-4 - Hypersensitivity Pneumonitis	B279
1371-4 - Allergen Latex (k82), IgE.....	B340
1372-2 - Allergen D. Farinae (Dust Mite) (d2), IgE	B173
1373-0 - Allergen R. Nigricans (m11), IgE	B465
1380-5 - Allergen Celery (f85), IgE	B134
1383-9 - Magnesium, RBC	B369
1389-6 - Allergens Pineapple (f210), IgE.....	B454
1390-4 - Allergen Tea (rf222), IgE	B526
1396-1 - Giardia Antigen (EIA).....	B243
1397-9 - Resp. Syncytial Virus.....	B470
1402-7 - Cytology, Respiratory System	B127
1404-3 - Allergen Hazelnut (f17), IgE	B280
1405-0 - Allergen Brazil Nut (f18), IgE.....	B91
1406-8 - Allergen Alpha Lactalbumin (f76), IgE	B44
1409-2 - Allergen Lettuce (f215), IgE.....	B341
1410-0 - Influenza A/B, RSV BY PCR.....	B213
1412-6 - Allergen Cashew (rf202), IgE.....	B569
1414-2 - Allergen Sunflower Seed (k84), IgE.....	B505
1422-5 - Allergen Black Pepper (f250), IgE.....	B86
1423-3 - Allergen Chili Pepper (rf279), IgE	B138
1424-1 - Allergen Paprika (f218),IgE.....	B443
1426-6 - Allergen Turkey Feathers (e89), IgE.....	B536
1436-5 - Allergen Cotton-Fiber (o1), IgE	B155
1438-1 - Lead, Blood (Adult)	B445
1461-3 - Allergen Basil (rf269), IgE	B81
1462-1 - Allergen Mushroom (rf212), IgE.....	B375
1463-9 - Allergen Anchovy (rf313), IgE	B47
1468-8 - Cyclospora ID, Modified Acid Fast	B127

1476-1 - Protein Total, Pleural Fluid B0

1478-7 - Cocaine, Urine, Screen (100) B113

1479-5 - Allergen Walnuts (f256), IgE..... B557

1480-3 - Allergen Eucalyptus (t18), IgE..... B200

1481-1 - Allergen Marsh Elder (w16), IgE B365

1484-5 - Allergen Pumpkin Seed (rf226), IgE..... B463

1485-2 - Allergen Halibut (rf303), IgE..... B280

1486-0 - Allergen Blackberry (rf211), IgE..... B86

1487-8 - Allergen Blueberry (rf228), IgE B90

1488-6 - Allergen Pumpkin (f225), IgE B463

1497-7 - WBC B555

1498-5 - RBC..... B466

1499-3 - Hemoglobin B281

1500-8 - Urethral Smear..... B543

1576-8 - Cannabinoids (Marijuana), Urine, Screen (20)..... B513

1583-4 - Allergen Cherries (f242), IgE B569

1584-2 - Allergen Plum (f255), IgE..... B456

1589-1 - Vancomycin, Trough, Serum B548

1592-5 - Allergen Pecan Nuts (f201), IgE..... B449

1602-2 - Rickettsia Antibody Panel with Reflex to Titers B0

1605-5 - Allergen Storage Mite (Acarus) (d70), IgE B503

1606-3 - Hepatitis B Surface Antibody, Quant. B262

1608-9 - Hepatitis B Surface Antigen, Quant. B262

1615-4 - Lyme Ab., IgM B344

1620-4 - Allergen Whey (f236), IgE B559

1627-9 - Allergen Rough Careless Pigweed (w14), IgE B479

1629-5 - Allergen Buckwheat (f11), IgE..... B94

1630-3 - Allergen Malt (f90), IgE B362

1631-1 - Allergen House Dust (Greer) (h1), IgE B291

1632-9 - Panel A (Inhalants)..... B439

1639-4 - Lipoprotein-a (Lp-a)..... B343

1644-4 - Immunotyping, Urine B82

1648-5 - Insulin, Fasting B305

1649-3 - Insulin, 1st HR. Specimen B304

1659-2 - Potassium, Urine, Timed Specimen B424

1660-0 - Calcium, Urine, Timed Specimen..... B130

1661-8 - Creatinine Urine, Timed or Random B543

1662-6 - Protein,Urine,Timed or Random B463

1663-4 - Uric Acid, Urine, Timed Spec.....	B540
1664-2 - Allergen Cucumber (f244), IgE.....	B159
1679-0 - BUN, Post Dialysis.....	B76
1687-3 - Post Dialysis Creatinine	B119
1688-1 - Post Dialysis Potassium	B423
1689-9 - PTH, Intact w/ Calcium	B438
1698-0 - Urea Clearance, 24hr Urine.....	B538
1715-2 - Amantadine, Serum/Plasma	B0
1741-8 - Aspergillus Antibodies, Immunodiffusion, Serum	B32
1744-2 - Metanephrines Fractionated, Urine, Random.....	B0
1754-1 - Beta-2-Microglobulin, Urine, Random	B55
1763-2 - Immune Deficiency Panel CD4/CD8 w/CBC	B508
1765-7 - H. Pylori,IgG (For Panel #0953)	B255
1766-5 - H. Pylori,IgA (For Panel #0953).....	B255
1793-9 - Sodium, Urine, Random.....	B384
1794-7 - Chloride, Urine, Random	B145
1795-4 - AAMI Water Analysis	B555
1814-3 - Allergen Brussel Sprouts (rf217), IgE.....	B93
1815-0 - Allergen Cauliflower (rf291), IgE	B133
1817-6 - Aluminum, Water, Dialysis	B20
1823-4 - Allergen Mustard (f89), IgE	B375
1824-2 - Allergen Johnson Grass (g10), IgE.....	B317
1825-9 - Allergen Bahia Grass (g17), IgE	B79
1827-5 - Culture, Seminal Fluid	B498
1829-1 - U.R.R. PANEL.....	B542
1841-6 - Grass Mix (gx1), IgE.....	B251
1843-2 - Rast/Seafood Mix(fx2) IgE +	B497
1846-5 - Procainamide, Ser.(+napa #0525)	B84
1871-3 - Tacrolimus (FK506)	B460
1894-5 - Allergen Helminthosporium (m8), IgE.....	B281
1900-0 - Infant Tot.Bili (tot/dir=#0046)	B377
1901-8 - Infant Dir.Bili (tot/dir=#0046).....	B377
1930-7 - GENTAMICIN, TROUGH	B229
1936-4 - Meprobamate	B367
1937-2 - Allergen Fig (f328), IgE.....	B219
1938-0 - Allergen Pecan/Hickory Tree (t22), IgE.....	B450
1939-8 - Allergen Cocklebur (w13), IgE.....	B148
1941-4 - Allergen Watermelon (f329), IgE	B558

1942-2 - Allergen Papaya (f293), IgE	B442
1945-5 - Allergen Avocado (f96), IgE	B54
1948-9 - Allergen Sheep Epithelium (e81), IgE	B499
1962-0 - PAP, Liquid-Based	B406
1976-0 - Glucose.....	B245
2063-6 - PDF Creatinine (Timed Test code)	B120
2073-5 - Anemia III Profile.....	B24
2075-0 - Hepatitis B Recovery	B264
2083-4 - Urea, Urine, Timed	B539
2087-5 - Bioavailable Testosterone	B556
2088-3 - PSA Free + Total	B435
2100-6 - Cytology, Urine With Interpretation	B128
2103-0 - Allergen Isocyanate TD1 (k75), IgE	B315
2104-8 - Allergen Isocyanate MD1 (k76), IgE	B314
2106-3 - Allergen Formaldehyde (k80), IgE.....	B222
2115-4 - CF/Fragile X/SMA.....	B105
2122-0 - Amniotic fluid AFP.....	B12
2126-1 - Hepatitis A Ab.(non-reflex)	B261
2130-3 - CA 15-3, Serum.....	B98
2132-9 - Tobramycin, Trough	B517
2136-0 - Chromosome Analysis,Amnio.FLD(IG)	B142
2161-8 - Hepatitis C Genotype	B263
2163-4 - Troponin I.....	B522
2167-5 - Allergen Penicillium Notatum Mold (m1), IgE	B450
2177-4 - Weed Mix (wx1), IgE.....	B558
2178-2 - Tree Mix (tx1), IgE	B533
2179-0 - Allergen Mosquito (i71), IgE.....	B372
2180-8 - Animal Mix (ex1), IgE	B48
2181-6 - Mold Mix (mx1), IgE	B370
2182-4 - Dust Mix (hx2), IgE	B176
2183-2 - Allergen Nut Mix (fx1), IgE	B387
2185-7 - Cereal Mix (fx3), IgE	B134
2186-5 - Pediatric Mix (fx5), IgE.....	B450
2187-3 - Northern Allerpanel (3 mix).....	B382
2188-1 - Perennial Allerpanel (3 mix)	B418
2189-9 - Food Allerpanel (4 mix)	B214
2192-3 - Myeloperoxidase Ab.....	B354
2194-9 - LDL Direct	B328

2237-6 - PT + PTT	B436
2256-6 - Post Dialysis Albumin	B14
2263-2 - Allergen Goat Milk (rf300), IgE.....	B247
2265-7 - Allergen Eggplant (rf262), IgE	B194
2276-4 - Citric Acid, Urine with Creatinine.....	B0
2297-0 - Lyme Disease AB RFLX IGM TO WB G/M	B337
2302-8 - Hepatitis Comprehensive II	B143
2308-5 - Allergen Cardamon (rf267), IgE.....	B132
2317-6 - Estrogen Panel	B189
2318-4 - Gonadotropin I.....	B235
2319-2 - Pituitary I	B420
2321-8 - Hypercoagulation	B497
2322-6 - Thyroiditis	B516
2326-7 - Ova and Parasite (O&P) Panel, Stool W/ Trichrome(3 specimens)	B388
2344-0 - Allergen Apricot (f237), IgE	B51
2345-7 - Allergen Cinnamon (rf220), IgE	B145
2346-5 - Allergen Ginger (rf270), IgE.....	B244
2347-3 - Allergen Vanilla (rf234), IgE.....	B552
2350-7 - Allergen Clove (rf268), IgE.....	B147
2351-5 - Allergen Oregano (rf283), IgE.....	B398
2352-3 - Allergen Nutmeg (rf282), IgE	B387
2356-4 - Allergen Curry (rf281), IgE.....	B161
2361-4 - Allergen Coriander (rf317), IgE.....	B144
2376-2 - Toxoplasma IgG W/Reflex.....	B518
2377-0 - CMV Ab.,IGG *POS. Reflex To IGM	B111
2379-6 - HCG Quant.,Progesterone.....	B258
2382-0 - Lyme Ab.,Total (No Reflex)	B337
2383-8 - Fragile X.....	B215
2392-9 - Atenolol Level	B33
2394-5 - Allergen Brome Grass (G11), IgE.....	B93
2395-2 - Allergen Common Reed (G7) IgE+	B151
2406-7 - Allergen Cultivated Rye (G12), IgE	B159
2407-5 - Allergen False Ragweed (w4), IgE.....	B216
2420-8 - Allergen Chestnut (t206), IgE.....	B18
2439-8 - Mercury, Urine, Random	B0
2440-6 - Zinc, Urine, Random.....	B568
2441-4 - H. Pylori Antigen, Stool	B254
2446-3 - S.PNEUMONIAE IGG AB.PAN.....	B0

2469-5 - Allergen Profile, Insects (6)	B553
2475-2 - LAL (Limulus Amebocyte Lysate)	B327
2476-0 - Allergen Curvularia Lunata (m16), IgE	B162
2477-8 - Allergen Phoma Betae (m13), IgE	B453
2509-8 - Arsenic, Urine, Random	B0
2523-9 - Mycoplasma/Ureaplasma Culture.....	B361
2534-6 - Ova and Parasite (O&P), Urine.....	B388
2543-7 - Allergen Scallop (f338), IgE	B497
2555-1 - Basic Metabolic Profile (BMP)	B70
2565-0 - Measles,Mumps,Rubella.....	B353
2570-0 - Allergen Chicken Feathers (e85), IgE.....	B137
2596-5 - Drug Confirmation, Opiates, Urine (GC/MS).....	B390
2603-9 - Urine Culture, Catheterized	B133
2642-7 - Allergen Melaleuca Tree (t212), IgE.....	B366
2663-3 - Proteinase-3 AutoAbs	B424
2670-8 - Isospora Belli ID, Modified Acid Fast.....	B308
2673-2 - Microsporidium ID, Stool.....	B350
2675-7 - Beta-2-Glycoprotein I (G,A,M).....	B58
2696-3 - Hereditary Hemochromatosis DNA Mutation Analysis (NY only).....	B269
2721-9 - Toxoplasma Antibodies, (IgG, IgM).....	B518
2765-6 - Prenatal + HIV.....	B425
2797-9 - Food Panel II	B214
2861-3 - Vitamin Panel	B550
2922-3 - Creatine Kinase-MB Fraction(CKMB).....	B109
2948-8 - THYROID AB.....	B0
3001-5 - Tissue Path.(POC).....	B516
3110-4 - Glucose, Post Dialysis.....	B232
3116-1 - Estriol, Unconjugated.....	B198
3121-1 - Allergen Oyster (f290), IgE	B400
3122-9 - Allergen Kidney Bean (rf287), IgE	B325
3134-4 - Phosphorus, Urine, Random.....	B419
3135-1 - Magnesium, Urine, Random	B345
3138-5 - Hepatitis B Core W/Reflex To IgM	B262
3147-6 - Nefazodone Level	B378
3151-8 - Rubella Ab, IgG, Qual, Immune Status.....	B470
3157-5 - Glucose, 2HR. PP W/Glucola.....	B234
3158-3 - Alpha-Fetoprotein (AFP 4) includes HCG, UE3, Inhibin-A.....	B464
3171-6 - Microalbumin, Urine, 24Hr.....	B350

3174-0 - Ecstasy (MDMA), Urine, Screen, Qual. (300).....	B347
3188-0 - Decalcification (Nails/Bone)	B166
3221-9 - Allergen Cultivated Wheat (g15), IgE	B159
3233-4 - Benzodiazepines, Urine, Screen (100)	B58
3239-1 - Methaqualone, Urine, Screen (300)	B464
3245-8 - PBG Deaminase, WB	B410
3247-4 - Mycoplasma Culture, Blood	B360
3267-2 - Ashkenazi Jewish DNA Carrier	B31
3275-5 - Legionella Pneumophila Antibody, IgM by IFA	B0
3283-9 - Acute Hepatitis Panel	B11
3284-7 - Renal Function (Medicare).....	B467
3320-9 - C-Reactive Protein (hs-CRP), High Sensitivity	B95
3325-8 - Vitamin B2.....	B478
3329-0 - Protein/Creatinine w/Ratio	B434
3371-2 - Microsatellite Instability (MSI) by PCR.....	B356
3374-6 - Fecal Lactoferrin	B208
3376-1 - Hepatitis C RNA, Quant., PCR.....	B259
3389-4 - Hepatitis B Virus, DNA, Quantitative, RT-PCR	B282
3399-3 - Allergen White Pine (t16), IgE	B560
3400-9 - Culture, Genital (Beta strep. Group B).....	B243
3412-4 - HPV High/Low Risk by ISH	B274
3420-7 - Hereditary Hemochromatosis DNA Mutation Analysis (Non-NY).....	B269
3422-3 - Hepatic Function Panel	B264
3423-1 - LAMBDA LT.CHAIN,QUANT.24HR.URN	B0
3427-2 - Comprehensive Metabolic Panel.....	B115
3442-1 - Mycoplasma pneumoniae IgG/IgM Antibody, IFA (CSF) (NON-NY).....	B361
3456-1 - Allergen Mango (f91), IgE.....	B363
3465-2 - Keppra Level (Levetiracetam).....	B319
3472-8 - Transglutaminase, IgA Abs.	B520
3473-6 - Transglutaminase IGG/A,AB, Evaluation	B524
3476-9 - Citalopram, Serum/Plasma	B0
3480-1 - Allergen Passion Fruit (rf294), IgE.....	B445
3481-9 - Sirolimus, Blood	B471
3508-9 - Allergen Red Raspberry (rf343), IgE.....	B466
3514-7 - Beta-Galactosidase, Leukocytes	B59
3520-4 - Allergen Mandarin (rf302), IgE	B362
3521-2 - Allergen Lentil (f235), IgE	B341
3522-0 - Allergen Lime (RF306), IgE	B43

3553-5 - B-Type Natriuretic Peptide (BNP)	B70
3556-8 - Allergen Horsefly (i204), IgE.....	B290
3566-7 - Allergen Guinea Pig (e6), IgE	B253
3573-3 - Anthrax Culture Screen.....	B25
3576-6 - Fetal Fibronectin	B209
3582-4 - Prenatal Rapid Aneuploidy FISH	B304
3586-5 - LDL SUBFRACTIONS.....	B0
3605-3 - Allergen Dill (f277), IgE	B174
3617-8 - Absolute Neutrophil Count	B10
3618-6 - HIV-1 PCR QUANT, HIV GENOTYPE BY SEQUENCING	B0
3638-4 - Allergen Pignolia (Pine) Nuts (t213), IgE.....	B454
3642-6 - Allergen Bamboo Shoot (f51), IgE.....	B79
3643-4 - Allergen Trout (f204), IgE	B536
3644-2 - Allergen Macadamia Nuts (RF345), IgE	B362
3704-4 - Allergen Poplar White IgE	B457
3707-7 - Anti-Mullerian Hormone.....	B27
3709-3 - Mesoridazine	B367
3745-7 - Chlamydia Pneumoniae Antibodies (IgG, IgA, IgM).....	B0
3747-3 - CCP (CYCLIC CITRULLINATED PEPTIDE) IGG, SERUM	B102
3767-1 - Pyruvic Acid (Pyruvate), Blood	B464
3773-9 - Culture, Anal (Beta strep. Group B).....	B160
3788-7 - Enterovirus RNA, Qualitative.....	B183
3792-9 - Chloramphenicol	B107
3796-0 - Allergen Cranberry (rf341), IgE.....	B157
3802-6 - Allergen Western Ragweed (w2), IgE.....	B559
3803-4 - Allergen Beetroot (f319), IgE	B82
3813-3 - Allergen Glycyphagus Domesticus (d73), IgE.....	B247
3822-4 - Gliadin IgG, ABS.....	B232
3825-7 - dsDNA Antibody IgG, Crithidia IFA	B570
3828-1 - Allergen StemphyliumBotryosum(m10),IgE.....	B493
3829-9 - Allergen Trichoderma Viride (m15), IgE	B502
3835-6 - Allergen White Mulberry (T70), IgE.....	B18
3837-2 - Allergen Australian Pine, IgE (t73)	B53
3851-3 - High Risk STI CT/GC	B123
3852-1 - Chlamydia NAAT	B138
3853-9 - Gonorrhea NAAT.....	B249
3874-5 - Legionella Culture.....	B329
3879-4 - Allergen Marguerite/Ox-eye Daisy (w7), IgE.....	B364

3882-8 - Alpha-2-Macroglobulin.....	B19
3889-3 - Chlamydia Pneumoniae Antibody (IGM)	B0
3891-9 - Modafinil/Armodafinil, Blood.....	B0
3893-5 - Free Kappa + Lambda Light Chain (serum).....	B319
3898-4 - Giardia Antigen EIA (2ND SPECIMEN)	B230
3899-2 - Giardia Antigen EIA (3RD SPECIMEN).....	B230
3904-0 - Beta-2 Glycoprotein 1 Antibody, IgG.....	B26
3905-7 - Beta-2 Glycoprotein 1 Antibody, IgA.....	B25
3906-5 - Beta-2 Glycoprotein 1 Antibody, IgM.....	B26
3911-5 - 6-Monoacetyl Morphine (6-MAM), Urine, Screen (10).....	B353
3935-4 - Mitotane, Serum/Plasma.....	B0
3946-1 - Protein Electrophoresis, CSF.....	B135
3956-0 - Buprenorphine (Suboxone), Urine, Screen (5)	B77
3968-5 - Allergen Ispaghula (Psyllium) (k72), IgE	B315
3971-9 - Cannabinoids (Marijuana), Urine, Screen (50).....	B513
3997-4 - Interferon-ALPHA, EIA.....	B305
4000-6 - Biopsy, general-@ (AREA).....	B61
4002-2 - Tissue Pathology, Gross Only.....	B516
4003-0 - Biopsy Gen 3@AREA.....	B62
4004-8 - Biopsy Gen 4-@ (AREA)	B62
4005-5 - Biopsy Gen 5-@ (AREA)	B63
4006-3 - Biopsy Gen 6-@ (AREA)	B63
4007-1 - Biopsy Gen 2-@ (AREA)	B62
4009-7 - Nasal Smear: Complete WBC w/ Diff	B385
4030-3 - BUN + Creatinine with Ratio.....	B77
4033-7 - Lyme AB West-Blot IgG/IgM	B337
4291-1 - Lupus Antibody Panel (SLE).....	B336
4300-0 - Cytology, Thyroid By FNA.....	B127
4447-9 - PET II (Dialysis Fluid)	B418
4488-3 - VASCULITIS PROFILE.....	B0
4514-6 - AFP/HCG Tumor Marker	B13
4534-4 - Creatinine, Arterial	B120
4535-1 - Creatinine, Systemic.....	B121
4633-4 - Basic Food Allergy Panel.....	B56
4742-3 - Heavy Metals Profile, Blood (3).....	B0
5000-5 - Biopsy Pod -@	B61
5002-1 - Biopsy Pod 2-@	B60
5003-9 - Biopsy Pod 3-@	B61

5009-6 - Liver Biopsy	B67
5011-2 - Congo Red Amyloid Stain	B117
5024-5 - AML - M2 (AML/ETO) (RUNX1/RUNX1T1) by FISH.....	B37
5025-2 - AML - M4Eo by FISH	B38
5026-0 - Lymphoma - Mantle Cell (MCL), BCL1/IGH by FISH.....	B573
5027-8 - Lymphoma - Burkitt, MYC/IGH by FISH	B359
5028-6 - Lymphoma - Diffuse large B-Cell, BCL6 by FISH.....	B169
5030-2 - Automated UroVysion By FISH	B34
5031-0 - Lymphoma - T-Cell, TCR Gamma by PCR.....	B512
5032-8 - Lymphoma - B-Cell, IGH by FISH.....	B293
5034-4 - AML - M2 (AML/ETO) by RT-PCR.....	B38
5038-5 - TEL/AML1: t(12;21) by FISH.....	B512
5042-7 - Glucose, Urine, Random	B234
5066-6 - X/Y Bone Marrow Transplant Monitoring	B72
5093-0 - Vaginitis/Vaginosis, DNA - BD Affirm	B55
5097-1 - Prostate Carcinoma IHC Panel (CK903, p63, P504S).....	B461
5099-7 - Iron, Liver Studies, Quant.....	B307
5101-1 - Hodgkin Lymphoma IHC Panel.....	B381
5102-9 - Cutaneous Immunofluorescence	B0
5103-7 - Reactive Hyperplasia vs. Lymphoma (NHL) IHC Panel	B467
5104-5 - Prostate Cancer (Stain Only)	B460
5108-6 - GP/AE1/AE3 (Stain Only	B12
5111-0 - Pathology Slide Consultation	B409
5112-8 - Cerebrospinal Fluid (CSF) Cytology	B121
5116-9 - Acute Leukemia/MDS IHC Panel.....	B37
5119-3 - GP/Bile Duct Brushing.....	B235
5122-7 - Prostate Cancer/ Lymph Nodes.....	B428
5126-8 - Metastatic Melanoma - Lymph Node.....	B347
5127-6 - ER / PR / HER2 (Tech Only)w/digital imaging.....	B92
5136-7 - IHC Single Antibody With Interpretation	B299
5137-5 - IHC Antibody (2) With Interpretation	B297
5138-3 - IHC Antibody (3) With Interpretation	B298
5139-1 - IHC Antibody (4) With Interpretation	B298
5141-7 - IHC Antibody (5) With Interpretation	B298
5144-1 - IHC 1 Antibody/Stain Only	B296
5145-8 - IHC 2 Antibody/Stain Only	B296
5146-6 - IHC 3 Antibody/Stain Only	B296
5147-4 - IHC 4 Antibody/Stain Only	B297

5148-2 - IHC 5 Antibody/Stain Only	B297
5151-6 - Cathepsin D by IHC.....	B101
5152-4 - Ki-67 by IHC	B325
5153-2 - p53 by IHC w/digital imaging.....	B572
5155-7 - Minimal Residual Disease for CLL/SLL.....	B111
5158-1 - GP/Level III Histo. TC (1 Slide)	B237
5159-9 - GP/Level III Histo (2 Slides)	B236
5160-7 - ER/PR /DNA/ Ki-67 - w/Interp w/digital imaging	B186
5161-5 - ER/PR/HER2 - w/Interpretation w/digital imaging.....	B186
5163-1 - ER/PR/ Ki-67/ HER2 - w/Interp w/digital imaging.....	B187
5164-9 - Level III Histo. (3 Slides).....	B331
5170-6 - Breast Carcinoma/Lymph Nodes/MM Metastasis	B75
5171-4 - HER2 by IHC w/digital imaging.....	B266
5172-2 - Level III Histo. (4 Slides).....	B331
5173-0 - Level III Histo.(5 Slides).....	B332
5174-8 - Prostate Core Biopsy (20 Spec.).....	B429
5175-5 - GIST Profile (IHC).....	B231
5176-3 - Breast Carcinoma, Rule Out Microinvasion, with Interp.....	B101
5177-1 - Small Round Cell Tumor IHC Panel	B488
5179-7 - C-KIT (D816V).....	B507
5181-3 - Adenocarcinoma vs. Mesothelioma IHC Panel.....	B11
5182-1 - Chronic Eosinophilic Leukemia: FIP1L1-PDGFRΑ by FISH	B210
5190-4 - Tumor of Unknown Primary IHC Panel	B524
5191-2 - Spindle Cell Neoplasm IHC Panel	B491
5193-8 - Hepatocellular vs. Cholangiocarcinoma vs. Metastasis IHC Pnl	B265
5194-6 - Germ Cell Tumor IHC Panel	B230
5195-3 - Neuroendocrine Neoplasm IHC Panel.....	B378
5198-7 - ER/PR - w/Interpretation w/digital imaging.....	B185
5199-5 - Bone Marrow Morphology	B70
5200-1 - Bone Marrow Morphology - Core Only	B71
5202-7 - Bone Marrow Analysis - Clot Only	B71
5207-6 - Bone Marrow Slide Consult.....	B71
5211-8 - Bone Marrow Analysis - Smears	B72
5218-3 - Cytopathology Urine-TC (Thinprep)	B128
5221-7 - T-NK CD57 Flow Cytometry Analysis.....	B102
5223-3 - CLL Prognosis: IgVH Mutation Analysis by PCR	B310
5224-1 - SURGICAL PATHOLOGY LEVEL 2.....	B495
5233-2 - Biopsy Level II, PC ONLY	B59

5234-0 - Biopsy Level III, PC ONLY	B59
5235-7 - Biopsy Level IV, PC ONLY	B60
5236-5 - Biopsy Level V, PC ONLY	B60
5237-3 - Biopsy, Level VI Global	B67
5240-7 - PAP W/O Image Reflex HR HPV.....	B405
5242-3 - Fine Needle Aspirate	B213
5243-1 - Cytology, Alcohol-Based FNA	B213
5244-9 - Cytology, Thinprep Fluids	B386
5245-6 - Cytology, Alcohol-Based Fluids	B386
5247-2 - HER-2/neu,FISH (Manual)	B266
5249-8 - UroVysion (Tech Only).....	B545
5250-6 - Oncology Chromosome Analysis- BM/PB	B163
5253-0 - Comprehensive Urine Pathology (Technical Only).....	B545
5254-8 - Cytopathology, Urine (Thinprep)	B541
5256-3 - Comprehensive Slide Consult	B116
5259-7 - Automated HER2 FISH TC Only	B268
5260-5 - AML - M3 (PML/RARA) by FISH	B39
5261-3 - AML - M3 (PML/RARA) by RT-PCR	B38
5262-1 - HER2 by FISH with Interpretation.....	B267
5265-4 - BCR/ABL: t(9;22) by FISH	B453
5270-4 - Lymphoma - Follicular, BCL2/IGH by FISH	B574
5271-2 - Lymphoma - Follicular, BCL2/IGH by PCR	B574
5273-8 - Lymphoma - Non-Hodgkin, BCL1, BCL2, BCL6 (FISH Panel)	B211
5278-7 - Lymphoma - B-Cell, IGH by PCR.....	B309
5280-3 - Chronic Lymphocytic Leuk.(CLL): FISH Pnl(+12, 13q, 17p, 11q).....	B402
5281-1 - MDS FISH Panel.....	B375
5282-9 - Multiple Myeloma FISH Panel	B211
5283-7 - Path Cons On Ref.BLK/SLD GLOB.....	B409
5285-2 - Bone Marrow Smear Interpretation	B236
5286-0 - ISH Single Antibody, W/Interp.....	B236
5289-4 - ISH Single Antibody Stain Only	B54
5308-2 - Imaged IHC (1) W/Interp.	B299
5314-0 - Cytopatho. Urine-PC (Alcohol)	B128
5324-9 - Comprehensive Urine Pathology.....	B545
5394-2 - Chromosome Analysis w/Reflex Microarray.....	B108
5405-6 - HER2 by IHC (Tech Only)w/digital imaging.....	B265
5406-4 - ER/PR/HER2w/digital imaging- Rflex to HER2 by FISH	B186
5407-2 - ER/ PR/ KI-67/ HER2 - Reflex to HER2 by FISH.....	B187

5408-0 - ER/ PR/ DNA/ KI-67/ HER2 - Reflex to HER2 by FISH	B185
5410-6 - ZAP-70 BY IHC	B564
5416-3 - Histology Processing	B271
5425-4 - ER/PR (Tech Only)	B185
5427-0 - AML - M4 + M5 by FISH	B352
5428-8 - HER2 by IHC - reflex +2 to HER2 by FISH	B265
5433-8 - ER / PR / Ki-67/ HER2 (Tech Only)	B184
5434-6 - CK-903 for Prostate Cancer	B109
5500-4 - OneCheck Hematopathology	B392
5564-0 - PNH by Flow Cytometry	B422
5579-8 - Plasma Cell Dyscrasia/Multiple Myeloma IHC Panel	B421
5602-8 - Skin Biopsy, (1 Jar).....	B371
5603-6 - Skin Biopsy, (2 Jars).....	B486
5604-4 - Skin Biopsy, (3 Jars).....	B486
5605-1 - Skin Biopsy, (4 Jars).....	B487
5606-9 - Skin Biopsy, (5 Jars).....	B487
5607-7 - Skin Biopsy, (6 Jars).....	B487
5611-9 - Prostate Core Biopsy (1 Spec.).....	B428
5612-7 - Prostate Core Biopsy (2 Spec.).....	B430
5613-5 - Prostate Core Biopsy (3 Spec.).....	B429
5614-3 - Prostate Core Biopsy (4 Spec.).....	B429
5615-0 - Prostate Core Biopsy (5 Spec.).....	B430
5616-8 - Prostate Core Biopsy (6 Spec.).....	B431
5617-6 - Prostate Core Biopsy (7 Spec.).....	B430
5618-4 - Prostate Core Biopsy (8 Spec.).....	B428
5619-2 - Surgical Pathology,Level IV.....	B495
5622-6 - Prostate Core Biopsy (12 Spec.).....	B431
5630-9 - Esophageal Biospy X1	B189
5631-7 - Biopsy, GI (1 Jar)	B63
5632-5 - Biopsy, GI (2 Jars)	B64
5633-3 - Biopsy, GI (3 Jars)	B64
5634-1 - Biopsy, GI (4 Jars)	B65
5635-8 - Biopsy, GI (5 Jars)	B65
5636-6 - Biopsy, GI (6 Jars)	B65
5640-8 - Gastric Biopsy (1 Specimen).....	B225
5643-2 - Gastric Washings/Brushings	B225
5645-7 - H.Pylori (Stain Ident.).....	B255
5652-3 - Colon Biopsy (2 Specimens)	B114

5660-6 - Cervical Biopsy	B135
5662-2 - Endo-Cervical Biopsy.....	B182
5664-8 - Endometrial Biopsy	B182
5672-1 - Cone Biopsy/LEEP	B116
5695-2 - Level 3 Surgical Pathology	B329
5696-0 - Level 4 Surgical Pathology	B330
5697-8 - Level 5 Surgical Pathology	B330
5704-2 - Activated Protein C Resistance	B432
5705-9 - Activated Protein C Resistance Modified.....	B431
5706-7 - Anti-Annexin V. Antibodies, IgG, IgM	B25
5707-5 - Factor X Chromogenic.....	B205
5714-1 - Antithrombin III Activity (AT III)	B28
5715-8 - Antithrombin III Antigen (ATIII)	B29
5722-4 - Phosphatidyl Ethanolamine	B419
5723-2 - Antithrombin III Anti-Xa Assay	B33
5724-0 - Factor II Activity	B201
5725-7 - Factor V Activity.....	B202
5726-5 - Factor V Mutation (Leiden)	B204
5728-1 - Factor VII Activity.....	B202
5729-9 - Factor VIII Activity.....	B203
5730-7 - DRVVT Complete	B174
5735-6 - Factor X Activity	B203
5736-4 - Factor X Antigen	B205
5737-2 - Factor XI Activity.....	B204
5738-0 - Factor XII Activity.....	B204
5739-8 - Factor XIII Activity.....	B205
5743-0 - D-Dimer	B165
5744-8 - Collagen Binding Assay	B114
5746-3 - Heparin Platelet Antibody	B456
5749-7 - Lupus-Hexagonal Phase	B326
5757-0 - Inhibitor Titer	B303
5761-2 - Heparin Level (ANTI-XA) - LMWH	B263
5769-5 - Factor IX Activity.....	B201
5771-1 - Part.Thromb.Time(PTT)-LA	B439
5775-2 - Plasminogen Quant.....	B421
5776-0 - Platelet Antibody Identification	B422
5784-4 - Platelet Neutralizing Procedure.....	B422
5785-1 - Protein C Activity.....	B462

5786-9 - Protein C Antigen	B432
5787-7 - Protein S Activity	B433
5788-5 - Protein S Antigen, Free	B434
5789-3 - Protein S Antigen Total	B433
5790-1 - Antiphosphatidylserine Antibodies, IgG, IgM, Special Coagulation.....	B28
5795-0 - Prothrombin G20210A Mutation.....	B436
5797-6 - Prothrombin Fragment 1+2	B292
5798-4 - COAG/Prothrombin Time	B435
5799-2 - Reptilase Time	B467
5800-8 - General Inhibitor Screen.....	B227
5805-7 - Ristocetin Cofactor ACT.	B554
5807-3 - Thrombin Time.....	B514
5808-1 - Thrombin/Antithrombin	B514
5811-5 - Activation Markers.....	B36
5815-6 - Von Willebrand Factor, Antigen	B575
5817-2 - Lupus Anticoagulant Screen	B344
5818-0 - Specific Factor Inhibitor	B491
5820-6 - Prothrombin ABS IGG	B435
5825-5 - Von Willebrand Screen- Basic	B550
5830-5 - Mixing Studies-APTT	B351
5834-7 - Imaged IHC, 4 Slides TC Only	B301
5837-0 - Imaged IHC, 3 Slides TC Only	B301
5838-8 - Imaged IHC, 2 Slides TC Only	B300
5839-6 - DIC Assessment Panel Level-1	B168
5840-4 - Mixing Studies-PT	B351
5842-0 - Imaged IHC, 1 Slide TC Only	B300
5843-8 - Imaged IHC, 4 Slides W/Interp.	B302
5844-6 - Imaged IHC, 3 Slides W/Interp.	B301
5845-3 - Imaged IHC, 2 Slides W/Interp.	B300
5848-7 - RET-HE	B476
5850-3 - Bleeding Diatheses Evaluation Profile I.....	B67
5851-1 - Abnormal APTT Assessment	B35
5853-7 - Abnormal PT and APTT Assessment	B35
5858-6 - BCR/ABL Quantification by RT-PCR	B57
5860-2 - Fibrinolytic Evaluation - Standard	B210
5868-5 - Gliadin IGA, ABS.	B231
5869-3 - Transglutaminase, IGG ABS.....	B520
5871-9 - HIT Profile-Level 1	B271

5872-7 - DIC Level 2	B168
5874-3 - AtheroSclerotic Risk (Coag)	B113
5875-0 - Venous Prothrom. Risk Factor	B548
5886-7 - Vulvar Biopsy	B551
5889-1 - Vaginal Biopsy	B547
5917-0 - Tramadol, Urine, Screen, Qual. (200).....	B520
5924-6 - Lorazepam (Ativan)), Urine	B53
5989-9 - Venous Prothrombic Risk	B548
5990-7 - Expanded Anti-Phos/Lupus	B191
5992-3 - DIC Assessment	B168
5993-1 - Hyperactive Coag. Panel	B279
5994-9 - Basic Antiphospholipid/Lupus	B56
5997-2 - Bleeding Diathesis Evaluation Profile II	B68
6108-5 - Enterococcus/VRE Screen-@	B183
6116-8 - High Risk-STI CT/GC/Trich-TP/SP	B139
6119-2 - MRSA Screen	B355
6155-6 - Lp-Pla2 (PLAC)	B343
6167-1 - Lipase, Peritoneal Fluid.....	B335
6169-7 - CEA, Peritoneal Fluid	B103
6180-4 - 5-Flucytosine, S	B8
6182-0 - Occult Blood, Stool, Immunochem (FIT).....	B218
6193-7 - Vitamin C Dilution	B553
6194-5 - Fentanyl, Urine, Screen, Qual. (2)	B208
6195-2 - Pap (Reflex to HPV DNA Genotyping 16,18 when ASC-US).....	B440
6196-0 - Pap (Reflex to HPV DNA Genotyping 16, 18 when ASC-US or >).....	B573
6203-4 - Cortisone, 24HR Urine	B0
6206-7 - Allergen Marjoram (f274), IgE	B365
6218-2 - Amylase,Pleural Fluid	B22
6219-0 - Silicate, Water	B485
6223-2 - Naproxen (Aleve,Anaprox,Naprosyn)	B384
6235-6 - Phenosense Entry (incl.Fuzeon).....	B418
6236-4 - H.Pylori Urea Breath Test	B538
6238-0 - Peritoneal KT/V, Measured.....	B324
6239-8 - Drug Confirmation, Ethanol, Urine (GC/HS).....	B190
6240-6 - Lambda Free Light Chain Serum	B327
6241-4 - Kappa Free Light Chain Serum	B319
6246-3 - Combined First (NT, Total BHCG, PAPP-A).....	B220
6258-8 - Opiates, Urine, Screen (200).....	B390

6261-2 - Warfarin, Pharmacogenomic Assay	B156
6266-1 - Allergen Sardine (Pilchard) (f308), IgE	B496
6267-9 - Allergen Herring (f205), IgE.....	B18
6272-9 - Lyme Blot, IgG W/Bands	B338
6273-7 - Lyme Blot, IgM W/Bands.....	B338
6274-5 - GP/Non-Automated Urovysion BY FISH	B87
6284-4 - Alcohol Biomarkers (ETG), Qual. (500).....	B14
6285-1 - Dihydropyrimidine Dehydrogenase (DPD) Gene Mutation Analysis.....	B170
6286-9 - EDDP (Methadone Metabolite), Urine, Screen, Qual. (300)	B177
6289-3 - Zinc Protoporphyrin (ZPP),Whole Blood, Adult	B567
6293-5 - Tysabri Antibodies	B525
6305-7 - Allergen Cedar Red (t212), IgE.....	B17
6308-1 - Allergen Nettle (w20), IgE.....	B378
6310-7 - Acylglycines, Quant. Urine.....	B39
6334-7 - Allergen Mountain Juniper (T6), IgE.....	B372
6337-0 - Chlamydia Trachomatis by Multiplex PCR	B107
6338-8 - GC By Multiplex PCR.....	B385
6339-6 - Ureaplasma By Multiplex PCR.....	B539
6340-4 - Mycoplasma Genitalium By Multiplex PCR	B360
6341-2 - Mycoplasma Hominis By Multiplex PCR (NON-NY).....	B360
6342-0 - Trichomonas By Multiplex PCR.....	B534
6351-1 - Herpes Simplex I by Multiplex (Non NY).....	B286
6352-9 - Herpes Simplex II by Multiplex (Non NY).....	B286
6353-7 - Cytomegalovirus (CMV) by Multiplex PCR (NON-NY).....	B111
6354-5 - T.Pallidum By Multiplex PCR (NON-NY)	B509
6356-0 - Lymphogranuloma Venereum (LGV) (Non-NY)	B334
6364-4 - HIV-1,RNA Viral Reflex Genotype.....	B272
6368-5 - Chlamydia Trachomatis, Urine, rRNA	B97
6369-3 - N. Gonorrhoea, Urine, rRNA.....	B376
6375-0 - Integrated Part 1 (NT, PAPP-A).....	B312
6376-8 - Integrated Part 2 (AFP, UE3, HCG, INHIBIN-A)	B312
6380-0 - Spinal Muscular Atrophy (SMA) Carrier (SMN1 Dosage)	B488
6431-1 - PAP/CT/GC (Reflex HPV PLUS when ASC-US)	B407
6626-6 - Anal Cytopathology, ThinPrep	B129
6629-0 - PAP (RFLX HPV HIGH ASC-US)	B403
6630-8 - PAP (RFLX HPV HIGH ASC-US or >)	B403
6794-2 - GTT,GEST,3HR.(ACOG,CARP/COUS).....	B239
6800-7 - General Thrombophilia Assessment.....	B228

6801-5 - Preg.+Throm.Assessment.....	B425
6802-3 - Atherosclerotic Risk Expanded	B33
6820-5 - Region 13 Allergy Profile (S. CA Coast).....	B474
6821-3 - Region 14 Allergy Profile (Central CA).....	B474
6822-1 - Region 17 Allergy Profile (North West)	B474
6823-9 - Region 1 Allergy Profile (North East).....	B473
6824-7 - Region 2 Allergy Profile (Mid Atlantic).....	B474
6825-4 - Region 3 Allergy Profile (South East).....	B475
6826-2 - Region 4 Allergy Profile (S.Florida).....	B475
6827-0 - Region 10 Allergy Profile (TX,OK).....	B473
6829-6 - Region 5 Allergy Profile (IN, OH, TN, WV, KY).....	B475
7301-5 - Testosterone Free + Total.....	B512
7736-2 - H. PYLORI IgM.....	B254
8017-6 - Electrolytes (24 hr urine).....	B181
8431-9 - Perimenopausal Hormonal Profile.....	B451
8545-6 - Meperidine, Urine, Screen, Qual. (200).....	B172
8573-8 - Vit. D, 25-Hydroxy Tot.,D2+D3	B554
8601-7 - Sequential Screen Part 2 (AFP, UE3, HCG, Inhibin-A).....	B483
8602-5 - Sequential Screen Part 1 (NT, Total BHCG, PAPP-A)	B483
8622-3 - Combined First Plus (NT, Total BHCG, Inhibin-A, PAPP-A)	B220
8689-2 - HPV High by ISH Stain ONLY	B274
8696-7 - HPV High/Low Risk by ISH (Stain Only)	B275
8714-8 - HPV High Risk by ISH	B274
8860-9 - AML - CEBPA Mutational Analysis.....	B103
8963-1 - Allergen Tilapia (RF414), IgE	B529
8979-7 - Glucose Tolerance 2 Hour	B233
9081-1 - Gentamicin, Peak + Trough	B229
9102-5 - PAP + HPV HI	B404
9149-6 - General Health Panel	B241
9150-4 - Comp. Metabolic/GGTP.....	B116
9235-3 - PAP + HPV HI + CT/GC	B404
9241-1 - PAP + CT/GC.....	B403
9891-3 - Comprehensive Thyroid-Expanded.....	B152
9892-1 - Arthritis Profile - Expanded	B30
9893-9 - Anemia Profile, Expanded	B24
9985-3 - PAP/CT/GC(RFLX HPV PLUS ASC-US or >).....	B442
9990-3 - CALL RESULT TO DOCTOR	B100
9991-1 - Send Copy Of Results To:.....	B483

9993-7 - Fax Result To Doctor	B207
A044-0 - OnkoSight FLT3	B212
A111-7 - GC, Oral Aptima	B226
A112-5 - GC, Anal Aptima	B226
A113-3 - Chlamydia, Oral Aptima.....	B123
A114-1 - Chlamydia, Anal Aptima	B122
A130-7 - Lung - Adenocarcinoma vs. Squamous cell carcinoma.....	B492
A205-7 - Trofile DNA	B522
A228-9 - Level V Histo,TC	B333
A229-7 - Level II Histo,TC.....	B330
A230-5 - Level VI Histo,TC	B334
A231-3 - Level IV Histo, TC	B332
A233-9 - ANA with Reflex to IFA	B23
A240-4 - EGFR by FISH	B179
A241-2 - EML4-ALK by FISH (FDA): Lung.....	B16
A250-3 - Pap and HPV mRNA Genotyping 16, 18/45.....	B441
A251-1 - Pap (Reflex to HPV mRNA Genotyping 16, 18/45 when ASC-US).....	B441
A252-9 - PAP/CT/GC/ HPV mRNA,RFX 16/18/45 GENOTYPE.....	B407
A253-7 - PAP RFLX HPV mRNA ASCUS,RFX 16/18/45 w/CT/GC	B404
A254-5 - PAP Reflex to HPV mRNA when ASCUS or greater,RFX 16,18/45	B405
A255-2 - PAP RFX HPV Plus mRNA ASCUS>,RFX 16/18/45 w/CT/GC.....	B407
A265-1 - Culture, Genital - Beta strep. Group B (Penicillin Allergic).....	B253
A269-3 - PAI 4G/5G Polymorphism.....	B402
A278-4 - ER,PR,KI-67,P53,HER2 FISH, TC	B184
A282-6 - Ashkenazi Jewish DNA Panel (NON NY).....	B14
A292-5 - Cannabinoid, synthetic (K2,Spice), Urine, Screen, group 1, Qual.(20).....	B318
A293-3 - Biopsy, GI (7 JARS).....	B66
A294-1 - Biopsy, GI (8 JARS).....	B66
A295-8 - Biopsy, GI (9 JARS).....	B66
A296-6 - Biopsy, GI (10 JARS).....	B64
A301-4 - HER-2, IHC.....	B266
A304-8 - Interleukin 28B POLYMORPHISM (NON-NY).....	B305
A311-3 - 1P/19Q,D(1;19)DELETION BY FISH	B5
A313-9 - MLH-1 Promo. Methylation.....	B352
A314-7 - P53 by FISH.....	B5
A315-4 - Epstein-Barr Virus By PCR	B184
A321-2 - Allergen Privet Tree (T210), IgE.....	B459
A333-7 - Allergen Hops (RF324), IgE	B290

A334-5 - Serum Integrated Part 2 (AFP, UE3, HCG, Inhibin-A)	B484
A335-2 - Serum Integrated Part 1 (PAPP-A)	B484
A346-9 - KI-67 Imaged By IHC-Tech Only	B320
A427-7 - HER2 FISH Global Gastric+GEJ	B267
A428-5 - HER2 FISH TC-Gastric+GEJ	B284
A473-1 - Alcohol, Methyl	B42
A489-7 - PAP(RFX HPV Plus)	B405
A502-7 - Culture, Fungal (Yeast)	B527
A518-3 - TSH W/ Reflex To Free T4	B523
A531-6 - HPV DNA 16, 18, NON 16/18	B273
A583-7 - Prenatal Targeted Chromosomal Microarray	B0
A587-8 - Karyotype CVS	B141
A591-0 - Prenatal Whole Genome Chromosomal Microarray	B0
A592-8 - Reflex to Prenatal Targeted Chromosomal Microarray from Chromosomes	B0
A595-1 - Oxycodone, Urine, Screen (300)	B400
A612-4 - MET by IHC	B349
A614-0 - Vancomycin (Random)	B547
A615-7 - Tobramycin (Random)	B517
A616-5 - Gentamicin (Random)	B229
A643-9 - Familial Hyperinsulinism	B206
A644-7 - Dihydropyridine Dehydrogenase	B169
A645-4 - Joubert Syndrome	B316
A646-2 - Nemaline Myopathy 2	B385
A647-0 - Usher Syndrome Type 3	B546
A648-8 - Walker Warburg Syndrome	B555
A649-6 - Usher Syndrome Type 1F	B542
A650-4 - Maple Syrup Urine Disease(Type 1A & 1B)	B346
A651-2 - Glycogen Storage, Type IA (GSDIA)	B235
A652-0 - Tay-Sachs Disease Carrier Test	B511
A653-8 - Niemann-Pick Type A+B Carrier Test	B380
A654-6 - Mucopolidosis Type IV	B357
A655-3 - Gaucher Disease Carrier Screening	B240
A656-1 - Fanconi Anemia, Type C	B206
A657-9 - Canavan Disease Test	B100
A658-7 - Bloom Syndrome Carrier Test	B68
A659-5 - Inherigen Plus	B303
A662-9 - Inherigen	B303
A673-6 - Allergen Trichophyton Rubrum (m205), IgE	B535

A697-5 - Familial Dysautonomia	B206
A744-5 - BRL FUNGUS CULTURE+STAIN,MISC	B0
A780-9 - Reflex to GenCerv	B477
A861-7 - Trichomonas Vaginalis, Probe rRNA	B521
A868-2 - -7/7Q- BY FISH-TC ONLY	B2
A869-0 - -5/5Q- BY FISH-TC ONLY	B1
A870-8 - Trisomy 8 By FISH -TC Only	B522
A871-6 - 20Q12 BY FISH-TC ONLY	B6
A873-2 - T(11;14)+11 BY FISH-TC Only	B508
A874-0 - 17P13 BY FISH-TC ONLY	B572
A875-7 - T(4;14) BY FISH-TC Only	B508
A876-5 - 1Q21 BY FISH-TC ONLY	B110
A877-3 - Trisomy 12 by FISH-TC Only	B521
A878-1 - 13Q14.3 BY FISH-TC ONLY	B165
A879-9 - 11Q22.3 BY FISH-TC ONLY	B34
A881-5 - AML - APL PML/RARA t(15;17) - Tech Only	B0
A888-0 - MAG ANTIBODY, IGM	B0
A889-8 - Myelin Associated Glycoprotein (MAG)-SGPG Antibody (IgM)	B0
A890-6 - GANGLIOSIDE GD1B,AB,IGG	B0
A891-4 - GANGLIOSIDE GD1B,AB,IGM	B0
A892-2 - GANGLIOSIDE GD1A,AB,IGG	B0
A893-0 - GANGLIOSIDE GD1A AB,IGM	B0
A894-8 - GANGLIOSIDE ASIALO-GM-1,AB,IGG	B0
A895-5 - GANGLIOSIDE ASIALO-GM1 AB,IGM	B0
A943-3 - Mismatch Repair(MMR) by IHC w/Interp	B355
A944-1 - Mismatch Repair(MMR) by IHC, Tech Only	B356
A951-6 - ER/PRKI-67/HER2/P53 BY IHC,TC	B188
A952-4 - Breast Carcinoma, R/O Microinvasion BY IHC-TC	B75
A954-0 - ER/PRKI-67/HER2/P53 BY IHC with Interp. w/digital imaging	B187
A972-2 - HER2 BY IHC GI/ESO Carcinoma	B268
A985-4 - PTEN BY IHC With Interpretation	B437
B012-6 - LEVEL II-PC SURGICAL PATHOLOGY	B331
B013-4 - LEVEL III-PC SURGICAL PATH	B332
B014-2 - LEVEL IV-PC SURGICAL PATHOLOGY	B333
B015-9 - LEVEL V-PC SURGICAL PATHOLOGY	B333
B016-7 - LEVEL VI-PC SURGICAL PATHOLOGY	B334
B024-1 - IHC 1 Antibody Interp.	B295
B124-9 - Neutrophil Function Panel	B379

B125-6 - Hepatitis C Antibody W/ Reflex RT PCR	B283
B137-1 - Clostridium Difficile Toxins A+B Genes By DNA Amplification.....	B102
B139-7 - Smith-Lemli-Opitz Synd.Prenat (Non-NY)	B489
B148-8 - Abetalipoproteinemia	B9
B149-6 - Tyrosinemia, Type 1	B525
B150-4 - Hermansky-Pudlak Syndrome 3.....	B285
B151-2 - Familial Hypercholesterolemia LDLR.....	B217
B155-3 - Galactosemia Carrier Screen	B225
B158-7 - Reflex to Prenatal Whole Genome Chromosomal Microarray from Chromosomes.....	B0
B164-5 - Wilson Disease Carrier Screen.....	B556
B222-1 - Reflex to Microarray from nl FISH (chromosomes if FISH abnl)	B0
B232-0 - Hold Cells For Further Testing	B273
B235-3 - Allergen Blomia Tropicalis (d201), IgE	B87
B236-1 - ALK BY FISH (TC).....	B16
B269-2 - HER2 IHC Gastric/GEJ, TC Only w/digital imaging	B268
B275-9 - Comprehensive Common Cancer Panel	B154
B331-0 - Pediatric Blood Culture #1	B416
B332-8 - Pediatric Blood Culture #2	B416
B333-6 - Pediatric Blood Culture #3.....	B417
B334-4 - ROS1 BY FISH Global	B469
B335-1 - MET BY FISH-Global.....	B348
B336-9 - MET BY FISH-Tech Only.....	B349
B337-7 - ROS-1 BY FISH-TECH Only	B469
B342-7 - HPV mRNA screen, if pos, Reflex to 16, 18/45.....	B275
B343-5 - Pancreatic Cancer Panel.....	B284
B352-6 - StormPath PAP, Liquid based	B406
B362-5 - BRCA1/2 Sequencing and Del/Dup Analysis.....	B75
B363-3 - Rest of Comprehensive Cancer Panel	B472
B370-8 - FAMILY MEMBER TESTING.....	B0
B394-8 - Renal Cancer Panel	B0
B424-3 - Allergen Sardine, IgG (NON-NY)	B482
B506-7 - D7S486/XCE 7 FOR -7/7q- BY FISH	B2
B507-5 - D5S2064 FOR -5/5q- BY FISH	B1
B508-3 - D20C108 FOR DEL (20q) BY FISH	B6
B510-9 - D13S319 (13q14.3) BY FISH.....	B4
B512-5 - ATM FOR 11q22.3 BY FISH.....	B3
B514-1 - XCE 8 BY FISH	B1
B518-2 - Urinalysis w/Reflex to Culture	B537

B522-4 - Lynch/Colorectal High Risk Panel.....	B0
B534-9 - Chromosome Analysis- Breakage, Fanconi Anemia.....	B108
B542-2 - JCV Ab w/ Index w/ RFX.....	B316
B544-8 - MDM2 Gene Amplification by FISH.....	B346
B545-5 - FISH, IGH/MAFB, t(14;20).....	B294
B556-2 - EWSR1(22Q12) Gene Rearrang. FISH.....	B200
B557-0 - FISH, FGFR1, 8p11-12.....	B211
B565-3 - ALK (EML4/ALK) BY IHC.....	B15
B608-1 - Liver Cytosol (LC-1) Autoantibodies.....	B328
B609-9 - ACTIN (SMOOTH MUSCLE) ANTIBODY (IgG).....	B0
B611-5 - t(1;19) (q23.3: p13.3) FISH.....	B574
B656-0 - Culture, Wound, Superficial.....	B562
B664-4 - PROCESS AND HOLD FOR FISH.....	B426
B665-1 - FOXO1 (FKHR) (13q14) GENE BY FISH.....	B215
B669-3 - Alpha-Globin mut. Fetus.....	B20
B670-1 - Histo Chemical Micro Organism TCx3.....	B271
B671-9 - Histo Chemical Micro Organism TCx2.....	B269
B672-7 - Histo Chemical Micro Organism TC.....	B270
B673-5 - Histo Chemical Micro Organism x3 W/IN.....	B270
B674-3 - Histo Chemical Micro Organism x2 W/IN.....	B270
B675-0 - Histo Chemical Micro Organism W/INTERP.....	B271
B684-2 - ALK-2P23: Anaplastic Lymphoma.....	B17
B685-9 - TAO Brush Endometrial Cytology.....	B511
B688-3 - HIV AG/AB 4th Generation.....	B272
B700-6 - HOLD BLOCK FOR ADDT'L TESTING.....	B273
B735-2 - InheriGenTx.....	B311
B749-3 - OncoGene Dx Custom Panel.....	B0
B756-8 - Special Stain Enzyme Constituents - PC.....	B490
B757-6 - Special Stain Enzyme Constituents - Global.....	B490
B758-4 - Special Stain Enzyme Constituents -TC.....	B490
B760-0 - FISH, IGH/MAF, t(14;16)- TC.....	B294
B761-8 - Beta-Hexosaminidase B Enzyme Analysis (Sandhoff Disease).....	B496
B770-9 - Allergen Mucor Racemosus IgG.....	B358
B784-0 - Drug Confirmation, Ethanol, Blood (GC/HS).....	B199
B796-4 - 15q Del Characterization, FISH.....	B424
B812-9 - HER2 FISH-Reflex Equivocal to HER2 IHC.....	B267
B818-6 - OnkoSight MPN Panel Sequencing.....	B395
B819-4 - OnkoSight Lung Cancer Panel.....	B394

B820-2 - OnkoSight Melanoma Panel	B396
B821-0 - OnkoSight Solid Tumor Panel	B397
B822-8 - OnkoSight Colorectal Cancer Panel	B393
B823-6 - OnkoSight Myeloid Disorder Panel	B396
B832-7 - Accutype Ribavirin	B308
B843-4 - Bone Marrow Morph. Clot Only- Tech Only	B69
B844-2 - Bone Marrow Morph. Smear Only- Tech Only	B69
B845-9 - Bone Marrow Morph. Core Only- Tech Only	B69
B846-7 - Bone Marrow Morph.- Tech Only	B68
B854-1 - PTEN BY FISH	B437
B889-7 - AML-M4 and M5 by Fish- TC ONLY	B352
B891-3 - AML-M2 (AML/ETO) by Fish- TC ONLY	B21
B893-9 - OnkoSight MDS Panel Sequencing	B395
B894-7 - OnkoSight AML Panel Sequencing	B393
B906-9 - Hypersensitivity Pneumonitis Extended Panel	B0
B960-6 - ISH Single Antibody -PC Only	B307
B964-8 - Prostate Cancer Triple Stain By IHC (PC)	B427
B975-4 - Pap + Pap Dependent HPV	B439
B976-2 - Vulvovaginal Candidiasis Panel by PCR w/ reflex to species	B101
B978-8 - ExtraTube(s) Drawn	B191
B984-6 - Rest of common cancer panel after Lynch Syndrome Panel B985	B471
B993-7 - PD-L1 IHC 22C3 pharmDx, Keytruda W/INTERP	B412
B994-5 - PD-1 BY IHC	B410
F021-1 - InheriGen Plus w/o Fragile X	B311
F100-3 - Chronic Fatigue Syndrome	B143
F101-1 - Gastric Distress	B239
F238-1 - Female High Risk Sexual Health Panel	B556
F239-9 - Male High Risk Sexual Health Panel	B348
F245-6 - Sinus Culture W/ Gram Stain	B486
F306-6 - Comprehensive Heart Health Baseline	B115
F307-4 - Comprehensive Heart Health Monitoring	B115
F368-6 - PID/Infertility W/O PAP Expanded	B420
F371-0 - Vaginitis/Vaginosis W/O PAP W/O HPV Expanded	B546
H001-1 - Chlamydia, Gonorrhea, and Trichomonas rRNA,TMA Aptima Urine	B125
H125-8 - OnkoSight Lung Cancer Panel, ALK and ROS1 by FISH	B389
H191-0 - Milk Component Panel IgE	B351
H474-0 - Iodine, Urine Panel	B306
H510-1 - Fibrotest/Actitest	B219

H575-4 - Bacterial Vag. W/RFX TP-NY.	B78
H578-8 - KIT and PDGFRA Sequencing- Non Ny	B321
H579-6 - EGFR and KRAS, if both neg, reflex ALK, if neg reflex ROS1 Non NY	B180
H580-4 - KRAS by NGS, if negative reflex to BRAF; + MSI by PCR, MSI-High, reflex to MLH1	B323
H581-2 - EGFR by NGS, ALK, ROS1, and MET by FISH- Non NY.....	B180
H583-8 - IBD Differentiation Panel	B292
H760-2 - Herpes Virus I/II DNA.....	B277
H852-7 - HPV PRIMARY SCREENING >=25 YRS, WHEN + RFX TO PAP	B275
H870-9 - BV PANEL BY RT-PCR	B77
H906-1 - Estrogens, Fractionated	B199
H940-0 - S. Cerevisiae Antibodies (IGG/IGA).....	B481
H959-0 - Glycemic Control Profile.....	B246
H991-3 - Total IgG and subclasses IgG 1, 2, 3 & 4	B531
H992-1 - Egg Components IgE.....	B192
H993-9 - Peanut Components IgE.....	B448
J002-7 - Ova and Parasite (O&P), Sputum.....	B388
J010-0 - Culture, Catheter Tip	B125
J011-8 - Culture, Actinomyces	B160
J013-4 - Cystatin C	B162
J015-9 - ALK BY FISH, IF NEG REFLEX TO ROS1 BY FISH	B16
J018-3 - Cadmium, Blood	B99
J019-1 - Arsenic, Blood.....	B29
J020-9 - KIT Sequencing (Solid Tumors) Non NY	B320
J021-7 - BRAF Sequencing- Non NY.....	B73
J022-5 - EGFR Sequencing- Non NY.....	B510
J023-3 - NRAS Sequencing- Non NY	B382
J024-1 - KRAS Sequencing- Non NY.....	B322
J025-8 - HPV DNA, High Risk Anal - TP	B276
J032-4 - Human Epididymal Protein 4	B260
J044-9 - TP53 Sequencing.....	B519
J045-6 - PDGFRA Sequencing- Non NY	B414
J047-2 - PTH, Intact, FNA	B438
J048-0 - Pneumococcal Antibody Panel (23 Serotype)	B495
J050-6 - C-PEPTIDE, 24-HOUR URINE	B0
J054-8 - Reflex to Remainder of Common Cancer Panel (after J055).....	B473
J055-5 - Breast Cancer Management Panel.....	B92
J057-1 - Anti CMV Avidity, IgG	B26
J076-1 - Iodine, Serum/Plasma	B306

J090-2 - Cannabinoid, synthetic (K2,Spice), Urine, Screen, group 2, Qual.(10)	B318
J096-9 - Abacavir Genotype HLA B57	B35
J097-7 - HLA Typing for Celiac Disease	B273
J104-1 - KRAS, if neg reflex to BRAF- Non NY	B322
J105-8 - EGFR, if neg reflex to ALK- Non Ny	B177
J106-6 - EGFR, if neg reflex to ALK, if neg reflex to ROS1- Non NY	B179
J107-4 - BRAF, if neg reflex to MLH1 Hypermethylation- Non NY	B73
J114-0 - GLYCOMARK	B2
J130-6 - OnkoSight MYD88 Sequencing (Non-NY)	B396
J132-2 - Onkosight NGS ABL1 Sequencing	B9
J137-1 - Onkosight NGS NPM1 Sequencing	B382
J138-9 - Onkosight NGS KIT Sequencing (Non-NY)	B320
J141-3 - MYC-BA, 8q24 Rearrangement FISH	B359
J143-9 - Gleason Score 3+4, Reflex To PTEN FISH	B231
J144-7 - MYC-BA, 8q24 Rearrangement FISH- TECH ONLY	B358
J146-2 - SOX-10 BY IHC TECH Only	B489
J147-0 - SOX-10 BY IHC W/Interp	B489
J148-8 - 4K Prostate Score (Serum)	B7
J150-4 - Pregnenolone, Serum	B425
J154-6 - Onkosight NGS FLT3 IF NEG. RFLX TO CEBPA+NPM1 NGS	B212
J155-3 - Onkosight NGS FLT3 AND NPM1 Sequencing	B212
J157-9 - XL RARA (17q21) BY FISH	B562
J158-7 - H. Pyloril BY IHC W/ Interp	B256
J159-5 - H. Pylori BY IHC Tech Only	B253
J160-3 - XL RARA (17q21) BY FISH - TC (NON-NY)	B563
J164-5 - HERPES SIMPLEX VIRUS I DNA	B277
J165-2 - HERPES SIMPLEX VIRUS II DNA	B278
J185-0 - HIV-1 Genotype	B0
J202-3 - EGFR, if neg reflex to ALK - TC only Non NY	B178
J211-4 - G. vaginalis by RT-PCR	B224
J212-2 - A. Vaginalis by RT-PCR	B8
J213-0 - Lactobacillus (Species) by PCR	B326
J214-8 - Megasphaera Type 1 by RT PCR	B347
J215-5 - BVAB2 by RT-PCR	B78
J216-3 - B. fragilis by RT PCR	B54
J217-1 - Mobiluncus (Species) by RT PCR	B353
J218-9 - Hemoglobin A1c with reflex to Glycomark	B281
J220-5 - Leptin (ELISA)	B329

J221-3 - PD-L1 IHC 28-8 pharmDx (OPDIVO) W/INTERP.....	B413
J223-9 - Interleukin 5.....	B306
J257-7 - PD-L1 (Clone E1L3N) IHC w/Interp.	B412
J258-5 - Estradiol (E2), Ultrasensitive, LC/MS, Serum.....	B189
J259-3 - Estrone (E1), Serum.....	B190
J264-3 - 4K PSA rflx >3.0	B7
J265-0 - TSH (Pregnancy).....	B523
J275-9 - Syphilis Reverse Algorithm	B342
J280-9 - Allergen Ovomuroid (f233), IgE.....	B399
J281-7 - Allergen, Ovalbumin (F232), IgE.....	B399
J283-3 - Peanut Component Ara (h2), IgE	B447
J284-1 - Allergen Peanut Component Ara (h1), IgE.....	B446
J285-8 - Peanut Component Ara (h3), IgE	B447
J286-6 - Peanut Component Ara (h9), IgE	B448
J287-4 - Peanut Component Ara (h8), IgE	B448
J291-6 - Tick Identification.....	B528
J292-4 - H. Pylori Urea Breath Test Pediatric.....	B0
J293-2 - Ashkenazi Jewish Expanded Profile.....	B13
J294-0 - Ehlers-Danlos Syndrome, Type VIIC.....	B194
J295-7 - Osteopetrosis, Autosomal Recessive 1	B511
J296-5 - Multiple Sulfatase Deficiency	B374
J297-3 - Arthrogyrosis, Mental Retardation, and Seizures.....	B52
J298-1 - Spastic Tetraplegia, Thin Corpus Callosum and Progressive Microcephaly.....	B501
J299-9 - Dyskeratosis Congenita, Autosomal Recessive 5	B176
J301-3 - Polycystic Kidney Disease, Autosomal Recessive (ARPKD).....	B457
J302-1 - Phosphoglycerate dehydrogenase deficiency	B453
J303-9 - Peroxisome Biogenesis Disorder 5A (Zellweger).....	B451
J304-7 - Phenylketonuria (PKU) Carrier Screening.....	B452
J305-4 - Mitochondrial Complex I Deficiency	B370
J306-2 - Congenital Amegakaryocytic Thrombocytopenia (CAMT).....	B154
J307-0 - Familial Mediterranean Fever	B217
J308-8 - Deafness, Autosomal Recessive 1A (GJB2)	B171
J310-4 - Factor XI carrier screening (Hemophilia C).....	B216
J311-2 - Ciliary Dyskinesia, Primary 1	B459
J313-8 - Retinitis Pigmentosa 59	B477
J314-6 - Smith-Lemli-Opitz Syndrome	B500
J315-3 - Carnitine Palmitoyltransferase Def. Type 2 Carrier Screening.....	B119
J316-1 - Alport Syndrome, Autosomal Recessive.....	B45

J317-9 - Bardet-Biedl Syndrome 2.....	B80
J322-9 - Coenzyme Q10	B149
J331-0 - Synovial Sarcoma 18q11.2 FISH.....	B492
J355-9 - Solid Tumor Gene Fusion Assay	B240
J356-7 - p16 by IHC.....	B571
J361-7 - Reflex to Noonan from Array	B0
J362-5 - Reflex to Noonan from Chromosome	B0
J371-6 - Bladder Washing, FISH	B546
J405-2 - Zika Serum RNA PCR w/reflex to IgM.....	B565
J406-0 - Zika Serum RNA PCR	B565
J407-8 - Zika Urine RNA PCR	B566
J434-2 - Expanded Cystic Fibrosis Panel	B201
J435-9 - Cystic Fibrosis 40- NGS.....	B162
J440-9 - SPECIAL STAIN ONLY, VON KOSSA- Tech Only	B0
J441-7 - CoEnzyme Q10 Deficiency, Primary, 7	B147
J451-6 - Total Testosterone by LC/MS/MS	B526
J452-4 - KRAS/NGS if neg rflx to BRAF/ MSI-PCR, if MSI-High, rflx MLH(NY).....	B323
J453-2 - KIT Sequencing (Solid Tumors) NY Only	B321
J454-0 - EGFR Sequencing NY Only	B180
J455-7 - NRAS Sequencing NY Only.....	B383
J456-5 - BRAF Sequencing NY Only	B74
J458-1 - KRAS Sequencing NY Only	B323
J459-9 - PDGFRA Sequencing.....	B414
J460-7 - EGFR, if neg reflex to ALK, if neg reflex to ROS1 NY Only	B179
J461-5 - EGFR, if neg reflex to ALK NY Only	B178
J462-3 - KRAS, if neg reflex to BRAF NY Only	B322
J463-1 - BRAF, if neg reflex to MLH1 Hypermethylation NY Only.....	B73
J464-9 - KIT and PDGFRA Sequencing NY Only	B321
J465-6 - EGFR and KRAS, if both neg, reflex ALK, if neg reflex ROS 1 NY Only.....	B181
J466-4 - EGFR, if neg reflex to ALK-TC only- NY Only	B178
J496-1 - Minimal Residual Disease for CLL/SLL- Tech Only	B355
J497-9 - Anti-Zika Virus Elisa IgM Serum	B566
J510-9 - Thyroid Stimulating Immunoglobulin	B524
J520-8 - Cystic Fibrosis (CFTR) Intron 9 Poly-T Analysis (5T/7T/9T).....	B106
J529-9 - Cytology, Bladder Washing (50% Alcohol).....	B163
J534-9 - Cytology, Bladder Washing (ThinPrep)	B163
J546-3 - p57 by IHC.....	B572
J571-1 - Potassium, RBC	B458

J632-1 - JAK2 (including V617F and exon 12), MPL, and CALR if negative, reflex to MPN Panel	B316
J633-9 - MPN: JAK2/MPL/CALR by NGS	B354
J634-7 - OnkoSight Glioma Panel	B394
J635-4 - Glioma: IDH1/IDH2 by NGS	B244
J655-2 - Process and Hold for Cytogenetics	B460
J661-0 - Citrated Platelet Count.....	B363
J664-4 - Rest of Comprehensive Cancer Panel (after Prostate Panel)	B472
J665-1 - Hereditary Prostate Cancer Panel	B461
J745-1 - Cryptosporidium ID, Modified Acid Fast.....	B158
J761-8 - MSI: if MSI-H reflex to BRAF, if BRAF neg. reflex to MLH1	B356
J785-7 - Interleukin 6.....	B299
J791-5 - Anal GC on ThinPrep (Non NY).....	B47
J792-3 - Anal CT on ThinPrep (Non NY)	B47
J795-6 - Rest of Demyelinating CMT Panel (after PMP22 Del/Dup).....	B476
J834-3 - GC MultiPlex PCR - swab	B226
J835-0 - Chlamydia Trachomatis MultiPlex PCR - swab	B139
J836-8 - Trichomonas by MultiPlex PCR - swab.....	B535
J837-6 - Mycoplasma Genitalium by MultiPlex PCR- swab	B345
J839-2 - Varicella-Zoster Virus.....	B551
J847-5 - HSV 1 MultiPlex PCR	B0
J848-3 - HSV 2 MultiPlex PCR	B0
J857-4 - Histology Processing 6 H+E	B288
J858-2 - Histology Processing 8 H+E	B287
J859-0 - Histology Processing 7 H+E	B288
J860-8 - Histology Processing 9 H+E	B289
J861-6 - Histology Processing 10 H+E	B288
J862-4 - Histology Processing (GI)	B287
J865-7 - Opiates, Serum Screen.....	B397
J876-4 - MGMT Promoter Region Methylation	B350
J892-1 - Dihydrotestosterone	B167
J897-0 - Malaria/Babesia/Other Blood Parasites	B88
J900-2 - RUBELLA IGM.....	B479
J913-5 - Chromosome Analysis (Oncology)- Tissue	B131
J916-8 - IHC Stain Only Ref To ARUP additional Antibody	B0
J920-0 - Fecal Calprotectin	B207
J949-9 - Non-Gyn Other (Slides) - Interp Only	B386
J950-7 - RUBELLA VIRUS RNA, QUALITATIVE, RT-PCR (NON-NY).....	B0

J964-8 - Lactobacillus reuteri/rhamnosus (Non NY)..... B340

J965-5 - Group A Strep (Non NY)..... B237

J966-3 - Group B Strep (Non NY)..... B238

J967-1 - Enterococcus faecalis (Non NY)..... B177

J968-9 - Escherichia coli (Non NY)..... B176

J970-5 - Staphylococcus aureus (Non NY)..... B481

K002-6 - Celiac Disease Eval + HLA Genetics..... B104

K135-4 - 4K PSA Total + Free Reflex >3.0 to 4Kscore B7

K335-0 - EGFR by NGS, ALK, ROS1, and MET by FISH (Ny Only)..... B180

K361-6 - Food Allergy Profile w/ rfx to Milk and Peanut..... B221

K426-7 - Food Allergy Profile B221

K559-5 - Zika Molecular Serum+Urine B565

K583-5 - Zika Molecular Serum w/rfx and Urine B566

K630-4 - Testosterone, Total/Free (LC/MS)..... B532

K726-0 - Childhood Allergy Profile..... B138

K835-9 - BV Panel by RT-PCR B78

K853-2 - Vulvovaginal Candidiasis Panel by PCR..... B554

K908-4 - Region 11 Allergy Profile (AZ, ID, NM, WY, CO, UT)..... B473

K909-2 - Culture, Tracheal Airway..... B161

K924-1 - OnkoSight Lung Comprehensive..... B394

K935-7 - COA/BMP B112

K943-1 - HCV Viral Load reflex to NS3/NS5 (Non-NY) B259

K987-8 - Anal GC/CT ThinPrep (Non-NY)..... B23

L008-2 - HSV 1+ 2/ VZV Assay B277

L009-0 - HSV 1+2 Reflex to VZV B276

L010-8 - HSV 1+2 (NON NY) B279

L012-4 - ROMA..... B468

L047-0 - High Risk STI-CT/GC/Trich-Aptima Swab..... B505

L048-8 - High Risk STI- CT/GC/Trich/M.Gen-Aptima Swab..... B506

L049-6 - Aptima Swab Cervicitis (NON-NY)..... B506

L050-4 - Vaginitis/Vaginosis by RT-PCR, Aptima Tube..... B52

L062-9 - High Risk STI- CT/GC- Aptima Swab..... B505

L063-7 - Serum Drug Screen..... B498

L068-6 - Vaginitis/Vaginosis by PCR, Aptima Tube B0

L079-3 - StormPath PAP, liquid based + HPV DNA Genotyping 16/18 B503

L080-1 - StormPath PAP, liquid based + HPV DNA + CT/GC..... B504

L113-0 - OnkoSight Lung Panel, if EGFR- neg. rfx to (J015) ALK by FISH, if neg. rfx to ROS1 by FISH B395

L202-1 - PANEL L202.....	B402
L204-7 - Drug Confirmation, Clonazepam, Urine (LC/MS).....	B147
L218-7 - DLBCL FISH PANEL: MYC BA, BCL2, BCL6	B169
L282-3 - H. Pylori by IHC. If+ rflx to H. pylori PCR with rflx to H. pylori Clarithromycin resistance detecti	B256
L325-0 - Kidney Profile	B325
L326-8 - Pap (Reflex to HPV DNA Genotyping 16,18 when ASC-US or UNSAT)	B440
L327-6 - PAP/CT/GC (Reflex HPV PLUS when ASC-US or UNSAT)	B407
L344-1 - Chlamydia/N. Gonorrhoea CT/GC PCR Urine	B124
L393-8 - PDL1 (22C3 Keytruda) IHC + EGFR NGS + ALK FISH + ROS1 FISH	B415
L399-5 - Parvovirus B19 Ab (IgG/IgM).....	B444
L416-7 - Candida Species by M. PCR - Swab.....	B131
L570-1 - BV If Negative Reflex to Aerobic Vaginitis (Non NY).....	B78
L571-9 - Bacterial Vaginosis and Aerobic Vaginitis (Non NY)	B77
L580-0 - Tick-borne Panel by RT-PCR (Non-NY).....	B0
L583-4 - GC/CT SWAB BY PCR.....	B227
L592-5 - URINE FREE KAPPA+LAMBDA	B541
L626-1 - Rheumatoid Factor IgG, IgA, IgM	B478
L767-3 - AV Panel w/o Lacto Reuteri and Rhamnosis (Non NY)	B34
M063-6 - Lipid Screen W/Reflex to Direct LDL	B342
M064-4 - Iron Reflex to TIBC.....	B313
M065-1 - Peanut w/Reflex to Component Panel	B449
M066-9 - Egg White w/Reflex to Component Panel	B193
M067-7 - PSA Reflex to Free PSA.....	B436
M123-8 - COVID-19, FLU A and FLU B	B0
M160-0 - COVID-19 Ab. Qual/Quant.....	B0
M162-6 - COVID-19 Ab. Qual/Quant w/ Reflex to TH99	B0
M175-8 - Region 8 ImmunoCap Respiratory Allergy Profile.....	B0
M176-6 - Region 7 ImmunoCap Respiratory Allergy Profile.....	B0
M237-6 - Drug Confirmation, Fentanyl Urine (LC/MS)	B0
M255-8 - HSV 1+2 PANEL - TMA.....	B0
M264-0 - MALE/FEMALE STI PROFILE	B0
M280-6 - PAP + HPV MRNA (NO GENO)	B0
M281-4 - PAP + HPV MRNA (NO GENO) + CT/GL	B0
M289-7 - CORD BLOOD/INFANT ABO Rh TYPE w/Direct Coombs Test	B0
M322-6 - COVID-19, Flu A/B, RSV by PCR	B0
P006-0 - High Risk STI CT/GC/Trich/M.Gen.....	B125
P010-2 - Chlamydia + Gonorrhoea ,Oral, Swab, Aptima	B139

P011-0 - Chlamydia + Gonorrhoeae, Anal, Swab, Aptima	B124
P079-7 - Pap and HPV DNA Genotyping 16,18.....	B441
P114-2 - Urogenital Mycoplasmas	B541
P150-6 - Expanded CF, Fragile X, SMA.....	B190
P210-8 - Bone/Tissue Culture.....	B72
P216-5 - STI Serology	B493
P223-1 - Ashkenazi Jewish DNA (NY patients only)	B31
P244-7 - Menopause/Perimenopause Risk Management Profile	B348
P258-7 - Polycyst.Ovary Synd. Conf.Panel	B423
P259-5 - Polycyst.Ovary Synd. Treat. Panel	B423
P260-3 - ALL Prognosis Panel by FISH	B17
P261-1 - AML Diagnostic Panel by FISH	B21
P264-5 - MSI by PCR if MSI-H Reflex to MLH1	B370
P267-8 - Allergy Outdoor Panel.....	B19
P268-6 - Indoor Allergy Panel	B302
P310-6 - Male STI Urine Profile (GP)	B346
P372-6 - PAP + HPV Plus + CT + GC.....	B406
P374-2 - PAP+CT+GC(RFX HPV PLUS)	B406
P383-3 - PAP+CT+GC(RF HPV+ LO ASC OR >)	B405
P576-2 - Vaginitis/Vaginosis by DNA probe, Thin Prep.....	B57
P658-8 - Ashkenazi Jewish 11	B30
P659-6 - Ashkenazi Jewish 18	B30
P699-2 - PAP/HPVHI/GC RFX>GENCERV-P+HPV	B408
P701-6 - PAP/HPV+/GC RFX>GENCERV-PAP+HP (NON-NY).....	B408
P734-7 - Pap + Pap Dependent HPV + CT/GC.....	B440
P764-4 - HSV/LGV Lesion Profile (NON-NY).....	B278
P782-6 - GTT,GEST.,2HR.(ADA) PANEL.....	B238
P791-7 - Fractional Excretion Of Sodium	B215
P863-4 - Female STI Urine Profile	B218
P960-8 - MDS FISH Panel- TC Only.....	B376
P961-6 - Multiple Myeloma FISH Panel - TC Only	B345
P962-4 - CLL FISH Panel - TC Only.....	B110
Q032-5 - Bean Profile	B57
Q069-7 - Protein Electrophoresis, 24Hr Urine	B538
Q161-2 - Ashkenazi Jewish 25.....	B31
Q265-1 - Cervicitis W/O PAP W/O HPV PLUS.....	B105
Q301-4 - STI Lesion (NON-NY)	B493
Q471-5 - Heart Health Baseline	B260

Q472-3 - Heart Health Monitoring..... B260

Q605-8 - G.T.T. 2 hr (includ. 30 min)..... B224

Q723-9 - Stormpath PAP + HPV HR W/Rfx to HPV 16 and 18/45..... B494

Q724-7 - Stormpath PAP+CT/GC+HPV HR w/RFX to HPV 16 and 18/45. B494

R163-7 - ALDOSTERONE, Urine, Random B0

R732-9 - Aerobic Vaginitis Profile (Non NY) B41

R733-7 - POST TREATMENT SYPHILIS (RPR)..... B469

R734-5 - Lymphoproliferative Disorders Analysis (TC Only) and StormPath BMM
ReportTemplate B143

R735-2 - Acute Leukemia + Myeloid Disorders Analysis (Technical Only) and StormPath BMM
Report Template B37

R736-0 - Plasma Cell Analysis by Flow Cytometry (Tech Only) and StormPath BMM Report
Template B421

T007-4 - Rubeola/Measles IgG W/Reflex..... B366

T008-2 - Varicella Zoster Ab IgG W/Reflex B553

T009-0 - Mumps Virus Ab IgG W/Reflex B374

T385-4 - Rubella Ab IgG W/Reflex B479

T397-9 - Reverse Triiodothyronine Serum/Plasma (rT3) B573

T407-6 - Chlamydia (CT) PCR Urine..... B122

T408-4 - N. Gonorrhea/NG (GC) PCR Urine B379

T409-2 - Special Stain, Halls Bile - TC..... B84

T411-8 - ERG by IHC (EP111) B188

T412-6 - ERG by IHC (EP111) - Tech Only B188

T413-4 - PD-L1 (SP142) , TECENTRIQ. B411

T420-9 - LabFinder Report B339

T421-7 - B-CrossLaps/Serum B150

T447-2 - Ostase - Alkaline Phosphatase Bone Specific B91

T449-8 - Candida krusei B96

T536-2 - Candida parapsilosis B96

T745-9 - Candida tropicalis by PCR..... B97

T781-4 - PDGFRB by FISH (Global)- Non NY..... B414

T782-2 - PDGFRB by FISH (Tech Only) Non NY B415

T783-0 - Candida glabrata by PCR B95

T786-3 - Candida dubliniensis by PCR..... B96

T790-5 - Parvovirus B19 Ab, IgG..... B444

T791-3 - Parvovirus B19 Ab, IgM B444

T814-3 - QFT - TB GOLD PLUS (QUANTIFERON)..... B295

T824-2 - Pregnancy, Serum (Qualitative) B465

T830-9 - Hereditary MDS/Leukemia Panel	B284
T832-5 - Chlamydia Swab by PCR	B241
T833-3 - Gonorrhea Swab by PCR.....	B242
T835-8 - Intrinsic Factor AB.....	B313
T841-6 - Histone Ab, IgG.....	B289
T850-7 - Storage Sample Clinical Trials	B0
T969-5 - Procalcitonin.....	B410
T974-5 - Anti-Jo-1, IgG	B49
T979-4 - OnkoSight BRAF Sequencing - BM/PB.....	B74
T980-2 - OnkoSight NGS TP53 Sequencing for Solid Tumors	B519
T982-8 - Rule Out Mosaicism Chromosome Analysis, Blood	B354
T983-6 - Feulgen Stain - Tech Only	B218
T988-5 - PD-L1 IHC, Keytruda for Cervical Cancer W/INTERP.....	B411
TA08-0 - Rheumatoid Factor IgA	B468
TA09-8 - Rheumatoid Factor IgM	B468
TA14-8 - Babesia Microti by PCR (Non NY)	B0
TA15-5 - Borrelia mayonii by RT-PCR (Non NY).....	B0
TA16-3 - Borrelia burgdorferi sensu stricto by RT-PCR (Non NY).....	B0
TA17-1 - Breast Basal Layer Triple Stain by IHC Tech Only	B93
TA18-9 - Breast Basal Layer Triple Stain by IHC Global	B91
TA31-2 - PD-L1 IHC , Keytruda for Urothelial Carcinoma W/INTERP.....	B413
TA34-6 - Chromogranin A	B141
TB33-6 - StormPath PAP(RFLX to HPV mRNA GENOTYPING 16,18/45 w/ASC-US)	B503
TB34-4 - BRAF/EGFR NGS if EGFR neg rfx to ALK FISH if neg rfx ROS1	B74
TB35-1 - PD-L1 (Clone E1L3N) IHC Tech Only.....	B413
TB61-7 - IGH/MAFB T(14;20) - TC ONLY.....	B295
TF69-2 - Storage Sample Clinical Trials	B0
TF79-1 - Inhibin-A	B311
TF85-8 - Varicella Zoster Virus by IHC	B551
TF87-4 - Human Telomerase Reverse Transcriptase - Tech Only	B571
TF88-2 - Varicella Zoster Virus by IHC Tech Only.....	B552
TF92-4 - MUC4 by IHC Tech Only	B357
TF93-2 - MUC4 by IHC.....	B357
TG28-6 - Urine Cytology, if Atypical, reflex to hTert	B544
TG35-1 - Peanut Component Ara (h6), IgE.....	B447
TG39-3 - Medical Implants/Devices Culture	B160
TG46-8 - FGFR1 by FISH (Tech Only)	B209
TH18-5 - ClariTest Core Non-Invasive Prenatal Screen (NON-NY)	B381

TH19-3 - ClariTest Core Non-Invasive Prenatal Screen with 22q11.2 Microdeletion (DiGeorge Syndrome)(NON-NY)..... B381

TH31-8 - Scarlet Eligible Order B0

TH39-1 - NGS MPN Cascade Sequencing Report..... B0

TH48-2 - OnkoSight Advanced Prostate NGS Panel B41

TH49-0 - OnkoSight Advanced Lung Cancer NGS Panel B40

TH50-8 - OnkoSight Advanced Colorectal Cancer NGS Panel..... B119

TH53-2 - OnkoSight Advanced GYN Tumor NGS B40

TH55-7 - OnkoSight Advanced Chronic Lymphoid Neoplasm NGS Panel..... B110

TH57-3 - OnkoSight Advanced CNS NGS Panel..... B39

TH58-1 - OnkoSight Advanced Gastrointestinal Stromal Tumor NGS Panel B40

TH59-9 - OnkoSight Advanced Melanoma NGS Panel B41

TH61-5 - Acute 10 Color Flow Cytometry Leukemia/Lymphoma Panel (Global)..... B0

TH62-3 - Keytruda PD-L1 Esophageal Squamous Cell Carcinoma by IHC..... B415

TH63-1 - Keytruda PD-L1 Head & Neck Squamous Carcinoma (HNSCC) B412

TH64-9 - Acute 10 Color Flow Cytometry Leukemia/Lymphoma Panel (TC Only) B0

TH65-6 - Chronic 10 Color Flow Cytometry Lymphoproliferative Panel (Tech Only)..... B0

TH66-4 - Chronic 10 Color Flow Cytometry Lymphoproliferative Panel (Global)..... B0

TH67-2 - RCCA TACROLIMUS B0

TH68-0 - Novel Coronavirus COVID-19 Nasal/Nasopharynx..... B118

TH69-8 - Novel Coronavirus COVID-19 Oropharynx B118

TH71-4 - Novel Coronavirus COVID-19 Pooled N/NP/OP B118

TH76-3 - Urine Culture, Prenatal, W/GBS B544

TJ06-6 - Throat Culture, R/O Group A Strep. ONLY B527

TJ12-4 - CD71 IHC B0

TJ13-2 - CD71 IHC- TC..... B0

TJ16-5 - OnkoSight Advanced Solid Tumor NGS Panel B41

TJ45-4 - OnkoSight NGS JAK2 V617F Reflex to MPL & CALR Sequencing..... B0

TJ93-4 - OnkoSight Advanced Chronic Lymphoid NGS Panel B0

TK69-1 - 10-Color Screening Flow Cytometry (Short) Panel (Global) B0

TK78-2 - OnkoSight Advanced Comprehensive Lung B0

TK81-6 - 10-Color Screening (Short) Flow Cytometry Panel (Tech Only) B0

TK84-0 - OnkoSight Advanced Pancreatic and Biliary Tract Tumor Panel B0

TK85-7 - PD-L1 22C3 for Keytruda Other Tumor Types..... B0

TK96-4 - Plasma 10-Color Flow Cytometry Myeloma/Monoclonal Gammopathy Panel (Tech Only)..... B0

TK97-2 - Plasma 10-Color Flow Cytometry Myeloma/Monoclonal Gammopathy Panel (Global) ... B0

TL02-0 - PD-L1 22C3 for Keytruda (Triple-Neg Breast Cancer) B0

TL40-0 - Novel Coronavirus COVID-19 for International Travel..... B0

TL42-6 - ClariTest Core Non-Invasive Prenatal Screen (NY) B0

TL43-4 - ClariTest Core Non-Invasive Prenatal Screen with 22q11.2 Microdeletion (NY)..... B0

TL62-4 - Urine Culture, Special B0

TL67-3 - OnkoSight Advanced Plasma Cell Myeloma NGS Panel B0

TL68-1 - OnkoSight Advanced NGS Multiple Myeloma Panel B0

TL69-9 - SOX-11, Stain Only B0

TL70-7 - SOX-11, IHC, Global..... B0

TL74-9 - CORD BLOOD ABO and Rh Blood Typing B0

TL75-6 - Mycoplasma genitalium Aptima Tube B0

TL76-4 - Trichomonas - Aptima Tube B0

TL87-1 - HSV 1- Aptima Tube B0

TL88-9 - HSV 2- Aptima Tube B0

TL94-7 - OnkoSight Advanced NGS Rapid AML Panel..... B0

TL95-4 - OnkoSight Advanced NGS Myeloid Panel..... B0

TL96-2 - OnkoSight Advanced NGS MPN Panel B0

TL97-0 - OnkoSight Advanced NGS MPN Cascade Panel B0

TL98-8 - OnkoSight Advanced NGS MDS Panel B0

TL99-6 - OnkoSight Advanced NGS JAK2, MPL and CALR Panel B0

TM07-7 - OnkoSight Advanced NGS JAK2, MPL and CALR Reflex to OnkoSight Advanced NGS Myeloid Panel B0

TM70-5 - CLL/MCL Prognosis: IgVH Mutation Analysis by NGS B0

TM71-3 - Chlamydia trachomatis DNA..... B0

TM72-1 - Neisseria gonorrhoea DNA..... B0

TM73-9 - Trichomonas vaginalis DNA B0

TM74-7 - Mycoplasma genitalium DNA..... B0

TM75-4 - Prenatal 46, XY Disorders of Sex Development Panel B0

TM76-2 - Prenatal Tuberous Sclerosis B0

TM79-6 - GeneStratNGS Genomic Test..... B0

TM80-4 - HIV-1 Qualitative TMA..... B0

TN00-0 - Anaplasma phagocytophilum DNA, Qualitative Real-Time PCR..... B0

TN01-8 - Allergen Aspergillus flavus IgE..... B0

TN02-6 - Early Sjogren's Syndrome Profile B0

TN03-4 - Helicobacter pylori Culture w/Antimicrobial Suscept B0

TN04-2 - Cardio IQÂ® Insulin Resistance Panel with Score B0

TN05-9 - Sjogren's Profile Plus (RDL)..... B0

TN06-7 - Omega-3 and -6 Fatty Acids..... B0

TN07-5 - OmegaCheck..... B0

TN08-3 - Selenium, Blood	B0
TN09-1 - Titanium, Blood	B0
TN10-9 - TMAO (Trimethylamine N-Oxide).....	B0
TN11-7 - TSH Antibody	B0
TN12-5 - TSH with HAMA Treatment.....	B0
TN13-3 - Rheumatoid Factor (IgG)	B0
TN14-1 - Apixaban.....	B0
TN16-6 - Copeptin	B0
TN17-4 - Ku Autoantibodies	B0
TN19-0 - Muscle Fibers, Fecal.....	B0
TN20-8 - Natural Killer Cells, Functional	B0
TN21-6 - Neopterin, Serum.....	B0
TN22-4 - Digoxin, Free, Serum	B0
TN23-2 - LRP4 Autoantibody Test.....	B0
TN24-0 - Hepatitis C Virus (HCV) FibroSure	B0
TN26-5 - Anti-MDA-5 Ab (CADM-140)	B0
TN27-3 - Anti-Synthetase Profile.....	B0
TN28-1 - Crohn's Disease Prognostic Profile.....	B0
TN29-9 - ADmark® Phospho-Tau/Total-Tau/A Beta 42, Interp, CSF	B0
TN30-7 - Allergen Almond Food IgG4 (Non-NY).....	B0
TN31-5 - Allergen Alpha Lactalbumin IgG4 (Non-NY).....	B0
TN32-3 - Allergen Beta Lactoglobulin IgG4 (Non-NY).....	B0
TN33-1 - Allergen Casein IgG4 (Non-NY).....	B0
TN34-9 - CU Index®* Panel.....	B0
TN36-4 - Allergen Egg White IgG4 (Non-NY)	B0
TN38-0 - Allergen Hazelnut/Filbert Food IgG4 (Non-NY)	B0
TN39-8 - Allergen Helminthosporium sativum/Drechslera IgE	B0
TN40-6 - Immuknow Immune Cell Function	B0
TN41-4 - Allergen Milk Cow IgG4 (Non-NY).....	B0
TN42-2 - Allergen Ovalbumin IgG4 (Non-NY)	B0
TN43-0 - Allergen Ovomuroid IgG4 (Non-NY)	B0
TN44-8 - Allergen Peanut IgG4 (Non-NY)	B0
TN46-3 - Allergen Rabbit Urine Proteins IgE	B0
TN47-1 - Allergen Mite Tyrophagus putrescentiae IgE	B0
TN48-9 - Oxidized Low-density Lipoprotein (OxLDL)	B0
TN49-7 - Ri Antibody Screen with Reflex to Titer and Western Blot	B0
TN51-3 - IgH Gene Rearrangement by NGS	B0
TN53-9 - Myelin Oligodendrocyte Glycoprotein (MOG) Ab with Reflex, Ser	B0

TN55-4 - Cardio IQ[®] Lipoprotein Fractionation, Ion Mobility (Non-NY)B0

TN56-2 - SC5b-9 Level Terminal Complement ComplexB0

TN57-0 - Advanced Lipid Panel, Cardio IQB0

TN59-6 - Eastern Equine Encephalitis Virus IgM Antibody, IFA (Serum).....B0

TN60-4 - Myelin Oligodendrocyte Glycoprotein (MOG) Ab with Reflex, CSFB0

TN77-8 - OVA1plusB0

TN92-7 - Monkeypox PCRB0

TN93-5 - OnkoSight Advanced NGS CEBPAB0

TN95-0 - OnkoSight Advanced NGS MYD88B0

TN96-8 - OnkoSight Advanced NGS BRAF Sequencing.....B0

TP00-6 - OnkoSight Advanced NGS MYD88 and CXCR4 Panel.....B0

TP14-7 - Influenza A by MolecularB0

TP15-4 - Influenza B by MolecularB0

TP16-2 - Resp. Syncytial Virus by MolecularB0